The Spike2 script language

Version 5

Copyright © Cambridge Electronic Design Limited 1988-2007

Neither the whole nor any part of the information contained in, or the product described in, this guide may be adapted or reproduced in any material form except with the prior written approval of Cambridge Electronic Design Limited.

Version 5.00	March 2003
Version 5.01	April 2003
Version 5.02	May 2003
Version 5.03	June 2003
Version 5.04	October 2003
Version 5.05	January 2004
Version 5.06	April 2004
Version 5.07	June 2004
Version 5.08	September 2004
Version 5.09	November 2004
Version 5.10	January 2005
Version 5.12	May 2005
Version 5.13	October 2005
Version 5.15	May 2006
Version 5.17	March 2007

Published by:

Cambridge Electronic Design Limited Science Park Milton Road Cambridge CB4 0FE UK Telephone: Cambridge (01223) 420186 International +44 1223 420186 Fax: Cambridge (01223) 420488 International +44 1223 420488 Email: info@ced.co.uk Home page: www.ced.co.uk

Curve fitting functions are based on routines in Numerical Recipes: The Art of Scientific Computing, published by Cambridge University Press and are used by permission.

Trademarks and Tradenames used in this guide are acknowledged to be the Trademarks and Tradenames of their respective Companies and Corporations.

Introduc	ction	1-1
	Hello world	
	Views and view handles	
	Writing scripts by example	
	Using recorded actions	
	Differences between systems	
	Notation conventions	
	Sources of script information	
The scri	pt window and debugging	2-1
	Syntax colouring	
	Editing features for scripts	
	Debug overview	
	Preparing to debug	
	Inspecting variables.	2-4
	Call stack	
Scrint la	anguage syntax	3-1
Jonpula	Keywords and names	3_1
	Data types	
	Arrays	
	Result views as arrays	3_/
	Statement types	3_/
	Comments	2 /
	Variable deelerations	
	Constant delerations	
	Every Every and operators	
	Expressions and operators	
	Flow of control statements	/
	Channel specifiers	
	Channel specificis	
Commar	nds by function	4-1
	Windows and views	
	Channel commands	
	Time views	
	Result views	
	XY view	
	Vertical cursors	
	Horizontal cursors	
	Editing operations	
	Text files	
	Binary files	
	File system	
	User interaction commands	
	Analysis	
	Mathematical functions	
	Array and matrix arithmetic	
	String functions	
	Curve fitting	4-8
	Digital filtering	4-9
	Sampling configuration	4-9
	Runtime sampling	<u>4_</u> C
	Arbitrary waveform output	ر-ب. 4_1(
	Discriminator control	۲۰۱۲ + ۲۰ ۱∟۱۲
	Signal conditioner control	
	Debugging operations	+-۱۰ ۲۰ ۱۱ /
	Debugging operations	

	Environment	4-11
	Spike shape window	
	Spike Monitor	
	Multimedia files	
	Serial line control	
	Multiple monitor support	
Alphabet	tical command reference	
Translati	ing DOS scripts	
	What will translate	
	What will not translate	
	General translation issues	
	After translating	
XY views	8	7-1
	Creating an XY view	
	Overdrawing data	
Curve fit	ting from a script	
	Testing the fit	
	Fitting routines	

What is a script? For many users, the interactive nature of Spike2 may provide all the functionality required. This is often the case where the basic problem to be addressed is to present the data in a variety of formats for visual inspection with, perhaps, a few numbers to extract by cursor analysis and a picture to cut and paste into another application. However, some users need analysis of the form:

- 1. Find the peak after the second stimulus pulse from now.
- 2. Find the trough after that.
- 3. Compute the time difference between these two points and the slope of the line.
- 4. Print the results.
- 5. If not at the end, go back to step 1.

This could all be done manually, but it would be very tedious. A script can automate this process, however it requires more effort initially to write it. A script is a list of instructions (which can include loops, branches and calls to user-defined and built-in functions) that control the Spike2 environment. You can create scripts by example, or type them in by hand.

Hello world Traditionally, the first thing written in every computer language prints "Hello world" to the output device. Here is our version that does it twice! To run this, use the File menu New command to create a new, empty script window. Type the following:

Message("Hello world");
PrintLog("Hello world");

Click the Bun button to check and run your first script. If you have made any typing errors, Spike2 will tell you and you must correct them before the script will run. The first line displays "Hello world" in a box and you must click on OK to close it. The second line writes the text to the Log view. Open the Log view to see the result (if the Log window is hidden you should use the Window menu Show command to make it visible).

So, how does this work? Spike2 recognises names followed by round brackets as a request to perform an operation (called a *function* in computer-speak). Spike2 has around 400 built-in functions, and you can add more with the script language. You pass the function extra information inside the round brackets. The additional information passed to a function is called the function *arguments*.

Spike2 interprets the first line as a request to use the function Message() with the argument "Hello world". The message is enclosed in double quotation marks to flag that it is to be interpreted as text, and not as a function name or a variable name.

An argument containing text is called a *string*. A string is one of the three basic data types in Spike2. The other two are *integer* numbers (like 2, -6 and 0) and real numbers (like 3.14159, -27.6 and 3.0e+8). These data types can be stored in *variables*.

Spike2 runs your script in much the same way as you would read it. Operations are performed in reading order (left to right, top to bottom). There are also special commands you can insert in the script to make the it run round loops or do one operation rather than another. These are described in the script language syntax chapter.

Spike2 can give you a lot of help when writing a script. Move the text caret to the middle of the word Message and press the F1 key. Help for the Message() command appears, and at the bottom of the help entry you will find a list of related commands that you might also find useful.

handles

Views and view The most basic concept in a script is that of a view and the view handle that identifies it. A view is a window in Spike2 that the script language can manipulate. The running script is hidden from most commands, however you can obtain its handle using App (3) so you can show and hide it.

> There is always a *current view*. Even if you close all windows the Log view, used for text output by the PrintLog() command, remains. Whenever you use a built-in function that creates a new view, the function returns a view handle. The view handle is an integer number that identifies the view. It is used with the View() and FrontView() functions to specify the current view and the view that should be on top of all windows.

> Whenever a script creates a new view, the view becomes the current view. However, views are created invisibly so that they can be configured before appearing. You can use WindowVisible(1) to display a new window.

example

Writing scripts by To help you write scripts Spike2 can monitor your actions and write the equivalent script. This is a great way to get going writing scripts, but it has limitations. Scripts generated this way only repeat actions that you have already made. The good point of recording your actions is that Spike2 shows you the correct function to use for each activity.

> For example, let us suppose you use the Script menu Turn Recording On option, open a data file, select interval histogram analysis mode of channel 3, process all data in the file and end with the Script menu Stop recording command. Spike2 opens a new window holding the equivalent script (we have tidied the output up a little and added comments):

```
'declare two integer variables (% means integer)
var v6%,v7%;
v6% := FileOpen("demo.smr", 0, 3);
                                       'Open file, save view handle
Window(10, 10, 80, 50);
                                       'Set a window position
v7% := SetInth(3, 100, 0.01, 0);
                                       'Create invisible INTH view
WindowVisible(1);
                                       'make the INTH visible
Process(0, View(-1).Maxtime(), 0, 1); 'Add data to INTH view
```

Some of this is fairly straightforward. You can find the FileOpen(), SetInth(), and Process () functions described in this manual and they seem to map onto the actions that you performed. However, there is extra scaffolding holding up the structure.

In the first line, the var keyword creates two integer variables, v6% and v7%. These variables hold view handles returned by FileOpen() and SetInth(). Spike2 generates the names from the internal view numbers (so your script may not be exactly the same). The result of the functions is copied to the variable with the := operator. In English, the second line of the script could be read as: Copy the result of the FileOpen command on file demo.smr to variable v6%.

The SetInth() command makes a new window to hold an interval histogram of the data in channel 3 with 100 bins of 0.01 seconds width and with the first bin starting at an interval of 0 seconds. The WindowVisible(1) command is present because the new window created by SetInth() is hidden. Spike2 creates invisible windows so that you can size and position them before display to prevent excessive screen repainting.

The View (-1). syntax accesses data belonging to a view other than the current view. The current view when the Process () command is used is the result view and we want to access the maximum time in the original time view. The negative argument tells the script system that we want to change the current view to the time view associated with this result view. The dot after the command means that the swap is temporary, only for the duration of the Maxtime() function.

actions

Using recorded You can now run this script with the control buttons. The script runs and generates a new result view, repeating your actions. Now suppose we want to run this for several files, each one selected by the user. You must edit the script a bit more and add in some looping control. The following script suggests a solution. Notice that we have now changed the view handle variables to names that are a little easier to remember.

```
var fileH%, inthH%;
                                      'view handle variables
fileH% := FileOpen("", 0, 1);
                                      'blank for dialog, single window
'FileOpen returns -ve if no file
while fileH\% > 0 do
   inthH% := SetInth(3, 100, .01, 0);WindowVisible(1);
   Process(0, View(-1).Maxtime(), 0, 1);
   Draw();
                                      'Update the INTH display
   fileH% := FileOpen("", 0, 1); 'ask for the next file
wend;
```

This time, Spike2 prompts you for the file to open. The file identifier is negative if anything goes wrong opening the file, or if you press the Cancel button. We have also included a Draw() statement to force Spike2 to draw the data after it calculates the interval histogram. There is a problem with this script if you open a file that does not contain a channel 3 that holds events or markers of some kind. We will deal with this a little later.

However, you will find that the screen gets rather cluttered up with windows. We do not want the original window once we have calculated the histogram, so the next step is to delete it. We also have added a line to close down all the windows at the start, to reduce the clutter when the script starts.

```
var fileH%, inthH%;
FileClose(-1);
                                  'close all windows except script
fileH% := FileOpen("", 0, 1);
                                  'use a blank name to open dialog
while fileH\% > 0 do
                                  'FileOpen returns -ve if no file
   inthH% := SetInth(3, 100, .01, 0);WindowVisible(1);
   Process(0, View(-1).Maxtime(), 0, 1);
   View(fileH%).FileClose();
                                  'Shut the old window
                                  'Update the INTH display
   Draw();
   fileH% := FileOpen("", 0, 1); 'ask for the next file
wend:
```

This seems somewhat better, but we still have the problem that there will be an error if the file does not hold a channel 3, or it is of the wrong type. The solution to this is to ask the user to choose a channel using a dialog. We will have a dialog with a single field that asks us to select a suitable channel:

```
var fileH%, inthH%, chan%;
                                    'Add a new variable for channel
FileClose(-1);
                                    'close all windows to tidy up
fileH% := FileOpen("", 0, 1);
                                    'use a blank name to open dialog
while fileH\% > 0 do
                                    'FileOpen returns -ve if no file
   DlgCreate("Channel selection"); 'Start a dialog
   DlgChan(1, "Choose channel for INTH", 126); 'all but waveform
   if (DlgShow(chan%) > 0) and
                                    'User pressed OK and..
                                    '...selected a channel?
      (chan \ge 0) then
      inthH% := SetInth(chan%, 100, .01, 0);WindowVisible(1);
Process(0, View(-1).Maxtime(), 0, 1);
                                   'Shut the old window
      View(fileH%).FileClose();
                                    'Update the display
      Draw();
   endif
   fileH% := FileOpen("", 0, 1); 'ask for the next file
wend;
```

The DlgCreate() function has started the definition of a dialog with one field that the user can control. The DlgChan() function sets a prompt for the field, and declares it to be a channel list from which we must select a channel (or we can select the No channel entry). The DlgShow() function opens the dialog and waits for you to select a channel and press OK or Cancel. The if statement checks that all is well before making the histogram.

Differences between systems	If you need to write code that can run on different systems, you should try to avoid system dependent features, or if this is impossible, make use of the System() function to find out where you are. The Macintosh versions of Spike2 stop at version 2, so if you use any version 3, 4 or 5 features, your script will not run on a Macintosh.						
Notation conventions	Throughout this manual we use the font of this sentence to signify descriptive tex Function declarations, variables and code examples print in a monospaced font, for example $a := \text{View}(0)$. We show optional keywords and arguments to functions is curly braces:						
	<pre>Func Example(needed1, needed2 {,opt1 {,opt2}});</pre>						
	In this example, the first two arguments are always required; the last two are optional. Any of the following would be acceptable uses of the function:						
	<pre>a := Example(1,2); 'Call omitting the optional arguments a := Example(1,2,3); 'Call omitting one argument a := Example(1,2,3,4); 'Call using all arguments</pre>						
	A vertical bar between arguments means that there is a choice of argument type:						
	<pre>Func Choice(i% r str\$);</pre>						
	In this case, the function takes a single argument that could be an integer, a real or string. The function will detect the type that you have passed and may perform a different action depending upon the type.						
	Three dots (\ldots) stand for a list of further, similar items.						
Sources of script information	You will find that the rest of this manual is a reference to the script language and to the built-in script functions. Once you are familiar with the scripting system it will be your most useful documentation. There are example and utility script provided with Spike2. These are copied to the scripts folder within the folder that contains Spike2.						
	There is a separate manual provided with Spike2 that has been used for our Spike2 training courses, held at CED and around the world. This manual contains many annotated examples and tutorials. Some of the scripts in this manual are useful in their own right; others provide skeletons upon which you can build your own applications.						
	Our web site at www.ced.co.uk has example scripts and script updates that you can download.						

Script window

You use the script window when you write and debug a script. Once you are satisfied that your script runs correctly you would normally run a script from the script menu without displaying the source code. You can have several scripts loaded at a time and select one to run with the Script menu Run Script command.

The script window is a text window with a few extra controls including a "splitter" control so that you can view two parts of the

🔛 Script1	
ОК.	Functions 💌 📑 🖑 🐀 ?
PrintLog("Hello world"); 's	say hello!
Margin	Splitter control

script at the same time. To use the splitter control, drag it down the window. To cancel it, drag it to the top or bottom of the window. There is no script language control of the splitter.

To the left of the text area is a margin where you can set break points (one is shown already set), bookmarks and where the current line of the script is indicated during debugging. Above the text area is a message bar and several controls. The controls have the following functions:

Function

This control is a quick way to find any func or proc in your script. Click on this to display a list, in alphabetical order, of the names of all user-defined routines. Select one, and the window will scroll to it. To be located, the keywords func and proc must be at the start of a line and the name of the routine must be on the same line.

Compile

The script compiler checks the syntax of the script and if no errors are found it creates the compiled version, ready to run. If the script has not been changed since the last compile and no other script has been compiled, the button is disabled as there is no need to compile again. Spike2 can have one compiled script in memory at a time.



If the script has not been compiled it is compiled first. If no errors are found, Spike2 runs the compiled version, starting from the beginning. Spike2 skips over proc ... end; and func ... end; statements, so the initial code can come before, between or after any user-defined procedures and functions. This button is disabled once the script has started to run.

Set break point

This button sets a break point on the line containing the text caret, or clears a break if one is already set. A break point stops a running script when it reaches the start of the line containing the break point. You can also set and clear break points by moving the mouse pointer over the margin on the left of the script and double clicking.

Not all statements can have break points set on them. Some statements, such as var, const, func and proc compile to entries in a symbol table; they generate no code. If you set a break point on one of them the break point will appear at the next statement that is "breakable". If you set break points before you compile your script, you may find that some break points move to the next breakable line when you compile.

Clear all break points

oints 🛞 This button is enabled if there are break points set in the script. Click this button to remove all break points from the script. Break points can be set and cleared at any time, even before your script has been compiled.



This button provides help on the script language. It takes you to an alphabetic list of all the built-in script functions. If you scroll to the bottom of this list you can also find links to the script language syntax and to the script language commands grouped by function. Within a script, you can get help on keywords and built in commands by clicking on the keyword or command and pressing the F1 key.

Syntax colouring	Spike2 supports syntax colouring for both the script language and also for the output sequencer editor. You can customise the colouring (or disable it) from the Editor settings section of the Edit menu Preferences. The language keywords have the standard colour blue, quoted text strings have the standard colour red, and comments have the standard colour green. You can also set the colour for normal text (standard colour black) and for the text background (standard colour white).
	The syntax colouring options are saved in the Windows registry. If several users share the same computer, they can each have their own colouring preferences as long as they log on as different users.
Editing features for scripts	There are some extra editing features that can help you when writing scripts. If you like to indent your text, you can indent and outdent selected blocks by selecting one or more characters, then use the Tab key to indent and the Shift-Tab combination to outdent. The text moves by one tab stop. Tab stop sizes are set in the Edit menu Preferences.
Debug overview	Despite all our efforts to make writing a script easy, and all your efforts to get things right, sooner or later (usually sooner), a script will refuse to do what you expect. Rather than admit to human error, programmers attribute such failures to "bugs", hence the act of removing such effects is "debugging". The term dates back to times when an insect in the high voltage supply to a thermionic valve really could cause hardware problems.
	To make bug extermination a relatively simple task, Spike2 has a "debugger" built into the script system. With the debugger you can:
	 Step one statement at a time Step into or over procedures and functions Step out of a procedure or function Step to a particular line
	 Enter the debugger on a script error to view variable values View local and global variables Edit variable values See how you reached a particular function or procedure
	• Set and clear break points With these tools at your disposal, most bugs are easy to track down.
Preparing to debug	A Spike2 script does not need any special preparation for debugging except that you must set a break point or include the Debug() command at the point in your script where you want to enter the debugger. Alternatively, you can enter the debugger "on demand" by pressing the Esc key; you may need to hold it down for a second or two, depending on what the script is doing. If the Toolbar() or Interact() commands are active, hold down Esc and click on a button. This is a very useful if your script is running round in a loop with no exit! You can disable this with the Debug(0) command, but we suggest that this feature is added after the script has been tested! Once you have disabled debugging, there is no way out of a loop.
	You can also choose to enter debug on a script error by checking the Enter debug on script error box in the General tab of the Edit menu Preferences option. Depending on the error, this may let you check the values of variables to help you fix the problem.

When your script enters the debugger, the debug toolbar opens if it was not visible. The picture shows the toolbar as a floating window,



but you can dock it to any side of the Spike2 window by dragging over the window edge.



Stop running the script. There is no check that you really meant to do this, as we assume that if you know enough to get into the debugger, you know what you are doing! You can use the Debug() command to disable the debugger.



Display the current line in the script. If the script window is not visible, this will make it visible, then bring it to the top and scroll the text to the current line.



If the current statement contains a call to a user-defined Proc or a Func, step into it, otherwise just step. This does not work with the Toolbar() command which is not user-defined, but which can cause user-defined routines to be called. To step into a user-defined Func that is linked to a Toolbar() command, set a break point in the

Func.



Step over this statement to the next statement. If you have more than one statement on a line you will have to click this button once for each statement, not once per line.

If you are in a procedure or function, step until you return from it. This does not {}• work if you are in a function run from the Toolbar() command as there is nowhere to return to. In this case, the button behaves as if you had pressed the run button.



Run until the script reaches the start of the line containing the text caret. This is slightly quicker than setting a break point, running to it, then clearing it (which is what this does).



Run the script. This disables the buttons on the debug toolbar and the script runs until in reaches a break point or the end of the script.



Show the local variables for the current user-defined func or proc in a window. If there is no current routine, the window is empty. You can edit a value by double clicking on the variable. Elements of arrays are displayed for the width of the text window. If an array is longer than the space in the window, the text display for the array ends with ... to show that there is more data.



Show the global variable values in a window. You can edit a global variable by double clicking on it. The very first entry in this window lists the current view by handle, type and window name.

Display the call stack (list of calls to user-defined functions with their arguments) on the way to the current line in a window. If the Toolbar() function has been used, the arguments for it appear, but the function name is blank.

The debug toolbar and the locals, globals and the call window close at the end of a script. The buttons in the debug toolbar are disabled if they cannot be used. If you forget what a particular button does, move the mouse pointer over the button. A "Tool tip" window will open next to the button with a short description; if the Status bar is visible, a longer description can be seen there.

Enter debug on error There is an option in the Edit menu Preferences dialog General tab that allows you to enter the debugger if an error occurs. After an error you can inspect the values of local and global variables and the contents of the call stack, but you are not allowed to continue running the script.



The example above shows a script that prompts the user for a data file, opens it and prints a summary of the data channels to the log window. The user set two break points, ran the script then clicked on the Step button five times and clicked the Globals button.

Inspecting variables

If the locals or globals windows are open, they display a list of variables. If there are more variables than can fit in the window you can scroll the list up and down to show them all. Simple variables are followed by their values. If you double click on one a new window opens in which you can edit the value of the variable.

If you double click on an array, a new window opens that lists the values of the elements of the array. You must double click on an element to edit the value. There is a limit of 32000 to the number of array elements that can be displayed and edited. This should not be a problem for most users. Function and procedure arguments show the name and cannot be edited.

The call stack can sometimes be Call stack useful to figure out how your script arrived at a position in your code. This is particularly true if your script makes recursive use of functions. A function is recursive when it calls itself, either directly, or indirectly through other functions. This example calculates factorials using recursion. We have set a break point and then displayed the call stack so you can see all the calls, and the arguments for each

🎥 Fact.s2s	💱 Calls for Fact.s2s 📃 🗖 🗙
Executing script Functions	Fact% (10)
<pre>var n%; n% := Fact%(10); halt;</pre>	Fact% (8) Fact% (7) Fact% (6) Fact% (5) Fact% (4)
<pre>func Fact%(n%) fif n% <= 1 then return 1; else return Fact%(n%-1); endif; end;</pre>	

call. A common fault with scripts is to have mutually recursive user options. This leads to users burrowing deeper and deeper into the call stack until they run out of memory. The call stack can help to detect such problems.

Script format A script consists of lines of text. Each line can be up to 240 characters long, however we suggest a maximum line length of 78 characters as experience shows that this makes printing and transfer of scripts to other systems simple.

The script compiler treats consecutive white space as a single space except within a literal string. White space characters are end of line, carriage return, space and tab. The compiler treats comments as white space.

The maximum size of a script is limited by the number of instructions that it compiles into. This number is displayed in the status bar of the script window when you compile. The limit is currently 1,000,000 instructions, which is a very large script, probably around 160,000 lines of typical script code.

Keywords and names All keywords, user-defined functions and variable names in the script language start with one of the letters a to z followed by the characters a to z and 0 to 9. Keywords and names are not case sensitive, however users are encouraged to be consistent in their use of case as it makes scripts easier to read. Variables and user-defined functions use the characters % and \$ at the end of the name to indicate integer and string type.

User-defined names can extend up to a line in length. Most users will restrict themselves to a maximum of 20 or so characters.

The following keywords are reserved and cannot be used for variables or function names:

and	band	bor	bxor	case
const	diag	do	docase	else
end	endcase	endif	for	func
halt	if	mod	next	not
or	proc	repeat	resize	return
step	then	to	trans	until
var	view	wend	while	xor

Further, names used by Spike2 built-in functions cannot be redefined as user functions or global variables. They can be redefined as local variables (not recommended).

Data types There are three basic data types in the script language: real, integer and string. The real and integer types store numbers; the string type stores characters. Integer numbers have no fractional part, and are useful for indexing arrays or for describing objects for which fractions have no meaning. Integers have a limited (but large) range of allowed values.

Real numbers span a very large range of number and can have fractional parts. They are often used to describe real-world quantities, for example the weight of an object.

Strings hold text and automatically grow and shrink in length to suit the number of text characters stored within them.

Real data type This type is a double precision floating point number. Numbers are stored to an accuracy of at least 16 decimal digits and can have a magnitude in the range 10^{-308} to 10^{308} . Variables of this type have no special character to identify them. Real constants have a decimal point or the letter e or E to differentiate from integers. White space is not allowed in a sequence of characters that define a real number. Real number constants have one of the following formats where digit is a decimal digit in the range 0 to 9:

```
{-}{digit(s)}digit.{digit(s)}{e|E{+|-}digit(s)}
{-}{digit(s)}.digit{digit(s)}{e|E{+|-}digit(s)}
{-}{digit(s)}digitE|e{+|-}digit(s)
```

A number must fit on a line, but apart from this, there is no limit on the number of digits. The following are legal real numbers:

1.2345 -3.14159 .1 1. 1e6 23e-6 -43e+03

E or e followed by a power of 10 introduces exponential format. The last three numbers above are: 1000000 0.000023 -43000.0. The following are not real constants:

1 e6	White space is not allowed	1E3.5	Fractional powers are not allowed
2.0E	Missing exponent digits	1e500	The number is too large

Integer data type The integer type is identified by a % at the end of the variable name and stores 32-bit signed integer (whole) numbers in the range -2,147,483,648 to 2,147,483,647. There is no decimal point in an integer number. An integer number has the following formats (where digit is a decimal digit 0 to 9, and hexadecimal-digit is 0 to 9 or a to f or A to F, with a standing for decimal 10 to f standing for decimal 15):

> {-}{digit(s)}digit {-}0x|X{hexadecimal-digit(s)}hexadecimal-digit

You may assign real numbers to an integer, but it is an error to assign numbers beyond the integer range. Non-integral real numbers are truncated (towards zero) to the next integral number before assignment. Integer numbers are written as a list of decimal digits with no intervening spaces or decimal points. They can optionally be preceded by a minus sign. The following are examples of integers:

1 -1 -2147483647 0 0x6789abcd 0X100 -0xced

Integers use less storage space than real numbers and are slightly faster to work with. If you do not need fractional numbers or huge numeric ranges, use integers.

String data type Strings are lists of characters. String variable names end in a \$. String variables can hold strings up to 65534 characters long. Literal strings in the body of a program are enclosed in double quotation marks, for example:

"This is a string"

A string literal may not extend over more than one line. Consecutive strings with only white space between them are concatenated, so the following:

"This string starts on one lin" "e and ends on another"

is interpreted as "This string starts on one line and ends on another". Strings can hold special characters, introduced by the escape character backslash:

- \" The double quote character (this would normally terminate the string)
- The Backslash character itself (beware DOS paths) $\backslash \backslash$
- The Tab character \t
- The New Line character (or characters, depending on the system) \n
- The Carriage Return character (ASCII code 13) ١r

Conversion between data You can assign integer numbers to real variables and real numbers to integer variables types (unless the real number is out of the integer range when a run-time error will occur). When a real number is converted to an integer, it is truncated. The Asc(), Chr\$(), Str\$() and Val() functions convert between strings and numbers.

Arrays of data The three basic types (integers, reals and strings) can be made into one-dimensional and two-dimensional arrays. We call the one-dimensional array a vector, and the two-dimensional array a matrix. Declare arrays with the var statement:

var v[20], M[10][1000];

This declares a vector with 10 elements and a matrix with 10 rows and 1000 columns. To reference array elements, enclose the element number in square brackets (the first element in each dimension is number 0):

v[6] := 1; x := M[8][997];

Vector subsets Use v[start:size] to pass a vector or a subset of a vector v to a function. start is the index of the first element to pass, size is the number of elements. Both start and size are optional. If you omit start, 0 is used. If you omit size, the sub-set extends to the end of the vector. To pass the entire vector use v[0:], v[:], v[] or just v.

For example, consider the vector of real numbers declared as var data[15]. This has 15 elements numbered 0 to 14. To pass it to a function as a vector, you could specify:

data or data[]This is the entire vector. This is the same as data[:] or data[0:15].data[3:9]This is a vector of length 9, being elements 3 to 11.data[:8]This is a vector of length 8, being elements 0 to 7.

()	1	2	3	4	5	6	7	8	9	10	11	12	13	14	var data[15]
[data[3:9]

Matrix subsets With a matrix you have more options. You can pass a single element, a vector sub-set, or a matrix sub-set. Consider the matrix of real numbers defined as var sq[6][7]. You can pass this as a vector or a matrix to a function as (a, b, c and d are integer numbers):

sq[a][b:c] sq[a][] sq[a:b][c] sq[][c]	a vector of length c a vector of length 7 a vector of length b a vector of length 6
sq[a:b][c:d]	a matrix of size [b] [d]
sq or sq[][]	a matrix of size [6] [7]

This diagram shows how sub-sets are constructed. sq[1:4][0] is a 4 element vector. This could be passed to a function that expects a vector as an argument. sq[5][1:6] is a vector with 6 elements. sq[2:2][2:4] is a matrix, of size [2][4].



Transpose of a vector or
matrixYou can pass the transpose of a vector or matrix to a function with the trans() operator,
or by adding ` (back-quote) after the array or matrix name. The transpose of a matrix
swaps the rows and columns. To be consistent with normal matrix mathematics, a one-
dimensional array is treated as a column vector and is transposed into a matrix with 1
row. That is given var data[15], trans(data) is a matrix of size [1][15].

var M[5][3], v[5], W[5][5];	
<pre>PrintLog(M, M`);</pre>	'Print M and its transpose
<pre>PrintLog(M[][], trans(M[][]));</pre>	'Exactly the same as last line
<pre>MatMul(W, M, M[][]`);</pre>	'set W to M times its transpose
MatMul(W , v, v`);	'set W to v times its transpose

Diagonal of a matrix You can pass the diagonal of a matrix to a function using the diag() operator. This expects a matrix as an argument and produces a vector whose length is the smaller of the dimensions of the matrix. Given a matrix M[10][10], diag(M) is a 10 element vector.

Result views as arrays	The script language treats a result view as vectors of real numbers, one vector per channel. To access a vector element use View(v%, ch%).[index] where v% is the view, ch% is the channel and index is the bin number, starting from 0. You can pass a channel as an array to a function using View(v%, ch%).[], or View(v%, ch%).[a:b] to pass a vector subset starting at element a of length b. You can omit ch%, in which case channel 1 is used. You can also omit View(v%, ch%), in which case channel 1 in the current view is used. See the View() command for more information.
	If you change a visible result view, the modified area is marked as invalid and will update at the next opportunity.
Statement types	 The script language is composed of statements. Statements are separated by a semicolon. Semicolons are not required before else, endif, case, endcase, until, next, end and wend, or after end, endif, endcase, next and wend. White space is allowed between items in statements, and statements can be spread over several lines. Statements may include comments. Statements are of one of the following types: A variable or constant declaration An assignment statement of the form: variable := expression; Set the variable to the value of the expression
	 variable += expression; Add the expression value to the variable variable -= expression; Subtract the expression value from the variable variable *= expression; Multiply the variable by the expression value to the variable variable /= expression; Divide the variable by the expression value to the variable with strings (a\$+=b\$ is the same as a\$:=a\$+b\$, but is more efficient). A flow of control statement, described below A procedure call or a function with the result ignored, for example View (vh%);
Comments in a script	A comment is introduced by a single quotation mark. All text after the quotation mark is ignored up to the end of the line.
	View(vh%); 'This is a comment, and extends to the end of the line
Variable declarations	Variables are created by the var keyword. This is followed by a list of variable names. You must declare all variable names before you can use them. Arrays are declared as variables with the size of each dimension in square brackets. The first item in an array is at index 0. If an array is declared as being of size n, the last element is indexed $n-1$.
	<pre>var myInt%,myReal,myString\$; 'an integer, a real and a string var aInt%[20],arl[100],aStr\$[3] 'integer, real and string vectors var a2d[10][4]; '10 rows of 4 columns of reals var square\$[3][3]; '3 rows of 3 columns of strings</pre>
	You can define variables define in the main program, or in user-defined functions. Those defined in the main program are global and can be accessed from anywhere in the script after their definition. Variables defined in user-defined functions exist from the point of definition to the end of the function and are deleted when the function ends. If you have a recursive function, each time you enter the function you get a fresh set of variables.
	The dimensions of global arrays must be constant expressions. The dimensions of local arrays can be set by variables or calculated expressions. Simple variables (not arrays) can be initialised to constants when they are declared. The initialising expression may not

include variables or function calls. Uninitialised numeric variables are set to 0, uninitialised strings are empty.

var Sam%:=3+2, jim := 2.3214, sally\$:= "My name is \"Sally\"";

Constant declarations Constants are created by the const keyword. A constant can be of any of the three basic data types, and must be initialised as part of the constant declaration. Constants cannot be arrays. The syntax and use of constants is the same as for variables, except that you cannot assign to them or pass them to a function or procedure as a reference parameter.

const Sam%:=3+2, jim := 2.3214, sally\$:= "My name is \"Sally\"";

operators

Expressions and Anywhere in the script where a numeric value can be used, so can a numeric expression. Anywhere a string constant can be used, so can a string expression. Expressions are formed from functions, variables, constants, brackets and operators. In numerical expressions, the following operators are allowed, listed in order of precedence:

Numeric operators

	Operators	Names
Highest	`, [], ()	Matrix transpose, subscript, round brackets
	-, not	Unary minus, logical not
	*, /, mod	Multiply, divide and modulus (remainder)
	+, -	Add and subtract
	<, <=, >, >=	Less, less or equal, greater, greater or equal
	=, 🗇	Equal and not equal
	and, band	Logical and, bitwise and
Lowest	or, xor, bor, bxor	Logical or, exclusive or and bitwise versions

The order of precedence determines the order in which operators are applied within an expression. Without rules on the order of precedence, 4+2*3 could be interpreted as 18 or 10 depending on whether the add or the multiply was done first. Our rules say that multiply has a higher precedence, so the result is 10. If in doubt, use round brackets, as in 4+(2*3) to make your meaning clear. Extra brackets do not slow down the script.

The divide operator returns an integer result if both the divisor and the dividend are integers. If either is a real value, the result is a real. So 1/3 evaluates to 0, while 1.0/3, 1/3.0 and 1.0/3.0 all evaluate to 0.333333...

The minus sign occurs twice in the list because minus is used in two distinct ways: to form the difference of two values (as in fred:=23-jim) and to negate a single value (fred :=-jim). Operators that work with two values are called *binary*, operators that work with a single value are called *unary*. There are four unary operators, [], (), – and not, the remainder are binary.

There is no explicit TRUE or FALSE keyword in the language. The value zero is treated as false, and any non-zero value is treated as true. Logical comparisons have the value 1 for true. So not 0 has the value 1, and the not of any other value is 0. If you use a real number for a logical test, remember that the only way to guarantee that a real number is zero is by assigning zero to it. For example, the following loop may never end:

```
var add:=1.0;
repeat
  add := add - 1.0/3.0; ' beware, 1/3 would have the value 0!
                         ' beware, add may never be exactly 0
until add = 0.0;
```

Even changing the final test to add<=0.0 leads to a loop that could cycle 3 or 4 times depending on the implementation of floating point numbers.

The result of the comparison operators is integer 0 if the comparison is false and integer 1 if the comparison is true. The result of the binary arithmetic operators is integer if both operands are integers, otherwise the result is a real number. The result of the logical operators is integer 0 or 1. The result of the exclusive or operator is true if one operand is true and the other is false.

The bitwise operators band, bor and bxor treat their operands as integers, and produce an integer result on a bit by bit basis. They are not allowed with real number operands.

String operators

	Operators	Names
Highest	+	Concatenate
	<, <=, >, >=	Less, less or equal, greater, greater or equal
Lowest	=, <>	Equal and not equal

The comparison operators can be applied to strings. Strings are compared character by character, from left to right. The comparison is case sensitive. To be equal, two strings must be identical. You can also use the + operator with strings to concatenate them (join them together). The character order for comparisons (lowest to highest) is:

> space !"#\$%&'()*+,-./0123456789:;<=>?@ ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^ abcdefghijklmnopqrstuvwxyz{|}~

Do not confuse assignment := with the equality comparison operator, =. They are entirely different. The result of an assignment does not have a value, so you cannot write statements like a:=b:=c; (which was allowed in the MS-DOS version of Spike2).

Examples of expressions The following (meaningless) code gives examples of expressions.

var jim,fred,sam,sue%,pip%,alf\$,jane\$; jim := Sin(0.25) + Cos(0.25);fred := 2 + 3 * 4; 'Result is 14.0 as * has higher precedence fred := (2 + 3) * 4;'Result is 20.0 fred += 1;'Add 1 to fred sue% := 49.734; 'Result is 49 sue% := -49.734; 'Result is -49 pip% := 1 + fred > 9; 'Result is 1 as 21.0 is greater than 9 jane\$:= "Jane"; alf\$:= "alf"; sam := jane\$ > alf\$; 'Result is 0.0 (a is greater than J) sam := UCase\$(jane\$)>UCase\$(alf\$); 'Result is 1.0 (A < J)</pre> sam := "same" > "sam"; 'Result is 1.0 pip% := 23 mod 7; 'Result is 2 'Result is 3.5 jim := 23 mod 6.5; 'Result is -2.0 jim := -32 mod 6; sue% := jim and not sam; 'Result is 0 (jim = -2.0 and sam = 1.0) pip% := 1 and 0 or 2>1; 'Result is 1 sue% := 9 band 8; 'Result is 8 (9=1001 in binary, 8=1000) sue% := 9 bxor 8; 'Result is 1 sue% := 9 bor 8; 'Result is 9

Mathematical constants We don't provide maths constants as built-in symbols, but the two most common ones, e and π are easily generated within a script; *e* is Exp(1.0) and π is 4.0*ATan(1.0).

Flow of control statements If scripts were simply a list of commands to be executed in order, their usefulness would be severely limited. The flow of control statements let scripts loop and branch. It is considered good practice to keep flow of control statements on separate lines, but the script syntax does not require this. There are two branching statements, if...endif and docase...endcase, and three looping statements, repeat...until, while...wend and for...next. You can also use user-defined functions and procedures with the Func and Proc statements.

if...endif The if statement can be used in two ways. When used without an else, a single section of code can be executed conditionally. When used with an else, one of two sections of code is executed. If you need more than two alternative branches, the docase statement is usually more compact than nesting many if statements.

```
if expression then 'The simple form of an if
  zero or more statements;
endif;
if expression then 'Using an else
  zero or more statements;
else
  zero or more statements;
endif;
```

If the expression is non-zero, the statements after the then are executed. If the expression is zero, only the statements after the else are executed. The following code adds 1 or 2 to a number, depending on it being odd or even:

docase...endcase These keywords enclose a list of case statements forming a multiway branch. Each case is scanned until one is found with a non-zero expression, or the else is found. If the else is omitted, control passes to the statement after the endcase if no case expression is non-zero. Only the first non-zero case is executed (or the else if no case is non-zero).

```
docase
  case exp1 then
    statement list;
  case exp2 then
    statement list;
    ...
  else
    statement list;
endcase;
```

The following example sets a string value depending on the value of a number:

```
var base%:=8,msg$;
docase
case base%=2 then msg$ := "Binary";
case base%=8 then msg$ := "Octal";
case base%=10 then msg$ := "Decimal";
case base%=16 then msg$ := "Hexadecimal";
else msg$ := "Pardon?";
endcase;
```

repeat...until The statements between repeat and until are repeated until the expression after the until keyword evaluates to non-zero. The body of a repeat loop is always executed at least once. If you need the possibility of zero executions, use a while loop. The syntax of the statement is:

```
repeat
   zero or more statements;
until expression;
```

For example, the following code prints the times of all data items on channel 3 (plus an extra -1 at the end):

```
var time := -1; 'start time of search
repeat
time := NextTime(3, time); 'find next item time
PrintLog("%f\n", time); 'display the time to the log
until time<0; 'until no data found</pre>
```

while...wend The statements between the keywords are repeated while an expression is not zero. If the expression is zero the first time, the statements in between are not executed. The while loop can be executed zero times, unlike the repeat loop, which is always executed at least once.

```
while expression do
    zero or more statements;
wend;
```

The following code fragment, finds the first number that is a power of two that is greater than or equal to some number:

for...next A for loop executes a group of statements a number of times with a variable changed by a fixed amount on each iteration. The loop can be executed zero times. The syntax is:

```
for v := exp1 to exp2 {step exp3} do
    zero or more statements;
next;
```

- v This is the loop variable and may be a real number, or an integer. It must be a simple variable, not an array element.
- exp1 This expression sets the initial variable value before the looping begins.
- exp2 This expression is evaluated once, before the loop starts, and is used to test for the end of the loop. If step is positive or omitted (when it is 1), the loop stops when the variable is greater than exp2. If step is negative, the loop stops when the variable is less than exp2.
- exp3 This expression is evaluated once, before the loop starts, and sets the increment added to the variable when the next statement is reached. If there is no step exp3 in the code, the increment is 1. The value of exp3 can be positive or negative.

The following example prints the squares of all the integers between 10 and 1:

```
var num%;
for num% := 10 to 1 step -1 do
    PrintLog("%d squared is %d\n", num%, num% * num%);
next;
```

If you want a for loop where the end value and/or the step size are evaluated each time round the loop you should use a while...wend or repeat...until construction.

Halt The Halt keyword terminates a script. A script also terminates if the path of execution reaches the end of the script. When a script halts, any open external files associated with the Read() or Print() functions are closed and any windows with invalid regions are updated. Control then returns to the user.

Functions and A user-defined function is a named block of code. It can access global variables and procedures create its own local variables. Information is passed into user-defined functions through arguments. Information can be returned by giving the function a value, by altering the values of arguments passed by reference or by changing global variables.

> User-defined functions that return a value are introduced by the func keyword, those that do not are introduced by the proc keyword. The end keyword marks the end of a function. The return keyword returns from a function. If return is omitted, the function returns to the caller when it reaches the end statement. Arguments can be passed to functions by enclosing them in brackets after the function. Functions that return a value or a string have names that specify the type of the returned value. A function is defined as:

```
func name({argument list})
                                 or
                                        proc name({argument list})
                                        {var local-variable-list;}
{var local-variable-list;}
statements including return x;
                                        statements including return;
end;
                                        end;
```

There is no semicolon at the end of the argument list because the argument list is not the end of the func or proc statement; the end keyword terminates the statement. Functions may not be nested within each other.

- **Argument lists** The argument list is a list of variable names separated by commas. There are two ways to pass arguments to a function: by value and by reference:
 - Value Arguments passed by value are local variables in the function. Their initial values are passed from the calling context. Changes made in the function to a variable passed by value do not affect the calling context.
 - Reference Arguments passed by reference are the same variables (by a different name) as the variables passed from the calling context. Changes made to arguments passed by reference do affect the calling context. Because of this, reference arguments must be passed as variables (not expressions or constants) and the variable must match the type of the argument (except we allow you to pass a real variable where an integer variable is expected).

Simple (non-array) variables are passed by value. Simple variables can be passed by reference by placing the & character before the name in the argument list declaration.

Arrays and sub-arrays are always passed by reference. Array arguments have empty square brackets in the function declaration, for example one [] for a vector and two [][] for a matrix. The number of dimensions of the passed array must match the declaration. The array passed in sets the size of each dimension. You can find the size of an array with the Len () function. An individual array element is a treated as a simple variable.

If you use the trans () or diag () operators to pass the transpose or diagonal of an array to a function, the array is still passed by reference and changes made in the function to array elements will change the corresponding elements in the original data.

return The return keyword is used in a user-defined function to return control to the calling context. In a proc, the return must not be followed by a value. In a func, the return should be followed by a value of the correct type to match the function name. If no return value is specified, a func that returns a real or integer value returns 0, and a func that returns a string value returns a string of zero length.

Examples of user-defined functions	<pre>proc PrintInfo() PrintLog(ChanTitle\$(1)); PrintLog(ChanComment\$(1)); return; end;</pre>	<pre>'no return value, no arguments 'Show the channel title 'and the comment 'return is optional in this case as 'end forces a return for a proc</pre>
	<pre>func sumsq(a, b) return a*a + b*b; end;</pre>	'sum the square of the arguments
	<pre>func removeExt\$(name\$) var n := 0, k := 1; repeat</pre>	'remove text after last . in a string
	<pre>k:=InStr(name\$,".",k); if (k > 0) then</pre>	'find position of next dot 'if found a new dot 'remember where
	until k=0;	'until all found
	return name\$; else	'no extension
	<pre>return Left\$(name\$,n-1); end;</pre>	;
	<pre>proc sumdiff(&arg1, &arg2) arg1 := arg1 + arg2; arg2 := arg1 - 2*arg2;</pre>	'returns sum and difference of args 'sum of arguments 'original arg1-arg2
	return; end;	'results returned via arguments
	<pre>func sumArr(data[]) var sum:=0.0; var i%; function 0 to Log(data[]) 1</pre>	<pre>'sum all elements of a vector 'initialise the total 'index ''</pre>
	<pre>for 1%:=0 to Len(data[])-1 sum := sum + data[i%]; next; return sum;</pre>	do 'of course, ArrSum() is much faster!
	end;	
	<pre>Func SumArr2(data[][]) var rows%,cols%,i%,sum; rows% := Len(data[0][]) cols% := Len(data[1][0]);</pre>	'Sum of all matrix elements 'sizes, index and sum, all set to 0 'get sizes of dimensions
	<pre>if rows%>cols% then for i%:=0 to cols%-1 do st else</pre>	<pre>'choose more efficient method um += ArrSum(data[i%][]) next;</pre>
	<pre>for i%:=0 to rows%-1 do su endif; return sum; end;</pre>	um += ArrSum(data[][i%]) next;
	Variables declared within a functi- cannot be used from elsewhere. Y functions. Each variable is separate it calls itself), each time vou enter th	on exist only within the body of the function. They ou can use the same name for variables in different . In addition, if you call a function recursively (that is he function, you have a fresh set of variables.
Scope of user-defined functions	Unlike global variables, which are are declared onwards, and local function only, user-defined function define two functions that call each o	only visible from the point in the script in which they variables, which are visible within a user-defined ns are visible from all points in the script. You may other, if you wish.

Functions as arguments The script language allows a function or procedure to be passed as an argument. The function declaration includes the declaration of the function type to be passed. Functions and procedures can occur before or after the line in which they are used as an argument.

```
proc Sam(a,b$,c%)
. . .
end;
func Calc(va)
return 3*va*va-2.0*va;
end;
func PassFunc(x, func ff(realarg))
return ff(x);
end;
func PassProc(proc jim(realarg, strArg$, son%))
jim(1.0,"hello",3);
end;
val := PassFunc(1.0, Calc);
                              'pass function
PassProc( Sam );
                               'pass procedure
```

The declaration of the procedure or function argument is exactly the same as for declaring a user-defined function or procedure. When passing the function or procedure as an argument, just give the name of the function or procedure; no brackets or arguments are required. The compiler checks that the argument types of a function passed as an argument match those declared in the function header. See the ToolbarSet() function for an example.

Although user-defined functions and built-in functions are very similar, you are not allowed to pass a built-in function as an argument to a built-in function. Further, you cannot pass a built-in function to a user-defined function if it has one or more arguments that can be of more than one type. For example, the built-in sin() function can accept a real argument, or a real array argument, and so cannot be passed.

- **Channel specifiers** Several built-in script commands use a channel specifier to define a list of 1 or more channels. This argument is always called cSpc in the documentation. This argument stands for a string, an integer array or an integer.
 - CSpc\$ This holds a list of channel numbers and channel ranges, separated by commas. A channel range is a start channel number followed by an end channel number separated by two dots. The end channel number can be less that the start channel number. For example "1, 3, 5...7, 12...9" is a list of channels 1, 3, 5, 6, 7, 9, 10, 11 and 12.
 - cSpc%[] The first element of the vector is the number of channels. The remaining elements are channel numbers. It is an error for the vector passed in to be shorter than the number of channels + 1.
 - cSpc% This is either a channel number, or -1 for all channels, -2 for visible channels or -3 for selected channels.

In a result view, channel numbers start at 1, but we accept 0 (meaning 1) to avoid breaking scripts written for older versions of Spike2. Some commands expect channels of specific types; channels that do not meet the type requirements are removed from the list. It is usually an error for a channel specification to generate an empty list.

Functional	command
	aroups

This section of the manual lists commands by function. The next section lists the command alphabetically with a description of the command arguments and operation.

Windows and views These commands open, close, manipulate, position, display, size, colour and create windows (views). These commands apply to windows in general. See the sections on Time, Result, XY and Text windows for more specific commands.

aqA	Get the application view handle and version and special views
ChanNumbers	Hide and show channel numbers
Colour	Get or set the palette entry associated with a screen item
Draw	Draw invalid regions of the view (and set x axis range)
DrawAll	Update all invalid regions in all views
EditCopy	Copy window to clipboard
FileClose	Closes a window or windows
FileCommentS	Gets and sets the file comment for time views
FileConvertS	Convert a foreign file to a Spike? file and open it
FileNameŚ	Gets the file name associated with a window
FileNew	Opens an output file or a new text or data window
FileOpen	Opens an existing file (in a window)
FilePrint	Prints a range of data from the current view
FilePrintVisible	Prints the current view
FilePrintScreen	Prints all text-based time and result views
FileOuit	Closes the application
FileSave	Save a view
FileSaveAs	Save a view in variety of formats
FontGet	Read back information about the font
FontSet	Set the font for the current window
FocusHandle	Get handle of script-controllable window with the focus
FrontView	Get or set the front window on screen
Grid	Get or set the visibility of the axis grid
LogHandle	Gets the view handle of the log window
PaletteGet	Get the RGB colour of a palette entry
PaletteSet	Set the RGB colour of a palette entry
View	Change or override current view and get view handle
ViewColour	Override application colours for a view
ViewFind	Get a view handle from a view title
ViewKind	Get the type of the current view
ViewList	Form a list of handles to views that meet a specification
ViewStandard	Returns a window to a standard state
ViewUseColour	Get and set monochrome/colour use for view
Window	Sets the window size and position
WindowGetPos	Get window position
WindowSize	Changes the window size
WindowTitle\$	Gets or changes the window title
WindowVisible	Sets or gets the visibility of the window (hide/show)
XAxis	Get or set visibility of the x axis
XAxisMode	Control how the x axis is drawn
XAxisStyle	Control x axis tick spacing and time view hh:mm:ss mode
XScroller	Get or set visibility of x axis scroll bar
XUnits\$	Get units of the x axis
YAxis	Get or set the visibility of the y axis
YAxisMode	Control how the y axis is drawn
YAxisStyle	Control y axis tick spacing

Channel commands These commands manipulate data channels.

 1	
BinError	Get or set error bar information
Binsize	The sample interval for waveforms, otherwise time resolution
BurstMake	Extract bursts to a memory channel from an event channel
BurstRevise	Modify a list of bursts in a memory channel
BurstStats	Collect statistics on bursts in preparation for BurstRevise
Chan\$	Get the symbolic name for a channel
ChanCalibrate	Calibrate waveform or WaveMark channels from known values
ChanColour	Override drawing mode based channel colours
ChanComment\$	Get or set the channel comment
ChanData	Fill an array with a waveform or event times
ChanDelete	Delete a channel from a time view
ChanDuplicate	Duplicate channels in a time view
ChanHide	Make a channel invisible
ChanKind	Get the type of a channel
ChanList	Get a list of channels meeting a specification
ChanMeasure	Make the same measurements as the Cursor Regions dialog
ChanNew	Create a new channel for ChanWriteWave() to write to
ChanOffset	Set waveform and WaveMark value for input of zero
ChanOrder	Modify the channel order and y axis grouping
ChanProcessAdd	Add new processing option to a waveform or RealWave channel
ChanProcessArg	Read or modify an argument of a channel process
ChanProcessClear	Delete a process, all processes for a channel, or all processes
ChanProcessInfo	Get information about the processes attached to a channel
ChanPort	Get the physical sampling port of a time view data channel
ChanScale	Set waveform and WaveMark channel scale factor
ChanSave	Copy time view channels between views with optional time shift
ChanSearch	Search a channel for a feature (peak, level crossing, slope)
ChanSelect	Select and report on selected state of channels
ChanShow	Make a channel visible
ChanTitle\$	Get or set the channel title string
ChanUnits\$	Get or set the channel units
ChanValue	Get channel data at a particular time or x axis position
ChanVisible	Get the visibility state of a channel
ChanWeight	Change the relative vertical space of a channel
ChanWriteWave	Write data to a waveform or RealWave channel
Count	Count items in a time range, sum bins in result view
DrawMode	Get or set display mode for a channel
DupChan	Get information about duplicate channels
EventToWaveform	Convert an event channel into a waveform channel
FitLine	Fit a straight line to waveform or result channel
LastTime	Find the previous item in a channel (and return values)
MarkEdit	Edit marker codes and other information
MarkInfo	Get information on extended marker types
MarkMask	Set the marker filter for a channel
MarkSet	Set the marker codes for data in a time range
Maxtime	Time of last item on the channel
MemChan	Create a channel in memory
MemDeleteItem	Delete one or more items from a memory channel
MemDeleteTime	Delete one or more items based on time in a memory channel
MemGetItem	Get information on item in memory channel
MemImport	Import items into a memory channel
MemSave	Write the contents of a memory channel to the data file
MemSetItem	Edit or add an item in a memory channel
Minmax	Find minimum and maximum values (and positions)
NextTime	Find the next item in a channel (and return values)

-

ime views	These command	ls manipu	late time views.
-----------	---------------	-----------	------------------

	BinToX	Convert time units to x axis units
	DrawMode	Get or set display mode for a channel
	Dup	Get the view handle of duplicates of the current view
	EventToWaveform	Create a waveform channel from events
	ExportChanFormat	Set channel format for export
	ExportChanList	Set list of channels for export
	ExportRectFormat	Set rectangular time view export
	ExportTextFormat	Set format for text output of channels
	FileApplyResource	Apply a resource file to the current time view
	FileGlobalResource	Set a global resource file
	FileSaveResource	Create a resource file matching the current time view
	Maxtime	Maximum x axis value for any channel in the view
	Optimise	Set reasonable y range for channels with axes
	ReRun	Get or set the rerunning state of a channel
	ViewTrigger	Control on-line triggered displays and off-line display stepping
	WindowDuplicate	Duplicate a time view
	XLow	Time of the start of the displayed area in seconds
	XHigh	Time of the end of the displayed area in seconds
	XRange	Set the x axis range for next draw
	XScroller	Show or hide the x axis scroll bar and controls
	XTitle\$	Get x axis title
	XToBin	Convert x axis units to time units
	YAxisLock	Get and set the locking state of grouped channels
	YRange	Set y axis range for a channel
	YHigh	Get upper limit of y axis for a channel
	YLow	Get lower limit of y axis for a channel
Beault viewe	These commands manit	pulate result views. Desult views can also be treated as arrays: so
Result views	all the array arithmetic	pommands are available
		Cat ar act area has information
	BinError	Width of each bin in the w direction
	Binsize	with of each off in the x direction
	BINTOX DucuMada	Convent on number to x axis value
	DrawMode	Number of hims in the result view
	Maxtime	Find activities and accessing and accessing the second accessing and the second accessing and the second accessing a
	Minmax	Find minimum and maximum values (and positions)

DrawMode	Get or set display mode for a channel
Maxtime	Number of bins in the result view
Minmax	Find minimum and maximum values (and positions)
Optimise	Set reasonable y range for display
RasterAux	Deprecated, use RasterSort() and RasterSymbol()
RasterGet	Read back raster information
RasterSet	Set raster information
RasterSort	Get and set the sorting values
RasterSymbol	Get and set the time values for raster symbols
Sweeps	Gets and sets the number of items added to result view
XHigh	Bin number of the end of the displayed area
XLow	Bin number of the start of the displayed area
XRange	Set the x axis range for next draw
XScroller	Show or hide the x axis scroll bar and controls
XTitle\$	Get and set x axis title
XToBin	Convert x axis value to bin number
YAxisLock	Get and set the locking state of grouped channels
YHigh	Get upper limit of y axis
YLow	Get lower limit of y axis
YRange	Set y axis range

XY views	XY views hold from 1 to 256 channels of (x, y) co-ordinate pairs that can be displayed in
	a variety of styles. Most functions that work for views in general will work for an XY
	view. The following are the XY view specific commands.

MeasureToXY	Create an XY view and associated measurement process
XYAddData	Add data points to a channel of an XY view
XYColour	Set the colour of a channel
XYCount	Return the number of XY data points in a channel
XYDelete	Delete one or more data points from a channel
XYDrawMode	Control how a channel is drawn
XYGetData	Read data back from an XY channel
XYInCircle	Count the number of XY points within a circle
XYInRect	Count the number of XY points inside a rectangle
XYJoin	Get or set the point joining method
ХҮКеу	Control the display of the XY view channel key
XYRange	Get rectangle containing one or more channels
XYSetChan	Create or modify an XY channel
XYSize	Get or set maximum size of an XY channel
XYSort	Change the sort (and draw) order of a channel

Vertical cursors The following commands control the vertical cursors.

-	
Cursor	Set or get the position of a cursor
CursorActive	Set or get the active cursor mode
CursorActiveGet	Read back active cursor parameters
CursorDelete	Delete a designated cursor
CursorExists	Test if a cursor exists
CursorLabel	Set or get the cursor label style
CursorLabelPos	Set or get the cursor label position
CursorNew	Add a new cursor (at a given position)
CursorRenumber	Renumber the cursors in ascending position order
CursorSearch	Trigger active cursor searches from the script
CursorSet	Set the number (and position) of vertical cursors
CursorValid	Test if an active cursor search succeeded
CursorVisible	Get or set cursor visibility

Horizontal cursors Horizontal cursors are controlled from a script with the following commands.

HCursor	Set or get the position of a horizontal cursor
HCursorChan	Gets the channel that a horizontal cursor belongs to
HCursorDelete	Delete a designated horizontal cursor
HCursorExists	Test if a cursor exists
HCursorLabel	Gets or sets the horizontal cursor style
HCursorLabelPos	Gets or sets the horizontal cursor position
HCursorNew	Add a new horizontal cursor on a channel (at a given position)
HCursorRenumber	Renumbers the cursors from bottom to the top of the view

Editing operations These functions mimic the Edit menu commands and provide additional functionality.

EditClear	Delete text from a text window at the caret
EditCopy	Copy the current selection to the clipboard
EditCut	Delete the current selection to the clipboard
EditPaste	Paste the clipboard into the current text field
EditSelectAll	Select the entire text or cursor window contents
MoveBy	Move relative to current position
MoveTo	Move to a particular place
Selection\$	This function returns the text that is currently selected

Text files Spike2 can create, read and write text files. You can also open a text file into a window. FileNew Open a new text file in a window Open a text file in a window or for reading and writing FileOpen Save a view in variety of formats, including text FileSaveAs Write formatted output to a file or log window Print Read Extract data from a text file Set separators and delimiters for Read and ReadStr ReadSetup

Binary files Spike2 can read and write binary files. Binary files provide links to other programs and are generally more efficient than text for transferring large quantities of data.

FileClose	Close a file opened in binary mode
FileOpen	Open an external file in binary mode
BRead	Extract 32-bit integer, 64-bit real and string data from a file
BReadSize	Extract 8 and 16-bit integer and 32-bit real data from a file
BRWEndian	Change byte order for read/write of numbers in binary files
BSeek	Change the current file position for next read or write
BWrite	Write 32-bit integer and 64-bit real data to a file
BWriteSize	Write 8 and 16-bit integer, 32-bit real and string data to a file

File system Spike2 can read file information, change the current directory, delete and copy files and convert foreign files into Spike2 format. You can also execute external files.

FileCo	nvert\$ Co	nvert a foreign file to a Spike2 file and open it
FileCo	ру Сс	py one file to a new location
FileDa	te\$ Re	trieve date of creation of Spike2 data file
FileDe	lete De	lete one or more files
FileLi	st Ge	t a list of files or directories
FilePa	th\$ Ge	t the current directory or directory for new data files
FilePa	thSet Ch	ange the current directory or directory for new data files
FileTi	me\$ Re	trieve time of creation of Spike2 data file
FileTi	meDate Re	trieve time and date of creation of Spike2 data file as numbers

commands

User interaction These commands exchange information with the user or let the user interact with the data.

Inkey	Return key code of last key user pressed
Input	Prompt user for a number in a defined range
Input\$	Prompt user for a string with a list of acceptable characters
Interact	Allow user to interact with data
Keypress	Detect if Inkey function would read a character
Message	Display a message in a box, wait for OK
Print	Formatted text output to a file or window
PrintLog	Formatted text output to the Log window
Print\$	Formatted text output to a string
Query	Ask a user a question, wait for response
Sound	Play a tone or a .wav file
Speak	Convert text to speech on systems that support this
Yield	Give idle time to the system; delay for a time
YieldSystem	Surrender current time slice and sleep Spike2

You can build simple dialogs, with a set of fields stacked vertically or you can build freeformat dialogs (but with more work to define the positions of all the dialog fields):

DlgAllow	Set allowed actions and change and idle functions
DlgButton	Add buttons to the dialog and modify existing buttons
DlgChan	Define a dialog entry as prompt and channel selection
DlgCheck	Define a dialog item as a check box
DlgCreate	Start a dialog definition
DlgEnable	Enable and disable dialog items in change or idle functions
DlgGroup	Position a group box in the dialog
DlgInteger	Define a dialog entry as prompt and integer number input
DlgLabel	Define a dialog entry as prompt only
DlgList	Define a dialog entry as prompt and selection from a list
DlgReal	Define a dialog entry as prompt and real number input
DlgShow	Display the dialog, get values of fields
DlgString	Define a dialog entry as prompt and string input
DlgText	Define a fixed text string for the dialog
DlgValue	Gives access to dialog fields in change and idle functions
DlgVisible	Show and hide dialog items in change or idle functions
DlgXValue	Define a dialog entry to collect an x axis value

These commands control the various toolbars and link script functions to the toolbar.

Арр	Get the view handle of the toolbars
Toolbar	Let the user interact with the toolbar
ToolbarClear	Remove all defined buttons from the toolbar
ToolbarEnable	Get or set the enables state of toolbar buttons
ToolbarSet	Add a button (and associate a function with it)
ToolbarText	Display a message using the toolbar
ToolbarVisible	Get or set the visibility of the toolbar
SampleBar	Controls the sample toolbar buttons
ScriptBar	Controls the script toolbar buttons

Analysis These functions create new result views and analyse data into the result views.

BinError	Get or set error bar information
MeasureToXY	Create an XY view and associated measurement process
MeasureX	Set the x part of a measurement
MeasureY	Set the y part of a measurement
RasterAux	Set additional data for sorted rasters (deprecated)
RasterGet	Read back raster information
RasterSet	Set raster information
RasterSort	Get and set the sorting values
RasterSymbol	Get and set the time values for raster symbols
SetEvtCrl	Create result view for event correlation
SetEvtCrlShift	Set shift for shuffled event correlation
SetINTH	Create result view for interval histogram
SetPhase	Create result view for phase histogram
SetPower	Create result view for power spectrum
SetPSTH	Create result view for post/peri-stimulus histogram
SetWaveCrl	Create result view for waveform correlation
SetWaveCrlDC	Set DC use/ignore for SetWaveCrl()
SetAverage	Create result view for waveform copy/average/sum
SetResult	Create result view for user-defined data
Sweeps	Number of items added to result view
Process	Process data into result view
ProcessAll	Process all result views attached to a time view
ProcessAuto	Process data automatically during sampling
ProcessTriggered	Process triggered by data during sampling

Mathematical The following mathematical functions are built into Spike2. You can apply many of these functions functions to real arrays by passing an array to the function.

Abs	Absolute value of a number or array
ATan	Arc tangent of number or array
Cos	Cosine of a number or array
Exp	Exponential function of a number or array
Frac	Remove integral part of a number or array
GammaP	Incomplete Gamma function, used to calculate probabilities
Ln	Natural logarithm of a number or array
LnGamma	Natural logarithm of the Gamma function (use for big factorials)
Log	Logarithm to base 10 of a number or array
Max	Finds maximum of array or variables
Min	Finds minimum of array or variables
Pow	Raise a number or an array to a power
Rand	Generate pseudo-random numbers with uniform density
RandExp	Generate random numbers with exponential density
RandNorm	Generate random numbers with normal density
Round	Round a real number to the nearest integral value
Sin	Sine of a number or array
Sqrt	Square root of a number or an array
Tan	Tangent of a number or array
Trunc	Remove fractional part of number or array

Array and matrix These functions are used with arrays and result views to speed up data manipulation. You **arithmetic** can also use the built-in mathematical functions directly on an array.

ArrAdd	Adds an array or constant to an array
ArrConst	Copies an array, or sets an array to a constant value
ArrDiff	Replaces an array with an array of simple differences
ArrDiv	Divides an array by another array or a constant
ArrDivR	Divides array into another array or constant
ArrDot	Forms the dot product (sum of products) of two arrays
ArrFFT	Fourier transforms and related operations
ArrFilt	Applies a FIR filter to an array
ArrIntgl	Integrates array; inverse of ArrDiff()
ArrMul	Multiples an array by another array or constant
ArrSort	Sort an array and optionally order others in the same way
ArrSpline	Interpolate one array to another using cubic splines
ArrSub	Subtract constant from array, or difference of two arrays
ArrSubR	Subtract array from constant, or reversed difference of arrays
ArrSum	Sum, mean and standard deviation of an array
Len	Returns the length of a string or array
MATDet	Calculate the determinant of a matrix
MATInv	Invert a matrix
MATMul	Multiply matrices and vectors.
MATSolve	Solve a set of linear equations
MATTrans	Transpose a matrix (also see the trans () operator)
PCA	Principal component analysis (singular value decomposition)

String functions These functions manipulate strings and convert between strings and other variable types.

Asc	ASCII code of first character of a string
Chr\$	Converts a code to a one character string
DelStr\$	Returns a string minus a sub-string
InStr	Searches for a string in another string
LCase\$	Returns lower case version of a string
Left\$	Returns the leftmost characters of a string
Len	Returns the length of a string or array
Mid\$	Returns a sub-string of a string
Print\$	Produce formatted string from variables
ReadStr	Extract variables from a string
ReadSetup	Set separators and delimiters for ReadStr and Read
Right\$	Returns the rightmost characters of a string
Str\$	Converts a number to a string
UCase\$	Returns upper case version of a string
Val	Converts a string to number

ChanFitFit function to time, result or XY chaFitExpFits multiple exponentials		Fit function to time, result or XY channels and display fitted line
		Fits multiple exponentials
	FitGauss	Fit multiple Gaussians (normal distributions)
	FitLine	Fit straight line to time or result view data
	FitLinear	Fit to linear combination of user-defined functions
	FitNLUser	Fit to non-linear user function
	FitPoly	Fit polynomial to the data
	FitSin	Fit multiple sinusoids

Digital filtering These functions create and apply digital filters and manipulate the filter bank.

_		
	ArrFilt	Array arithmetic routine to apply FIR coefficients to an array
	FiltApply	Apply a set of coefficients or a filter bank filter to a waveform
	FiltAtten	Set the desired attenuation of a filter in the filter bank
	FiltCalc	Force coefficient calculation of a filter in the filter bank
	FiltComment\$	Get or set comment for a filter in the filter bank
	FiltCreate	Create a new filter definition in the filter bank
	FiltInfo	Retrieve information about a filter in the filter bank
	FiltName\$	Get or set the name of a filter in the filter bank
	FiltRange	Get the useful sampling rate range for a filter bank filter
	FIRMake	Generate FIR filter coefficients in an array
	FIRQuick	Generate FIR filter coefficients with desired attenuation
	FIRResponse	Calculate frequency response of array of coefficients
	IIRBp	Create and/or apply an IIR bandpass filter
	IIRBs	Create and/or apply an IIR bandstop filter
	IIRHp	Create and/or apply an IIR highpass filter
	IIRLp	Create and/or apply an IIR lowpass filter
	IIRNotch	Create and/or apply an IIR notch filter
	IIRReson	Create and/or apply an IIR resonator filter

Sampling configuration

Sampling These commands set the sampling configuration to use when you create a new data file.

configuration	SampleAutoComment	Control automatic prompt for file comment after sampling
•	SampleAutoCommit	Controls how frequently a data file is flushed to disk
	SampleAutoFile	Set if file automatically written to disk at end of sampling
	SampleAutoName\$	Set file name template for automatic file saving
	SampleCalibrate	Set calibration of waveform or WaveMark channel
	SampleComment\$	Set or get the channel comment
	SampleClear	Sets the Sampling configuration to a known state
	SampleEvent	Adds a channel to sample event data
	SampleDigMark	Adds a channel to sample digital marker data
	SampleLimitSize	Set or clear the size limit on a file
	SampleLimitTime	Set or clear the time limit for sampling
	SampleMode	Set the sampling mode (Continuous, Timed, Triggered)
	SampleOptimise	Set methods for optimising waveform and WaveMark rates
	SampleRepeats	Set the number of files to sample in sequence
	SampleSeqCtrl	Set or get the allowed sequencer jump methods
	SampleSequencer	Set the name of the sequencer file
	SampleSequencer\$	Get the name of the sequencer file
	SampleStartTrigger	Set triggered or non-triggered sample start
	SampleTextMark	Adds a channel for text notes
	SampleTimePerAdc	Set the base ADC conversion interval
	SampleTitle\$	Set and get the title of a channel
	SampleTrigger	Enable peri-trigger sampling of groups of channels
	SampleUsPerTime	Set and get the basic timing interval
	SampleWaveform	Adds a channel to sample waveform data
	SampleWaveMark	Adds a channel to sample WaveMark data

Runtime	samp	lina	The

These commands control data sampling. There is only one new data file at a time, and these commands refer to it, regardless of the current view.

SampleAbort	Stop sampling and throw data away
SampleHandle	Gets the view handle of sampling windows and controls
SampleKey	Adds to the keyboard marker channel, controls output sequencer

	SampleReset SampleSeqStep SampleSeqTable SampleSeqVar SampleStart SampleStatus SampleStop SampleText SampleWrite	Clear the current data file and restart sampling Get the current sequencer step Get and set output sequencer table data Set an output sequencer variable Start sampling after creating a new time view Get the sampling state Stop sampling and keep the data Adds a string to the text marker channel Control writing data to sampling file
Arbitrary waveform output	These commands conta Waveforms are identified PlayOffline PlayWaveAdd PlayWaveChans PlayWaveCopy PlayWaveCopy PlayWaveCycles PlayWaveDelete PlayWaveEnable PlayWaveEnable PlayWaveEnable PlayWaveLink\$ PlayWaveLink\$ PlayWaveRate PlayWaveStatus\$ PlayWaveStop PlayWaveStop PlayWaveTrigger	rol the output of waveforms during data sampling and offline. ed by a key code and up to 10 can be defined. Play offline via 1401 or sound card Add waveform to the on-line play list Get or change the DAC channels of an area before sampling Update output waveforms in 1401 after sampling starts Get or change the repeats of a waveform Delete one or more waveforms from the play wave list Controls which waveforms in the list are available for playing Get waveform information including list of all waveform keys Get or set the label associated with each waveform Get or set the links between waveforms Get or set the desired base wave replay rate Scale factor for the wave replay rate - can be used on-line Get information about waveform output during sampling Stop the current output immediately or set one more cycle Get or set the trigger state for a waveform
Discriminator control	These functions control DiscrimChanGet DiscrimChanSet DiscrimClear DiscrimLevel DiscrimMode DiscrimMonitor DiscrimTimeOut	the CED 1401-18 discriminator card. Get the settings of a discriminator channel Change discriminator channel settings Set all channels to a known state Get or set discriminator thresholds Get or set discriminator modes Get or set the waveform monitor channel Get or set the time-out for modes 7 and 8
Signal conditioner control	These functions control CondFilter CondFilterList CondGain CondGainList CondGet CondOffset CondOffsetLimit CondRevision\$ CondSet CondSourceList CondType	serial line controlled signal conditioners. Get or set the conditioner low-pass or high-pass filter Get a list of possible low-pass or high-pass filter settings Get or set the conditioner gain Get a list of the possible gains for the conditioner Get all the settings for one channel of the conditioner Get or set the conditioner offset for a channel Get or set the conditioner offset range for a channel Get or set the conditioner offset for a channel Single call to set all channel parameters Get names of the signal sources available on the conditioner Get the type of signal conditioner

Debugging operations Th

These functions can be used when debugging a script.			
Debug	Set a permanent break point and disable/enable debugging		
DebugList	List internal Spike2 script objects		
DebugOpts	Gets and optionally sets system level debugging options		
Eval	Convert the argument to text and display		

Environment	These functions don't fit well into any of the other categories!	
	Арр	Get the program serial number
	Date\$	Get system date in a string in a variety of formats
	Error\$	Convert a runtime error code to a message string
	Profile	Access to the system registry and to Spike2 Preferences
	ProgKill	Terminate a process started by ProgRun
	ProgRun	Run a program, optionally position the window
	ProgStatus	Test if a program started by ProgRun is still active
	ScriptRun	Set a script to run when the current script ends
	Seconds	Get or set current relative time in seconds
	Sound	Play a tone or a .wav file
	Speak	Convert text to speech on systems that support this
	System	Get system revision as number
	System\$	Get system name as a string, access environment variables
	Time\$	Get system time in a string in a variety of formats
	TimeDate	Get system date and time as numbers

Spike shape window These functions can be used with a spike shape window.

	· ·
Cursor	Moving cursor 0 sets the time to read spikes from
FileClose	Close the spike shape window
HCursor	Set the horizontal cursor positions
Optimise	Optimise the y range and horizontal cursors
SSButton	Set and read button states in the window
SSChan	Set and get channel information, save channel configuration
SSClassify	Classify and create WaveMark data
SSOpen	Open spike shape windows and get handle of open window
SSParam	Equivalent to the spike shape parameters dialog
SSRun	Get and set the run state and speed
SSTempDelete	Delete templates based on index and code
SSTempGet	Read back template shape
SSTempInfo	Get and set template information (width, locking, counts)
SSTempSet	Create and modify templates
SSTempSizeGet	Get template size and offset in display and displayed points
SSTempSizeSet	Set template size and position and display size in points
ViewLink	Get the associated time view handle
Window	Set the position of the spike shape window
WindowGetPos	Get the window position
WindowVisible	Show and hide the spike shape window
XHigh	Get the start of the time view range to process
XLow	Get the end of the time view range to process
XRange	Set time range in associated time view to process
YHigh	Read back the y axis range
YLOW	Read back the y axis range
YRange	Set displayed y range

Spike Monitor	These function control the Spike monitor window.		
	SMControl SMOpen ViewLink Window WindowGetPos WindowVisible	Get and set the states of the control bar items Open and or get the handle of a spike monitor window Get the associated time view handle Set the position of the spike shape window Get the window position Show and hide the spike shape window	
Multimedia files	These function control a MMAudio MMImage MMOpen MMPosition	avi files associated with data fles. Get information on the audio content Get a video image in a variety of formats Open and or get the handle of a multimedia file Set the play position and play state of a multimedia file	
	MMRate MMVideo	Set the desired frames per second during recording Get information about the video format	
Serial line control	These functions let the computer. This feature of SerialOpen SerialWrite SerialRead SerialCount SerialClose	e script writer read to and write from serial line ports on their would normally be used to control equipment during data capture. Open a serial port and configure it (set Baud rate, parity etc.) Write characters to the serial port Read characters from the serial port Count the number of data items available to read Release a previously opened serial port	
Multiple monitor support	These functions give th Spike2 windows in a mu System Window WindowVisible	e script writer contol over positioning the Spike2 application and altiple monitor environment. Get monitor count, positions and identify the primary monitor. Position a script-controllable window relative to monitors. Maximise the application over the entire desktop.	


reference

Alphabetical command This section of the manual lists commands alphabetically. If you are not sure which command you require, look in the previous chapter, Commands by function. You might also find the index useful as it cross-references commands and common keywords.

Abs()

This evaluates the absolute value of an expression as a real number. This can also form the absolute value of a real or integer array.

Func Abs(x|x[]|x[]]);

A real number, or a real or integer vector or matrix Х

Returns If x is an array, this returns 0 if all was well, or a negative error code if integer overflow was detected. Otherwise it returns \times if \times is positive, otherwise -x.

See also:ATan(), Cos(), Exp(), Frac(), Ln(), Log(), Max(), Min(), Pow(), Rand(), Round(),Sin(),Sqrt(),Tan(),Trunc()

App()

This returns view handles for system specific windows and program information.

Func App({type%});

- type% -3 The program serial number.
 - -2 The highest type% that returns a value.
 - -1 The program revision multiplied by 100.

The remaining values return view handles. If type% is omitted, 0 is used.

- 4 Edit toolbar 0 Application 8 Sample control panel 1 Main toolbar 5 Play wave bar 9 Sequence control panel
- 2 Status bar 6 Script bar

3 Running script

- 7 Sample bar
- Returns The requested information or a handle for the selected window. If the requested window does not exist, the return value is 0.

See also:View(), Dup(), LogHandle(), System(), Window()

ArrAdd

This function adds a constant or an array to a real or integer array.

Func ArrAdd(dest[]{[]}, source[]{[]}|value);

- The destination array (vector or matrix). dest
- source An array of reals or integers with the same number of dimensions as dest. If the arrays have different sizes, the smaller size is used for each dimension.
- A value to be added to all elements of the destination array. value
- Returns The function returns 0 if all was well, or a negative error code for integer overflow. Overflow is detected when adding a real array to an integer array and the result is set to the nearest valid integer.

In the following examples we assume that the current view is a result view:

```
var fred[100], jim%[200], two[3][30], add[3][30];
ArrAdd(fred[],1.0);
                             'Add 1.0 to all elements of fred
ArrAdd(fred[],jim%[]);
                             'Add elements 0-99 of jim% to fred
ArrAdd([],fred[10]);
                             'Add fred[10] to result channel 1
                             'add corresponding elements
ArrAdd(two[][], add[][]);
See also:ArrConst(), ArrDiff(), ArrDiv(), ArrDivR(), ArrDot(), ArrFFT(),
       ArrFilt(), ArrIntgl(), ArrMul(), ArrSub(), ArrSum(), Len()
```

ArrConst()	This function sets an array or result view to a array or result view to another array or result	a constant value, or copies the elements of an view.	
	<pre>Func ArrConst(dest[]{[]} dest%[]{</pre>	<pre>[]}, source[]{[]} value);</pre>	
	dest The destination real or integer array	(vector or matrix).	
	source An array of reals or integers with th arrays have different sizes, the smal	e same number of dimensions as dest. If the ler size is used for each dimension.	
	value A value to be copied to all elements of the destination array.		
	 Returns The function returns 0, or a negative error code. If an integer overflow element is set to the nearest integer value to the result. The function performs the operations listed below. The indices <i>j</i> and <i>i</i> mean reperoperation for all values of the indices. Both ald and bld are one dimensional a ald and bld are two dimensional arrays. The arrays and value may be integer or reference. 		
	<pre>Function ArrConst(ald[], value); ArrConst(ald[], bld[]); ArrConst(a2d[][], value); ArrConst(a2d[][], b2d[][]);</pre>	Operation ald[<i>i</i>] := value ald[<i>i</i>] := bld[<i>i</i>] a2d[<i>j</i>][<i>i</i>] := value a2d[<i>j</i>][<i>i</i>] := b2d[<i>j</i>][<i>i</i>]	
	See also:ArrAdd(),ArrDiff(),ArrDiv() ArrIntgl(),ArrMul(),ArrSub(,ArrDivR(),ArrDot(),ArrFilt(),),ArrSubR(),ArrSum(),Len()	
ArrDiff()	This function replaces an array or result vie this as a crude form of differentiation, however	w with an array of differences. You can use ver ArrFilt() provides a better method.	
	<pre>Proc ArrDiff(dest[]);</pre>		

dest[] A real or integer vector that is replaced its differences. The first element of dest is not changed.

The effect of the ArrDiff() function can be undone by ArrIntgl(). The following block of code performs the same function as ArrDiff(work[]):

```
var work[100],i%;
for i%:=99 to 1 step -1 do work[i%] -= work[i%-1]; next;
See also:ArrAdd(),ArrConst(),ArrDiv(),ArrDivR(),ArrDot(),ArrFFT(),
```

```
ArrDiv()
```

This function divides a real or integer array by an array or a constant. Use ArrDivR() to form the reciprocal of an array. Division by zero and integer overflow are detected.

ArrFilt(),ArrIntgl(),ArrMul(),ArrSub(),ArrSubR(),ArrSum()

Func ArrDiv(dest[]{[]}, source[]{[]}|value);

- dest An array (vector or matrix) of real or integer values.
- source An array of reals or integers with the same number of dimensions as dest, used as the denominator of the division. If the arrays have different sizes, the smaller size is used for each dimension.
- value A value used as the denominator of the division.
- Returns 0 or a negative error code if integer overflow or division by zero occurs.

If there was integer overflow when assigning the result to an integer array the result is set to the nearest allowed integer value. If division by zero occurs, the associated destination element is not changed. The function performs the operations listed below. The indices j and i mean repeat the operation for all values of the indices. Both ald and bld are vectors, ald and bld are matrices. The arrays and value may be integer or real.

Function	Operation
<pre>ArrDiv(ald[], value);</pre>	ald[i] := ald[i] / value
<pre>ArrDiv(ald[], bld[]);</pre>	ald[i] := ald[i] / bld[i]
<pre>ArrDiv(a2d[][], value);</pre>	a2d[j][i] := a2d[j][i] / value
<pre>ArrDiv(a2d[][], b2d[][]);</pre>	a2d[j][i] := a2d[j][i] / b2d[j][i]
See also:ArrAdd(), ArrConst(), ArrDiff ArrFilt(), ArrIntgl(), ArrMul	(),ArrDivR(),ArrDot(),ArrFFT(), (),ArrSub(),ArrSubR(),ArrSum()

ArrDivR()

This function divides a real or integer array into an array or a constant.

Func ArrDivR(dest[]{[]}, source[]{[]}|value);

- dest An array (vector or matrix) of reals or integers used as the denominator of the division and for storage of the result.
- source A real or integer array with the same number of dimensions as dest used as the numerator of the division. If the arrays have different sizes, the smaller size is used for each dimension.
- value A value used as the numerator of the division.

Returns 0 or a negative error code if integer overflow or division by zero occurs.

If there was integer overflow when assigning the result to an integer array the result is set to the nearest allowed integer value. If division by zero occurs, the associated destination element is not changed.

The function performs the operations listed below. The indices j and i mean repeat the operation for all values of the indices. Both ald and bld are vectors, a2d and b2d are matrices. The arrays and value may be integer or real.

Function	Operation
<pre>ArrDivR(ald[], value);</pre>	ald[i] := value / ald[i]
<pre>ArrDivR(ald[], bld[]);</pre>	ald[i] := bld[i] / ald[i]
<pre>ArrDivR(a2d[][], value);</pre>	a2d[j][i] := value / a2d[j][i]
<pre>ArrDivR(a2d[][], b2d[][]);</pre>	a2d[j][i] := a2d[j][i] / b2d[j][i]
See also:ArrAdd(),ArrConst(),ArrDi	ff(),ArrDiv(),ArrDot(),ArrFilt(),
ArrIntgl(),ArrMul(),ArrSu	b(),ArrSubR(),ArrSum(),Len()

ArrDot()

This function multiplies a vector by another and returns the sum of the products (sometimes called the dot product). The vectors are not changed.

Func ArrDot(arr1[], arr2[]);

- arr1 A real or integer vector.
- arr2 A real or integer vector.
- Returns The function returns the sum of the products of the corresponding elements of the two vectors.

See also:ArrAdd(), ArrConst(), ArrDiff(), ArrDiv(), ArrDivR(), ArrFFT(), ArrIntgl(), ArrMul(), ArrSub(), ArrSubR(), ArrSum(), Len()

ArrFFT()	This command performs spectral analysis on a result view, or on an array of data. Variants of this command produce log amplitude, linear amplitude, power and relative phase as well as an option to window the original data. The command has the syntax:				
	<pre>Func ArrFFT(dest[], mode%);</pre>				
	dest A real vector to process. It should be a power of two points long, from 8 points upwards; the upper size is limited by available memory. If the number of points is not a power of two, the size is reduced to the next lower power of two points.				
	mode $%$ The mode of the command, in the range 0 to 5. Modes are defined below.				
	Returns The function returns 0 or a negative error code.				
	This command uses real arithmetic to calculate the fast Fourier transform. This is more precise, but somewhat slower than the integer transform used by SetPower().				
	Modes 1 and 3-5 take an array of data that is a set of equally spaced samples in some unit (usually time). If this unit is xin, the output is equally spaced in units of $1/xin$. In the normal case of input equally spaced in seconds, the output is equally spaced in 1/seconds, or Hz. If there are n input points, and the interval between the input points is t, the spacing between the output points is $1/(n*t)$. The transform assumes that the sampled waveform is composed of sine and cosine waves of frequencies: $0, 1/(n*t), 2/(n*t), 3/(n*t)$ up to $(n/2)/(n*t)$ or $1/(2*t)$.				
Display of phase in result views	The phase information sits rather uncomfortably in a result view. When it is drawn, the x axis has the correct increment per bin, but starts at the wrong frequency. If you need to draw it, the simplest solution is to copy the phase information to bin 1 from bin $n/2+1$ and set bins 0 and $n/2$ to 0 (this destroys any amplitude information):				
	<pre>ArrConst([1:],[n/2+1:]); 'Copy phase to the start of the view [0]:=0; [n/2]:=0; 'Set phase of the DC and Nyquist points Draw(0,n/2+1); 'Display the phase</pre>				
Mode 0: Window data	This mode is used to apply a raised cosine window to the data array. See SetPower() for an explanation of windows. The selected data is multiplied by a raised cosine window of maximum amplitude 1.0, minimum amplitude 0.0. This window causes a power loss in the result of a factor of 3/8.				
	You can supply your own window to taper the data, using the array arithmetic commands. The raised cosine is supplied as a general purpose window.				
Mode 1: Forward FFT	This mode replaces the data with its forward Fast Fourier Transform. You would most likely use this to allow you to remove frequency components, then perform the inverse transform. The output of this mode is in two parts, representing the real and imaginary result of the transform (or the cosine and sine components). The first $n/2+1$ points of the result hold the amplitudes of the cosine components of the result. The remaining $n/2-1$ points hold the amplitudes of the sine components. In the case of an 8 point transform, the output has the format:				
	$ \begin{array}{cccccccccccccccccccccccccccccccccccc$				
	There is no sine amplitude at a frequency of $4/(n*t)$, the Nyquist frequency, as this sine wave would have amplitude 0 at all sampled points.				

Mode 2: Inverse FFT	This mode takes data in the format produced by the forward transform and converts in back into a time series. In theory, the result of mode 1 followed by mode 2, or mode 2 followed by mode 1, would be the original data. However, each transform adds some noise due to rounding effects in the arithmetic, so the transforms do not invert exactly One use of modes 1 and 2 is to filter data. For example, to remove high frequency noise use mode 1, set unwanted frequency bins to 0, and use mode 2 to reconstruct the data.			
Mode 3: dB and phase	This mode produces an output with the fir power spectrum in dB, and the second n/2 data. In the case of our 8 point transform t	st n/2+1 po 2-1 points l he output fo	ints holding holding the p ormat would	the log amplitude of the hase (in radians) of the be:
	$\begin{array}{ccc} \text{point} & \text{frequency} & \text{contents} \\ 0 & \text{DC} & \log \text{ amplitude in dB} \\ 1 & 1/(n^*t) & \log \text{ amplitude} \\ 2 & 2/(n^*t) & \log \text{ amplitude} \\ 3 & 3/(n^*t) & \log \text{ amplitude} \end{array}$	point 4 5 6 7	frequency 4/(n*t) 1/(n*t) 2/(n*t) 3/(n*t)	contents log amplitude in dB phase in radians phase in radians phase in radians
	There is no phase information for DC or for for both of these points is zero. If you w $(180^{\circ}/\pi)$. The log amplitude is calculated mode 1 above) and forming:	or the point ant the pha by taking t	at 4/(n*t). T ase in degree he result of a	his is because the phase es, multiply by 57.3968 a forward FFT (same as
	dB = 10.0 Log(power)	T	he <i>power</i> is c	alculated as for Mode 5
Mode 4: Amplitude and phase	This mode produces the same output form amplitude. The amplitude is calculated by mode 1 above), and forming:	at as mode y taking th	3, but with a e result of a	mplitude in place of log forward FFT (same as
	$amplitude = (\cos^2 + \sin^2)^{0.5} $ Th	ere is no si	n component	for the DC and Nyquist
Mode 5: Power and phase	$amplitude = (\cos^2 + \sin^2)^{0.5}$ TheThis mode produces the same output for power. The sum of the power component original data divided by the number of data result of a forward FFT (same as mode 1 at	ere is no sir rmat as mo nts is equa ta points. T above), and	des 2 and 3 l to the sum he power is a forming:	for the DC and Nyquist , but with the result as n of the squares of the calculated by taking the
Mode 5: Power and phase	<i>amplitude</i> = $(\cos^2 + \sin^2)^{0.5}$ The This mode produces the same output for power. The sum of the power component original data divided by the number of dat result of a forward FFT (same as mode 1 at <i>power</i> = $(\cos^2 + \sin^2) * 0.5$ <i>power</i> = DC ² or Nyquist ²	ere is no sin rmat as mo nts is equa ta points. T above), and for all cor	des 2 and 3 l to the sum he power is 6 forming: mponents exc for the DC a	for the DC and Nyquist , but with the result as n of the squares of the calculated by taking the rept the DC and Nyquist nd Nyquist components
Mode 5: Power and phase	$amplitude = (\cos^2 + \sin^2)^{0.5}$ TheThis mode produces the same output for power. The sum of the power component original data divided by the number of data result of a forward FFT (same as mode 1 at $power = (\cos^2 + \sin^2) * 0.5$ $power = DC^2$ or Nyquist ² You can compare the output of this mode waveform channel on channel 1 in view 1,	ere is no sin rmat as mo nts is equa ta points. T above), and for all cor e with the re , and do the	des 2 and 3 l to the sum he power is 6 forming: mponents exc for the DC a esult of SetH following:	for the DC and Nyquist , but with the result as n of the squares of the calculated by taking the ept the DC and Nyquist nd Nyquist components Power (). If you have a
Mode 5: Power and phase	<pre>amplitude = (cos²+sin²)^{0.5} Th This mode produces the same output for power. The sum of the power component original data divided by the number of dat result of a forward FFT (same as mode 1 at power = (cos²+sin²) * 0.5 power = DC² or Nyquist² You can compare the output of this mode waveform channel on channel 1 in view 1, var spView%, afView%; spView% := SetPower(1,1024); Process(0,1024*View(-1).Binsize WindowVisible(1); afView% := SetAverage(1,1024,0), Process(0,0);WindowVisible(1); View(afView%); ArrFFT([], 0);ArrFFT([], 5); Draw(0,500);Optimise(0);</pre>	ere is no sin rmat as monts is equal ta points. The above), and for all correct with the rected assume 'assume 'select (1)); 'prited 'make with 'make with 'make with 'make with 'bow the 'Apply with 'Show 50 'Look att 'Show sat	des 2 and 3 l to the sum he power is a forming: mponents exc for the DC a esult of SetH following: in a time power spe cocess fir ndow visi y 1024 poi he data the view yindow, ta 00 bins of c SetPower ame bins o	for the DC and Nyquist , but with the result as n of the squares of the calculated by taking the eept the DC and Nyquist nd Nyquist components eower(). If you have a view ctrum st 1024 points ble nts holding data ke power spectrum Power () result f power spectrum
Mode 5: Power and phase	<pre>amplitude = (cos²+sin²)^{0.5} Th This mode produces the same output for power. The sum of the power component original data divided by the number of dat result of a forward FFT (same as mode 1 and power = (cos²+sin²) * 0.5 power = DC² or Nyquist² You can compare the output of this mode waveform channel on channel 1 in view 1, var spView%, afView%; spView% := SetPower(1,1024); Process(0,1024*View(-1).Binsize WindowVisible(1); afView% := SetAverage(1,1024,0), Process(0,0);WindowVisible(1); View(afView%); ArrFFT([], 0);ArrFFT([], 5); Draw(0,500);Optimise(0); View(spView%); Draw(0,500);Optimise(0); The results are identical except that the and view generated by SetPower(). The read command compensates for the effect of th result by 8/3. To produce the same numering proceeding of the same numering of the same numering of the same numering set by 8/3. To produce the same numering of the same numering proceeding of the same numering of the same numerin</pre>	ere is no sin rmat as monts is equal ta points. The above), and for all correct with the rection assume 'select (1)); 'pri'make wid ;'To copy 'show the 'Move to 'Apply widdle 'Show 50 'Look at assume 'Show select 'Apply widdle assume 'show the 'show the 'show the 'show select 'assume 'show the 'show the 'show select 'show the 'show select 'show the 'show select 'show the 'show select 'show the 'show select 'show	des 2 and 3 l to the sum he power is a forming: mponents exc for the DC a esult of SetH following: in a time power spe cocess fir ndow visi 7 1024 poi he data the view window, ta 00 bins of c SetPower ame bins o view is 3/8 a difference i it uses interr altiply by 8/3	for the DC and Nyquist , but with the result as n of the squares of the calculated by taking the eept the DC and Nyquist nd Nyquist components Power (). If you have a view ctrum st 1024 points ble nts holding data ke power spectrum Power () result f power spectrum of the amplitude of the s that the SetPower () nally by multiplying the

ArrFilt()

This function applies a FIR (Finite Impulse Response) filter to a real array. You can use FiltCalc() and FIRMake() to generate filter coefficients for a wide range of filters.

Func ArrFilt(dest[], coef[]);

- dest[] A real vector holding the data to filter. It is replaced by the filtered data.
- coef[] A real vector of filter coefficients. This is usually an odd number of data points long so that the result is not phase shifted.

Returns The function returns 0 if there was no error, or a negative error code.



This diagram shows how the FIR filter works. Open circles represent filter coefficients, solid circles are the input and output waveforms. Each output point is generated by multiplying the waveform by the coefficients and summing the result. The coefficients are then moved one step to the right and the process repeats.

From this description, you can see that the filter coefficients (from right to left) are the *impulse response* of the filter. The impulse response is the output of a filter when the input signal is all zero except for one sample of unit amplitude. In the example above with seven coefficients, there is no time shift in the output. If the filter has an even number of coefficients, there is an output time shift of half a sample.

The filter operation is applied to every vector element. There is a problem at the start and end of the vector where some coefficients have no corresponding data element.



The simple solution is to take these missing points as copies of the first and last points. This is usually better that taking these points as 0. You should remember that the first and last (nc+1)/2 points are unreliable, where *nc* is the number of coefficients.

A simple use of this command is to produce three point smoothing of data, replacing each point by the mean of itself and the two points on either side:

```
var data[1000],coef[3]; 'Arrays of data and the coefficients
... 'Fill data[] with values
coef[0]:=0.33333; coef[1]:=0.33333; coef[2]:=0.33333;
ArrFilt(data[],coef[]); 'smooth the data.
```

If you use this command to apply a causal filter, that is, one with all coefficients that use data points ahead of the current point set to zero, you must still provide these coefficients. If you omit the trailing zero coefficients, the output will be time shifted backwards by half the number of coefficients you do supply.

```
See also:ArrDiff(), ArrDiv(), ArrDot(), ArrFFT(), ArrIntgl(), ArrMul(),
ArrSub(), ArrSubR(), ArrSum(), FiltCalc(), FIRMake(), Len()
```

ArrIntgl()

This function is the inverse of ArrDiff(), replacing each point by the sum of the points from the start of the array up to the element. The first element is unchanged.

Proc ArrIntgl(dest[]);

dest A vector of real or integer data.

The function is equivalent to the following:

for i%:=1 to Len(dest[])-1 do dest[i%] += dest[i%-1]; next;

See also:ArrAdd(), ArrConst(), ArrDiff(), ArrDiv(), ArrDivR(), ArrDot(), ArrMul(), ArrSub(), ArrSubR(), ArrSum(), Len()

ArrMul()

This command is used to form the product of a pair of arrays, or to scale an array by a constant. A less obvious use is to negate an array by multiplying by -1.

Func ArrMul(dest[]{[]}, source[]{[]}|value);

- dest A vector or matrix of real or integer numbers. If dest is integer, the multiplication is done as reals and truncated to integer.
- source A real or integer array with the same number of dimensions as dest. If the arrays have different sizes, the smaller size is used for each dimension.
- value A value to multiply the data in dest.

Returns The function returns 0 if all was well, or a negative error code.

The function performs the operations listed below. The indices j and i mean repeat the operation for all values of the indices. Both ald and bld are vectors, a2d and b2d are matrices. The arrays and value may be integer or real.

Operation

Function

<pre>ArrMul(a1d[], value);</pre>	ald[i] := ald[i] * value
<pre>ArrMul(a1d[], b1d[]);</pre>	ald[i] := ald[i] * bld[i]
<pre>ArrMul(a2d[][], value);</pre>	a2d[j][i] := a2d[j][i] * value
ArrMul(a2d[][], b2d[][]);	a2d[j][i] := a2d[j][i] * b2d[j][i]
See also:ArrAdd(), ArrConst(), ArrDiff((),ArrDiv(),ArrDivR(),ArrDot(),
<pre>ArrIntgl(),ArrSub(),ArrSubR(</pre>	(),ArrSum(),Len(),MATMul()

ArrSort()

This function sorts an array of any data type and optionally orders additional arrays to match the sorted array.

Proc ArrSort(sort[]{, opt%{, arr1[]{, arr2[]{, ...}}});

- sort[] An array of integer, real or string data to sort.
- opt% This optional argument holds the sorting options. If omitted, the value 0 is used. It is the sum of the following flag values:
 - 1 Sort in descending order. If omitted, the data is sorted in ascending order.
 - 2 Case sensitive sort when sort[] is an array of strings. String sorts are usually case insensitive. If omitted, the sort is case insensitive.
- arrn[] If present, these arrays are sorted in the same order as the sort[] array. The arrays can be of any type. You can sort up to 18 additional arrays.

ArrSpline()	This fund array sar has contr from poi are set to source po	ction interpolates an array of real or integer data sampled at one rate into another npled at a different rate using cubic splines. This assumes that the source data nuous first and second derivatives and that the second derivatives vary linearly nt to point. The second derivatives at the first and last point of the source data to zero. We interpolate from up to one input sampling interval before the first point to up to one sampling interval after the last source point.
	Proc Ai	<pre>crSpline(dest[], source[]{, ratio{, start}});</pre>
	dest	A real or integer vector that holds the interpolated result.
	source	A vector that is integer if dest is integer or real if dest is real.
	ratio	If present, this is the ratio of the sampling interval of dest to the sampling interval of source (or the source frequency divided by the dest frequency). If omitted, we assume that the first and last points of the two arrays are aligned. There are no restrictions on ratio, it can even be negative (in which case the output is backwards relative to the input).
	start	If present, this is the position of the first point of dest relative to the first point of source in units of the sample interval of the source array. If omitted, it is assumed to be 0 (the first points of both arrays fall at the same position).
	You will time rang falls by a supply m	get the best results if you can supply source data before and after the output ge. The effect of a source point on the interpolation of an interval n points away factor of approximately 4 each time n increases by 1. There is rarely the need to hore than 15 data points before and after the interpolation range.
An example	Suppose sampled at 5.335 0.1. The 0.01, the	we have a source vector of 100 points sampled at 100 Hz, with the first point at 5.02 seconds, and we want to generate an array sampled at 1000 Hz that starts and lasts 0.5 seconds. In this case, the value of ratio is 0.001/0.01, which is value of start is 31.5, which is the time difference (5.335 - 5.02) divided by sample interval of the source channel.
	See also:	ChanProcessAdd(),DrawMode()
ArrSub()	This fun- the destin	ction forms the difference of two arrays or subtracts a constant from an array. If nation is an integer array, overflow is detected when subtracting real values.
	Func Ai	<pre>crSub(dest[]{[]}, source[]{[]} value);</pre>
	dest	A real or integer vector or matrix that holds the result.
	source	A real or integer array with the same number of dimensions as dest. If the arrays have different sizes, the smaller size is used for each dimension.
	value	A real or integer value.
	Returns	0 if all is well or a negative error code if integer overflow is detected.
	The function operation matrices	tion performs the operations listed below. The indices j and i mean repeat the for all values of the indices. Both ald and bld are vectors; ald and bld are vectors; ald and bld are vectors; all value may be integer or real.
	Functior	Operation

FUNCTION	Operation
<pre>ArrSub(ald[], value);</pre>	ald[i] := ald[i] - value
<pre>ArrSub(ald[], bld[]);</pre>	ald[i] := ald[i] - bld[i]
<pre>ArrSub(a2d[][], value);</pre>	a2d[j][i] := a2d[j][i] - value
<pre>ArrSub(a2d[][], b2d[][]);</pre>	a2d[j][i] := a2d[j][i] - b2d[j][i]
See also:ArrAdd(), ArrConst(), ArrDiff	(),ArrDiv(),ArrDivR(),ArrDot(),
<pre>ArrIntgl(),ArrMul(),ArrSubR</pre>	(),ArrSum(),Len()

ArrSubR()

This function forms the difference of two arrays or subtracts an array from a constant. If the destination is an integer array, overflow is detected when subtracting real values.

Func ArrSubR(dest[]{[]}, source[]{[]}|value);

- dest A real or integer vector or matrix.
- source A real or integer array with the same number of dimensions as dest. If the arrays have different sizes, the smaller size is used for each dimension.
- value A real or integer value.

Returns 0 if all is well or a negative error code if integer overflow is detected.

The function performs the operations listed below. The indices j and i mean repeat the operation for all values of the indices. Both ald and bld are vectors, a2d and b2d are matrices. The arrays and value may be integer or real.

Function	Operation
<pre>ArrSubR(ald[], value);</pre>	ald[i] := value - ald[i]
<pre>ArrSubR(ald[], bld[]);</pre>	ald[i] := bld[i] - ald[i]
<pre>ArrSubR(a2d[][], value);</pre>	a2d[j][i] := value - a2d[j][i]
<pre>ArrSubR(a2d[][], b2d[][]);</pre>	a2d[j][i] := b2d[j][i] - a2d[j][i]

See also:ArrAdd(), ArrConst(), ArrDiff(), ArrDiv(), ArrDivR(), ArrDot(), ArrIntgl(), ArrMul(), ArrSub(), ArrSum(), Len()

ArrSum()

This function forms the sum of the values in a vector or matrix, and optionally forms the mean and standard deviation of the vector or of each matrix column. The standard deviation of m data points with a sum of squared difference from the mean of errSq is sqrt(errSq/(m-1)). There are two command variants:

```
Func ArrSum(arr[]|arr%[]{, &mean{, &stDev}});
Func ArrSum(arr[][]|arr%[][]{, mean[], stDev[]});
```

- arr A real or integer vector or matrix to process.
- mean If present it is returned holding the mean of the values in the array. The mean is the sum of the values divided by the number of array elements. If arr is an m by n matrix, mean must be a vector of at least n elements and is returned holding the mean of each column of arr.
- stDev If present, this returns the standard deviation of the array elements around the mean. If the array has only one element the result is 0. If arr is a m by n matrix, stDev must be a vector with at least n elements and is returned holding the standard deviation of each column of arr.

Returns The function returns the sum of all the array elements as a real number.

See also:ArrAdd(), ArrConst(), ArrDiff(), ArrDiv(), ArrDivR(), ArrDot(), ArrIntgl(), ArrMul(), ArrSub(), ArrSubR(), Len()

Asc()

This function returns the ASCII code of the first character in the string as an integer.

Func Asc(text\$);

text\$ The string to process.

```
See also:Chr$(), DelStr$(), LCase$(), Left$(), Len(), Mid$(), Print$(),
Right$(),Str$(),UCase$(),Val()
```

ATan()

This function returns the arc tangent of an expression, or the arc tangent of an array:

```
Func ATan(s|s[]|s[][] {,c});
```

- s If the only argument, the function uses this for the arc tangent calculation. s can also be a real array (in which case c must not be present).
- c If this is present, the function uses s/c for the calculation.
- Returns If s is an array, each element of s is replaced by its arc tangent in the range $-\pi/2$ to $\pi/2$ radians. The function returns 0 if all was well or a negative error code.

When s is not an array, if s is the only argument, the function returns the arc tangent of s in the range $-\pi/2$ to $\pi/2$. If c is present, the function calculates the result of ATan (s/c) and uses the signs of s and c to decide the quadrant of the result. With the second argument, the result is in the range $-\pi$ to π .

BinError()

This function is used in a result view with error bins enabled to access the error information. Error bins are enabled for a result view created with SetAverage() or with SetResult() with 4 added to flags%. If you are setting the error information you must set the sweep count with Sweeps() first as the sweep count is used to convert the standard deviation into the internal storage format. There are two command variants; the first transfers data for a single bin, the second for an array of bins:

Func BinError(chan%, bin% {,newSD});
Func BinError(chan%, bin%, sd[]{, set%});

- chan% The channel number in the result view.
- bin% The first bin number for which to get or set the error information.
- newSD If present, this sets the standard deviation for a single bin.
- sd[] An array used to hold standard deviation values. Values are transferred starting at bin bin% in the result view. If the array is too long, extra bins are ignored.
- set% If present and non-zero, the values in sd[] are copied to the result view. If omitted or zero, values are copied from the result view into sd[].
- Returns The first command variant returns the standard deviation at the time of the call. The second variant returns the number of bins copied. If there are 0 or 1 sweeps of data or errors are not enabled, the result is 0.

To illustrate how errors are calculated, we will assume that we are dealing with an average that is set to display the mean of the data in each bin. In terms of the script language, if the array s[] holds the contribution of each sweep to a particular bin, the mean, standard deviation and standard error of the mean are calculated as follows:

```
var mean, sd:=0, i%, diff, sem;
for i%:= 0 to Sweeps()-1 do mean += s[i%] next; 'form sum
mean /= Sweeps(); 'form mean data value
for i%:= 0 to Sweeps()-1 do
    diff := s[i%]-mean; 'difference from mean
    sd += diff*diff; 'sum squares of differences
    next;
sd := Sqrt(sd/(Sweeps()-1)); 'the standard deviation
sem := sd/Sqrt(Sweeps()); 'the standard error of the mean
```

We divide by Sweeps()-1 to form the standard deviation because we have lost one degree of freedom due to calculating the mean from the data.

See also:BinSize(), BinToX(), SetAverage(), SetResult(), Sweeps()

F	Rin	siz	e()

In a result view, this returns the x axis increment per bin. In a time view, the value returned depends on the channel type.

```
Func Binsize({chan%});
```

- chan% In a time view this is the channel from which to return information. If you omit the argument, the function returns the file time resolution in seconds. In a result view, chan% is ignored, and should be omitted.
- Returns In a time view, the sampling interval between points is returned for Waveform, WaveMark and RealWave channels or a negative number if the channel does not exist. Otherwise the underlying time resolution of the file in seconds is returned.

See also:BinToX(), XToBin()

BinToX()

This converts bin numbers to x axis units in the current result view. If the current view is a time view, it converts the underlying Spike2 time units into time in seconds.

Func BinToX(bin);

bin A bin number in the result view. You can give a non-integer bin number without error. If you give a bin number outside the result view, the bin number is limited to the range of the result view before it is converted to an x axis value.

The x axis range of a result view is BinToX(0) to BinToX(MaxTime()). Do not confuse this range with XLow() to XHigh(), which is the visible range of the x axis in the current view.

In a time view, this is in the underlying time units. If the value is beyond the x axis range, it is limited to the x axis range. The value need not be integral, but you should note that all data items in the time view have time stamps that correspond with integral values of bin. The returned value is in seconds.

Returns It returns the equivalent x axis position.

See also:BinSize(), MaxTime(), XHigh(), XLow(), XToBin()

BRead()

This reads data into variables and arrays from a binary file opened by FileOpen(). The function reads 32-bit integers, 64-bit IEEE real numbers and zero-terminated strings.

Func BRead(&arg1|arg1[]|&arg1%|arg1%[]|&arg1\$|arg1\$[] {,...});

- Arguments may be of any type. Spike2 reads a block of memory equal in size to the combined size of the arguments and copies it into the arguments. Strings or string arrays are read a byte at a time until a zero byte is read.
- Returns It returns the number of data items for which complete data was read. This will be less than the number of items in the list if end of file was reached. If an error occurred during the read, a negative code is returned.

See also:FileOpen(), BReadSize(), BRWEndian(), BSeek(), BWrite()

BReadSize()	This co The fun and 64- fixed-si current	nverts da iction rea bit IEEE ze region file posit	ta into varia ds 8, 16 and real number is in the file ion. The cur	bles and arrays from a binary file opened by FileOpen(). 32-bit integers and converts them to 32-bit integers, and 32 is and converts them to 64-bit reals. It also reads strings from (zero bytes are ignored during the read). The read is from the rent position after the read is the byte after the last byte read.
	Func E	ReadSiz	ze(size%,	&arg arg[] &arg% arg%[] &arg\$ arg\$[]{,});
	size%	The byt	tes to read fo	r each argument. Legal values depend on the argument type:
		Integer	1, 2 or 4 -1, -2	Read 1, 2 or 4 bytes and sign extend to 32-bit integer. Read 1 or 2 bytes and zero extend to 32-bit integer.
		Real	4 8	Read 4 bytes as 32-bit real, convert to 64-bit real. Read 8 bytes as 64-bit real.
		String	n	Read n bytes into a string. Null characters end the string.
	arg	The targ	get variable(s) to be filled with data. size% applies to all targets.
	Returns	It return be less occurre	ns the number than the num d during the	er of data items for which complete data was read. This will aber of items in the list if end of file was reached. If an error read, a negative code is returned.
	See also):FileOp	oen(),BRea	d(),BRWEndian(),BSeek(),BWrite()
BRWEndian()	This ge with Bl use little the Mac are little	ts and se Read(), e-endian cintosh, u e-endian	ets the "endi BReadSize data (least s ise big-endia by default.	anism" of binary data files. This affects numeric data used (), EWrite() and EWriteSize(). PC programs normally significant byte at lowest address). Some systems, including in data (most significant byte at lowest address). Binary files
	Most us	sers do no	ot need to us dian comput	se this routine. Only use it if you are writing binary files for ter or reading binary files that were generated with a big-

endian system.

Func BRWEndian({new%});

Omit or set -1 for no change. Set 0 for little-endian and 1 for big-endian. new%

Returns The current endianism as 0 for little, 1 for big or a negative error code.

See also:FileOpen(), BRead(), BReadSize(), BSeek(), BWrite(), BWriteSize()

BSeek(

This function moves and reports the current position in a file opened by FileOpen() with a type% code of 9. The next binary read or write operation to the file starts from the position returned by this function.

Func BSeek({pos% {, rel%}});

- The new file position. Positions are measured in terms of the byte offset in the pos% file from the start, the current position, or from the end. If a new position is not given, the position is not changed and the function returns the current position.
- rel% This determines to what the new position is relative:
 - 0 Relative to the start of the file (same as omitting the argument).
 - 1 Relative to the current position in the file.
 - 2 Relative to the end of the file.

Returns The new file position relative to the start of the file or a negative error code.

See also:FileOpen(), BReadSize(), BRead()

BurstMake()

This function extracts burst start and end times from an event channel and writes them to a memory channel. Three parameters control the formation of bursts: the interval between the first two events in a burst, the interval between the last event in a burst and any following event, and the minimum number of events in a burst.



The diagram illustrates the method of forming bursts. In the first case, with a short maximum interval between events, the algorithm finds three bursts. In the second case with a longer period, the algorithm detects one burst. The command for this function is:

Func BurstMake(mChan%,eChan%,sTime,eTime,maxT{,firstT{,minE%}});

- mChan% The channel number of an event or marker channel created by MemChan() for the output. If this is a marker channel, the start of each burst is given a first marker code of 00 and the end has a first marker code of 01.
- eChan% The channel number of an event channel to search for bursts.
- sTime The start of the time range to search for bursts.
- eTime The end of the time range to search for bursts.
- maxT The maximum time between two events (after the first pair) for the events to lie in the same burst.
- firstT The maximum time between the first two events in a burst. If omitted, maxT sets the time interval for the first pair of events.
- minE% The minimum number of events that can make up a burst. The default is 2.
- Returns The function returns the number of bursts found, or a negative error code.

See also:BurstRevise(),BurstStats(),MemChan()

BurstRevise()

This command modifies a list of times, indicating bursts and produces a new list of bursts that have inter-burst intervals and burst durations greater than specified minimum times (See BurstStats() for more information).

Func BurstRevise(mChan%, eChan%, sTime, eTime, minI, minD);

- mChan% A memory channel created by MemChan() holding event data to which the output of this command is added. This can be the same channel as eChan%, but if it is, the output is generated by deleting unwanted events from the channel.
- eChan% The event or marker channel to read burst times from.
- sTime The start of the time range to revise.
- eTime The end of the time range to revise.
- minI The minimum interval between the end of one burst and the start of the next. Bursts with shorter intervals are amalgamated.
- minD The minimum duration of a burst. Shorter bursts are deleted.

Returns The number of bursts in the output or a negative error code.

See also:BurstMake(),BurstStats(),MemChan()

BurstStats()

This command returns statistics on bursts (possibly made by BurstMake()). Additional rules can be applied to the bursts before the statistics are calculated to amalgamate bursts that are too close together and to delete bursts that are too short. The statistics are the mean and standard deviation of the duration of the bursts and the intervals between the start of one burst and the start of the next.



- eChan% The event channel containing burst data. The first event found in the time range is assumed to be the start of a burst, the second the end, and so on.
- sTime The start of the time range to process.
- eTime The end of the time range to process.
- mean I This is a real variable that is returned holding the mean interval between the starts of bursts (as long as at least 2 bursts were found).
- sdI This is a real variable that is returned holding the standard deviation of the mean interval (as long as at least 3 bursts were found).
- meanD This is a real variable that is returned holding the mean burst duration (as long as at least 1 burst was found).
- sdD This is a real variable that is returned holding the standard deviation of the burst durations (as long as at least 2 bursts were found).
- minI This optional value sets the minimum interval between the end of one burst and the start of the next. It takes the value 0.0 if omitted. Bursts that are closer than this are amalgamated for the purpose of forming statistics.
- minD This optional value (taken as 0.0 if omitted) sets the minimum burst duration. Bursts shorter than this (after amalgamations) are ignored in the statistics.

Returns It returns the number of bursts used for the statistics, or a negative error code.

This command is used in scripts that optimise the values of minI and minD for forming bursts. For example, when bursting is known to follow a cyclical pattern, one would like to find values of minI and minD that minimise the total coefficient of variance:

Total coefficient of variance = $(sdD/meanD)^2 + (sdI/meanI)^2$

Once suitable values are found, the BurstRevise() command can generate a memory channel holding bursts based on the optimised parameters.

See also:BurstMake(), BurstRevise(), MemChan()

BWrite()

This function writes binary data values and arrays into a file opened by FileOpen() with a type% code of 9. The function can write 32-bit integers, 64-bit IEEE real numbers and strings. The output is at the current position in the file. The current position after the write is the byte after the last byte written.

Func BWrite(arg1 {,arg2 {,...}});

Arguments may be of any type, including arrays. Spike2 fills a block of memory equal in size to the combined size of the arguments with the data held in the arguments and copies it to the file. An integer uses 4 bytes and a real uses 8 bytes. A string is written as the bytes in the string plus an extra zero byte to mark the end. Use BWriteSize() to write a fixed number of bytes. Returns It returns the number of arguments for which complete data was written. If an error occurred during the write, a negative code is returned.

See also:FileOpen(), BRead(), BReadSize(), BRWEndian(), BWriteSize()

BWriteSize()

This function writes variables or arrays as binary into a file opened by FileOpen() with a type% code of 9. It writes 8, 16 and 32-bit integers and 32 and 64-bit reals and strings. It allows you to write formats other than the 32-bit integer and 64-bit real used internally by Spike2 and to write variable-length strings into fixed-size fields in a binary file.

Func BWriteSize(size%, arg1 {,arg2 {,...}});

size% Bytes to write for each argument (or array element if the argument is an array). Legal values of size% depend on the argument type:

1, 2	Write least significant 1 or 2 bytes.
4	Write all 4 bytes of the integer.
4	Convert to 32-bit real and write 4 bytes.
8	Write 8 bytes as 64-bit real.
n	Write n bytes. Pad with zeros if the string is too short.
	1, 2 4 4 8 n

arg The target variable(s) to be filled with data. size% applies to all targets.

Returns It returns the number of data items for which complete data was written or a negative error code.

See also:FileOpen(), BRead(), BReadSize(), BRWEndian(), BWrite()

Chan\$()

This function converts a channel number or a list of channel numbers into a string. Memory, duplicate and virtual channels are listed as they appear on screen (not just as numbers). If a channel does not exist in the current view, it is represented as a number.

Func Chan\$(chan%|chan%[]);

chan% Either a channel number or an array of integers in the same format as a channel specification (the first element holds the number of items, followed by the channel numbers).

Returns A channel specification string, for example "1, 3, 5..8, m2, v1a".

See also:ChanList()

This fun and offs user-de	nction is equivalent to the Analysis menu Calibrate command. It changes the scale set of Waveform and WaveMark channels and rewrites RealWave channels so that fined data sections have user-defined values. The command is:
Func C	ChanCalibrate(cSpc,mode%,cfg%,t1,t2,units\$,v1{,v2{,t3,t4}});
cSpc	A channel specifier for the channels to calibrate or -1 for all, -2 for visible and -3 for selected channels. You cannot calibrate processed channels if the process changes the scale or offset or any processed RealWave channel.
mode%	 The calibration mode. The items in brackets are the required optional arguments. Mean levels of two time ranges (v2, t3, t4). Values at two time points (v2). Set offset from mean of time range. Set scale from mean of time range. Square wave, upper and lower level (v2). Servere wave, upper and lower level (v2).
	This fur and offs user-de: Func C cSpc mode%

- 6 Peak to peak amplitude and mean (v2).
- 7 RMS amplitude about mean (v2).
- 8 Area under curve, assume zero at end.
- 9 Areas under curve, two time ranges (v2, t3, t4).
- 10 Slope of best-fit line to the data, the offset not changed.
- cfg% If non-zero, the new calibration is saved in the sampling configuration.
- t1, t2 In all modes except mode 1, these are the start and end of the first time range. In mode 1 these times correspond to the two calibration values v1 and v2.
- units\$ The units to apply to the channel.
- v1, v2 The values in user units that correspond to the times or time ranges. v1 is always used; v2 is only used in modes 0, 1, 4, 6, 7 and 9.
- t3, t4 A second time range. These values are used in modes 0 and 9.

Returns The return value is the sum of the following values:

- 1 A channel in the list did not exist or was the wrong type.
- 2 A channel was processed OK.
- 4 Unknown or unimplemented calibration mode.
- 8 A time range had the end past the start or two time ranges overlapped or the two times in mode 1 were the same.
- 16 The v1 and/or v2 values were too big, too small or too similar.
- 32 Not enough data to process at least one channel in the list.
- 64 The data is unsuitable. For example, in mode 0 mean levels must differ by at least the standard deviation of the data around the mean.

See also:ChanOffset(), ChanScale(), ChanUnits\$()

ChanColour()

This returns and optionally sets the colour of a channel in a time or result view. This colour overrides the application colour set for the drawing mode of the channel.

Func ChanColour(chan%, item%{, col%});

- chan% A channel in the time or result view.
- item% The colour item to get or set; 0=background, 1=primary, 2=secondary colour.
- col% If present, the new colour index for the item. There are 40 colours in the palette, indexed 0 to 39. Use -1 to revert to the application colour for the drawing mode.
- Returns The palette colour index at the time of the call, -1 if no colour is set or a negative error code if the channel does not exist.

See also:Colour(), PaletteGet(), PaletteSet(), ViewColour(), XYColour()

ChanComment\$()

This returns or sets the comment string for a channel in a time view. It returns an empty string in a result view. The comment can be up to 71 characters long. It is an error to use this in any other view type.

Func ChanComment\$(chan%{, new\$});

- chan% A channel in the time view.
- new\$ An optional string with a new comment. If the string is too long, it is truncated.
- Returns It returns the comment string for the designated channel. If the channel does not exist, the function does nothing and returns an empty string. Setting a comment for a result view does nothing, but is not an error.

See also:ChanTitle\$(), FileComment\$()

ChanData()

Fills an array with waveform data from a waveform or RealWave channel or with event

times from all other channel types. Use NextTime() and LastTime() to get data

	attached	to WaveMark, RealMark and TextM	lark channels and marker codes.
	Func C	hanData(chan%, arr[] arr%[],	sTime, eTime {,&fTime});
	chan%	The data channel in the current time	view to read data from.
	arr	A real or integer array. Real arrays of times in seconds. Integer arrays coll- in the underlying time units (as return	collect waveform data in user units and event ect waveforms in ADC units and event times med by Binsize()).
	sTime	The first data value returned is at or	after this time in seconds.
	eTime	The last data point returned is before	e or at this time in seconds.
	fTime	This optional argument is a variable	that is set to the time of the first data point.
	Returns	The number of data values placed is is returned; gaps terminate a read.	n the array. Only contiguous waveform data
	See also	:Binsize(),ChanValue(),ChanW	riteWave(),NextTime(),ChanScale()
ChanDelete()	This de time vie memory duplicat	letes a channel from a time, result of ew channel deletion if the channel is v buffer channels are not confirme e channels. In an XY view you can	r XY view. You can make the user confirm stored in the file. Duplicated channels, and d. In a result view you may only delete not delete the last XY channel as XY views
	must alv an XY channel	ways have at least one channel. Chan view, so if you delete a channel, the s will change.	nnels are always numbered consecutively in e channel numbers of any higher numbered
	Func C	hanDelete(cSpc {,query%});	
	cSpc	A channel specifier for the channels selected. In an XY view only channel	to delete or -1 for all, -2 for visible or -3 for el numbers greater than 0 are allowed.
	query%	If present and non-zero, the user is a part of a time view. You cannot dele	sked to confirm the deletion if the channel is the channels that are being sampled.
	Returns	0 if the channel was deleted or a n operation or tried to delete the last X	legative error code if the user cancelled the CY channel or for other problems.
	See also	ChanDuplicate(),MemChan(),X	YDelete(),XYSetChan()
ChanDuplicate()	This du view to marker	plicates a time or result view channe delete all duplicates of the current c filter setting for each template. Use D	I. It can also be used in the Edit WaveMark hannel and generate duplicates with suitable upChan() to find duplicates of a channel.
	Func C	hanDuplicate({chan%});	
	chan%	For time or result views, this is cha exist. This argument must be omitted	nnel number to duplicate; this channel must d for Edit WaveMark views.
	Returns	For time or result views it returns the Edit WaveMark view, it returns the error stops the script.	he channel number of the duplicate. For an number of duplicates that were created. Any
	For time visible. followir	e and result views the new channel is In an Edit WaveMark view, chang example duplicates a time or result	not displayed. Use ChanShow() to make it nnels are made visible immediately. The view channel and makes it visible:
	var ch ChanSh	%;ch% := ChanDuplicate(1); ow(ch%);	'create a duplicate 'make visible

See also:ChanShow(), ChanDelete(), DupChan(), MemChan()

	This f	notion together with charteleness () and charteleness () incomparison the
ChanFit()	function initialise current	ality of the Analysis menu Fit Data dialog. ChanFitShow() incorporates the ality of the Analysis menu Fit Data dialog. ChanFit() has three variants that: e ready for a new fit, perform the fit and return information about the last fit. The window must be a time, result or XY view to use these functions.
Initialise fit information	This co- limits a existing	mmand associates a fit with a channel. The fit parameters and the coefficient re reset to their default values, the coefficient hold flags are cleared and any fit for this channel is removed.
	Func C	<pre>hanFit(chan%, type%, order%);</pre>
	chan%	The channel number to work on. Each channel in a time, result or XY view can have one fit associated with it.
	type%	0=Clear fit, 1=Exponential, 2=Polynomial, 3=Gaussian, 4=Sine, 5=Sigmoid.
	order%	The order of the fit. This is 1 or 2 for an exponential or Gaussian fit, 1 for a Sine or Sigmoid fit and 1 to 5 for a polynomial fit. If type% is 0 this should also be 0.
	Returns	0 if the command succeeded.
Porform the fit	This you	iant of the command does the fit set by the previous variant
r enorm the fit	Func C	<pre>hanFit(chan%, opt%, start start\$, end end\${, ref ref\${,</pre>
	chan%	A channel number in the current view that has had a fit initialised.
	opt%	This is the sum of:
		 Estimate the coefficients before fitting, else use current values. Draw the fit over the user-defined range, not the fit range.
	start	This is the start of the fit range as a value or as a dialog expression string that is to be evaluated.
	end	The end of the fit range in x axis units as a value or a string to evaluate.
	ref	The reference time as a value or a string to evaluate. If omitted start is used.
	err	If present, this optional variable is updated with the chi-squared or least-squares error between the fit and the data.
	max1%	If present, this sets the maximum number of iterations. If omitted, the current number set for the channel is used. The system default number is 100.
	iTer%	If present, this integer variable is updated with the count of iterations done.
	covar	An optional two dimensional array of size at least $[nCoef][nCoef]$ that is returned holding the covariance matrix when the fit is complete. It is changed if the return value is -1, 0 or 1. However, the values it contains are probably not useful unless the return value is 0.
	Returns	0 if the fit is complete, 1 if max iterations done, or a negative error code: -1=the fit is not making progress (results may be OK), -2=the fit failed due to a singular matrix, -5=the fit caused a floating point error, -6=too little data for the number of coefficients, -7=unknown fitting function, -8=ran out of memory during the fit (too many data points), -9=the fit was set up with bad parameters.
Get fit information	This var	iant of the command returns information about the current fit set for a channel.
	Func C	<pre>hanFit(chan%{, opt%});</pre>
	chan%	The channel number of the fit to return information about.
	opt%	This determines what information to return. If omitted, the default value is 0. Positive values return information about the fit that is set-up to be done next. Negative value return information about the last fit that was done and that can be displayed. The returned information for each value of opt% is:

		opt% 0 -1 -2 -3 -4 -5 -6 -7	Returns Fit type of next fit 1=a fit exists, 0=no fit exists Type of existing fit or 0 Order of existing fit Chi or least-squares error Fit probability (estimated) X axis value at fit start X axis value at fit end	opt% 1 -8 -9 -10 -11 -12 -13 -14	Returns Fit order of next fit Reference x value User-defined x draw start User-defined x draw end 1=chi-square, 0=least-square Last fit result code Number of fitted points Number of fit iterations used
	Returns	The in	formation requested by the opt%	argumer	it or 0 if opt% is out of range.
	See also	ChanE: FitPc	'itCoef(),ChanFitShow(),Ch ly()	anFitV	alue(),FitExp(),
ChanFitCoef()	This co result or and lim	mmand XY vi its and l	gives you access to the fit coeffi ew. You can return the values from hold values fixed for iterative fits.	cients for n any ty There a	or a channel in the current time, pe of fit and set the initial values re two command variants:
Set and get coefficients	This co coeffici	mmand ent valu	l variant lets you read back the les and limits for iterative fitting:	current	coefficient values and set the
	Func C	hanFi	tCoef(chan%{, num%{, new{,	, lower	<pre>c{, upper}}});</pre>
	chan%	The cl	nannel number of the fit to access.		
	num%	If this If pres presen no fit, iterativ	is omitted, the return value is the sent, it is a coefficient number. Th it, the return value is the coefficie the coefficient value that would we fit is returned.	number e first co nt value be used	of coefficients in the current fit. befficient is number 0. If num% is for the existing fit, or if there is as the starting point for the next
	new	If pres	ent, this sets the value of coefficiel.	ent num ⁹	For the next iterative fit on this
	lower	If pres on thi There	sent, this sets the lower limit for or s channel. There is currently no is also no check made that the lim	coefficie way to its are s	nt num% for the next iterative fit read back the coefficient limits. et to sensible values.
	upper	If pres on this	sent, this sets the upper limit for os channel.	coefficie	nt num% for the next iterative fit
	Returns	The nu	umber of coefficients or the value	of coeff	icient num%.
Get and set the hold flags	This co Data dia	mmand alog Co	variant sets the hold flags (equiverent efficients tab).	alent to	the Hold checkboxes in the Fit
	Func C	hanFi	tCoef(chan%, hold%[]);		
	chan%	The cl	nannel number of the fit to access.		
	hold%	An arr extra o Set ho less th corres	ray of integers to correspond with elements are ignored. If it is too s bld%[i%] to 1 to hold coefficien han 0, the hold state is not char ponding coefficient is held and to	the coest short, ex t i% and ged, bu 0 if it is	fficients. If the array is too long, tra coefficients are not affected. d to 0 to fit it. If hold%[i%] is t hold%[i%] is set to 1 if the not held.
	Returns	This a	lways returns 0.		
	<pre>See also:ChanFit(), ChanFitShow(), ChanFitValue(), FitExp(), FitPoly()</pre>				

ChanFitShow()	This controls the display of data fitted to a channel in the current time, result or XY view.		
	<pre>Func ChanFitShow(chan%{, opt%{, start start\${, end end\$}});</pre>		
	chan% The channel number of the fit to access.		
	opt%If present and positive, this is the sum of:1Display the fitted data.2Use the user-defined display range rather than the fitting range.		
	If opt % is omitted or positive, the return value is the current option value. Use negative values to return the user-defined display range: -1=return the start, -2=return the end.		
	start If present, this is an x axis value or a string holding a dialog expression to be interpreted as an x axis value that sets the start of the user-defined display range.		
	end If present, this is an x axis value or a string holding a Dialog expression to be interpreted as an x axis value that sets the end of the user-defined display range		
	Returns The current opt% value or the information requested by opt%. If there is no fit defined for the channel, the return value is 0.		
	See also:ChanFit(), ChanFitShow(), ChanFitCoef(), FitExp(), FitPoly()		
ChanFitValue()	This function returns the value at a particular x axis value of the fitted function to a channel in the current time, result or XY view.		
	<pre>Func ChanFitValue(chan%, x);</pre>		
	chan% The channel number of the fit to access.		
	x The x axis value at which to evaluate the current fit. You should be aware that some of the fitting fuctions can overflow the floating point range if you ask for x values beyond the fitted range of the function.		
	Returns The value of the fitted function at x. If the result is out of floating point range, the function may return a floating point infinity or a NaN (Not a Number) value or a 0. If there is no fit, the result is always 0.		
	See also:ChanFit(), ChanFitShow(), ChanFitCoef(), FitExp(), FitPoly()		
ChanHide()	Hide a channel, or a list of channels. Hiding a channel that doesn't exist has no effect.		
	<pre>Proc ChanHide(cSpc {,cSpc});</pre>		
	cSpc A channel specifier or -1 for all, -2 for visible, and -3 for selected channels.		
	See also:ChanShow(), ChanList()		
r			
ChanKind()	This returns the type of a channel in the current time, result or XY view.		
	Func chanking (chang)		
	meaning the first channel to be compatible with Spike2 version 3 scripts.		
	ReturnsA code for the channel type or -3 if this is not a time, XY or result view:0None/deleted3Event (Event+)6WaveMark9Real wave1Waveform4Level7RealMark120XY channel2Event (Event-)5Marker8TextMark127Result channel		

The result and XY codes changed at version 3.16. Use ViewKind() to detect result and XY views if the script must be compatible with all versions of Spike2.

See also:ChanList(), MemChan(), ViewKind()

ChanList()

This function generates a channel list in an array for a time view, result or an XY view. The channels can be filtered to show only a subset of the available channels. There are two command variants:

```
Func ChanList(list%[]{, types%});
Func ChanList(list%[], str${, types%});
```

- list% An integer array to fill with channel numbers. Element 0 is set to the number of channels returned. The remaining elements are channel numbers. If the array is too short, enough channels are returned to fill the array. It is unnecessary to list all the channel numbers for an XY view, since they are numbered contiguously.
- types% This specifies the channels to return. If omitted, all channels are returned. The values are the same as those defined for the DlgChan() function. XY views only allow exclude visible (1024) and exclude hidden (2048) channels.
- str\$ A channel specification such as "1..10,v1,m1" or "v1a". Only channels that exist in the current view are returned in list%. If types% is provided, only channels that match both the string and types% are returned in list%.
- Returns The number of channels that would be returned if the array were of unlimited length or 0 if the view is not the correct type.

See also:ChanKind(), ChanShow(), DlgChan()

	T1	
ChanMeasure()	I his per	forms the cursor regions measurements on a channel in a time or result window.
	Func C	<pre>hanMeasure(chan%, type%, sPos, ePos{, &data%{, kind%}});</pre>
	chan%	The number of the channel to measure. In a result view, channels start at 1, but we accept 0 as meaning the first channel to be compatible with version 3.
	type%	The type of measurement to take, see the documentation of the Cursor Regions window for details of these measurements. The possible values are:
		1 Area5 Area (scaled)9 Minimum13 Abs max.17 Mean in X2 Mean6 Curve area10 Peak to Peak14 Peak18 SD in X3 Slope7 Modulus11 RMS Amplitude15 Trough19 Mean of abs4 Sum8 Maximum12 Standard deviation16 RMS error
	sPos	The start position for the measurement in x axis units. If sPos is greater than or equal to ePos, the result is 0. sPos is converted to a bin number in result views; the conversion is equivalent to $Trunc(XTOBin(sPos))$.
	ePos	The end position in x axis units. In a result view, this is converted to a bin number which must be greater than the bin number obtained from $sPos$.
	data%	Optional variable returned as 1 if a valid result was obtained or as 0 if there is no result (equivalent to a blank cell in the Cursor Regions dialog).
	kind%	This optional variable forces a WaveMark channel to be treated as a waveform or as events for this measurement. If 0 or omitted, the channel is treated as a waveform if drawn in waveform. WaveMark or Cubic Spline mode, otherwise it

Returns The function returns the requested measurement value.

is treated as an event channel. 1 forces waveform and 2 forces event.

See also:ArrSum(), ChanData(), ChanValue(), XToBin()

ChanNew()	This function the created this to c	nction creates a new channe ted channel is permanent; a reate channels for use with	el in the current time any data written to it ChanWriteWave ()	view. Unlike a memory channel, occupies disk space. You can use and ChanSave().
	Func (ChanNew(chan%, type%{	, size%{, binsz{	<pre>, pre%{, trace%}}});</pre>
	chan%	The channel number to us file (usually 32) or 0 for the routine to overwrite an existing the second	se in the range 1 to the first unused disk-basisting channel; use Chan	ne number of channels in the data ased channel. You cannot use this nanDelete() to remove it first.
	type%	The type of channel to cre 1 Waveform 2 Event (Event-) 3 Event (Event+)	ate. Codes are: 4 Level (Event+-) 5 Marker (default) 6 WaveMark	7 RealMark 8 TextMark 9 Real wave
	size%	Used for TextMark, Real number of characters, real use this to set the disk l recommend that you use 0	Mark and WaveMar ls or waveform point buffer size in bytes for the default buffer	k channels to set the maximum s to attach to each item. You can for other channel types, but we r size.
	binsz	Used for waveform and W waveform points. This is resolution. If you set this (/aveMark data to spe rounded to the neares) or negative, the sma	cify the time interval between the st multiple of the underlying time illest bin size possible is set.
		If binsz is not a multiple for the file, versions of Sp	of the underlying tin ike2 before 4.03 canr	ne resolution times <i>time per ADC</i> not read it.
	pre%	This must be present for V	VaveMark data to set	the number of pre-trigger points.
	trace8	Optional, default 1, the nu	mber of interleaved t	races for WaveMark data.
	Returns The channel number if the channel was created, or a negative error code			
	Channe these w ChanOf wavefor	ls created in this way are g ith the ChanTitle\$(), Ch ffset() routines. The fol rm channel) in channel wTo	given default titles, u aanUnits\$(), Chan(llowing code creates %:	nits and comments. You can set Comment\$(), ChanScale() and a copy of channel wFrom% (a
	func C var er if Cha if Cha ChanNe if err Char Char Char Char Char char char if err if char if err if err i	CopyWave% (wFrom%, wTo cr%, buffer%[8192], n unKind (wFrom%) <>1 the unKind (wTo%) <>0 then cw (wTo%, 1,0,BinSize (c%= wTo% then Scale (wTo%, ChanScal Offset (wTo%, ChanOff Units\$ (wTo%, ChanUni Title\$ (wTo%, "Copy") Comment\$ (wTo%, "Copi eat c := ChanData (wFrom%, c n% > 0 then	<pre>%) 'Copy wavefor %, sTime := 0.0; n return -1 endi ChanDelete(wTo%) wFrom%)); e(wFrom%)); set(wFrom%)); ts\$(wFrom%)); ; ed from channel buffer%[], sTim</pre>	<pre>tm to a memory channel f; 'Not a waveform! endif; 'Check if used Create waveform channel Created OK? Copy scaleand offsetand units Set our own title "+Str\$(wFrom%)); he, MaxTime(), sTime); read ok?</pre>
	er if	n% := ChanWriteWave(dif; n% > 0 then sTime +	wTo%, buffer%[:n = n% * BinSize(w	u%], sTime) /To%) endif;
	Unti Char endif;	LL N& <= U; NShow(wTo%);	if anosted of	display new channel
	end;	i erra; Keturns U	II Created UK	
	See also):ChanData(),ChanDelet FileNew(),MemChan()	ce(),ChanSave(),	ChanWriteWave(),

ChanNumbers()	You can show and hide channel numbers in the current view and get the channel number state with this function. It is not an error to use this with data views that do not support channel number display, but the command has no effect.		
	<pre>Func ChanNumbers({show%});</pre>		
	show% If present, 0 hides the channel number, and 1 shows it. Other values are reserved (and currently have the same effect as 1).		
	Returns The channel number display state at the time of the call.		
	See also:YAxis(), YAxisMode()		
ChanOffset()	Waveform and WaveMark data is stored as 16-bit integers with a scale factor and offset to convert to user units. This function gets and/or sets the y axis value that corresponds to a 16-bit waveform data value of zero. The y axis user units for a channel are:		
	y axis value = (16-bit value)*scale/6553.6 + offset		
	<pre>Func ChanOffset(chan% {,offset});</pre>		
	chan% The channel number.		
	offset If present, this sets the new channel offset. There are no limits on the value.		
	Returns The channel offset at the time of the call if this is a waveform or WaveMark channel, or 0 if it is not.		
	See also:ChanCalibrate(),ChanScale(),Optimise()		
ChanOrder()	This command changes the order of channels in a time or result view, groups channels with a y axis so that they share a common y axis and sets the channel sorting order.		
	Func ChanOrder(dest%, pos%, cSpc); Func ChanOrder(dest%, opt%); Func ChanOrder(order%);		
	dest% The destination channel number or channel number used with opt%.		
	The drop position relative to the destination channel: -1=drop before, 0=drop on top, 1=drop after. If you drop between grouped channels, the dropped channels become members of the group (as long as they have a y axis).		
	cSpc A channel specifier for the channels to move.		
	opt% 0=returns the number of channels in the group that dest% belongs to or 0 if not grouped. 1-n returns the channel number of the nth channel in the group or 0 if no channel1 ungroups the group and returns the number of changed channels.		
	order [®] This form of the command sorts all the channels into numerical order. Set -1 for low, 1 for high numbered channels at the top and 0 to use the default channel		

ordering set by the Edit menu Preferences. Returns With cSpc it returns the number of moved channels. With order%, it returns -1 if low and 1 if high numbered channels were previously at the top.

See also:ChanWeight(), ViewStandard(), Yaxis(), YaxisLock()

ChanPort()	This ret	urns the physical hardware port that sampled data in a time view channel. ChanPort (chan%);
	chan%	The channel number in the current time view.
	Returns	The physical data port or -1 if the view is not a time view, or if the channel was not sampled. Both event and ADC ports are enumerated from 0.
	See also	<pre>D:SampleEvent(),SampleWaveform(),SampleWaveMark()</pre>
ChanProcessAdd()	This ac matchir	lds a channel process to a waveform or RealWave channel in a time view, and the effect of the Add button in the Channel Process dialog.
	Func (ChanProcessAdd(chan%, PType% {, arg1, arg2,});
	chan%	A waveform or RealWave channel number in the current time view.
	PType	 The process type. The following process types are currently defined: Rectify. Positive values are unchanged, negative values are negated. Data from waveform channels with a non-zero channel offset may be limited. Smooth. This has one argument, a time range in seconds. The output at time t is the average of the input data from times t-arg1 to t+arg1. DC Remove. There is one argument, a time range in seconds. The output at time t is the input minus the mean input from t-arg1 to t+arg1. Slope. There is one argument, a time range in seconds. The slope at time t is calculated from the points in the time range t-arg1 to t+arg1. If you apply this process the channel scale and y axis units change and the channel offset becomes 0. The output is in input units per second. Time shift. There is one argument, a time shift in seconds. A positive value shifts the trace right in the window, a negative value shifts it left. Down sample. There is one argument, the down sample ratio. Interpolate. Argument 1 is the channel to match. RMS amplitude. Argument 1 is the time range.
	arg#	Optional process arguments to match the Channel Process dialog arguments for each process type. Each process has default arguments. It is an error to provide too many arguments or for an argument to be incorrect for the process.
	Returns	The index of the added process in the list of processes for the channel.
	See also	D:ChanProcessArg(),ChanProcessClear(),ChanProcessInfo()
ChanProcessArg()	This gi Process	ves you access to the parameters of channel processes added by the Channel dialog or by the ChanProcessAdd() command.
	Func (ChanProcessArg(chan%, id% {,n% {,arg}});
	chan%	The channel number in the current time view.
	id%	The process index. The first one added is number 1, the second is 2, and so on.
	n%	An optional argument number. The first argument is 1, the second 2, and so on. If you omit n $%$, the function returns the number of arguments for this process. If you supply n $%$, the function returns the argument value at the time of the call.
	arg	An optional argument that changes process argument n %. Arguments that are out of range or illegal are ignored

Returns If n% is omitted this returns the number of arguments this process expects. If n% is present it returns the value of argument n% at the time of the call.

See also:ChanProcessAdd(), ChanProcessClear(), ChanProcessInfo()

ChanProcessClear()

This removes a one or all processes from a channel or all processes from all channels.

Func ChanProcessClear({chan% {, id%}});

- chan% A channel number in the current time view. If you omit this argument, or set it to -1, all processes for all channels are removed.
- id% The index of the process to delete. Set -1 to delete all processes for channel chan%. After deleting, any remaining processes are renumbered.

Returns 0 if at least one process was deleted, -1 if no process was deleted.

See also:ChanProcessAdd(), ChanProcessArg(), ChanProcessInfo()

ChanProcess	Info()	
-------------	--------	--

This returns information about processes attached to a channel.

Func ChanProcessInfo(chan% {,id% {, arg%}});

chan% The channel number in the current time view.

- id% Omit this to return the count of processes for the channel or set it to the index of the process to return information on. The first added process is number 1.
- arg% Omit this to return the type of process index id%; see ChanProcessAdd() for the type codes. If present, the command returns the type of argument arg%: 0=integer, 1=real, 2=string (unused), 4=channel. The first argument is number 1.

Returns The information requested or a negative error code.

See also:ChanProcessAdd(), ChanProcessArg(), ChanProcessClear()

ChanSave()

This function copies data with an optional time shift from source channels in the current time view to disk-based or memory buffer destination channels in any time view. If a destination channel exists, ChanSave() adds data to it, otherwise it creates a disk-based channel to match the source channel. You can set the number of disk-based channels in a new data file with the nChans% argument of FileNew().

The source data must be compatible with the destination. All event-based data types are mutually compatible. Destination data that is not present in the source is set to 0. Event times are set as close to the original times as the time bases of the source and destination views allow. If the conversion causes multiple events to fall at the same time, the converter adds 1 tick to separate the times. However, it will not do this for more consecutive events that the ratio of the two time bases.

If the source channel has a marker filter set, only events in the source channel that match the marker filter are copied to the destination channel. If the destination channel has a marker filter, the filter is unchanged for an existing channel and is set to accept all data for a new channel.

Adc and RealWave channels are compatible; conversion between them uses the channel scale and offset values. If the source and destination sample rates do not match, the data is interpolated using cubic splines. This is also done for non-matching WaveMark data.

If a destination disk-based channel already contains data, any added event-based data must be written after all data already contained in the channel. You can over-write existing waveform data with the same restrictions as for ChanWriteWave().

Func ChanSave(cSpc, dest%{, dh%{, sTime{, sTo{, dTime}}});

- cSpc A channel specifier for the source channels or -1 for all, -2 for visible and -3 for selected channels. If you set multiple channels, dest% must be set to 0 or -1.
- dest% This sets the destination channel or channels in the view identified by dh%. If cSpc refers to a single channel, this can be a channel number of a disk-based channel or an existing memory channel; otherwise dest% must be 0 or -1.

If dest% is 0, the lowest unused disk-based channels are used. If dest% is -1, the same channel number as the source is used. If the destination channel does not exist, it is created using the source channel settings.

- dh% The handle of the destination view. Omit or set 0 to use the current time view.
- sTime The start time to read source data. This is 0.0 if omitted.
- sTo Source data is read up to this time. If omitted, the end of the data is used.
- dTime The destination time. If this is omitted, sTime is used. The source data is time shifted by dTime-sTime before being copied to the destination channels.
- Returns With a single source channel, the return value is the destination channel number. With multiple source channels, the return value is the number of channels that copied without error. Negative return values mean: -1 = cSpc illegal, -2 = dest illegal, -3 = dh is not a time view handle, -4 = could not create a destination channel or an existing destination channel was not compatible with the source, -5 = a file system error occurred when copying data.

Append files example This command has many applications. This example creates a new file and adds the contents of two data files to it. We assume the files hold the same types of channels. In this case we create an output file with the same time resolution as the first input file.

var fh%,dh%; 'file handles of the source and destination files fh% := FileOpen("file1.smr", 0, 0); 'open file1 dh% := FileNew(7, 1, View(fh%).Binsize()*1e6, 1, 1); 'make output View(fh%);ChanSave(-1, -1, dh%);FileClose(0,-1); 'copy and close fh% := FileOpen("file2.smr", 0, 0); 'open file2 and append ChanSave(-1, -1, dh%, 0, MaxTime(), View(dh%).MaxTime()+1.0); FileClose(0,-1); 'close file2 View(dh%);ChanShow(-1);Draw(0, MaxTime()); 'Display the result

See also:ChanNew(), ChanOffset(), ChanScale(), ChanWriteWave(), FileNew()

ChanScale()

This function gets and/or sets the scaling between the 16-bit integer data used to store waveform and WaveMark data and the real units of the channel. The 16-bit waveform data has values between +32767 and -32768. The y axis user units for a channel are:

```
y axis value = (16-bit value)*scale/6553.6 + offset
```

With the standard ± 5 Volt 1401 ADC inputs, 6553.6 is equivalent to 1 Volt at the 1401 input. In this case, a scale of 1.0 and an offset of 0.0 gives a v axis calibrated in Volts.

```
Func ChanScale(chan% {,scale});
```

chan% The channel number.

scale If present, sets the channel scaling. We suggest you don't set 0.0 as a scale!

Returns The scale at the time of the call for waveform or WaveMark channels, or 0.

See also:ChanCalibrate(), ChanOffset(), Optimise()

ChanSearch()

This searches a time view channel for a user-defined feature. It is the same as an active cursor search, but does not use or move cursors. Searches are done in the current channel drawing mode. WaveMark data drawn as a waveform is searched as a waveform; to search it as events, set an event drawing mode. There is no need for the channel to be visible; you only need to set the display mode before searching.

```
Func ChanSearch(chan%, mode%, sT, eT{, sp1{, sp2{, width
                                             {, flag%{, lv2}}};
```

- chan% The number of a channel in the time view to search.
- This sets the search mode, as for the active cursors. Modes in italics (15 and 16) mode% cannot be used: mode 16 is not relevant and mode 15 can be emulated by mode 6, 7 or 8. See the cursor mode dialog documentation for details of each mode.
 - 8 Falling threshold 1 Maximum value 15 Repolarisation %
 - 2 Minimum value
 - 9 Steepest rising
 - 3 Max excursion 10 Steepest falling
 - 4 Peak find

6 Threshold

- 11 Steepest slope (+/-) 18 Slope% 12 Slope peak
- 5 Trough find 13 Slope trough
- 20 Inside thresholds

16 Expression

17 Turning point

19 Outside thresholds

- 7 Rising threshold 14 Data points
- sT, eT The search start and end times. If eT is less than sT, the search is backwards.
- This is the amplitude for peaks, threshold level for threshold crossings, baseline sp1 level for maximum excursion and number of points for mode 14. It is in the y axis units of the search channel (y axis units per second for slopes). If omitted, the value 0.0 is used. Set it to 0 if sp1 is not required for the mode.
- sp2 This is hysteresis for threshold crossings and percent for Slope%. If omitted, the value 0.0 is used. Set it to 0 if sp2 is not required for the mode.
- This is the width in seconds for slope measurements. It sets the minimum time width that the data must be above/below a level for threshold measures. It sets the maximum allowed width of a peak or trough. If omitted, 0 is used. Set 0 if width is not required for the mode or you do not want a time constraint.
- This is the sum of flag values and is 0 if omitted. At the moment, the only value flag% defined is 1=ignore gaps in waveform data. If this value is not set, a search of a waveform channel will stop at a gap in the data.
- This is the second level for modes 19 and 20. If omitted, 0 is used. 1 v 2

Returns A positive time if the search succeeds or -1 if it fails or a negative error code.

See also:ChanMeasure(), ChanValue(), CursorActive()

ChanSelect()	This function is used to report on the selected/unselected state of a channel in a time or
	result view, and to change the selected state of a channel.
	<pre>Func ChanSelect(chan% {,new%});</pre>
	chan% The channel number or -1 for all visible channels. Hidden channels and non- existent channels always appear to be unselected.
	new% If present it sets the state: 0 for unselected, not 0 for selected. If omitted, the state is unchanged. Attempts to change invisible channels are ignored.
	Returns The channel state at the time of the call, 0 for unselected, 1 for selected. If you set chan% to -1, the function returns the number of selected channels.
	See also:ChanHide(),ChanOrder(),ChanShow(),ChanWeight()
Ob as Ob as ()	Display a channel, or a list of channels in a time, result or XY view. Turning on a channel
ChanShow()	that is on has no effect. Turning on a channel that doesn't exist has no effect.
	<pre>Proc ChanShow(cSpc {,cSpc});</pre>
	cSpc A channel specification for the channels to show.
	See also:ChanHide(),ChanList(),ChanVisible()
ChanTitle\$()	This returns or sets the channel title string in a time, result or XY view. A title string is up to 9 characters long. In an XY view the channel titles are visible in the Key window.
	Setting a duplicate channel title does not change the original channel. Setting the original channel title sets the titles of duplicates unless the duplicates have their own title. If you set the title of a duplicate to "", the title reverts to the title of the channel it duplicates.
	<pre>Func ChanTitle\$(chan%{,new\$});</pre>
	chan% The channel number. In XY views you can also use 0 to change the y axis title.
	new\$ An optional string holding the new channel title.
	Returns The channel title string for the channel. If the channel does not exist, the function does nothing and returns an empty string.
	See also:ChanComment\$(), XYKey()
ChanUnits\$()	This gets or sets the units for waveform or WaveMark channels and result or XY views.
	From 5.03 you can read back and set the displayed units of time view event-based channels. Changes to event-based channels are lost when the display mode changes.
	<pre>Func ChanUnits\$(chan%{,new\$});</pre>
	<pre>chan% A channel in the time view. In a result view, channels start at 1, but we accept 0 as meaning the first channel to be compatible with version 3. In an XY view the channel number is ignored. We suggest you use 0 to match ChanTitle\$().</pre>
	new\$ An optional string holding the new y axis units for the channel (up to 5 characters). If the string is too long, it is truncated.
	Returns It returns the old units of the designated channel. If the channel does not exist or is not of a suitable type, the function does nothing and returns an empty string.
	See also:ChanScale(), ChanOffset(), ChanTitle\$()

ChanValue()

In time views this returns the value of a channel at a time. For a waveform channel it returns the waveform value, for other channel types, it returns a value in the y axis units of the channel display mode. If the channel has no y axis or is drawn in raster mode, the value is the time of the next event on the channel.

In a result view, this returns the value of the result corresponding to an x axis value. You can use the [bin] notation to access result views by bin number.

Func ChanValue(chan%,pos{,&data%{,mode%{,binsz{,trig%|edge%}}});

- chan% The channel number in the time or result view.
- pos In a time view, the time for which the value is needed. In a result view, this is the x axis value for which a result in needed.
- data% Returned as 1 if there is data at pos, 0 if not. For example, for a waveform if there was no data within Binsize (chan%) of pos, this would be set to 0.
- mode% If present, this sets the display mode in which to read the value from a time view. If mode% is inappropriate or absent, the current display mode is used. This parameter is ignored in a result view. The modes in a time view are:
 - 0 The current mode for the channel. Any additional arguments are ignored.
 - 1 Dots mode for events. Returns the event time at or after pos.
 - 2 Lines mode, result is the same as mode 1.
 - 3 Waveform mode.
 - 4 WaveMark mode, returns waveform values for WaveMark channels.
 - 5 Rate mode. The binSz argument sets the width of each bin.
 - 6,11 Mean frequency mode, binSz sets the time period, 11 is rate per minute.
 - 7,12 Instantaneous frequency, returns value at next event, 12 is per minute.
 - 8 Raster mode, trig% sets the trigger channel, result as for mode 1.
 - 13 Cubic spline for waveforms and WaveMark data.
 - 16 Skyline mode for waveform, RealWave and RealMark channels.
- binSz This sets the width of the rate histogram bins and the smoothing period for mean frequency mode when specifying your own mode.
- trig% The trigger channel for raster displays, we assume raster displays of level data are never required.
- edge% For level data event channels. This sets which edges of the signal to use for mean frequency, instantaneous frequency and rate modes: 0=both edges, 1=rising edges, 2=falling edges. If edge% is omitted, both edges are used.
- Returns It returns the value or 0 if no data is found. For waveform data, if there is no data within Binsize(chan%) of the time, the value is zero.

If data% is omitted any error stops the script. Errors include: no current window, current window not a time or result view, no data at pos, and pos beyond range of x axis. If data% is present, errors cause it to be set to 0.

See also:ChanData()

ChanVisible()

This returns the state of a channel in a time, result or XY view as 1 if the channel is visible and 0 if it is not. If you use a silly channel number, the result is 0 (not displayed).

Func ChanVisible(chan%);

chan% The channel to report on.

Returns 1 if the channel is displayed, 0 if it is not.

See also:ChanShow(), ChanHide()

ChanWeight()

This function sets the relative vertical space to give a channel or a list of channels. The standard vertical space corresponds to a weight of 1. When Spike2 allocates vertical space, channels are of two types: channels with a y axis and channels without a y axis. Spike2 calculates how much space to give each channel type assuming all channels have a weight of 1. Then the actual space allocated is proportional to the standard space multiplied by the weight factor. This means that if you increase the weight of one channel, all other channels get less space in proportion to their original space.

Func ChanWeight(cSpc{, new});

- cSpc The specification for the list of channels to process. See the *Script language syntax* chapter for a definition of channel specifiers.
- new If present, a value between 0.001 and 1000.0 that sets the weight for all the channels in the list. Values outside this range are limited to the range.

Returns The command returns the channel weight of the first channel in the list.

See also:ChanOrder(), ViewStandard()

ChanWriteWave()

This function writes real or integer data to a waveform (16-bit integer) or RealWave (32bit floating point) channel. The time gap between array points is the Binsize() value of the channel. You can overwrite existing data and add data to the end of the channel. You cannot fill in gaps in wave channels; values written into gaps in previously written data are ignored.

Func ChanWriteWave(chan%, arr[]|arr%[], sTime);

- chan% The waveform or RealWave channel in the current time view to write data to. This can be a duplicate channel, a disk channel or a memory channel. If you write to a duplicated channel, the original channel data is changed.
- A real or integer array to write to the channel. When writing real data to a waveform channel or integer data to a RealWave channel, the data is converted to match the channel format using the channel scale and offset. When writing to a waveform, output is limited to 16-bit integers in the range -32768 to 32767.
- sTime The first array point time. When overwriting, if the time does not align with existing data it is reduced by less than one sample interval to align it.
- Returns The number of points processed including points skipped due to gaps in existing channel data or a negative error code, for example if the file is read-only.

The function will cause a fatal script error if used on the wrong view type, the wrong channel type or if the system runs out of memory.

See also:Binsize(), ChanData(), ChanNew(), ChanOffset(), ChanScale()

Chr\$()

This function converts a code to a character and returns it as a single character string.

Func Chr\$(code%);

code% The code to convert. Codes that have no character representation will produce unpredictable results.

See also:Asc(), DelStr\$(), LCase\$(), Left\$(), Len(), Mid\$(), Print\$(), Right\$(), Str\$(), UCase\$(), Val()

Colour()	This fur palette, The co ChanCo	nction gets and/or sets the colours of items. Colours are set in terms of the colour not directly in terms of colours. XY channels are coloured using XYColour(). lours set for time and result view drawing modes can be overridden by lour().
	Func C	<pre>olour(item% {,col%});</pre>
	item%	This selects the item to be coloured. The available items are:
		1Time background11Rate fill29Controls (not used)2Waveform channel12Result background30Grid colour3Events as dots13Result lines31X and Y Axes4Events as lines14Result dots32XY background5Level events15Result skyline33Not saving to disk6Marker data16Result histogram34Result SEM and SD7Mean frequency17Result histogram fill35Fitted curves8Inst. Frequency18-26WaveMark codes 0-836Cluster background9Raster dots27Channel number37WaveMark background10Rate outline28CursorsSecuritySecurity
	col%	If present, this sets the index of the colour in the colour palette to be applied to the item. There are 40 colours in the palette, numbered 0 to 39.
	Returns	The index into the colour palette of the colour of the item at the time of the call.
	See also	ChanColour(), PaletteGet(), PaletteSet(), XYColour()
Conditioner commands	The Cor program Axon In	ad family of commands control signal conditioners. They support the CED 1902 mable signal conditioner, the Power 1401 programmable gain options and the struments CyberAmp. Other conditioners may be added in the future.
	These c of condi set the port% a attached	ommands do not select which serial port (if any) the conditioner uses or the type itioner supported. You choose the conditioner type when you install Spike2. You serial port in the Edit menu Preferences option. All these commands require a argument. This is the physical waveform input port number that the conditioner is I to. It is not a channel number in a time view.
	You can Configut lists of g	n access the built-in interactive support for the conditioner from the Sampling ration channel parameters dialog. This can be a useful short-cut to getting the gains and signal sources available on your conditioner(s).
	See also	<pre>:CondFilter(), CondFilterList(), CondFilterType(), CondGain(), CondGainList(), CondGet(), CondOffset(), CondOffsetLimit(), CondRevision\$(), CondSet(), CondSourceList(), CondType()</pre>
CondFeature()	This co	mmand gets and sets special signal conditioner features that are not general
	more de	tails of conditioner operation. There are four command variants:
Get feature count	Func C	ondFeature(port%);
	port%	The waveform port number that the conditioner is connected to.
	Returns	The number of special features supported by the signal conditioner.
Get feature information	The con as ["Noi such as call the	nmand supports two types of features: those that have a set of discrete values such ne", "Rectify"], and those that support a continuous range of floating point values 10.0 to 25.6, for example. To determine if the feature is continuous or discrete function with only the first 3 or 4 arguments.

	<pre>Func CondFeature(port%, feat%, &name\${,&flags%{, list\$[]}}); Func CondFeature(port%, feat%, &name\${,&flags%{, &low{, &high}}});</pre>				
	feat% The feature number, from 1 to the number of features available				
	name\$ Returned set to the name of the feature.				
	flags $%$ Returned set to the feature flags. Currently, none are defined, so this will be 0.				
	list\$ Returned set to an array of the possible settings (strings) for discrete type.				
	low Returned as the low limit for a feature with continuous values.				
	high Returned as the upper limit for a feature with continuous values				
	Returns The number of discrete feature values, or 0 if the feature supports continuate values over the range returned in low and high.				
Set feature value	<pre>Func CondFeature(port%, feat%{, val});</pre>				
	val If present, it sets the value for the feature set by feat% and port%. If this is a continuous feature, val sets the new value. If the feature has n discrete setting, val should be 0 to n-1 to select the feature corresponding to the feature description in list\$[val]. If val exceeds the allowed range, a continuous feature is set to the nearest allowed value and a discrete feature is unchanged.				
	Returns The feature value at the time of the call (before any change). This is an integer index for features with discrete values otherwise it is the feature value.				
	<pre>See also:CondFilterList(), CondGain(), CondGainList(), CondGet(), CondSet(), CondSourceList(), CondType()</pre>				

CondFilter()

This sets or gets the frequency of the low-pass or high-pass filter of the signal conditioner. See the CondSet () command for more details of conditioner operation.

Func CondFilter(port%, high% {,freq{, type%}});

- port% The waveform port number that the conditioner is connected to.
- high% This selects which filter to set or get: 0 for low-pass, 1 for high-pass.
- freq If present, this sets the desired corner frequency of the selected filter. See the CondSet() description for more information. Set 0 for no filtering. If omitted, the frequency is not changed. The high-pass frequency must be set lower than the frequency of the low-pass filter, if not the function returns a negative code.
- type% Optional, taken as 1 if omitted. The filter type to use when setting the filter.
- Returns The cut-off frequency of the selected filter at the time of call, or a negative error code. A return value of 0 means that there is no filtering of the type selected.
- See also:CondFilterList(), CondGain(), CondGainList(), CondGet(), CondOffset(), CondOffsetLimit(), CondRevision\$(), CondSet(), CondSourceList(), CondType()

CondFilterList()

This function gets a list of the possible filter frequencies of the conditioner. Conditioners that support continuous frequency ranges also supply a list of frequencies to match the list of frequencies shown in the conditioner control panel. See the CondSet() command for more details of conditioner operation.

Func CondFilterList(port%, high%, freq[]{, type%});

- port% The waveform port number that the conditioner is connected to.
- high% Selects which filter to get: 0 for low-pass, 1 for high-pass.

	Ereq[] An array of reals holding the cut-off frequencies of the selected filter. A value of 0 means no filtering, that is the filter can be switched out of circuit.		
	<pre>Cype% Optional, taken as 1 if omitted. The filter type in the range 1 to the number of types (as returned by CondFilterType()).</pre>		
	Returns The number of filtering frequencies of the conditioner or a negative error code.		
	<pre>See also:CondFilter(), CondFilterType(), CondGet(), CondOffset(), CondOffsetLimit(), CondSet(), CondSourceList(), CondType()</pre>		
CondFilterType()	This function returns information about the filter types supported by the conditioner for he low pass and high pass filters. For example, the CED 1902 mk IV supports a choice of filter types. There are three command variants:		
Get number of filter types	<pre>Func CondFilterType(port%, high%);</pre>		
	The waveform port number that the conditioner is connected to.		
	Selects which filter to return information for: 0 for low-pass, 1 for high-pass.		
	Returns The number of filter types.		
Get filter type in use	June CondFilterType(port%, high%, 0):		
Get miler type in use	Returns. The currently selected filter type, from 1 to the number of filter types		
	counts the currently selected inter type, non t to the number of inter types.		
Get filter type information	<pre>'unc CondFilterType(port%,high%,type%{,&name\${,&lower{,&upper}}});</pre>		
	The filter number, from 1 to the number of available filters.		
	name\$ If present, returned holding the name of the filter selected by type%.		
	ower If present, returned as the lowest filter frequency (excludes 0, meaning 'off').		
	upper If present, returned as the highest supported filter frequency.		
	Returns The number of frequency values the filter can be set to (including 0) or 0 if the filter corner frequency can be set to any value in the range lower to upper.		
	See also:CondFilter(),CondGain(),CondGainList(),CondGet(),CondSet(), CondSourceList(),CondType()		
CondGain()	This sets and gets the gain of the signal passing through the signal conditioner See the CondSet() command for more details of conditioner operation.		
	<pre>Func CondGain(port% {,gain});</pre>		
	The waveform port number that the conditioner is connected to.		
	Jain If present this sets the ratio of output signal to the input signal. If this argument is omitted, the current gain is returned. The conditioner will set the nearest gain it can to the requested value.		
	Returns The gain at the time of call, or a negative error code.		
	<pre>See also:CondFilter(), CondFilterList(), CondGainList(), CondGet(), CondOffset(), CondOffsetLimit(), CondSet()</pre>		

CondGainList()

This function gets a list of the possible gains of the conditioner for the selected signal source. See the CondSet() command for more details of conditioner operation.

Func CondGainList(port%, gain[]);

- port% The waveform port number that the conditioner is connected to.
- gain[] An array of reals holding the conditioner gains for the selected signal source. If a conditioner (for example, 1902) has a fixed set of gains, this is the set of gain values. If the conditioner supports continuously variable gain, the first two elements of this array hold the minimum and the maximum values of the gain.
- Returns The number of gain values if the conditioner has a fixed set of gains or 2 if the conditioner has continuously variable gain. In the case of an error, a negative error code is returned.

CondGet()

This function gets the input signal source of the signal conditioner, and the conditioner settings for gain, offset, filters and coupling. The settings are returned in arguments which must all be variables. See CondSet () for details of conditioner operation.

- port% The waveform port number that the conditioner is connected to.
- in% Returned as the zero-based index of the input signal source (see CondSet ()).
- gain Returned as the ratio of output to input signal amplitude (ignoring filtering).
- offs A value added to the input waveform to move it into a more useful range. Offset is specified in user units and is only meaningful when DC coupling is used.
- low Returned as the cut-off frequency of the low-pass filter. A value of 0 means that there is no low-pass filtering enabled on this channel.
- hi Returned as the cut-off frequency of the high-pass filter. A value of 0 means that there is no high-pass filtering enabled on this channel.
- notch% Returned as 0 if the mains notch filter is off, and 1 if it is on.
- ac% Returned as 1 for AC or 0 for DC coupling.
- typeL% Optional integer variable returned holding the low-pass filter type number as described for CondFilterType().
- typeH% Optional integer variable returned holding the high-pass filter type number.
- Returns 0 if all well or a negative error code.
- See also:CondFilter(), CondFilterType(), CondFilterList(), CondGain(), CondGainList(), CondOffset(), CondSet(), CondSourceList()

CondOffset()

This sets or gets the offset added to the input signal of the signal conditioner. See the CondSet() command for more details of conditioner operation.

Func CondOffset(port%{, offs});

- port% The waveform port number that the conditioner is connected to.
- offs The value to add to the input waveform of the conditioner to move it into a more useful range. If this argument is omitted, the current offset is returned. The conditioner will set the nearest value it can to the requested value.

Returns The offset at the time of call, or a negative error code.

See also:CondGain(), CondGainList(), CondGet(), CondOffsetLimit(), CondRevision\$(), CondSet(), CondSourceList(), CondType()

CondOffsetLimit()	This function gets the maximum and minimum values of the offset range of the conditioner for the currently selected signal source.			
	<pre>Func CondOffsetLimit(port%, offs[]);</pre>			
	port% The waveform port number that the conditioner is connected to.			
	offs[] This is an array of real numbers returned holding the minimum (offs[0]) and the maximum (offs[1]) values of the offset range of the conditioner for the currently selected signal source.			
	Returns 2 or a negative error code.			
	<pre>See also:CondGain(), CondGainList(), CondGet(), CondOffset(), CondRevision\$(), CondSet(), CondSourceList(), CondType()</pre>			
	This function returns the name and version of the signal conditioner as a string or a blank			
CondRevision\$()	string if there is no conditioner for the port.			
	<pre>Func CondRevision\$(port%);</pre>			
	port% The waveform port number that the conditioner is connected to.			
	Returns This returns a string describing the signal conditioner. The strings defined so far are: "1902ssh", where ss is the 1902 ROM software version number and h is the hardware revision level; and "CYBERAMP 3n0 REV x.y.z" where n is 2 or 8. If there is no conditioner attached to the port it returns an empty string.			
	See also:CondFeature(),CondFilter(),CondType()			
CondSet()	This sets the input signal source, gain, offset, filters and coupling of the conditioner. All values are requests. The actual values set will depend on the capabilities of the conditioner. In all cases, the command sets the nearest value to that requested. If it is important to know what has actually been set you should read back the values with CondGet() after setting them, or use the functions for reading specific values.			
	<pre>Func CondSet(port%, in%, gain, offs {,low, high, notch%, ac%</pre>			

- port% The waveform port number that the conditioner is connected to.
- in% A conditioner has one or more signal sources. For example, the CED 1902 supports Grounded, Single ended, Normal Diff, Inverted Diff, etc. Conditioners of the same type may have different sources. To select a source, set in% to its zero-based index in the list returned by CondSourceList().
- gain This is the desired ratio of output signal amplitude to the input signal amplitude (ignoring the effect of any filtering). The actual gain depends on the capabilities of the signal conditioner, see CondGainList(). The gain range may be altered by the choice of signal source. For example, the 1902 Isolated Amp input has a build-in gain of 100. This command sets the nearest gain to the requested value.
- offs This is the desired value in user units to add to the input waveform to move it into a more useful range. Offsets are only meaningful with DC coupling. Different conditioners have different offset ranges, and the offset range may be altered by the choice of signal source, see CondOffsetLimit(). The command will set the nearest offset it can to the desired value.
- 10W If present and greater than 0, it is the desired corner frequency of the low-pass filter. Low-pass filters are used to reduce the high frequency content of the signal, both to satisfy the sampling requirement, and in case where it is known that no useful information is to be found in the signal above a certain frequency.

If omitted, or 0, there is no low-pass filtering. The actual filter value set depends on the capability of the signal conditioner.

high If present and greater than 0, it is the high-pass filter corner frequency. Highpass filters reduce the low-frequency content of the signal. This must be set lower than the frequency of the low-pass filter; if not, the function returns a negative code. If omitted, or set to 0, there is no high-pass filtering.

Different signal conditioners have different ranges of frequency filtering. To find out the real filter frequency set, use CondFilter(). CondFilterList() returns the list of possible filter frequencies.

- notch% Some signal conditioners have a mains-frequency notch filter (usually 50 Hz or 60 Hz) used to reduce the effect of mains interference on low level signals. This filter will remove the fundamental 50 Hz or 60 Hz signal; it will not remove higher harmonics (for example 150 Hz). If notch% is present with a value greater than 0, the notch filter is on. If omitted, or 0, the notch filter is off.
- ac% The 1902 supports both AC and DC signal coupling. If you set AC coupling you should probably set the offset to zero. If ac% is greater than 0, the signal conditioner is AC coupled. If omitted or 0, the signal conditioner is DC coupled.
- typeL% Optional value, taken as 1 if omitted, that sets the low-pass filter type as described for CondFilterType() in the range 1 to the number of filter types.
- typeH% Optional value, taken as 1 if omitted, that sets the high-pass filter type.

Returns 0 if all is well or a negative error code.

See also:CondFilter(), CondFilterType(), CondGain(), CondGet(), CondOffset(), CondRevision\$(), CondType()

CondSourceList()	This function gets a list of the possible signal source names of the conditioner, or the specific signal source name with the given index number.			
	Func C	ondSourceList(port%, src\$[] src\$ {,in%});		
	port%	The waveform port number that the conditioner is connected to.		
	src\$	This is either a string variable or an array of strings that is returned holding the name(s) of signal sources. Only one name is returned per string.		
	in% This argument lets you select an individual source or all sources. If present a greater than or equal to 0, it is the zero-based index number of the signal sou to return. In this case, only one source is returned, even if src\$ is an array.			
		If omitted and src\$ is a string, the first source is returned in src\$. If src\$[] is an array of strings, as many sources as will fit in the string array are returned.		
	Returns If in% is greater than or equal to 0, it returns 1 or a negative error code. If in% is omitted, it returns the number of signal sources or a negative error code.			
	<pre>See also:CondFilterList(), CondGainList(), CondGet(), CondOffsetLimit(), CondRevision\$(), CondSet(), CondType()</pre>			
CondType()	This function returns the type of the signal conditioner.			
Func CondType(port%);		ondType(port%);		
	port%	The waveform port number that the conditioner is connected to.		
	Returns 0 for no conditioner or it is not the type set when installing, 1 for a CED 1902 for an Axon Instruments CyberAmp and 3 for Power1401 with gain controls.			
	See also	CondFilterList().CondGainList().CondRevision\$().		

CondSourceList(), CondGainList(), CondRevi
Cos()	This calculates the cosine of one or an array of angles in radians.
	Func Cos(x x[] x[][]);
	The angle, expressed in radians, or a real array of angles. The best accuracy of the result is obtained when the angle is in the range -2π to 2π .
	Returns When the argument is an array, the function replaces the array with the cosines of all the points and returns either a negative error code or 0 if all was well. When the argument is not an array the function returns the cosine of the angle.
	See also:Abs(), ATan(), Cosh(), Exp(), Ln(), Log(), Pow(), Sin(), Sqrt(), Tan()
Cosh()	This calculates the hyperbolic cosine of one value or an array of values.
	<pre>Func Cosh(x x[] x[][]);</pre>
	x The value, or a real array of values.
	Returns When the argument is an array, the function replaces each value with its hyperbolic cosine and returns 0. When the argument is not an array the function returns the cosh of the argument.
	See also:ATan(),Cos(),Exp(),Ln(),Log(),Pow(),Sinh(),Sqrt(),Tanh()
Count()	In a time view this counts events and forms the mean level of a waveform channel. It can also be used in a result view to sum bin contents between a start and end bin number.
	<pre>Func Count(chan%, start, finish);</pre>
	chan% The channel number in a time or result view. If the channel does not exist, -1 is returned.
	start The start time/bin. If start is greater than finish, the result is 0.
	finish The last time/bin. In a time view, if start equals finish, only items that fall exactly at the time count towards the result.
	Returns In a time view, for waveform channels it returns the mean waveform level in the time range (gaps in waveforms are ignored). If the time range falls entirely in a gap, the result is 0. For all other channels, it returns the number of events. In a result view it returns the sum of the bins in the range.
	See also:ArrSum(), ChanData(), ChanMeasure()
Cursor()	This returns a cursor position and optionally sets a new position. If you move cursor 0 in a time view, all active cursors 1-9 will search, equivalent to CursorSearch(1).
	<pre>Func Cursor(num% {,where});</pre>
	num% The cursor number to use in the range 1 to 9 (0 to 9 in a time view).
	where If present, the new cursor position. If this exceeds the x axis range, it is limited to the x axis. In a time view the position is in seconds. In a result view it is a bin number, use XTOBin() to convert an x axis value to a bin number. In an XY

Returns The old cursor position or -1 if the cursor doesn't exist.

view the position is in x axis units.

```
See also:BinToX(),XToBin(),CursorDelete(),CursorLabel(),
CursorLabelPos(),CursorNew(),CursorRenumber(),CursorSet(),
CursorVisible()
```

CursorActive()	This function retrieves the current active cursor mode and optionally sets a new mode and search parameters. This is valid for a time view only. The function is equivalent to the vertical cursor mode dialog. Once you have set a cursor mode you can command it to seek with CursorSearch() and tell if the search succeeded with CursorValid().	
	Func C [.] Func C [.]	ursorActive(num%); ursorActive(num%, mode%{, ch%, str\$ minSt, end\${, def\${, ampLev n%{, hyPer{, width{, ref\${, lev2}}}}});
	num%	This is the cursor number, from 0 to 9.
	mode%	If this argument is present, it sets the new cursor mode. Modes in <i>italics</i> cannot be used for cursor 0 and are converted to Static mode. See the documentation for the cursor mode dialog for details of each mode.
		0Static7Rising threshold14Data points1Maximum value8Falling threshold15Repolarisation %2Minimum value9Steepest rising16Expression3Maximum excursion10Steepest falling17Turning point4Peak find11Steepest slope (+/-)18Slope%5Trough find12Slope peak19Outside thresholds6Threshold13Slope trough20Inside thresholds
	ch%	The time view channel to search. Use 0 for modes that do not require a channel.
	str\$	This is an expression that sets the start time for the search when num% is not 0. In expression mode (16) this is the expression to evaluate.
	minSt	This sets the minimum step for cursor 0 in all modes except 0 and 16.
	end\$	This string expression sets the end limit of the search. This is ignored for cursor 0 operations when it should be an empty string.
	def\$	Optional. If a search fails, and this string evaluates to a valid time, the cursor is positioned at the time and the position is valid. This is ignored for cursor 0 operations when it should be an empty string.
	ampLev	A value or string that sets the amplitude for peaks, threshold level for threshold crossings and baseline level for maximum excursion. It is in the y axis units of the search channel (y axis units per second for slopes). If omitted, the value 0.0 is used. Set it to 0 or an empty string if ampLev is not required for the mode.
	n%	The data point count for mode 14. A value of 0 is treated as 1.
	hyPer	The hysteresis for threshold crossings and percent for modes 15 and 18. If omitted, 0 is used. Set it to 0 if hyper is not used by the mode.
	width	This is the width in seconds for all slope measurements. It sets the reference level measurement width in mode 15. It sets the minimum time that the data must be above/below a level for threshold measures. It sets the maximum allowed width of a peak or trough. If omitted, 0 is used. Set it to 0 if width is not used in the mode or you do not want a time constraint.
	ref\$	This string expression is used in mode 15 to set the time at which the 100% value is measured. The 0% value is measured at the start time.
	lev2	This sets the second threshold level in modes 19 and 20. If omitted, 0 is used. You can supply this as a value or as a string to evaluate.
	Returns	The active cursor mode at the time of the call.
	The arguto a time	uments str\$, end\$, def\$ and ref\$ are strings holding expressions that evaluate ne in seconds. They are typically of the form "Cursor(0)+1.3". They can

See also:CursorActiveGet(), CursorNew(), CursorSearch(), CursorValid(), MeasureChan(),MeasureX()

contain any expression that would be valid in the Cursor mode dialog.

CursorActiveGet()	CursorActiveGet() This gets active cursor parameters set by CursorActive() or the cursor mode dialog.		
	Func CursorActiveGet(num%, item% {,&val\$});		
	num% I his is the cursor number from 0 to 9.		
	1 ch% 3 end\$ 5 ampLevIn% 7 width 9 ref\$		
	2 str $4 def$ $6 hyPer$ $8 minSt$ $10 lev2$		
	val\$ This optional string variable is updated with the parameter value when item% is 2, 3, 4, 5, 9 or 10.		
	Returns The numerical value of items 1, 5, 6, 7, 8 and 10. Items 5 and 10 may return 0 if the item is an un-parseable string. Otherwise it returns the active cursor mode.		
	See also:CursorActive(),CursorSearch(),CursorValid(),MeasureX()		
CursorDelete()	Deletes a cursor. It is not an error to delete an unknown cursor (which has no effect). You cannot delete cursor 0; use CursorVisible(0,0) to hide cursor 0.		
	<pre>Func CursorDelete({num%});</pre>		
	num% The cursor number to delete. If omitted, the highest numbered cursor is deleted. Use -1 to delete all cursors 1-9.		
	Returns The value of nume if any cursors were deleted or 0 if no cursor was deleted.		
	See also:Cursor(), CursorNew(), CursorRenumber(), CursorSet()		
CursorExists()	Use this function to determine if a vertical cursor exists.		
	Func CursorExists(num%);		
	num [®] The cursor number in the range 0-9. Cursor 0 always exists in a time view and never exists in any other view.		
	Returns 0 if the cursor does not exist, 1 if it does.		
	See also:CursorDelete(),CursorNew(),CursorValid(),HCursorExists()		
	This command sets (or gets) the current view cursor label style. Cursors can be annotated		
CursorLabel()	with a position and/or the cursor number, or with a user-defined string:		
	<pre>Func CursorLabel({style%{, num%{, form\$}}});</pre>		
	style% Label styles are: 0=None, 1=Position, 2=Number, 3=Both, 4=User-defined. Unknown styles cause no change. Style 4 is used with a format string. Styles 0-3 set the styles of cursors selected by the num% field and the view style for new cursors. Style 4 is applied to cursors set by num%; it does not set the view style.		
	num% A value of -1 or omitting the argument selects all cursors and sets the view style for new cursors, 0-9 selects one cursor. In version 3, <=0 selected all cursors.		
	form\$ Cursor label string with replaceable fields %p, %n and %v(chan) for position, number and channel value; chan is the channel number whose value you require. %w.dp and %w.dv(chan) formats are allowed where w and d are numbers that set the field width and number of decimal places.		
	Returns The view cursor style as 0-3 before any change. If style% is omitted or not 0-3, the current view cursor style is not changed.		

See also:Cursor(), CursorLabelPos(), CursorNew(), CursorRenumber()

CursorLabelPos()	This lets you set and read the position of the cursor label.
	<pre>Func CursorLabelPos(num% {,pos});</pre>
	num% The cursor number. Setting a silly number does nothing and returns -1.
	pos If present, the command sets the label position as the percentage of the distance from the top of the cursor. Out-of-range values are set to the appropriate limit.
	Returns The cursor position before any change was made.
	See also:Cursor(),CursorLabel(),CursorNew(),CursorRenumber()
CursorNew()	This command adds a new cursor to the view at the designated position. You cannot use this to create cursor 0 in a time view as this cursor always exists. To show cursor 0 use CursorVisible(0,1). A new cursor is created in Static mode (not active).
	<pre>Func CursorNew({where{, num%}});</pre>
	where The cursor position. In a time view it is a time in seconds, in a result view, it is the bin number. Use XToBin() to convert x axis units to bin numbers. In an XY view it is in x axis units. The position is limited to the x axis range. If the position is omitted, the cursor is placed in the middle of the window.
	num% If this is omitted, or set to -1, the lowest-numbered free cursor is used. If this is a cursor number, that cursor is created. This must be a legal cursor number or -1.
	Returns It returns the cursor number as an integer, or 0 if all cursors are in use.
	See also:Cursor(), CursorActive(), CursorDelete(), CursorLabel(),
CursorRenumber()	This command renumbers the cursors from left to right in the view. It has no effect on cursor 0 . Active cursor and label information stays with the cursors, not the number
	Func CursorRenumber():
	Returns. The number of cursors that were renumbered (cursor 0 is not counted)
	See also:Cursor(), CursorActive(), CursorDelete(), CursorLabel(),
	CursorLabelPos(), CursorNew(), CursorSet()
CursorSearch()	This function causes active cursors in a time view to search according to the current cursor mode. You can cause all cursors to search, or a restricted range of cursor numbers. Moving cursor 0 with Cursor(0, new) also causes all cursors 1-9 to search.
	<pre>Func CursorSearch(num% {,stop%});</pre>
	num% This is the first cursor number to run the search defined by the active cursor mode. Set this to 0 to cause cursor 0 to search forwards and to -1 for cursor 0 to search backwards. CursorSearch(0) and CursorSearch(-1) are equivalent to the Ctrl+Shift+Right and Ctrl+Shift+Left key combinations.
	stop% This optional argument sets number of the last cursor to try to reposition. If you omit this argument, all cursors from num% upward will search according to their active mode. To reposition a single cursor set stop% the same as num%.
	Returns The time that cursor num% moved to or -1 if didn't move. You can also use CursorValid() to test if searches have succeeded.
	<pre>See also:Cursor(), CursorActive, CursorActiveGet(), CursorNew(), CursorValid(), MeasureChan(), MeasureX()</pre>

CursorSet()	This sets the number of vertical cursors in addition to cursor 0. It deletes cursors 1-9, then positions num% cursors equally spaced in the view, numbered in order from left to right. Finally, if any cursors positions are given, they are applied. The cursor labelling style is not changed. This destroys all active cursor information for the deleted cursors. Proc CursorSet(num% {,where1 {,where2 {,where3 {,where4}}}); num% The number of cursors to display in the range 0 to 9. It is a run-time error to ask for more than 9 or less than 0 cursors. If num% is 0, cursor 0 is hidden. whereN Optional positions of cursor N (1 to 9). Positions that are out of range are set to the nearest valid position. In a time view positions are in seconds. In a result view, they are the bin number; use XToBin() to convert x axis units to bin numbers. In an XY view the position is in x axis units. You cannot change the position of cursor 0 in a time view with this command.
	Examples:
	CursorSet(0); 'Delete 1-9, hide cursor 0 CursorSet(2,20,30); 'Delete 1-9, cursor 1 at 20, cursor 2 at 30
	<pre>See also:BinToX(), Cursor(), CursorNew(), CursorRenumber(), CursorVisible(), XToBin()</pre>
CursorValid()	Use this function to test if the last search of a cursor in a time view succeeded. Cursor positions are valid if a search succeeds or if the cursor is positioned manually or by a script command. The position of a newly created cursor is valid.
	Func CursorValid(num%);
	num% The cursor number to test for a valid search result in the range 0-9.
	Returns The result is 1 if the position of the nominated cursor is valid or 0 if it is invalid or the cursor does not exist.
	<pre>See also:CursorActive, CursorActiveGet(), CursorNew(), CursorSearch(), CursorVisible(), MeasureChan(), MeasureX()</pre>
CursorVisible()	Vertical cursors can be hidden without deleting them. Interactively you can hide cursor 0, but from a script you can show and hide any vertical cursor. Cursors are always made visible by the Ctrl+n key combination.
	<pre>Func CursorVisible(num% {,show%});</pre>
	num% The cursor number in the range 0-9 or -1 for all vertical cursors.
	show% If present set this to 0 to hide the cursor and non-zero to show it.
	Returns The state of the cursor at the time of the call (0=hidden, 1=visible) or -1 if the cursor does not exist. If num% is -1, the result is the number of vertical cursors. See also:CursorExists(), CursorNew(), CursorSearch(), CursorValid()

This function returns a string holding the date. Use TimeDate() to get the date as Date\$(numbers. For this description, we assume that today's date is Wednesday 1 April 1998, the system language is English and the system date separator is "/". Default argument values are shown **bold**. Func Date\$({dayF%, {monF%, {yearF%, {order%, {sep\$}}}}); dayF% This sets the format of the day field in the date. This can be written as a day of the week or the day number in the month, or both. The options are: Show day of week: "Wednesday". 1 2 Show the number of the day in the month with leading zeros: "01". 4 Show the day without leading zeros: "1". This overrides option 2. 8 Show abbreviated day of week: "Wed". 16 Show weekday name first, regardless of the order% field. Use 0 for no day field. Add the numbers for multiple options. For example, to return "Wed 01", use 11 (1+2+8) as the dayF% argument. If you add 8 or 16, 1 is added automatically. If you request both the weekday name and the number of the day, the name appears before the number. The format of the month field. This can be returned as either a name or a monF% number. If this argument is omitted, the value 3 is used. The options are: 0 No month field. 1 Show name of the month: "April". 2 Show number of month: "04" 3 Show an abbreviated name of month: "Apr" Show number of month with no leading zeros: "4" 4 yearF% The format of the year field. This can be returned as a two or four digit year. 0 No year is shown Year is shown in two digits: "98". 1 2 Year is shown in two digits with an apostrophe before it: "' 98". 3 Year is shown in four digits: "1998". order% The order that the day, month and year appear in the string. 0 Operating system settings 1 month/day/year 2 day/month/year 3 year/month/day sep\$ This string appears between the day, month and year fields as a separator. If this string is empty or omitted, Spike2 supplies a separator based on system settings. For example, Date\$(20, 1, 2, 1, " ") returns "Wednesday April 1 '98". As 20 is 16+4, we have the day first, even through the order% argument places the day in between the month and the year. Date\$() returns "01/Apr/98". See also:Seconds(), FileDate\$(), TimeDate(), Time\$()

Debug()

This command has two functions. It can open the debug window so you can step through your script, set breakpoints and display and edit variables. From version 3 it can be used to stop the user entering the debugger with the Esc key.

Proc Debug({msg\$}|{Esc%});

msg\$ When the command is used with no arguments, or with a string argument, the script stops as though the Esc key had been pressed and enters the debugger. If the debugging toolbar was hidden, it becomes visible. If the msg\$ string is present, the string is displayed in the bar at the top of the script window.

Esc[%] When the command is used with an integer argument, it enables and disables the ability of the user to break out of a running script. If Esc[%] is 0, the user cannot break out of a script into the debugger with the Esc key and must wait for it to finish. If Esc[%] is 1, the user can break out. Spike2 enables the Esc key each time a script starts, so make this the very first instruction of your script if you want to be certain that the user cannot break out.

This command was included for use in situations such as student use, where it is important that the user cannot break out of a script by accident. It is advisable to test your script carefully before using this option. Once set, you cannot stop a looping script except by forcing a fatal error. Make sure you save your script before setting this option.

See also:Eval()

De	bua	List()

This command is used for debugging problems in the system. It writes information to the Log view about the internal list of "objects" used to implement the script language.

Proc DebugList(list% {, opt%});

- list% This determines what to list. A value of 0 lists a summary of the options, 1 lists fixed objects (constants and operators), 2 lists permanent objects (constants, operators and built in commands). Values greater than 2 list information for the object with that number.
- opt% This optional argument (default value 0) sets the additional object information to list. 1= list the index number, 2= list the type, 3= list both index and type.

See also:Debug(), Eval(), DebugOpts()

DebugOpts()

This command is used for debugging problems in the system. It controls internal options used for debugging at the system level.

Func DebugOpts(opt% {,val%});

- opt% This selects the option to return (and optionally to change). A value of 0 prints a synopsis of available options to the Log view and the current value of each option. Values greater than 0 return the value of that option, and print the option information to the Log view. At the time of writing, only option 1, dump compiled script to the file default.cod is implemented.
- Val% If present, this sets the new value of the option.

See also:Debug(), Eval(), DebugList()

DelStr\$()	This function removes a sub-string from a string.
• 0	<pre>Func DelStr\$(text\$, index%, count%);</pre>
	text\$ The string to remove characters from. This string is not changed.
	index% The start point for the deletion. The first character is index 1. If this is greater than the length of the string, no characters are deleted.
	count% The number of characters to remove. If this would extend beyond the end of the string, the remainder of the string is removed.
	Returns The original string with the indicated section deleted.
	<pre>See also:Asc(), Chr\$(), InStr(), LCase\$(), Left\$(), Len(), Mid\$(), Print\$(), ReadStr(), Right\$(), Str\$(), UCase\$(), Val()</pre>

Discriminator (CED The Discrim... family of commands supports the 1401-18 event discriminator card, **1401-18)** support which is available for the standard 1401 and the 1401*plus* only. You can also control the 1401-18 card interactively; see the Sample menu.

> If you use these commands with the interactive discriminator dialog active, the dialog changes to show any changes made from the script. You cannot change the current channel the dialog displays (except with DiscrimClear()), so you will only see changes in the dialog if the channel is the same as the current channel in the dialog.

```
See also:DiscrimChanGet(), DiscrimChanSet(), DiscrimClear(),
       DiscrimLevel(), DiscrimMode(), DiscrimMonitor(),
       DiscrimTimeOut()
```

DiscrimChanGet()

This gets the input, spike2 event port, 1401 event input E1 or E3, mode, lower and higher trigger levels, and the time-out for modes 7 and 8 of a discriminator channel.

Func DiscrimChanGet (chan%, &in%, &out%, &mode%, &low, &hi, &tOut);

- chan% the discriminator channel number (0-7).
- in% Returned holding the source of the discriminator channel as:
 - 0 The discriminator is not used
 - 1 From the 1401 front panel digital input port
 - 2 From the front panel event input (channels 0-4) or ADC Ext (channel 5)
- out% Returned holding a code that describes how the discriminator output, event inputs and digital inputs are connected to the Spike2 event ports and the E1 (digital marker trigger) and E3 (start sampling trigger). The E1 and E3 values are only valid for channels 1 and 3.

	Spike2 event port	E1 (digital marker) or E3 (start sampling)
0	Disconnected.	Disconnected
1	1401 digital input	Disconnected
2	Discriminator output	Disconnected
4	Disconnected	To front panel E1 or E3
5	1401 digital input	To E1 or E3
6	Discriminator output	To E1 or E3
8	Disconnected	Discriminator output
9	1401 digital input	Discriminator output
10	Discriminator output	Discriminator output
	-	-

the mode of the discriminator channel is returned in this variable: mode%

- detect a level below a threshold 1
- 2 detect a level above a threshold
- pulse on rising signal through a threshold 3
- 4 pulse on falling signal through a threshold
- 5 detect a level inside a region set by two thresholds
- detect a level outside region set by two thresholds 6
- 7 pulse when signal returns below lower threshold, without going above the higher threshold, within a programmable time out period
- 8 pulse when signal returns above higher threshold, without going below the lower threshold, within a programmable time out period
- low Returned holding the lower threshold level in Volts.
- Returned holding the higher threshold level in Volts. high
- Returns the programmable time out period used in modes 7 and 8. tOut

Returns 0 if all well or a negative error code.

See also:DiscrimChanSet(), DiscrimClear(), DiscrimLevel(), DiscrimMode(), DiscrimMonitor(), DiscrimTimeOut()

DiscrimChanSet()

This sets the input, output, mode, lower and higher trigger levels, and the time-out (if the mode is 7 or 8) of the discriminator channel.

Func DiscrimChanSet(chan%, in%, out%{, mode%, low, high, tOut});

- chan% the discriminator channel number (0-7).
- in% Sets the source of the discriminator channel as:
 - 0 Disconnected
 - 1 From the 1401 front panel digital input port
 - 2 From the front panel event input (channels 0-4) or ADC Ext (channel 5). For discriminator channels 6 and 7, you cannot set in% to 2 as there is no 1401 event input available as the source.
- out% Sets how the discriminator output, event inputs and digital inputs connect to the Spike2 event ports and the E1 (digital marker) and E3 (start sampling) triggers. If in% is 0, out% is forced to 1. If both in% and out% are 1, out% is forced to 2. Possible out% values are:
 - 0 Spike2 event port disconnected
 - 1 Spike2 event port connected to 1401 digital input
 - 2 Spike2 event port connected to discriminator output

When the channel is 1 or 3 you can choose what the 1401 event input E1 (digital marker) or E3 (start sampling) trigger connects to by adding 0, 4 or 8 to out%. If in% is 0, adding 4 is forced. If in% is 2, adding 4 is treated as adding 8.

- 0 E1 or E3 is disconnected
- 4 E1 or E3 is connected to the front panel E1 or E3
- 8 E1 or E3 is connected to the discriminator output
- mode% The discriminator channel mode (1-8). If omitted, no change is made.
 - 1 detect a level below a threshold
 - 2 detect a level above a threshold
 - 3 pulse on rising signal through a threshold
 - 4 pulse on falling signal through a threshold
 - 5 detect a level inside a region set by two thresholds
 - 6 detect a level outside region set by two thresholds
 - 7 pulse when signal returns below lower threshold, without going above the higher threshold, within a programmable time out period
 - 8 pulse when signal returns above higher threshold, without going below the lower threshold, within a programmable time out period
- Low Sets the lower threshold level in volts between -5 Volts and the higher threshold level value. If omitted, no change is made. If the lower threshold level is given greater than the higher one, it will be set equal to the higher level.
- high Sets the higher threshold level in Volts between the lower threshold level value and +5 Volts. If omitted, no change is made.
- tout Sets the time out period in seconds (0.00002 0.65535). This is only meaningful in mode 7 or 8. If omitted, no change is made.

Spike2 stores and uses the values for the threshold levels and time out period in terms of the resolution of the 1401-18 card. To find out the real value of the threshold levels and time out, use DiscrimLevel() and DiscrimTimeOut(). The difference is usually very small.

Returns 0 if the operation completed without a problem, or a negative error code.

See also:DiscrimChanGet(), DiscrimClear(), DiscrimLevel(), DiscrimMode(), DiscrimMonitor(), DiscrimTimeOut()

DiscrimClear()

DiscrimLevel(

This sets the Discriminator configuration dialogue contents to a standard state. All the discriminator channels are disconnected. Spike2 event ports are connected to digital inputs. E1 (digital marker trigger) and E3 (start sampling trigger) are connected to the front panel E1 and E3. All the discriminator channels have mode 3, and the lower and higher trigger levels are set to 1.25 and 2.5 Volts. The ADC monitor channel is set to 15. If the discriminator dialog is open, channel 0 becomes the current channel.

Proc DiscrimClear();

```
See also:DiscrimChanGet(), DiscrimChanSet(), DiscrimLevel(),
DiscrimMode(), DiscrimMonitor(), DiscrimTimeOut()
```

This sets or gets the lower or higher threshold level values of the discriminator channel.

Func DiscrimLevel(chan%, which%, {level});

- chan% The discriminator channel number (0-7).
- which % Selects which level to set or get: 0 =lower, 1 =upper threshold level.
- level If present it sets the threshold value in Volts:

If which% is 0, it is the lower level. It should lie between -5 volts and the higher level. If it is set greater than the higher level, it is set equal to the higher one.

If which% is 1, it is the higher level. It should lie between the lower level and 5 Volts. If it is set less than the lower level, it is set equal to the lower one

- Returns the lower or higher threshold level of the discriminator channel at the time of call, or a negative error code.
- See also:DiscrimChanGet(), DiscrimChanSet(), DiscrimClear(), DiscrimMode(), DiscrimMonitor(), DiscrimTimeOut()

DiscrimMode()	This sets or gets the mode of the discriminator channel.	
	<pre>Func DiscrimMode(chan%, {mode%});</pre>	
	chan% The discriminator channel number (0-7).	
	The mode of the discriminator channel (1-8). If omitted, the current mode of the discriminator channel is returned.	he
	Returns The discriminator channel mode at the time of call, or a negative error code.	
	 detect a level below a threshold detect a level above a threshold 	
	3 pulse on rising signal through a threshold	
	 4 pulse on falling signal through a threshold 5 detect a level inside a region set by two thresholds 	
	6 detect a level outside region set by two thresholds	
	7 pulse when signal returns below lower threshold, without going above the	

higher threshold, within a programmable time out period
pulse when signal returns above higher threshold, without going below the
lower threshold, within a programmable time out period

See also:DiscrimChanGet(),DiscrimChanSet(),DiscrimClear(),

DiscrimLevel(), DiscrimMonitor(), DiscrimTimeOut()

DiscrimMonitor()	This sets or gets the ADC monitor channel in the Discriminator configuration dialog.
	<pre>chan% Sets the ADC monitor channel number (0-15). If this is omitted, the current ADC monitor channel number is returned. Returns the ADC monitor channel number at the time of call, or a negative error code. See also:DiscrimChanGet(), DiscrimChanSet(), DiscrimClear(), DiscrimLevel(), DiscrimMode(), DiscrimTimeOut()</pre>
DiscrimTimeOut()	This sets or gets the time out period of mode 7 or 8 for the discriminator channel.
	$abar \delta_{abar}$ the discriminator channel number (0.7)
	tout the programmable time out period in seconds (0.00002 - 0.65535). This is only meaningful in mode 7 or 8. If this is omitted, the current time out period is returned.
	Returns the time out period at the time of call, or a negative error code.
	<pre>See also:DiscrimChanGet(), DiscrimChanSet(), DiscrimClear(), DiscrimLevel(), DiscrimMode(), DiscrimMonitor()</pre>
Dialogs	You can define your own dialogs to get information from the user. You can define dialogs in a simple way, where each item of information has a prompt, and the dialog is laid out automatically, or you can build a dialog by specifying the position of every item. A simple dialog has the structure shown in the diagram:
	The dialog is arranged in terms of items. Unless you specifically request otherwise, the dialog items are stacked vertically above each other with buttons arranged at the bottom.
	The dialog has a title of your choosing at the top and OK and Cancel plus user-defined buttons at the bottom. When the dialog is used, pressing the Enter key is equivalent to clicking on OK.
	This form of dialog is very easy to program; there is no need to specify any position or size information, the system works it all out. Some users require more complicated dialogs, with control over the field positions. This is also possible, but harder to program. You are allowed up to 1000 fields in a dialog (was 256 in version 4).
	In more complex cases, you specify the position (and usually the width) of the box used for user input. This allows you to arrange data items in the dialog in any way you choose. It requires more work as you must calculate the positions of all the items.
Version 5 extensions	There are new script functions to add and manipulate buttons, to collect a time or x axis value and to add a group box. You can define script functions that are called in response to button presses and user changes to the dialog and an idle-time function that is called repeatedly whilst the dialog is waiting for user actions. All these functions can enable and disable, hide and show and modify dialog items. DlgChan() can now get a channel from an XY view. You can add tooltips for all the prompts in a dialog to give users a more detailed explanation of a field. Finally, there are extensions to the integer, real and string fields that allow you to define a drop-down list of selectable items to copy into the fields

and you can add a spin control to both integer and real numeric fields. All dialogs created in previous versions of Spike2 should work without any change.

Dialog units Positions within a dialog are set in *dialog units*. In the x (horizontal) direction, these are in multiples of the maximum width of the characters '0' to '9'. In the y (vertical) direction, these are in multiples of the line spacing used for simple dialogs. The maximum y position is 40. Unless you intend to produce complex dialogs with user-defined positions, you need not be concerned with dialog units at all.

Simple dialog example The simple example dialog shown above can be created by this code:

var ok%, item1%, item2%, item3, item4\$:= "more...", item5; DlgCreate ("Title for the dialog"); 'start new dialog DlgInteger(1,"Item &1 prompt",0,10,0,0,1); 'range 0-10,spinner DlgChan (2, "Item &2 prompt|Tip", 1); 'Waveform channel list DlgReal (3, "Item &3 prompt", 1.0, 5.0); 'real, range 1.0-5.0 DlgString (4, "&More prompts...|Tip", 6); 'string, any characters DlgCheck (5, "A &checkbox with a prompt");'a checkbox item DlgButton (2, "&Extra||Tooltip"); 'button 2, tooltip ok% := DlgShow(item1%, item2%, item3, item4\$, item5); 'show dialog

Prompts, & and tooltips In the functions that set an item with a prompt, if you precede a character in the prompt with an ampersand (&), the following character is underlined and is used by Windows as a short-cut key to move to the field or activate the button. All static dialog items except a group box allow you define a tooltip by appending a vertical bar followed by the tooltip text to the text\$ argument. For example:

DlgReal(3, "Rate|Enter the sample rate in Hz", 100, 500);

Buttons allow you to specify an additional activation key and an optional tooltip by adding a vertical bar followed by the key code and then another vertical bar followed by the tooltip. See the label\$ argument of ToolbarSet() for details of key codes.

More complex example This example shows how to respond to user actions within a dialog. In this case we use a checkbox to enable and disable a group of items and a button that displays the current values of dialog items. The numbered fields are:

- 1 An integer, range 0-10 with a spinner
- 2 A drop list of 4 items
- 3 A checkbox, used to enable items 4 and 5
- 4 A real number with a spinner
- 5 A string with a drop down list of items

We have added button 2 (buttons 0 and 1 are

Dialog with user p	lacement	×
Integer 0 to 10	0	-
List item	List 0	-
🔲 Checkbox enab	oling items	
Extra items Real 1 to 5	1	* *
String length 10	Text	-
Current	Cancel OK	

Cancel and OK) and a group box around items 4 and 5. To make room for the group box, the y positions of items 4 and 5 are set explicitly.

With DlgAllow() we have set Func Change%(item%) to be called whenever the user changes a selection or checkbox or when an editable field loses the input focus. The item% argument is set to the item number that changed or to 0 if the dialog is appearing for the first time. We are interested in item 3, the checkbox, and we use the state to enable or disable the group box and the items inside it.

Func Current%() is linked to the "Current" button and can be activated with the C and F1 keys and has a tooltip. The function displays a message box that lists the current values of items in the dialog.

DIgAllow()

```
var ok%, item1%, item2%, item3%, item4, item5$:= "Text", gp%;
DlgCreate("Dialog with user placement",0,0,40,7.5);
DlgInteger(1,"Integer 0 to 10",0,10,0,1,1);
                                                'Int with spinner
DlgList(2,"List item","List 0|List 1|List 2|List 3", 4, 0, 2);
DlgCheck(3, "Checkbox enabling items",0,3);
                                                'checkbox item
DlgReal(4, "Real 1 to 5",1.0,5.0,0,4.5,0.5);
                                                'Real with spinner
DlgString(5, "String length 10",10,"",0,5.5,
                                                'String item with
             "String 1|String 2|String 3");
                                                ' drop-down list
DlgButton(2,"&Current|0x70|Tooltip", Current%);'button F1+function
DlgAllow(0x3ff, 0, Change%); 'Allow all, no idle, change function
gp% := DlgGroup("Extra items",1,3.8,-1,2.9);
                                                'Group box
ok% := DlgShow(item1%,item2%,item3%,item4,item5$);
Halt;
Func Change%(item%)
var v%;
docase
  case ((item% = 3) or (item% = 0)) then '0 is initial setup
    v% := DlgValue(3);
                               'get checkbox state
    DlgEnable(v%, gp%, 4, 5); 'enable groupbox+items 4, 5
endcase;
                               'Return 1 to keep dialog running
return 1;
end;
Func Current%()
var v1%, v2%, v3%, v4, v5$;
v1% := DlgValue(1);
                               'Retrieve the current values
v2% := DlgValue(2); v3% := DlgValue(3);
v4 := DlgValue(4); v5$ := DlgValue$(5);
Message("Values are %d, %d, %d, %g and %s",v1%,v2%,v3%,v4,v5$);
return 1;
                               'Return 1 to keep the dialog running
end:
See also:DlgAllow(), DlgButton(), DlgChan(), DlgCheck(), DlgCreate(),
       DlgEnable(), DlgGroup(), DlgInteger(), DlgLabel(), DlgList(),
       DlgReal(), DlgShow(), DlgString(), DlgText(), DlgValue(),
       DlgVisible(), DlgXValue(), ToolbarSet()
```

DIgAllow()

Call this function after DlgCreate() and before DlgShow() to enable dialog idle time processing, advanced call-back features and dynamic access to the dialog fields. There are no restrictions on what call-back functions can do. However, it is not sensible to place time-consuming code in an idle call-back function or to do anything other than check dialog fields and possible display a warning message in a dialog-item-change function. Call-back functions use DlgValue(), DlgEnable() and DlgVisible() to manipulate the dialog fields. DlgAllow() is new in version 5.

Proc DlgAllow(allow% {,func id%() {, func ch%() }});

- allow% A number that specifies the actions that the user can and cannot take while interacting with Spike2. See Interact() for a full description.
- id%() This is an integer function with no arguments. Use the name with no brackets, for example DlgAllow(0,Idle%); where Func Idle%() is a script function. When DlgShow() executes, the function is called repeatedly in system idle time, as for the ToolbarSet() idle function.

If the function return value is greater than 0, the dialog remains open. A zero or negative return value closes the dialog and DlgShow() returns the same value.

If this argument is omitted or 0, there is no idle time function.

ch%() This is an integer function with one integer argument, for example Func Changed%(item%). You would use DlgAllow(0,0,Changed%); to link this

function to a dialog. Each time the user changes a dialog item, Spike2 calls the function with the argument set to the changed item number. There is an initial call with the argument set to 0 when the dialog is about to be displayed.

A field is deemed to change when the user clicks a checkbox or changes a selection in a list or moves the focus from an editable item after changing the text. For real and integer values, the new value must be in range.

If the change function returns greater than 0, the change is accepted. If the return value is zero, the change is resisted and the focus set back to the changed item. If the return value is negative, the dialog closes and DlgShow() returns this value and the arguments are not updated.

See also:DlgCreate(), DlgEnable(), DlgShow(), DlgValue(), DlgVisible(), Interact(), ToolbarSet()

DlgButton()

Dialogs created by DlgCreate() have Cancel and OK buttons. This function adds, deletes and changes buttons. You can link a script function to a button and use the function return value to decide if the dialog should close. Use this function after DlgCreate() and before DlgShow(). DlgButton() is new in version 5.

Proc DlgButton(but%, text\${, func ff% {, x, y}});

- but% The button number from 0 to 200. Button 0 is the cancel button, 1 is the OK button. Button numbers higher than 1 create new buttons.
- text\$ This sets the button label. Set an empty string to delete a button. You cannot delete button 1; the label is set back to OK if you try. The label text can be followed by an optional key code and an optional tooltip separated by vertical bars. See the label\$ argument of ToolbarSet() for details of the format.

If you set a key code, the button can be activated even when the dialog does not have the input focus as long as it is the topmost user dialog and you have not created a toolbar or interact bar from a function linked to the dialog. This allows you to drag cursors in a window, then use the key code without the need to click in the dialog to activate it first.

ff%() This is an integer function with no arguments that is called when the button is used. Set the argument to zero or omit it if you don't want a button function, in which case clicking the button closes the dialog, DlgShow() returns the button number and the DlgShow() arguments are updated for all buttons except 0.

If you supply a function, it is called each time the button is used and the function return value determines what happens next:

- <0 The button acts as Cancel. The dialog closes, DlgShow() returns this value and its arguments are not updated.</p>
- 0 The button acts as OK. The dialog closes and the DlgShow() return value is the button number and its arguments are updated.
- >0 The dialog continues to display.

The button function can use DlgEnable(), DlgValue() and DlgVisible().

x, y Set the button position, both or neither of these must be supplied. If the button position is not supplied it will be positioned at the bottom of the dialog.

See also:DlgEnable(), DlgShow(), DlgValue(), DlgVisible(), ToolbarSet()

DlgChan()

You often need to select a channel from a time, result or XY view. This function defines a dialog entry that lists channels that meet a specification. For simple dialogs, the wide, x and y arguments are not used. Channel lists are checked or created when the DlgShow() function runs. If the current view is not a time, result or XY view, the list is empty.

Proc DlgChan(item%, text\$|wide, mask%|list%[]{, x{, y}});

- item% This sets the item number in the dialog in the range 1 to the number of items.
- text\$ The prompt to display, optionally followed by a vertical bar and tooltip text.
- wide This is an alternative to the prompt. It sets the width in dialog units of the box used to select a channel. If wide is omitted the number entry box has a default width of the longest channel name in the list or 12, whichever is the smaller.
- mask% This determines the channels to display. Select channel types by adding the following codes, which are given as decimal and hexadecimal.
 - 1 0x1 Waveform or result view channel
 - 2 0x2 Event+ and Event- channels
 - 4 0x4 Event +- channels (level data)
 - 8 0x8 Marker channels
 - 16 0x10 WaveMark data
 - 32 0x20 TextMark data
 - 64 0x40 RealMark data
 - 128 0x80 Unused/deleted disk channels
 - 256 0x100 Deleted channels on disk
 - 512 0x200 Real wave channel

If none of the above values are used or this is an XY view, the list includes all channels. Add the following codes to exclude channels from the list:

1024	0x400	Exclude visible channels
2048	0x800	Exclude hidden channels
4096	0x1000	Exclude time view disk channels but not their duplicates
8192	0x2000	Exclude memory channels but not their duplicates
16384	0x4000	Exclude duplicated channels
32768	0x8000	Exclude selected channels
65536	0x10000	Exclude non-selected channels
2097152	0x200000	Exclude virtual channels but not their duplicates

Finally, the following codes add special entries to the list:

131072 0x20000 Add None as an entry in the list, returns 0

262144	0x40000	Add All	channels	s as an entry in th	e list, returns -1
524288	0x80000	$Add {\tt All}$	visible	channels as an	entry, returns -2
1048576	0x100000	Add Sele	ected as a	n entry in the list,	returns -3

- list% As an alternative to a mask, you can pass in a channel list (as constructed by ChanList()). This must be an array of channel numbers, with the first element of the array holding the number of channels in the list.
- x If omitted or zero, the selection box is right justified in the dialog box, otherwise this sets the position of the left end of the channel selection box in dialog units.
- Y If omitted or zero, this takes the value of item%. It is the position of the bottom of the channel selection box in dialog units.

The variable passed to DlgShow() for this field should be an integer. If the variable passed in holds a channel number in the list, the field shows that channel, otherwise it shows the first channel in the list (usually None). The result from this field in DlgShow() is a channel number, or 0 if None is selected, -1 if All channels is selected, -2 if All visible channels is selected or -3 if Selected is chosen.

See also:DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(), DlgList(), DlgReal(), DlgShow(), DlgString(), DlgText(), DlgXValue()

DlgCheck()	This defines a dialog item that is a check box (on the left) with a text string to its right. For simple dialogs, the x and y arguments are not used.		
	Proc I	<pre>DlgCheck(item%, text\${, x{, y}});</pre>	
	item%	This sets the item number in the dialog in the range 1 to the number of items.	
	text\$	The prompt to display, optionally followed by a vertical bar and tooltip text.	
	х,у	The position of the bottom left hand corner of the check box in dialog units. If omitted, x is set to 2 and y to $item$ %. When used without these fields, this behaves exactly like the simple dialog functions, and can be mixed with them.	
	The associated DlgShow() variable should be an integer. It sets the initial state (0 for unchecked, not 0 for checked) and returns the result as 0 (unchecked) or 1 (checked). This item does not have a prompt. If you use DlgValue(), a string value refers to the text and a numeric value refers to the check box state.		
	See also	D:DlgChan(),DlgCreate(),DlgInteger(),DlgLabel(),DlgList(), DlgReal(),DlgShow(),DlgString(),DlgText(),DlgXValue()	
	This Ca		
DIgCreate()	kills of	f any previous dialog that might be partially defined.	
	Func I	<pre>DlgCreate(title\${, x{, y{, wide{, high{, help% help\$}}}});</pre>	
	title	A string holding the title for the dialog.	
	х,у	Optional, taken as 0 if omitted. The position of the top left hand corner of the dialog as a percentage of the screen size. The value 0 means centre the dialog. Values out of the range 0 to 95 are limited to the range 0 to 95.	
	wide	The width of the dialog in dialog units. If this is omitted, or set to 0, Spike2 works out the width for itself, based on the items in the dialog.	
	high The height of the dialog in dialog units. If omitted, or set to 0, Spike2 out for itself, based on the dialog contents.		
	help	This is a string or numeric identifier that identifies the help page to be displayed if the user requests help when the dialog is displayed.	
	Returns	This function returns 0 if all was well, or a negative error code.	
	For simple use, only the first argument is needed. The remainder are for use with more complicated menus where precise control over menu items is required.		
	See also	D:DlgButton(),DlgChan(),DlgCheck(),DlgGroup(),DlgInteger(), DlgLabel(),DlgList(),DlgReal(),DlgShow(),DlgString(), DlgText(),DlgXValue()	
DlgEnable()	Use this argume enabled are also	s only from a dialog call-back function to enable or disable dialog items. With one nt, it returns the enabled state of an item; with two or more arguments it sets the l state of one or more items. Prompts and spin controls associated with the item o enabled or disabled. DlgEnable() is new in version 5.	
	Func I Func I	<pre>DlgEnable(en%, item% item%[]{, item% item%[]}); DlgEnable(item%);</pre>	
	en%	Set 0 to disable list items, 1 to enable them and 2 to enable and give the first item the input focus. Input focus changes should be used sparingly to avoid user confusion; they can cause button clicks to be missed.	
	item%	An item number or an array of item numbers of dialog elements. The item number is either the number you set, or the number returned by DlgText() or	

DlgGroup(), or -button, where button is the button number. You cannot access prompts separately from their items as this makes no sense.

Returns When called with a single argument it returns the enabled state of the item, otherwise it returns 0.

See also:DlgAllow(), DlgCreate(), DlgShow(), DlgValue(), DlgVisible()

DlgGroup()

This routine creates a group box, which is a rectangular frame with a text label at the top left corner. You can use this between calls to DlgCreate() and DlgShow(). There is nothing for the user to edit in this item, so you do not supply an item number and there is no matching argument in DlgShow(). However, the returned number is an item number (above the values used to match items to DlgShow() arguments) that you can use in callback functions to identify the group box. This function was added at version 5.

Func DlgGroup(text\$, x, y, width, height);

- text\$ The text to display at the top left of the group box. Any tooltip text is ignored.
- x, y The position of the top left corner of the group box.
- width If positive, the width of the group box. If negative, this is the offset of the right hand side of the group box from the right hand edge of the dialog.
- height The height of the group box.
- Returns The routine returns an item number so that you can refer to this in call-back functions to use DlgVisible() and DlgEnable().

See also:DlgCreate(), DlgEnable(), DlgShow(), DlgVisible()

Diginteger()

This function defines a dialog entry that edits an integer with an optional spin control or drop down list of selectable items. The numbers you enter may not contain a decimal point. For simple dialogs, the wide, x, y, sp% and li\$ arguments are not used. The sp and li\$ arguments are new in version 5.

Proc DlgInteger(item%, text\$|wide, lo%, hi%{, x{, y{, sp%|li\$}});

- item% This sets the item number in the dialog in the range 1 to the number of items.
- text\$ The prompt to display, optionally followed by a vertical bar and tooltip text.
- wide This is an alternative to the prompt. It sets the width in dialog units of the box in which the user types the integer. If the width is not given the number entry box has a default width of 11 digits (or width needed for the number range?).
- 10% The start of the range of acceptable numbers.
- hi% The end of the range of acceptable numbers.
- x If omitted or zero, the number entry box is right justified in the dialog box, otherwise this sets the position of the left end of the box in dialog units.
- y If omitted or zero, this takes the value of item%. It is the position of the bottom of the number entry box in dialog units.
- sp% If present and non-zero, this adds a spin box with a click increment of sp%.
- 11\$ If present, this argument is a list of items separated by vertical bars that can be selected into the integer field, for example "1|10|100".

The variable passed into DlgShow() should be an integer. The field starts with the value of the variable if it is in the range. Otherwise, it is limited to the nearer end of the range.

See also:DlgChan(), DlgCheck(), DlgCreate(), DlgLabel(), DlgList(), DlgReal(), DlgShow(), DlgString(), DlgText(), DlgXValue() DlgList()

DlgLabel()	This function sets an item with no editable part that is used as a label. For simple dialogs, the wide, x and y arguments are not used. You can add text to a dialog without using an item number with $DlgText()$.				
	<pre>Proc DlgLabel(item%, text\${, x{, y}});</pre>				
	item%	This sets the item number in the dialog in the range 1 to the number of items.			
	text\$	The prompt to display, optionally followed by a vertical bar and tooltip text.			
	х	If omitted, the text is left justified in the dialog box. Otherwise, this sets the position of the left end of the text in the dialog in dialog units.			
	У	If omitted or zero, this takes the value of item%. It is the position of the bottom of the text in the dialog in dialog units.			
	When y variable	rou call DlgShow(), you must provide a dummy variable for this field. The is not changed and can be of any type, but must be present.			
	See also	<pre>CDlgChan(),DlgCheck(),DlgCreate(),DlgInteger(),DlgList(), DlgReal(),DlgShow(),DlgString(),DlgText(),DlgXValue()</pre>			

This defines a dialog item for a one of n selection. Each of the possible items to select is identified by a string. For simple dialogs, the wide, x and y arguments are not used.

Proc DlgList(item%, text\$|wide, list\$|list\$[]{, n%{, x{, y}}});

- item% This sets the item number in the dialog in the range 1 to the number of items.
- text\$ The prompt to display, optionally followed by a vertical bar and tooltip text. The list box width is the smaller of the longest string in the list or 20.
- wide This is an alternative to the prompt. It sets the width of the box in which the user selects an item. If you need a prompt you can use DlgText().
- list\$ A string of list items separated by a "|" character or an array of strings, one per item. Long strings are truncated. The "|" method was new in version 3.
- n% The number of entries to display. If this is omitted, or if it is larger than the number of entries provided, then all the entries are displayed.
- x If omitted or zero, the selection box is right justified in the dialog, otherwise this sets the position of the left end of the selection box.
- y The bottom of the selection box position. If omitted, the value of item% is used.

The result obtained from this is the index into the list of the list element chosen. The first element is number 0. The variable passed to DlgShow() for this item should be an integer. If the value of the variable is in the range 0 to n-1, this sets the item to be displayed. Otherwise, the first item in the list is displayed.

The following example shows how to set a list:

```
var ok%,which%:=0; 'string list, test for OK, result
DlgCreate("List example"); 'Start the dialog
var list$[3]; 'these strings are the choices
list$[0] := "one"; list$[1] := "two"; list$[2] := "three";
DlgList(1,"Choose", list$[]); 'Add the list to the dialog
ok% := DlgShow(which%); 'Display dialog, wait for user
```

From version 3, you can replace the third, fourth and fifth lines with:

DlgList(1,"Choose","one|two|three"); 'version 3 onwards

See also:DlgChan(), DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(), DlgReal(), DlgShow(), DlgString(), DlgText(), DlgXValue()

DlgReal()

This function defines a dialog entry that edits a real number. For simple dialogs, the wide, x and y arguments are not used. The sp and li arguments are new in version 5.

Proc DlgReal(item%, text\$|wide, lo, hi{, x{, y{, sp|li\$}});

- item $\$ This sets the item number in the dialog in the range 1 to the number of items.
- text^{\$} The prompt to display, optionally followed by a vertical bar and tooltip text.
- wide This is an alternative to the prompt. It sets the width in dialog units of the box in which the user types a real number. If wide is not given the box has a default width of 12 digits.
- lo, hi The range of acceptable numbers.
- x If omitted or zero, the number edit box is right justified in the dialog, otherwise this sets the position in dialog units of the left end of the number entry box.
- y Bottom of the number edit box position. If omitted, the value of item% is used.
- sp If present and non-zero, this adds a spin box with a click increment of sp.
- 1i\$ If present, this argument is a list of items separated by vertical bars that can be selected into the editing field, for example "1.0|10.0|100.0".

The variable passed into DlgShow() should be a real number. The field will start with the value of the variable if it is in the range, otherwise the value is limited to lo or hi.

See also:DlgChan(), DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(), DlgList(), DlgShow(), DlgString(), DlgText(), DlgXValue()

DIgShow()

This function displays the dialog you have built and returns values from the fields identified by item numbers, or makes no changes if the dialog is cancelled. Once the dialog has closed, all information about it is lost. You must create a new dialog before you can use this function again.

Func DlgShow(&item1|item1[], &item2|item2[], &item3|item3[] ...);

Returns 0 if the user clicked on the Cancel button, or 1 if the user clicked on OK.

For each dialog item with an item number, you must provide a variable of a suitable type to hold the result. It is an error to use the wrong variable type, except an integer field can have a real or an integer variable. Items created with DlgLabel() must have a variable too, even though it is not changed.

The variables also set the initial values. If an initial value is out of range, the value is changed to the nearest legal value. In the case of a string, illegal characters are deleted before display. In addition to passing a simple variable, you can pass an array. An array with n elements matches n items in the dialog. The array type must match the items.

```
See also:DlgChan(), DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(),
DlgList(), DlgReal(), DlgString(), DlgText(), DlgXValue()
```

DlgString()

This defines a dialog entry that edits a text string. You can limit the characters that you will accept in the string. For simple dialogs, the wide, x and y arguments are not used.

Proc DlgString(item%, text\$|wide, max%{, legal\${, x{,y{,sel\$}}});

- item% This sets the item number in the dialog in the range 1 to the number of items.
- text\$ The prompt to display, optionally followed by a vertical bar and tooltip text.
- wide This is an alternative to the prompt. It sets the width in dialog units of the box in which the user types the string. If the width is not given the number entry box has a default width of max% or 60, whichever is the smaller.

max% The maximum number of characters allowed in the string.

- legal\$ A list of acceptable characters. See Input\$() for a full description. If this is omitted, or an empty string, all characters are allowed.
- x If omitted or zero, the string entry box is right justified in the dialog, otherwise this sets the position in dialog units of the left end of the string entry box.
- y If omitted or zero, this takes the value of item%. It is the position of the bottom of the string entry box in dialog units.
- sel\$ If this string is present, it should hold a list of items separated by vertical bars, for example "one|two|three". The field becomes an editable combo box with the items in the drop down list. This is new at version 5.

The result from this operation is a string of legal characters. The variable passed to DlgShow() should be a string. If the initial string set in DlgShow() contains illegal characters, they are deleted. If the initial string is too long, it is truncated.

```
See also:DlgChan(), DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(),
DlgList(), DlgReal(), DlgShow(), DlgText(), DlgXValue(), Input$()
```

DlgText()

This places non-editable text in the dialog box. This is different from DlgLabel() as you do not supply an item number and it does not require a variable in the DlgShow() function. It returns an item number (higher than item numbers for matching arguments in DlgShow()) that you can use to identify this field in call-back functions, for example DlgVisible(). There was no returned value in Spike2 version 4.

Func DlgText(text\$, x, y{, wide});

text\$ The prompt to display, optionally followed by a vertical bar and tooltip text.

- x, y The position of the bottom left hand corner of the first character in the string, in dialog units. Set x to 0 for the default label position (the same as DlgLabel()).
- wide Normally, the width of the field is set based on text\$. This optional argument sets the width in dialog units. This allows you to replace the text with a longer string from a call-back function.

Returns An item number to identify this field for call-back functions.

See also:DlgChan(), DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(), DlgList(), DlgReal(), DlgShow(), DlgString(), DlgXValue()

DlgValue() and DlgValue\$()

These functions can only be used from a dialog call-back function to get and optionally set the value of an item, item prompt or button text. They are new in version 5.

```
Func DlgValue(item%{, val});
Func DlgValue$(item%{, val$});
```

- item% This identifies the dialog item. For items with arguments in DlgShow(), use the item% value you set to create the field. For items created with DlgText() and DlgGroup(), use the returned item number. For buttons use minus the button number. To access the prompt for an item add 1000 to the item number.
- val This optional argument holds the new item value. Use val on numeric fields or to set a checkbox or an item number in a list. Use val\$ to set a prompt, button label or the text of an editable control or to select the first matching item in a list box. It is up to you to make sure the text is acceptable for editable items.
- Returns The returned value is the current value of the item. You can use DlgValue\$() on any item to get the current contents of the field, checkbox text, button or prompt as a text string. Use DlgValue() to collect numeric or checkbox values.

If there is a problem running the command, for example if the item does not exist, or an argument type is not appropriate for an item, the result is an empty string or the value 0.

```
See also:DlgCreate(), DlgEnable(), DlgShow(), DlgVisible()
```

DlgVisible()

This can only be used from a dialog call-back function to show or hide dialog items. There are two versions of this command. The version with a single argument returns the visible state of an item; the version with two or more arguments sets the visible state of one or more dialog items. When you show or hide an item, any prompt or spin control associated with the item is also shown or hidden. This function is new in version 5.

```
Func DlgVisible(show%, item%|item%[]{, item%|item%[]...});
Func DlgVisible(item%);
```

- show% Set this to 1 to show the items in the list and to 0 to hide them.
- item% An item number of an element of the dialog or an integer array containing a list of item numbers. The item numbers are either the number you set, or the number returned by DlgText() or DlgGroup(), or -button, where button is the button number. You cannot access prompts separately from their items as this makes no sense.
- Returns When called with a single argument it returns the visible state of the item, otherwise the return value is 0.

See also:DlgAllow(), DlgCreate(), DlgEnable(), DlgShow(), DlgValue()

DIgXValue()

This creates an editable combo box to collect an x axis value for the current time, result or XY view. The combo box drop-down list is populated with cursor positions and other window values when DlgShow() runs. If the current view is not suitable, the list is empty. This control accepts expressions, for example: (Cursor(1)+Cursor(2))/2. The matching DlgShow() argument is a real number to hold a time in seconds for a time view, or an x axis value for other views. This command is new in version 5.

Proc DlgXValue(item%, text\$|wide{, x{, y}});

- item% This sets the item number in the dialog in the range 1 to the number of items.
- text^{\$} The prompt to display, optionally followed by a vertical bar and tooltip text.
- wide This is an alternative to the prompt. It sets the width in dialog units of the combo box. If the width is not given the combo box has a default width of 18 numbers.
- x If omitted or zero, the string entry box is right justified in the dialog, otherwise this sets the position in dialog units of the left end of the string entry box.
- y If omitted or zero, this takes the value of item%. It is the position of the bottom of the string entry box in dialog units.

See also:DlgChan(), DlgCheck(), DlgCreate(), DlgInteger(), DlgLabel(), DlgList(), DlgReal(), DlgShow(), DlgString(), DlgText()

Draw()

This optionally positions the current view and allows invalid regions to update. Draw() with no arguments on a view that is up-to-date should make no change. The view is not brought to the front.

Proc Draw({from {, size}});

from The left hand edge of the window. For a time window, this is in seconds. For a result view, this is in bins. For an XY view, it is in x axis units.

size The width of the window in the same units as from. A negative size is ignored. With two arguments, the width is set (unless it is unchanged) and then it is drawn. With one argument, the view scrolls by an integral number of pixels such that from is in the first pixel. The following may not move the display if a pixel is more than a second wide:

Draw(XLow()+1.0); 'This may scroll by less than 1 second

Time views run from time 0 to the maximum time in the view. Result views have a fixed number of bins, set when they are created. XY view axes can be any positive length.

See also:DrawAll(), XRange(), XLow(), XHigh(), Maxtime()

DrawAll()

This routine updates all views with invalid regions. Nowadays, calling Yield() is a better solution as this also allows the system time to clean up unused resources.

Proc DrawAll();

See also:Draw(), Yield()

DrawMode()

This sets and reads the channel display mode in a time or result view. You can set the display mode for hidden channels. The first command variant returns information, the second is for time views, the third is for result views and the fourth is for sonograms.

```
Func DrawMode(chan%{,mode%});
Func DrawMode(cSpc, mode%{, dotSz%|binSz|flags%{, trig%|edge%}});
Func DrawMode(cSpc, mode%{, dotSz% {, opt%|err%{, sort}});
Func DrawMode(cSpc, 9{, fftSz%, wnd%, top, range, xInc%, skip%});
```

chan% The channel number used to read back information with a negative mode%.

- cSpc A channel specification or -1 for all, -2 for visible or-3 for selected channels. The first channel in a result view is 1; we allow 0 to mean 1 for backwards compatibility. Setting a draw mode for a bad channel number has no effect.
- mode% If mode% is omitted, the function returns the mode. If negative, see below for return values. Positive values set the display mode. If an inappropriate mode is requested, no change is made. Some modes require additional parameters (for example a bin size). If they are omitted, the last known value is used.

The mode values for setting modes in a time view are:

- 0 The standard drawing mode for the channel.
- 1 Dots mode for events or a waveform. dotSz% argument can be used.
- 2 Lines mode. Shows event data as vertical ticks on a horizontal line. If dotSz% is present and zero, the horizontal line is not drawn.
- 3 Waveform mode. Draws straight lines between waveform and WaveMark points. No WaveMark codes are displayed (faster than WaveMark mode).
- 4 WaveMark mode. Only use this is you wish to see WaveMark codes as it takes longer to draw than Waveform mode.
- 5 Rate mode. binSz sets the width of each bin in seconds.
- 6 Mean frequency mode. binSz sets the time period.
- 7 Instantaneous frequency mode. dotSz% can be used.
- 8 Raster mode. trig% sets the trigger channel, dotSz% is used.
- 9 Sonogram mode. In this mode the fftSz%, wnd%, top, range, xInc% and skip% arguments are allowed (or can be omitted from the right).
 - fftSz% The block size for the Fourier Transform used to generate the sonogram. Allowed values are 16, 32, 64, 256, 512, 1024, 2048 and 4096. Intermediate values set the next lower allowed value.

w

nd%	The window applied to the data. $0 = no$ window, $1 =$ Hanning, $2 =$
	Hamming. Values 3 to 9 set Kaiser windows with -30 dB to -90 dB
	sideband ripple in steps of 10 dB.

top The signal dB value relative to 1 bit at the ADC input to show as the maximum density output. Signal values above top dB are shown in maximum intensity. Values 0.0 to 100.0 are allowed.

- range The range in dB of output data below top that is mapped into the grey scale. Output below (top-range) is shown as minimum intensity. Values greater than 0.0 and up to 100.0 are allowed.
- xInc% The number of pixels (1 to 100) to step across the screen before calculating the next set of sonogram values. This is normally 1.
- skip% Normally 0. Set 1 to use one data block for each vertical strip of sonogram to reduce the calculation time for large files (but only shows samples of the sonogram). Values 0 and 1 are allowed.
- 10 WaveMark overdraw mode.
- 11 Same as mode 6, but display rate per minute rather than per second.
- 12 Same as mode 7, but display rate per minute rather than per second.
- 13 Cubic spline mode for waveform and WaveMark channels.
- 14 Text mode for Textmark channels. dotSz% is used. New in version 5.
- 15 State mode for marker channels. flags% is used. New in version 5.
- 16 Skyline mode for waveform, RealWave and RealMark channels. New in version 5.

For a result view channel the modes are listed below.

- 0 The standard drawing mode for the result view.
- 1 Draw as a histogram. err% can be used.
- 2 Draw as a line. err% can be used.
- 3 Draw as dots. err% and dotSz% can be used.
- 4 Draw as a skyline. err% can be used.
- 8 Draw raster as lines, dotSz% sets the line length, opt% is used.
- 9 Draw raster as dots, dotSz% sets the dot size, opt% is used.
- 13 Cubic spline mode. err% can be used. New in version 5.
- dotSz% Sets the dot or tick size to use on screen or the point size to use on a printer. 0 is the smallest size available. The maximum size is 10. Use -1 for no size change.
- binSz Sets rate histogram bin width and the mean frequency smoothing period.
- flags% For markers in State mode, this sets the additional information to display. 0=no extra information. Add 1 to draw the state number, 2 to display TextMark text.
- trig% This is the trigger channel for time view raster displays. We assume that you do not want to display a level event channel in raster mode.
- edge% Sets the edges of level event data to use for mean and instantaneous frequency and rate modes. 0 = both edges (default), 1 = rising edges and 2 = falling edges.
- opt% Used in result view raster channels. The default is 0. Add 1 for horizontal line in raster line mode. Add 2 for y axis as time of sweep or variable value in raster modes. Add 16 to show symbols. For backwards compatibility, if sort% is omitted, add 4 to select sort value 1 and 8 to select sort value 2. These values (4 and 8) are not returned when mode% is -12, use -15 to get the sort value.
- err% Result view error drawing style: 0=none (default), 1=1 SEM, 2=2 SEM, 3=SD.
- sort% Result view raster sort mode. 0=time (default), 1-4 = use sort variables 1 to 4.

Returns If a single channel is set it returns the previous mode%. For multiple channels or an invalid call, it returns -1. Negative mode% values return drawing parameters:

	i illi alla call, it	1000	ins it it gau		ouco fuido	s recarn aratting	parameters.
-1	Reserved	-4	trig%	-7	wnd%	-10 xInc%	-13 err%
-2	dotSz%	-5	edge%	-8	top	-11 skip%	-14 sort%
-3	binSz	-6	fftSz%	-9	range	-12 opt%	-15 flags%
See also:Draw(),SetEvtCrl(),SetPSTH(),ViewStandard(),XYDrawMode()							

Dup()	This gets the view handle of a duplicate of the current view or the number of duplicates. Duplicated views are numbered from 1 (1 is the original). If a duplicate is deleted, higher numbered duplicates are renumbered. See WindowDuplicate() for more information.					
	<pre>Junc Dup({num%});</pre>					
	num% The number of the duplicate view to find, starting at 1. You can also pass 0 (or omit num%) as an argument, to return the number of duplicates.					
	Returns If num% is greater than 0, this returns the view handle of the duplicate, or 0 if the duplicate does not exist. If num% is 0 or omitted, this returns the number of duplicates. The following code illustrates the use of Dup().					
	<pre>var maxDup%, i%, dvh%; 'declare variables maxDup% := Dup(0); 'Get maximum numbered duplicate for i% := 1 to maxDup% do 'loop round all possible duplicates dvh% := Dup(i%); 'get handle of this duplicate if (dvh% > 0) then 'does this duplicate exist? PrintLog(view(dvh%).WindowTitle\$()+"\n");'print window title endif; next;</pre>					
	See also:App(), View(), WindowDuplicate()					
DupChan()	This returns the number of duplicates of a channel in the current time or result view, or the channel number of the n th duplicate, or it identifies the channel on which a duplicate channel is based, or it identifies which duplicate of an original channel this is. See ChanDuplicate() for more information on duplicate channels.					
	Func DupChan(chan%{, num%});					
	chan% A channel number in the current view. The channel must exist. If this channel is a duplicate, the command behaves as though you passed the original channel on which the duplicates are based.					
	num% If greater than 0, this is the index of a duplicate channel in the range 1 to the number of duplicates. Use 0 to return the original channel that was duplicated. Use -1 as an argument (or omit num%), to return the number of duplicates of the channel. From version 5.10 onwards you can use -2 as an argument to return which duplicate a channel is: 0 means not a duplicate, 1 upwards is the duplicate index.					
	Returns If num% is greater than 0, this returns the channel number of the duplicate, or 0 if the duplicate does not exist. If num% is 0, this returns the channel number that chan% was duplicated from or chan% if it is the original. If num% is -1 or omitted, this returns the number of duplicates. If num% is -2, this returns the duplicate index of chan% (0 if not a duplicate, 1 for duplicate a and so on).					
	See also:ChanDuplicate()					
EditClear()	In a result view, this command zeros the histogram data. In a text view, this command deletes any selected text.					
	<pre>Func EditClear();</pre>					
	Returns The function returns 0 if nothing was deleted, otherwise it returns the number of items deleted or 1 if the number is not known, or a negative error code.					

See also:EditCut()

EditCopy()	This command copies data from the current view to the clipboard. Time views copy as a bitmap, as a scaleable image and as text in the format set by <code>ExportRectFormat()</code> , <code>ExportTextFormat()</code> , <code>ExportChanFormat()</code> and <code>ExportChanList()</code> . Result and XY views can copy as a bitmap, text and as a scaleable image. Multimedia views copy as a bitmap. Text windows copy as text.		
	 Func EditCopy({as%}); as% This sets how to copy data when several formats are possible. If omitted, al formats are used. Systems that cannot implement some or all of these feature ignore the argument. The format value is the sum of: Copy as a bitmap Copy as a scaleable image (Windows metafile) Copy as text Returns It returns 0 if nothing was placed on the clipboard, or the sum of the format 		
	<pre>Specifiers for each format used. See also:EditSelectAll(),ExportTextFormat(),ExportChanFormat(), ExportRectFormat(),ExportChanList(),EditCut(),EditPaste()</pre>		
EditCut()	<pre>This command cuts data from the current view to the clipboard and deletes the original. Func EditCut({as%}); as% This optional argument sets how data is cut and copied to the clipboard. Currently only text may be cut. If omitted, all allowed formats are used.</pre>		
EditPaste()	You may only paste into a text window when the clipboard contains text. The contents of the clipboard are inserted into the text at the current caret. If text is already selected, it is replaced by the clipboard contents. Func EditPaste(); Returns The function returns the number of characters inserted.		
EditSelectAll()	This function selects all items in the current view that can be copied to the clipboard. This is the same as the Edit menu Select All option. Func EditSelectAll(); Returns It returns the number of selected items that could be copied to the clipboard. See also:EditCopy(),EditClear(),EditCut()		

Error\$()	This function converts a negative error code returned by a function into a text string.					
	Func Error\$(code%);					
	code% A negative error code returned from a Spike2 function.					
	Returns It returns a string that describes the error.					
	See also:Debug(), Eval(), Print\$(), PrintLog()					
Eval()	Script window or Evaluate window message area, as appropriate, when the script ends.					
	Proc Eval (arg);					
	arg A real or integer number or a string.					
	If you use Eval() it will suppress any run-time error messages as it uses the same					
	mechanism as the error system. A common use of $Eval()$ in a script is to report an error condition during debugging, for example:					
	if val<0 then Eval("Negative value"); Halt; endif;					
	See also:Debug(), Error\$(), Print(), PrintLog()					
EventToWaveform()	This function creates a waveform channel from an event channel. The waveform channel contents depend on the frequency of the events. You can convert based on instantaneous					
	or smoothed frequency. VirtualChan () also converts events to a waveform.					
Instantaneous frequency	Spike2 calculates instantaneous frequency as the reciprocal of event intervals. Between					
	between the last two events. If it is less, the reciprocal of the last interval is used. If it is					
	more, the reciprocal of the interval from the last event is used.					
	$1/(t_2-t_1)$ $1/(t-t_2)$ $1/(t_3-t_2)$					
	<i>t1 t2 t t3</i>					
	The diagram shows the result for three events at times $t1$, $t2$ and $t3$. From $t1$ to $t2$, the					
	result depends on events before t_1 . From t_2 to $t_2+(t_2-t_1)$, the result is set by the previous interval. From this point onwards to t_3 , the output reduces until at t_3 it becomes $1/(t_3-t_2)$.					
Smoothed frequency	Spike2 calculates smoothed frequency by replacing (convoluting) each event with a					

waveform of unit area. The built-in waveforms are symmetrical about the event time. If the event is at time t, the waveforms extend from time *t*-smooth to *t*+smooth. You can

also provide your own waveforms, and these need not be symmetric.



The diagram shows the result of using the built-in raised cosine waveform. The area under the curve for each spike is 1. You would normally use a smoothing period that covered several events to obtain a smooth output.

The command Func EventToWaveform(smooth, eChan%, wChan%, wInt, sTime, eTime, maxF, meanF {,query% {,wave[]|wave%{, nPre}}});

- smooth The smoothing period in seconds. Set this 0 or negative for instantaneous frequency, otherwise an event at time t is spread over a time range t-smooth to t+smooth seconds for symmetric functions. If you omit the wave argument, the smoothing function is a raised cosine bell.
- eChan% The channel number to use as a source of event times. If the channel is a marker, and a marker filter is set, only filtered events are visible.
- wChan% The channel number to create as a waveform channel in the range 1 to the maximum allowed channels in the file.
- wInt The desired sampling interval for the new channel in seconds. Values outside the range 0.000001 to 1000 are errors. This interval is rounded to the nearest multiple of the microseconds per time unit set for the file. Use BinSize() to return the actual sampling interval. Before version 4.03, wInt was rounded to the nearest multiple of the base ADC convert interval for the file. Older Spike2 versions will not open the file if you set an interval that they cannot achieve.
- sTime The start time for output of the waveform data.
- eTime The last time for waveform output will be less than or equal to this time. The command processes events from sTime-smooth to eTime+smooth.
- maxF The output is stored as if it were a sampled signal. That is as a 16-bit integer in the range -32768 to 32767. maxF sets the frequency that corresponds to the 32767 value. It is the maximum frequency that can be represented in the result. Frequencies higher than this are limited to maxF.
- meanF This sets the frequency that corresponds to 0 in the 16-bit waveform data. Most users set this to 0. The lowest frequency that can appear in the result is maxF-2* (maxF-meanF). Frequencies lower than this are limited to this value.

If the result lies in the range freqtdelta, you get the best resolution by setting maxF to freq+delta and meanF to freq.

- guery% Set this to 0 or omit it to query overwriting an existing channel. If this is present and non-zero, an existing channel is overwritten with no comment.
- wave% This optional argument is used when smooth is greater than 0 and sets the smoothing function: 0=rectangular, 1=triangular, 2=raised cosine, 3=Gaussian (to 4 sigmas). If omitted, raised cosine smoothing is used. These functions are the same as are used for virtual channels.
- wave[] If you provide an array, it is scaled to span smooth seconds. The sum of the waveform values must be in the range 1.2e-38 to 3.4e+38. If nPre is omitted or negative, the smoothing function is symmetrical and you supply the right-hand

side only. You can set any array size; the built-in functions use 1000 points. This example matches the raised cosine generated with wave% set to 2.

var wave[1000], i%, nPre := -1.0; const pi := 3.14159; for i% := 0 to 999 do wave[i%]:=1 + Cos(i%*pi/1000) next;

nPre This sets the index in wave[] that matches the event time for asymmetric smoothing functions. The entire wave[] array is smooth seconds wide. If nPre is omitted or negative, a symmetric function is assumed. This example generates a typical asymmetric smoothing function in the wv[] array:

```
var wv[400],i%, nPre;
const a:=40, b:=80; 'The time constants in points
for i%:=0 to 399 do wv[i%]:=exp(-i%/a)*(1-exp(-i%/b)) next;
nPre := b*Ln(1+a/b); 'Calculate position of maximum
```

Returns The function returns 0 unless the user decided not to overwrite a channel, in which case the result is a negative error code. You can also get error -5799 if you run out of memory (unlikely) or disk space.

Uses of EventToWaveform This can be used to investigate frequency changes in an event train. If you have an event channel of heart beats at 70 beats per minute (1.17 Hz) and you are seeking a respiration-induced modulation at about 0.25 Hz, you can convert the heartbeats into a waveform with this command. The result would be a constant level of 1.17 Hz, plus a small oscillation around 0.25 Hz. The oscillation can be detected with SetPower().

For a 0.01 Hz resolution with SetPower(), you need 100 seconds of data per power spectrum block. With a 1024 point transform size, a wInt of 0.1 Hz gives a block size of 102.4 seconds. The SetPower() result will have 512 bins, spanning the range 0 to 5.00 Hz in steps of 0.01 Hz. A smooth period of 5-10 seconds would be appropriate.

See also:BinSize(), SetPower(), VirtualChan()

Exp()

This function calculates the exponential function for a single value, or replaces a real array by its exponential. If a value is too large, overflow will occur, causing the script to stop for single values, and a negative error code for arrays.

Func Exp(x|x[]|x[][]);

- x The argument for the exponential function or an array of real values.
- Returns With an array, the function returns 0 if all was well, or a negative error code if a problem was found (such as overflow).

With an expression, it returns the exponential of the number.

See also:Abs(), ATan(), Cos(), Frac(), Ln(), Log(), Max(), Min(), Pow(), Rand(), Round(), Sin(), Sqrt(), Tan(), Trunc()

ExportChanFormat()

This sets the channel text format used by ExportTextFormat(), FileSaveAs() and EditCopy(). It is equivalent to the Channel type parameters section of the File Dump Configuration dialog. See ExportTextFormat() for the command to reset the format.

Func ExportChanFormat(type%, synop%, data%{, as%});

type% The channel type to set the format for. Types 2, 3 and 4 share the same settings. Types 1 and 9 also share the same settings.

1 Waveform3 Event (Evt+)5 Marker7 RealMark9 RealWave2 Event (Evt-)4 Level (Evt+-)6 WaveMark8 TextMark

synop% Set this non-zero to enable synopsis output for this channel type.

- data% Set this non-zero to enable data output for this channel type.
- This sets the format for markers and descendent types. It is ignored and can be omitted for waveform and event channels. If omitted or set to 0, the channels are treated as their own type. The codes are the same as for the type% field.

Returns 0 if all was OK or a negative error code.

See also:FileSaveAs(), ExportChanList(), ExportTextFormat()

ExportChanList()

This command sets a time range and a channel list to export for use by FileSaveAs() and EditCopy() in a time view. The current view must be a time view. When used with ExportTextFormat(), repeated calls to this routine are additive. When used with ExportRectFormat(), only the first set of channels is used. Call it with one or no arguments to clear the list.

Proc ExportChanList(sTime, eTime, cSpc {,cSpc...});
Proc ExportChanList({flags%});

- sTime The start time of the range of data to save.
- eTime The end time of the range of data to save.
- cSpc A channel specifier for the channels to export.
- flags% This affects data written by FileSaveAs() with type% of 0 and is the sum of:
 1 Time shift data so that the earliest sTime appears at time 0 in the output file.
 2 White Particular Statement of the output file.
 - 2 Write RealWave channels as waveform data for backwards compatibility.

ExportRectFormat()

This command sets the text export format for use by FileSaveAs() and EditCopy() in a time view. This is an alternative to using ExportTextFormat() and generates equally spaced output in time for all channels. If used with no parameters, this will reset the rectangular text output configuration to the default of 100 Hz output, linear interpolation, WaveMark as event, RealMark as waveform, double-quotes delimiter and tab for the separator. Otherwise at least 3 arguments must be supplied.

Proc ExportRectFormat({freq, inter%{, delim\${, sep\$}}})

freq This sets the number of output points per channel generated per second.

- inter% This sets the interpolation method:
 - 0 Use the data point nearest to the time.
 - 1 Linear interpolation between the points either side.
 - 2 Cubic spline interpolation (waveform data only).
 - 3 Cubic spline interpolation if drawn splined, else linear (waveforms only).
- delim\$ This sets the character used to delimit the start and end of text strings. If omitted, text strings are surrounded by quote marks, for example, "Volts".

sep\$ The character to separate multiple data items on a line. If omitted, a Tab is used.

The following example exports channels 1, 2 and 3 at 200 Hz from times 10 seconds to 30 seconds as a text file. The times are shifted so that the first data item is at time 0. Waveform channels are mapped to 200 Hz using cubic spline interpolation.

ExportChanList(1); 'Clear export list, sets zero shift ExportChanList(10,30,1,2,3); 'Choose channels 1, 2 and 3 ExportRectFormat(200, 2, 0); 'Treat extended markers as events

```
FileSaveAs("rect.txt", 1); 'Export as text file, and...
EditCopy(); '...also copy to the clipboard
See also:FileSaveAs(), ExportChanList(), EditCopy()
```

ExportTextFormat()

This command sets the text export format for use by FileSaveAs() and EditCopy() in a time view. It is equivalent to the top section of the File Dump Configuration dialog. The form with no arguments resets the dialog to enable all check boxes, sets the columns to 1, sets the string delimiter to a double quote mark, sets the separator to a tab character, and sets all channel types to be treated as themselves.

```
Proc ExportTextFormat(head%, sum%, cols%{, delim${, sep$}});
Proc ExportTextFormat();
```

- head% If this is non-zero, Spike2 will output the header.
- sum% If this is non-zero, Spike2 outputs the channel summary information.
- cols% The number of columns to use when writing waveforms and event times.
- delim\$ This sets the character used to delimit the start and end of text strings. If omitted, text strings are surrounded by quote marks, for example, "Volts".
- sep\$ The character to separate multiple data items on a line. If omitted, a Tab is used.

The following example exports channels 1 and 2 from 90.0 to 100.0 seconds and from 110 to 120 seconds as a file 30 seconds long called fred.smr.

```
ExportChanList(1); 'Clear export list, sets zero shift
ExportChanList(90, 100, 1, 2);'set channels 1 and 2
ExportChanList(110,120, 1, 2);'add 10 more seconds after 10s gap
ExportTextFormat(); 'reset the file dump dialog
ExportChanFormat(1,0,1); 'turn off synopsis for waveform data
ExportTextFormat(0, 0, 1); 'No header or summary
FileSaveAs("fred.smr", 0); 'export to new data file
```

See also:FileSaveAs(), ExportChanList(), ExportChanFormat(), EditCopy()

FileApplyResource()

This applies a resource file to the current time window. If the window is duplicated, all other duplicates are deleted, then the resource file is applied, which may create duplicate windows. The current time window handle is not changed. Handles of duplicate windows are not preserved, even if the resource file creates the same number of duplicates.

Func FileApplyResource(name\$);

- name\$ The resource file to apply. You must include the .s2r file extension and required path. Remember that \ in a string must be entered as \\ or use /. If the name is "" or contains a "*" or a "?", the user is prompted for a file.
- Returns The number of modified windows (usually 1 unless the resource file generates duplicates), 0 if the resource file did not contain suitable information, -1 if the user cancelled the file dialog, -2 if the file could not be found or -3 if there is any other problem.

See also:FileGlobalResource(), FileOpen(), FileSaveResource()

FileClose()

This is used to close the current window or external file. You can supply an argument to close all windows associated with the current time or result view or to close all the windows belonging to the application.

Func FileClose({all% {,query%});

- all% This argument determines the scope of the file closing. Possible values are:
 - -1 Close all windows except loaded scripts and debug windows.
 - 0 Close the current view. This is the same as omitting all%.
 - 1 Close all windows associated with the current view.
- query% This determines what happens if a window holds unsaved data:
 - -1 Don't save the data or query the user
 - 0 Query the user about each window that needs saving. If the user chooses Cancel, the operation stops, leaving all unclosed windows behind. This is the same as omitting query%.
- Returns The number of views that have not been closed. This can occur if a view needs saving and the user requests Cancel.

If the current view is a time view and you have been sampling data, you must call SampleStop() before using this command to set a file name. If you do not, this is equivalent to aborting sampling and your data file will not be saved. Use FileSaveAs() just before calling FileClose() to set a file name.

Do not use View(fh%).FileClose() if fh% is the current view; you will get an error when Spike2 tries to restore the current view. Use View(fh%);FileClose() instead.

See also:FileOpen(), FileSave(), FileSaveAs(), FileNew(), SampleStop()

FileComment\$()

This function accesses the file comments in the file associated with the current time view. Files have five comment strings. Each string is up to 79 characters in length.

Func FileComment\$(n% {,new\$));

- n % The number of the file comment line in the range 1 to 5
- new\$ If present, the command replaces the existing comment with new\$.

Returns The comment at the time of the call or a blank string if n% is out of range.

See also:ChanComment\$()

FileConvert\$()

This function converts a data file from a "foreign" format into a Spike2 data file. The range of foreign formats supported depends on the number of import filters in the Spike5\import folder.

Func FileConvert\$(src\${, dest\${, flag%{ ,&err%}}});

- src\$ This is the name of the file to convert. The file extension is used to determine the file type (unless flag% bit 0 is set). Known file extensions include: abf, cfs, cnt, cut, dat, eeg, ewb, ibw, son and uff. We expect to add more.
- dest\$ If this is present, it sets the destination file. If this is not a full path name, the name is relative to the current directory. If you do not supply a file extension then Spike2 appends ".smr". If you set any other file extension, Spike2 cannot open the file as a Spike2 data file. If you do not supply this argument, the converted file will be written to the same folder as the source file using the original file name with the file extension changed to .smr.
- flag% This argument is the sum of the flag values: 1=Ignore the file extension of the source file and try all possible file converters, 2=Allow user interaction if required (otherwise sensible, non-destructive defaults are used for all decisions).
- err% Optional integer variable that is returned as 0 if the file was converted, otherwise it is returned holding a negative error code.

Returns The full path name of the created file, or an empty string if the file was not converted.

See also:FileOpen(), FilePath\$(), FilePathSet(), FileList()

FileCopy()

This function copies a source file to a destination file. File names can be specified in full or relative to the current directory. Wildcards cannot be used.

Func FileCopy(src\$, dest\${, over%});

- src\$ The source file to copy to the destination. This file is not changed.
- dest\$ The destination file. If this file exists you must set over% to overwrite it.
- over% If this optional argument is 0 or omitted, the copy will not overwrite an existing destination file. Set to 1 to overwrite.
- Returns The routine returns 1 if the file was copied, 0 if it was not. Reasons for failure include: no source file, no destination path, insufficient disk space, destination exists and insufficient rights.

See also:BRead(), BWrite(), FileDelete(), FileOpen(), ProgRun()

FileDate\$()

Since version 4.03, the time and date at which sampling started is saved in the data file. This function returns a string holding the date. Use FileTimeDate() to get the date as numbers. The current view must be a time view. The arguments are exactly the same as for the Date\$() command. If there is no date stored, the result is an empty string.

Func FileDate\$({dayF%, {monF%, {yearF%, {order%, {sep\$}}}}); See also:Date\$(),FileTimeDate(),FileTime\$()

FileDelete()

This deletes one or more files. File names can be specified in full or relative to the current directory. Windows file names look like $x:\folderl\folder2\folder2\folderl\folder2\folderl\file.ext$ or $\mbox{machine}\folder1\folder2\foldern\file.ext}$ across a network. If a name does not start with a $\$ or with $x:\$ (where x is a drive letter), the path is relative to the current directory. Beware that $\$ must be written $\$ in a string passed to the compiler.

Func FileDelete(name\$[]|name\${, opt%});

- name\$ This is either a string variable or an array of strings that holds the names of the files to delete. Only one name per string and no wildcard characters are allowed. If the names do not include a path they refer to files in the current directory.
- opt% If this is present and non-zero, the user is asked before each file in the list is deleted. You cannot delete protected or hidden or system files.
- Returns The number of files deleted or a negative error code.

See also:FilePath\$(), FilePathSet(), FileList()

FileGlobalResource()

This function is equivalent to the Global Resources dialog. Please see the documentation for the dialog for a full explanation. The function has two forms: the first is for setting values, the second is to read them back:

Func FileGlobalResource(flags%, loc%, name\${, path\$});
Func FileGlobalResource(&loc%, &name\$, &path\$);

flags% This is the sum of the following values:

- 1 Enable the use of global resources. If unset, no global resources are used.
- 2 Only use global resources if a data file has no associated resource file.
- 4 Only use if the data file is within the path set by path\$.
- 100% The location of the global resource file. 0=the folder that Spike2 was run from, 1=search the data file folder, then the Spike2 folder, 2=the data file folder only.
- name\$ The name of the global resource file excluding the path and .s2r file extension.
- path\$ The file path within which to use the global resources. If you omit this argument when setting values, the current path does not change.
- Returns When setting values, the return value is 0 or a negative error code. When reading back values, the return value is the flags% argument.

See also:FileApplyResource(), FileOpen(), FileSaveResource()

This function gets lists of files and sub-directories (folders) in the current directory and can also return the path to the parent directory of the current directory. This function can be used to process all files of a particular type in a particular directory.

Func FileList(names\$[]|&name\$, type%{, mask\$});

- name\$ This is either a string variable or an array of strings that is returned holding the name(s) of files or directories. Only one name is returned per string.
- type% The objects in the current directory to list. The Parent directory is returned with the full path, the others return the name in the current directory. Values are:
 - -3 Parent directory

All files

-2

-1

0

1

- 2 Output sequence files (*.pls)
- Sub-directories
- 3 Spike2 script files (*.s2s)
 4 Spike2 result view files (*.s2s)
- 4 Spike2 result view files (*.srf)
 6 Spike2 configuration files (*.s2c)
- Spike2 data files (*.smr) Text files (*.txt)
- 12 XY view data file (*.sxy)
- mask\$ This optional string limits the names returned to those that match it; * and ? in the mask are wildcards. ? matches any character and * matches any 0 or more characters. Matching is case insensitive and from left to right.
- Returns The number of names that met the specification or a negative error code. This can be used to set the size of the string array required to hold all the results.

See also:FilePath\$(), FilePathSet(), FileDelete(), FileName\$()

FileName\$()

This returns the name of the data file associated with the current view (if any). You can recall the entire file name, or any part of it. If there is no file the result is an empty string.

Func FileName\$({mode%});

mode% If present, determines what to return, if omitted taken as 0.

- 0 Or omitted, returns the full file name including the path
- 1 The disk drive/volume name
- 2 The path section, excluding the volume/drive and the name of the file
- 3 The file name up to but not including the last . excluding any trailing
 - number.
- 4 Any trailing numbers from 3.
- 5 The end of the file name from the last dot.

Returns A string holding the requested name, or a blank string if there is no file.

See also:FileList(), FilePath\$(), FilePathSet()

FileNew()

This is equivalent to the File menu New command. It creates a new window and returns the handle. You can create visible or invisible windows. Creating an invisible window lets you set the window position and properties before you draw it. The new window is the current view and if visible, the front view. Use FileSaveAs() to name created files.

Func FileNew(type%{, mode%{, upt, tpa%, maxT{, nChan%}});

type% The type of file to create:

- 0 A Spike2 data file based on the sampling configuration, ready for sampling. This may open several windows, including the floating command window and the sequencer control panel. Use SampleStart() to begin sampling and SampleStop() to stop sampling before calling FileSaveAs() to give the new file a name. During sampling, type 0 data files are saved in the folder set by the Edit menu preferences or the FilePathSet() command. Use FileSaveAs() command after SampleClose() to move the data file to its final position on disk.
- 1 A text file in a window
- 2 An output sequence file in a window
- 3 A Spike2 script file
- 7 An empty Spike2 data file (not for sampling). upt%, tpa% and maxT must be supplied. You can use ChanNew(), ChanWriteWave(), ChanSave(), MemSave() or to add data. Use FileSaveAs() to name the new file.
- 12 An XY view with one (empty) data channel. Use XYAddData() to add more data and XYSetChan() to create new channels.
- mode[%] This optional argument determines how the new window is opened. The value is the sum of these flags. If the argument is omitted, its value is 0. The flags are:
 - 1 Make the new window(s) visible immediately. If this flag is not set the window is created, but is invisible.
 - 2 For data files, if the sampling configuration holds information for creating additional windows, use it. If this flag is not set, data files extract enough information from the sampling configuration to set the sampling parameters for the data channels.
 - 4 Show the spike shape setup dialog if there are WaveMark channels in the sampling configuration for a data file. If the spike shape dialog appears, this function does not return until the user closes the dialog.
- upt Only for type = 7. The microseconds per unit time for the new file. This allows values in the range 0.001 to 32767. If upt is non-integral values, Spike2 versions prior to 4.02 will not open the file. This sets the time resolution of the new file and the maximum length of the file, in time. There are a maximum of 2^{31} -1 time units in a file, so a 1 µs resolution limits a file to 30 minutes in length.
- tpa% Only for type% = 7. The time units per ADC conversion for the new file. You can set values in the range 1 to 32767 (for compatibility with Spike2 sampling use 2 to 32767). The available sampling intervals for waveform data in the new file are n * upt% * tpa% microseconds where n is an integer.
- maxT Only for type% = 7. This sets the initial file size in time, which sets the value returned by MaxTime() until data written to the file exceeds this time.
- nChan% Only for type% = 7. This sets the number of channels in the file in the range 32 to 256. If you omit this argument, 32 channels are used. Spike2 version 4 can read files with 32-100 channels. Version 3 can only read files with 32 channels.
- Returns It returns the view handle (or the handle of the lowest numbered duplicate for a data file with duplicate windows) or a negative error.

See also:ChanNew(), ChanSave(), ChanWriteWave(), FileOpen(), FileSave(), FileSaveAs(), FileClose(), FilePathSet(), MaxTime(), MemSave(), SampleStart(), SampleStop(), XYAddDAta(), XYSetChan()

FileOpen()

This is the equivalent of the File Open... menu command. It opens an existing Spike2 data file or a text file in a window, or an external text or binary file. If the file is already opened, a handle for the existing view is returned. The window becomes the new current view. You can create windows as visible or invisible. It is often more convenient to create an invisible window so you can position it before making it visible.

Func FileOpen(name\$, type% {,mode% {,text\$}});

name\$ The name of the file to open. This can include a path. The file name is operating system dependent, see FileDelete(). If the name is blank or contains * or ? (Windows only), the file dialog opens for the user to select a file. For file types 8 and 9 only, name\$ can be of the form:

"Type 1 (*.fl1) |*.fl1||Type 2 (*.fl2;*.fl3) |*.fl2;*.fl3||"

This produces two file types, one with two extensions. There is one vertical bar between the description and the template and two after every template.

- type% The type of the file to open. The types currently defined are:
 - 0 Open a Spike2 data file in a window
 - 1 Open a text file in a window
 - 2 Open an output sequence file in a window
 - 3 Open a Spike2 script file in a window
 - 4 Open a result view file in a window
 - 6 Load configuration file or read configuration from data file
 - 8 An external text file without a window for use by Read() or Print()
 - 9 An external binary file without a window for use by BRead(), BWrite(), BSeek() and other binary routines.
 - 12 Open an XY view file in a window.
- mode[®] This optional argument determines how the window or file opens. If the argument is omitted, its value is 0. For file types 0 to 4 the value is the sum of:
 - 1 Make the new window(s) visible immediately. If this flag is not set the window is created, but is invisible.
 - 2 Read resource information associated with the file. This may create more than one window, depending on the file type. For data files, it restores the file to the state as it was closed. If the flag is unset, resources are ignored.
 - 4 Return an error if the file is already open in Spike2. If this flag is not set and the file is already in use, it is brought to the front and its handle is returned.

When used with file types 8 and 9 the following values of mode% are used. The file pointer (which sets the next output or input operation position) is set to the start of the file in modes 0 and 1 and to the end in modes 2 and 3.

- 0 Open an existing file for reading only
- 1 Open a new file (or replace an existing file) for writing (and reading)
- 2 Open an existing file for writing (and reading)
- 3 Open a file for writing (and reading). If the file doesn't exist, create it.
- text\$ An optional prompt displayed as part of the file dialog (not supported for type%=6 on the Macintosh 68k version).
- Returns If a file opens without any problem, the return value is the view handle for the file (if multiple views open, it is the handle for the first time view created). For configuration files (type% of 6), the return value is 0 if no error occurs. If the file could not be opened, or the user pressed Cancel in the file open dialog, the returned value is a negative error code.

If multiple windows are created for a data file, you can get a list of the associated view handles using ViewList(list%[],64).

See also:FileDelete(), FileNew(), FileSave(), FileSaveAs(), FileClose(), BRead(), BReadSize(), BSeek(), BWrite(), BWriteSize(), ViewList()

FilePath\$()	This fur file save	nction gets the "current directory". This is the place on disk where file open and e dialogs start. It can also get the path for created data files and auto-saving.			
	<pre>Func FilePath\$({opt%});</pre>				
	opt%	If 0 or omitted, this gets the current directory, 1 returns the path for Spike2 data files created by FileNew() as set by Edit menu Preferences, 2 gets the path to the Spike2 application, 3 returns the path for automatic file saving as set by the Sampling Configuration Automation tab.			
	Returns	A string holding the path or an empty string if an error is detected.			
	See also	:FileList(),FileName\$(),FilePathSet()			
FilePathSet()	This fur FileNe	action sets the "current directory/folder", and where Spike2 data files created by $w()$ are stored until they are sent to their final resting place by FileSaveAs().			
	Func F	<pre>ilePathSet(path\${,opt%{, make% text\$}});</pre>			
	path\$	A string holding the new path to the directory. The path must conform to the rules for path names on the host system and be less than 255 characters long. If the path is empty a dialog opens for the user to select an existing directory/folder.			
	opt%	If omitted or zero, this sets the current directory. 1 sets the path for Spike2 data files created by FileNew(), 2 is ignored, 3 sets the automatic file naming path.			
	make%	If this is present and non-zero, the command will create the directory/folder if all elements of the path exist except the last. You cannot use this option if path\$ is blank.			
	text\$	Optional prompt for use with the dialog.			
	Returns	Zero if the path was set, or a negative error code.			
	See also	:FileList(),FileName\$(),FilePath\$()			
FilePrint()	This fur current	nction is equivalent to the File menu Print command. It prints some or all of the view to the printer that is currently set for Spike2. If no printer has been set, the			
	current axis sca currentl	system printer is used. In a time or result view, it prints a range of data with the x ling set by the display. In a text or log view, it prints a range of text lines. There is y no script mechanism to choose a printer; you must do it interactively.			
	Func F	<pre>ilePrint({from{, to{, flags%}}});</pre>			
	from	The start point of the print. This is in seconds in a time view, in bins in a result view and in lines in a text view. If omitted, this is taken as the start of the view.			
	to	The end point in the same units as from. If omitted, the end of the view is used.			
	flags%	0=portrait, 1=landscape, 2=current setting. If omitted, the current value is used.			
	Returns The function returns 0 if all went well; otherwise it returns a negative error.				

The format of the printed output is based on the screen format of the current view. Beware that for time and result views the output could be many (very many) pages long.

All the FilePrintXXXX() routines allow you to choose the orientation of the printed output. If you do not set the orientation, the last print orientation you used in the current Spike2 session is set. If you have not printed, the orientation depends on the current orientation set for the system printer.

See also:FilePrintScreen(), FilePrintVisible()
FilePrintScreen()	This function is equivalent to the File menu Print Screen command. It prints all visible time, result, XY and text-based views to the current printer on one page. The page positions are proportional to the view positions in the Spike2 application window.
	<pre>Func FilePrintScreen({pTtl\${, vTtl%{, box%{,scTxt%{,flags%}}}});</pre>
	pTtl\$ This sets the page title string to print at the top of a page. If omitted or an empty string, there is no page title.
	vTtl% Set 1 or higher to print a title above each view, omitted or 0 for no title.
	box% Set 1 or higher for a box around each view. If omitted, or 0, no box is drawn.
	scTxt% Set 1 or higher to scale text differently in the x and y directions to match the original. If omitted or 0 scale both directions by the same scale factor.
	flags% 0=portrait, 1=landscape, 2=current setting. If omitted, the current value is used. In version 4, the default was portrait mode.
	Returns The function returns 0 if it all went well, or a negative error code.
	See also:FilePrint(), FilePrintVisible()
FilePrintVisible()	This function prints the current view as it appears on the computer screen to the current printer. In a text window, this prints the lines in the current selection. If there is no selection, it prints the line containing the cursor.
	<pre>Func FilePrintVisible({flags%});</pre>
	flags% 0=portrait, 1=landscape, 2=current setting. If omitted, the current value is used.
	Returns The function returns 0 if all went well; otherwise it returns a negative error.
	See also:FilePrint(), FilePrintScreen()
FileQuit()	This is equivalent to the File menu Exit command. You are asked if you wish to save any unsaved data before the application closes. If the user cancels the operation (because there were files that needed saving), the script terminates, but the Spike2 application is left running. Use FileClose $(-1, -1)$ before FileQuit() to guarantee to exit.
	<pre>Proc FileQuit();</pre>
	See also FileClose ()
FileSave()	This function is equivalent to the File menu Save command and saves the current view as
Theodve()	a file on disk. If the view has not been saved previously the File menu Save As dialog opens to request a file name. It cannot be used with external text or binary files. It cannot be used with a time view unless it has been sampled and has not been saved.
	<pre>Func FileSave();</pre>
	Returns The function returns 0 if the operation was a success, or a negative error code.
Use with a time view	You cannot use this command for a time view unless it has just been sampled; use FileSaveAs() instead. When used with a time view that has just been sampled, if an automatic file name is set, there is no prompt for a file name and the next automatic name is used. However, if automatic filing is enabled, the file is saved as soon as sampling stops, so this function is not needed.
	See also:FileOpen(), FileSaveAs(), FileSaveResource(), FileClose().

SampleAutoFile(), SampleAutoName\$()

FileSaveAs()	This fur picture. external FileNe	action saves the current view as a file on disk in its native format, as text or as a It is equivalent to the File menu Save As command. This cannot be used for text or binary files. This command can also name data files created with $w(7,)$ and after sampling, but read the remarks about time views below.
	Func F	<pre>ileSaveAs(name\${, type%{, yes%{, text\${, nChan% flag%}}});</pre>
	name\$	The output file name including the file extension. If the string is empty or it holds wild card characters * or ? or text\$ is not empty, then the File menu Save As dialog opens and name\$ sets the initial list of files.
	type%	The type to save the file as (if omitted, type -1 is used):
		-2 Valid for XY and result views only. Opens a File Save As dialog and allows the user to choose the file type to save as.
		-1 Save in the native format for the view.
		You can also use this to save and name a Spike2 data file immediately after it has been sampled. This will also save any multimedia data associated with the file. The correct sequence is:
		 Use FileNew(0, mode%) to create a data file Use SampleStart() to start data capture When capture is over use SampleStop() to stop capture, or wait for the SampleStatus() to indicate that sampling has finished. Use FileSaveAs(name\$,-1) to set the file name Use FileClose() to close the view
		Once a data file has been named you must use type 0 to save a selection of channels to a new file; you cannot use type -1 again in a time view.
		0 Save sections set by ExportChanList() to a new Spike2 data file. This will not save associated multimedia files.
		Save the contents of the current view as a Text file. If the current view is a time window, Spike2 saves the channels set by ExportChanList() in the text format set by ExportRectFormat() or ExportTextFormat() and ExportChanFormat(). This also can be used for result and XY views.
		2 Save a text view as an output sequence file.
		3 Save a text view as a Script file.
		4 Save as a result view (result views only).
		5 Save the screen image of a time, result or XY view as a metafile. The file extension sets the format; use WMF for a Windows metafile and EMF for an enhanced metafile. You cannot save a view as a bitmap file from the script.
		6 Save sampling configuration in a configuration file.
		12 Save as XY view (XY views only).
		13 Copy the screen image as bitmap (time, result and XY views only). If the view is obscured by another window or not visible, the result is undefined.
	yes%	If this operation would overwrite an existing file you are asked if you wish to do this unless yes% is present and non-zero.
	text\$	An optional prompt displayed as part of the file dialog to prompt the user. The File Save As dialog will appear if text\$ is not empty.
	nChan%	Used when $type$ is 0 to set the number of channels in the exported file in the range 32-256. If omitted, the file has the same number of channels as the source file. Spike2 version 4 can read files with up to 100 channels.
	flag%	This is used for result and XY views when t_{ype} is 1 (output as text) to override the normal behaviour, which is to save all channels and all data points. This optional argument has a default value of 0 and is the sum of:
		1 O the text of the strength

- 1 Output only visible channels
- 2 Output only selected channels, or visible channels if no channel is selected

	4 Output only data that is in the visible range (use 5 for visible data) Returns The function returns 0 if the operation was a success, or a negative error code.
Use with a newly sampled time view	This command can save a time view that has just been sampled as described above for $type$ % set to -1. However, if automatic filing is set with SampleAutoFile() or by the sampling configuration, the file is saved automatically as soon as sampling finishes and using a $type$ % of -1 will generate an error.
	When automatic file naming is enabled with SampleAutoName\$() or with the sampling configuration, you can use this command to override the next sequential name, as long as the file has not already been saved and has finished sampling.
	<pre>See also:EditCopy(), ExportChanFormat(), ExportChanList(),</pre>
FileSaveResource()	This saves a resource file for the current time window. If the window is duplicated, all the duplicates are saved.
	<pre>Func FileSaveResource({name\$})</pre>
	name\$ The resource file to save. You must include the .s2r file extension and required path. Remember that \ in a string must be entered as \\ or use /. If the name is "" or contains a "*" or a "?", the user is prompted for a file.
	If this argument is omitted, the current resource file that is in use is updated. If there is no current resource file, one is created to match the file name.
	Returns 0 if all was OK, -1 if the user cancelled the File Save dialog, -2 if the file could not be saved for any other reason.
	See also:FileGlobalResource(),FileApplyResource(),FileSave()
FileTime\$()	Since version 4.03, the time and date at which sampling started is saved in the data file. This function returns a string holding the time. Use FileTimeDate() to get the date as numbers. The current view must be a time view. The arguments are exactly the same as for the Time\$() command. If there is no time stored, the result is an empty string.
	<pre>Func FileTime\$({tBase%, {show%, {amPm%, {sep\$}}});</pre>
	<pre>See also:Time\$(), FileDate\$(), FileTimeDate()</pre>
FileTimeDate()	Since version 4.03, the time and date at which sampling started is saved in the data file. This function returns the time and date as numbers. Use FileTime\$() and FileDate\$() to get the result as strings. The current view must be a time view. The arguments are exactly the same as for the TimeDate() command. If there is no date stored, the returned values are all zero.
	<pre>Proc FileTimeDate(&s%, {&m%, {&h%, {&d%, {&mon%, {&y%, {&wDay%}}}}}); Proc FileTimeDate(now%[])</pre>
	<pre>See also:Date\$(), MaxTime(), Seconds(), Time\$(), TimeDate()</pre>

FiltApply()

Applies a set of filter coefficients or a filter in the filter bank to a source waveform channel in the current time view and places the resulting waveform in a destination channel. If there is a large amount of data in the source channel, you should filter directly to a disk based channel, rather than a memory channel to avoid running out of memory.

Each output point is generated from the same number of input points as there are filter coefficients. Half these points are before the output point, and half are after. Where more data is needed than exists in the source file (for example at the start and end of a file and where there are gaps), extra points are made by duplicating the nearest valid point.

Func FiltApply(n%|coef[], dest%, srce%, sTime, eTime {,scale%});

- n% Index of the filter in the filter bank to apply in the range 0-11, or
- coef[] An array holding a set of FIR filter coefficients to apply to the waveform. The size of the passed array sets the number of coefficients.
- dest% The channel to hold the filtered waveform: either an unused disk channel, a memory channel with the same sampling frequency as srce% or 0 to create a compatible memory channel and place the filtered waveform in it. When a new channel is created, the channel settings are copied from the old channel.
- srce% The source waveform channel. There must be at least half the number of sampling coefficients worth of data points before sTime if the output is to start at sTime. Similarly, the channel must extend for the same number of data points beyond eTime if the output is to extend to eTime.
- sTime Time to start the output of filtered data. There is no output for areas where there is no input data. If the filter has an even number of coefficients, the output is shifted by half a sample relative to the input.
- eTime The end of the time range for filtered data.
- scale% If present and non-zero, the output scale and offset values are optimised to give the best possible representation of the filtered data, but filtering takes twice as long. If there is data in the channel already, this data may lose some of its accuracy through this process. If you do not re-scale, the channel's scale and offset are unchanged and you run a slight risk of the waveform going outside the 16-bit range allowed for a waveform channel.
- Returns The channel number that the output was written to or a negative error code. A negative error code is also returned if the user clicks Cancel from the progress bar that may appear during a long filtering operation or if dest% is a disk channel that is in use. Delete an existing channel with ChanDelete(dest%).

See also:ChanDelete(), FiltAtten(), FiltCalc(), FiltComment\$(), FiltCreate(), FiltInfo(), FiltName\$(), FiltRange()

FiltAtten()

This set the desired attenuation for a filter in the filter bank. When FiltApply() or FiltCalc() is used, the number of coefficients needed to achieve this attenuation will be generated. A value of zero sets the attenuation back to the default (-65 dB).

Func FiltAtten(index%{, dB});

index% Index of the filter in the filter bank to use in the range 0-11.

dB If present and negative, this is the desired attenuation for stop bands in the filter.

Returns The desired attenuation for a filter at the time of the call.

FiltCalc()	ne calculation of filter of lculation of a filter for you do not force the owever, the coefficient hich may not be desirab	coefficients can take an appreciable time. This routine forces the a particular sampling frequency if it has not already been done. calculation, you can still use FiltApply() to apply a filter. calculation will then be done at the time of filter application, de if the filtering operation is time critical.
	nc FiltCalc(index	<pre>%, sInt{, coeff[]{, &dBGot {, nCoef%}}});</pre>
	ndex% Index of the filte	r in the filter bank to use in the range 0-11.
	Int% The sample inte returned by Bin	rval of the waveform you are about to filter. This is the value Size() for a waveform channel.
	beff An array to be small, as many e	filled with the coefficients used for filtering. If the array is too elements as will fit are set. The maximum size needed is 511.
	Got If present, set to	the attenuation in negative dB attained by the filter coefficients.
	Coef% If present, sets number for a ful	he number of coefficients used in the calculation (use an even l differentiator and an odd number for all other filter types).
	to FiltApply (coefficients generated by the filter. If you are passing the result), use coef[:n], not coef[] to pass n coefficients.
An example	uppose the first filter in Hz. If we know that we ter, we may want to cal	the bank (index 0) is a low pass filter with the pass band edge at ve will need to filter a channel 4 (sampled at 200 Hz) with this culate the coefficients needed in advance:
	ltCalc(0, BinSize	(4));
	his will calculate a filt equency of 200 Hz wit tenuation value for this	er corresponding to the specification of filter 0 for a sampling h an attenuation in the stop band of at least the current desired filter.
Constraints on filters	ne calculation of coeffi pating point rounding e eful result there is a lir e sampling frequency.	cients is a complex process and can produce silly results due to errors in some situations. To ensure that you will always get a nit to how small and how big a transition gap can be, relative to There is a similar limit on the width of a pass or stop band:
	The transition gap and of the sampling freque	d the width of a pass or stop band cannot be smaller than 0.0025 ency.
	The transition gap car	not be larger than 0.12 of the sampling frequency.
	nis function always cale order to do it (these ch	culates a set of coefficients, but may alter the filter specification anges are temporary, see later). This can happen in two cases:
	If the sampling freque pass and stop band we before calculating the with respect to a 12 $12000*0.0025 = 30$ H	ency is such that to produce the filter, the transition gap and/or dths are outside their limits, then the widths are set to the limits filter. In our 50 Hz low pass filter example, if we calculate it kHz sampling frequency, the minimum pass band width is z. So, the filter would be changed to a 60 Hz low pass filter.
	If half the sampling f or stop band, certain example, if we tried t that the Nyquist frequ a waveform sampled a filter. The filter specifi	requency (the Nyquist frequency) is less than an edge of a pass attributes of the filter are lost. In our 50 Hz low pass filter o calculate with a sampling frequency of 80 Hz, we would see ency is 40 Hz. No frequency above 40 Hz can be represented in at 80 Hz, so a 50 Hz low pass filter is equivalent to an "All pass" ication will be altered to reflect this before calculating.
	ny changes made to a ade with reference to the reference to the range of the range	filter specification to accommodate a particular calculation are ne original specification, not the specification that was last used
	e also:FiltApply(), FiltInfo(),F	<pre>FiltAtten(),FiltComment\$(),FiltCreate(), iltName\$(),FiltRange()</pre>

FiltComment\$()	This func	tion gets and sets the comr	nent ass	ociated with	n a filter in t	he filter bar	ık.
	Func Fi	ltComment\$(index% {,	new\$]);			
	indoxe	Index of the filter in the filt	ter hank	to use in th	e range 0-11		
	Index 6	If present sets the new een					
	newş .	The present, sets the new con					
	Returns	The previous comment for	the filte	r at the inde	Х.		
	See also:	FiltApply(), FiltAtte	n(),Fi o()	ltCalc(),	FiltCrea	te(),Filt	:Info(),
			C14 1			: (; ;	1
FiltCreate()	standard	name and comment.	inter ba	nk to the su	ppned spec	incation an	d gives it a
	Func Fi	ltCreate(index%, typ	pe%{, t	rW{, edge	el{, edge	2{,}	}});
	index%	Index of the new filter in t any existing filter at this in	the filter dex.	bank in the	e range 0-1	1. This action	on replaces
	type%	The type of the filter desire	ed (see ta	able).			
	trW	The transition width of the of a stop band and the edge	filter. T e of the a	his is the fro adjacent pas	equency inte s band.	erval betwee	en the edge
	edgen This is a list of edges of pass bands in Hz. See the table.						
	Returns	0 if there was no problem of	or a nega	ative error co	ode if the fil	lter was not	created.
	This tabl correspor is more th	e shows the relationship b ading arguments. The number an 1. An empty space in the	between bers in t he table	different fill brackets indi- means that	lter types an icate the nth the argument	nd the mean pass band nt is not requ	ning of the when there uired.
	type%	Name	trW	edge1	edge2	edge3	edge4
	0	All stop					
	1	All pass					
	2	Low pass	Yes	High			
	3	High pass	Yes	Low	11: 1		
	4	Band pass	Yes	LOW High(1)	Hign		
	6	Low pass differentiator	Ves	High	LOW(2)		
	7	Differentiator	105	mgn			
	8	1.5 Band Low pass	Yes	High(1)	Low(2)	High(2)	
	9	1.5 Band High pass	Yes	Low(1)	High(1)	Low(2)	
	10	2 Band pass	Yes	Low(1)	High(1)	Low(2)	High(2)
	11	2 Band stop	Yes	High(1)	Low(2)	High(2)	Low(3)

The values entered correspond to the text fields shown in the Filter edit dialog box.

See also:FiltApply(), FiltAtten(), FiltCalc(), FiltComment\$(), FiltInfo(), FiltName\$(), FiltRange()

FiltInfo()

Retrieves information about a filter in the bank.

Func FiltInfo(index%{, what%});

index% Index of the filter in the filter bank to use in the range 0-11.

what% Which bit of information about the filter to return:

- -2 Maximum what% number allowed
 - -1 Desired attenuation
 - 0 type (if you supply no value, 0 is assumed)

	<pre>1 Transition width 2-5 edge1-edge4 given in FiltCreate()</pre>
	Returns The information requested as real number.
	<pre>See also:FiltApply(), FiltAtten(), FiltCalc(), FiltComment\$(), FiltCreate(), FiltName\$(), FiltRange()</pre>
FiltName\$()	This function gets and/or sets the name of a filter in the filter bank.
	index% Index of the filter in the filter bank to use in the range 0-11.
	new\$ If present, sets the new name. Returns. The previous name of the filter at that index
	<pre>See also:FiltApply(), FiltAtten(), FiltCalc(), FiltComment\$(), FiltInfo(), FiltCreate(), FiltRange()</pre>
FiltRange()	Retrieves the minimum and maximum sampling rates that this filter can be applied to without the specification being altered. See the FiltCalc() command, <i>Constraints on filters</i> for more information.
	<pre>Proc FiltRange(index%, &minFr, &maxFr);</pre>
	index% Index of the filter in the filter bank to use in the range 0-11.
	minFr Returns the minimum sampling frequency you can calculate the filter with respect to, so that no transition width is greater than the maximum allowed and that no attributes of the filter are lost.
	maxFr Returns the maximum sampling frequency you can calculate the filter with respect to without the transition (or band) widths being smaller than allowed.
	It is possible to create a filter that cannot be applied to any sampling frequency without being changed. This will be apparent because minFr will be larger than maxFr.
	<pre>See also:FiltApply(), FiltAtten(), FiltCalc(), FiltComment\$(), FiltInfo(), FiltCreate(), FiltName\$()</pre>

FIRMake()

This function creates FIR filter coefficients and places them in an array ready for use by ArrFilt(). This command is very similar in operation to the DOS program FIRMake and has similar input requirements. Unless you need precise control over all aspects of filter generation, you may find it easier to use FiltCalc() or FIRQuick(). You will need to read the detailed information about FIR filters in the description of the Digital Filter dialog to get the best results from this command.

Proc FIRMake(type%, param[][], coef[]{, nGrid{, extFr[]}});

- The type of filter file to produce: 1=Multiband filter, 2=Differentiator, 3=Hilbert type% transformer, 4=Multiband pink noise (Multiband with 3 dB per octave roll-off).
- This is a 2-dimensional array. The size of the first dimension must be 4 or 5. The param size of the second dimension (n) should be the number of bands in your filter. You pass in 4 values for each band (indices 0 to 3) to describe your filter:

Indices 0 and 1 are the start and end frequency of each band. All frequencies are given as fraction of a sampling frequency and so are in the range 0 to 0.5.

Index 2 is the function of the band. For all filter types except a differentiator, this is the gain of the filter in the band in the range 0 to 1; the most common values are 0 for a stop band and 1 for a pass band. For a differentiator, this is the slope of the filter in the band, normally not more than 2. The gain at any frequency f in the band is given by f^* function.

Index 3 is the relative weight to give the band. The weight sets the relative importance of the band in multiband filters. The program divides each band into frequency points and optimises the filter such that the maximum ripple times the weight in each band is the same for all bands. The weight is independent of frequency, except in the case of the differentiator, where the weight used is weight/frequency.

If there is an index 4 (the size of the first dimension was 5), this index is filled in by the function with the ripple in the band in dB.

- coef An array into which the FIR filter coefficients are placed. The size of this array determines the number of filter coefficients that are calculated. It is important, therefore, to make sure this array is exactly the size that you need. The maximum number of coefficients is 511.
- nGrid The grid density for the calculation. If omitted or set to 0, the default density of 16 is used. This sets the density of test points in internal tables used to search for points of maximum deviation from the filter specification. The larger the value, the longer it takes to compute the filter. There is seldom any point changing this value unless you suspect that the program is missing the peak points.
- extFr An array to hold the list of extremal frequencies (the list of frequencies within the bands which have the largest deviation from the desired filter). If there are n coefficients, there are (n%+1)/2 extremal frequencies.

The parameters passed in must be correct or a fatal error results. Errors include: overlapping band edges, band edges outside the range 0 to 0.5, too many coefficients, differentiator slope less than 0. If the filter is not a differentiator the band function must lie between 0 and 1, the band weight must be greater than 0.

For example, to create a low pass filter with a pass band from 0 to 0.3 and a stop band from 0.35 to 0.5, and no return of the ripple, you would set up param as follows:

'No return of ripple, 2 bands var param[4][2] para[0][0] := 0; 'Starting frequency of pass band para[1][0] := 0.3; 'Ending frequency of pass band para[2][0] := 1; 'Desired gain (unity) 'Give this band a weighting of 1 para[3][0] := 1; para[0][1] := 0.35; 'Starting frequency of stop band 'Ending frequency of stop band para[1][1] := 0.5; para[2][1] := 0; 'Desired gain of 0 (stop band) para[3][1] := 10; 'Give this band a weighting of 10 See also:ArrFilt(), FiltApply(), FiltCalc(), FIRQuick(), FIRResponse()

FIRQuick()

This function creates a set of filter coefficients in the same way the FIRMake() does, but many of the parameters are optional, allowing the most common filters to be created with a minimal specification.

Func FIRQuick(coef[], type%, freq {, width {, atten}});

coef An array into which the FIR filter coefficients are placed. The size of this array should be 511 (was 255 in versions before 4.01). This is the maximum number of coefficients that can be created and this function reserves the right to return as many as it feels necessary up to that value to create a good filter.

type%	This sets the type of filter to create. 0=Low pass, 1=High pass, 2=Band pass,
	3=Band stop and 4=Differentiator.

freq This is a fraction of the sampling rate in the range 0 to 0.5 and means different things depending on the type of filter:

For Low pass, High pass and Differentiator types, this represents the cut-off frequency. This is the frequency of the higher edge of the first frequency band.

For Band pass and Band stop filters, this is the midpoint of the middle frequency band: the pass band in a Band pass filter, the stop band in a Band stop filter.

width For Low pass, High pass and Differentiator filters, this is the width of the transition gap between the stop band and the pass band. The default value is 0.02 and there is an upper limit of 0.1 on this argument.

For a Band pass, or Band stop filter, width is the width of the middle band. E.g. if you ask for a Pass band filter with the freq parameter to be 0.25 and the width to be 0.05, the middle pass band will be from 0.2 to 0.3. For these types of filter, you still need a positive transition width. This transition width is 0.02 and cannot be changed by the user.

- atten The desired attenuation in the stop band in dB. The default is 50 dB. This is analogous to the desired attenuation in the FiltAtten() command.
- Returns The number of coefficients calculated. If the array is not large enough the coefficient list is truncated (and the result is useless). If you are passing the result to FiltApply(), use coef[:n], not coef[] to pass n coefficients.

See also:ArrFilt(), FiltApply(), FiltCalc(), FIRMake(), FIRResponse()

FIRResponse()

This function retrieves the frequency response of a given filter as amplitude or in dB.

Func FIRResponse(resp[], coef[]{, as%{, type%}});

- resp The array to hold the frequency response. This array will be filled regardless of its size. The first element is the amplitude response at 0 Hz and the last is the amplitude response at the Nyquist frequency. The remaining elements are set to the response at a frequency proportional to the element position in the array.
- coef The coefficient array calculated by FIRMake(), FIRQuick() or FiltCalc().
- as% If this is 0 or omitted, the response is in dB (0 dB is unchanged amplitude), otherwise as linear amplitude (1.0 is unchanged).
- type% If present, informs the command of the filter type. The types are the same as those supplied for FIRQuick(): 0=Low pass, 1=High pass, 2=Band pass, 3=Band stop and 4=Differentiator. If a type is given, the time to calculate the response is halved. If you are not sure what type of filter you have, or you have type not covered by the FIRQuick() types, then do not supply a type to this command.

See also:ArrFilt(), FiltCalc(), FIRMake(), FIRQuick()

FitCoef()

This command gives you access to the fit coefficients for the next FitData() fit. You can return the values from any type of fit and set the initial values and limits and hold values fixed for iterative fits. There are two command variants:

Set and get coefficients

This command variant lets you read back the current coefficient values and set the coefficient values and limits for iterative fitting:

Func FitCoef({num%{, new{, lower{, upper}}});

	num%	If this is omitted, the return value is the number of coefficients in the current fi If present, it is a coefficient number. The first coefficient is number 0. If num% is present, the return value is the coefficient value for the existing fit, or if there is no fit, the coefficient value that would be used as the starting point for the nex- iterative fit is returned.
	new	If present, this sets the value of coefficient num% for the next iterative fit.
	lower	If present, this sets the lower limit for coefficient num% for the next iterative fr There is currently no way to read back the coefficient limits. There is also n check made that the limits are set to sensible values.
	upper	If present, this sets the upper limit for coefficient num% for the next iterative fit.
	Returns	The number of coefficients or the value of coefficient num%.
Get and set the hold flags	This control the next	nmand variant allows you to hold some coefficients at their current values durin fit.
	Func F	<pre>itCoef(hold%[]);</pre>
	hold%	An array of integers to correspond with the coefficients. If the array is too long extra elements are ignored. If it is too short, extra coefficients are not affected Set hold%[i%] to 1 to hold coefficient i% and to 0 to fit it. If hold%[i%] is less than 0, the hold state is not changed, but hold%[i%] is set to 1 if the corresponding coefficient is held and to 0 if it is not held.
	Returns	This always returns 0.
	See also	:FitData(),FitValue(),FitExp(),ChanFitCoef()
FitData()	This fu fitting f x and y point y	nction, together with FitCoef() and FitValue(), lets you apply the sam unctions that are available for channel data to data in arrays. You supply arrays of data points and an optional array holding the standard deviation of the input dat values. There are three command variants:
Initialise fit information	The firs fitting c are clea a guess	t variant sets the type of fit. If you select an iterative fit, the initial values of th oefficients are reset to standard values and any "hold" flags set by FitCoef(red. You can set your own initial values with the FitCoef() command or mak at the initial values when performing the fit.
	Func F	<pre>itData(type%, order%);</pre>
	type%	The fit type. 0=Clear any fit, 1=Exponential, 2=Polynomial, 3=Gaussian 4=Sine, 5=Sigmoid.
	order%	If positive, this is the order of the fit, if negative it is minus the number of fitte coefficients. See the information about each fit for the allowed values for each fit type. If type% is 0 this argument is ignored and should be 0.
	Returns	The number of fit coefficients for the fit or a negative error code.
Exponential fit	This fits	multiple exponentials by an iterative method. The data is fitted to the equation:
	$y = \mathbf{a}_0 \mathbf{e}_1$ $y = \mathbf{a}_0 \mathbf{e}_2$	$ \begin{array}{l} \exp(-x/\mathbf{a}_1) + \mathbf{a}_2 \exp(-x/\mathbf{a}_3) \dots \\ \exp(-x/\mathbf{a}_1) + \mathbf{a}_2 \exp(-x/\mathbf{a}_3) \dots + \mathbf{a}_n \end{array} \end{array} $ For an even number of coefficient For an odd number of coefficient
	You can coeffici Coeffici	a set up to 10 coefficients or orders 1 to 5. If you use a fit order, the number of ents is the order times 2. See the $FitExp()$ command for more information ent estimates are effective for orders 1 and 2.

- **Polynomial fit** This fits $y = \mathbf{a}_0 + \mathbf{a}_1 x + \mathbf{a}_2 x^2 + \mathbf{a}_3 x^3 \dots$ to a set of (x,y) data points. The fitting is by a direct method; there is no iteration. The fit order is the highest power of x to fit in the range 1 to 10. The number of coefficients is the fit order plus 1.
 - *Gaussian fit* This fits multiple Gaussians by an iterative method. The data is fitted to the equation:

 $y = \mathbf{a}_0 \exp(-\frac{1}{2}(x-\mathbf{a}_1)^2/\mathbf{a}_2^2) + \mathbf{a}_3 \exp(-\frac{1}{2}(x-\mathbf{a}_4)^2/\mathbf{a}_5^2) + \dots$

The fitted parameters (coefficients) are the \mathbf{a}_i . You can fit up to 3 Gaussians (order 1 to 3). The number of coefficients is given by the fit order times 3. Coefficient estimates become less useful as the order increases.

Sine fit This fits multiple sinusoids by an iterative method. The data is fitted to the equation:

 $y = \mathbf{a}_0 \sin(\mathbf{a}_1 x + \mathbf{a}_2) + \mathbf{a}_3 \sin(\mathbf{a}_4 x + \mathbf{a}_5) + \dots \{+ \mathbf{a}_{3n}\}$

The fitted coefficients are the \mathbf{a}_i . Angles are evaluated in radians. You can fit up to 3 sinusoids (order 3) and an optional offset. Although the function is given in terms of sine functions, you can easily convert to cosines by subtracting $\pi/2$ from the phase angle (\mathbf{a}_2 , \mathbf{a}_5 , \mathbf{a}_8) after the fit. The coefficient count can be set to 3, 4, 6, 7, 9 or 10. If you use orders, the number of coefficients is order times 3 and you cannot set an offset. A useful coefficient estimate is made for a single sinusoid fit.

Sigmoid fit This fits a single Boltzmann sigmoid by an iterative method. The data is fitted to the equation:

 $y = \mathbf{a}_0 + (\mathbf{a}_1 - \mathbf{a}_0)/(1 + \exp((\mathbf{a}_2 - x)/\mathbf{a}_3))$

In terms of the fitted result, \mathbf{a}_0 and \mathbf{a}_1 are the low and high fitting limits, \mathbf{a}_2 is the X50 point and \mathbf{a}_3 is a measure of the slope at the X50 point. You can set order 1 only, or 4 coefficients.

Perform the fit This variant of the command does the fit set by the previous variant. Use the FitCoef() command to preset fit coefficients and to read back the result of the fit.

- opt% 1=Estimate the coefficients before fitting, 0=use current values. Note that the estimates are usually only useful for a small number of coefficients.
- y[] An array of y values to be fitted.
- x[] A corresponding array of x values.
- s[] | s A corresponding array of standard deviations for the data points defined by y[] and x, or a single value, being the standard deviation of each point. If this value is omitted or set to 1.0, the result is a least squares fit. If standard deviations are supplied, the result is a chi-squared fit.
- err If present, this optional variable is updated with the chi-squared or least-squares error between the fit and the data.
- max1% If present, this changes the maximum number of iterations from 100.
- iTer% If present, this integer variable is updated with the count of iterations done.
- covar An optional two dimensional array of size at least [nCoef][nCoef] that is returned holding the covariance matrix when the fit is complete. It is changed if the return value is -1, 0 or 1. However, the values it contains are probably not useful unless the return value is 0.
- Returns 0 if the fit is complete, 1 if max iterations done, or a negative error code: -1=the fit is not making progress (results may be OK), -2=the fit failed due to a singular matrix, -5=the fit caused a floating point error, -6=too little data for the number

of coefficients, -7=unknown fitting function, -8=ran out of memory during the fit (too many data points), -9=the fit was set up with bad parameters.

Get fit information This variant of the command returns information about the current fit.

Func FitData({opt%});

opt% This determines what to return. opt% values match ChanFit(), where possible. The returned information for each value of opt% is:

opt%	Returns	opt%	Returns
0	Fit type of next fit	1	Fit order of next fit
-1	1=a fit exists, 0=no fit exists	-8	Not used
-2	Type of last fit or 0	-9	Not used
-3	Number of coefficients	-10	Not used
-4	Chi or least-squares error	-11	1=chi-square, 0=least-square
-5	Fit probability (estimated)	-12	Last fit result code
-6	Lowest x value fitted	-13	Number of fitted points
-7	Highest x value fitted	-14	Number of fit iterations used

Returns The information requested by the opt% argument or 0 if opt% is out of range.

See also:ChanFit(), FitCoef(), FitValue()

FitExp()

This command fits multiple exponentials to arrays of x,y data points, with an optional weight for each point. Fitting is by an iterative method. The data is fitted to the equation:

 $y = \mathbf{a}_0 \exp(-x/\mathbf{a}_1) + \mathbf{a}_2 \exp(-x/\mathbf{a}_3)...$ For an e y = $\mathbf{a}_0 \exp(-x/\mathbf{a}_1) + \mathbf{a}_2 \exp(-x/\mathbf{a}_3)... + \mathbf{a}_n$ For an

For an even number of coefficients For an odd number of coefficients

The fitted coefficients are the \mathbf{a}_i ; odd numbered \mathbf{a}_i are assumed to be positive. You can fit up to 5 exponentials or 4 exponentials and an offset. However, experience shows that trying to fit more than two exponentials requires care. The fit from even two exponentials should be viewed with caution, especially if the odd coefficients are similar. The commands to implement this are:

Set up the problem The first command sets the number of exponentials to fit and the data set to be fitted. You must call this function before you call any of the others.

func FitExp(nCoef%, y[], x[]{, s[]|s});

- nCoef[®] The number of coefficients to fit in the range 2 to 10. If this is even, the first form of the function above is used. If it is odd, the final coefficient is an offset.
- y[] An array of y data values. The length of the array must be at least nCoef%.
- x[] An array of x data values. The length of the array must be at least nCoef%.
- An optional argument which is either an array with one value for each y data point or a single value for all data points. If the value is the standard deviation (σ), then the error value returned when you iterate to find the best fit is the χ^2 value and the fit is a chi-squared fit.

If this value is proportional to the error in the y values, then the fit is still a chisquared fit, and the error returned is proportional to the chi-squared value. If you omit this array, the fit is a least-squares error fit, and the error value returned is the sum of the squares of the errors in the y values.

Returns The function returns 0. There is no other return value as all errors stop the script.

The number of data points is set by the smallest of the sizes of the y[], x[] or s[] (if present) arrays. The number of data points must be at least the number of coefficients.

Set coefficient values and ranges	This van the rang before th	riant of the function sets the initial value of each coefficient and optionally sets e of allowed values. You can call this function at any time after the setup call and he iterate call (below) has returned 0 indicating that the fit is completed.		
	It is imp you fit is constant overlapp wander damage convolu bit. Let	bortant that you give reasonable initial values for the coefficients, especially when more than one exponent. You should limit the odd coefficient values (the time is) so that they cannot be zero and make sure that multiple exponents do not have bing ranges. If two exponents get similar values, the fit is degenerate and will around forever without getting anywhere. However, setting too rigid a range may the fitting process as sometimes the minimisation process has to follow a ted n-dimensional path to reach the goal, and the path may need to wander quite a experience be your guide.		
	If you d values. does, we	o not give starting values, the command will make a simplistic guess at the fitting As we expect that you know more about the "right" answer than the command e suggest that you set the values you want.		
	func F	<pre>itExp(coef%, val{, lo, hi});</pre>		
	coef%	The coefficient to set. The first coefficient is number 0, the last is nCoef%-1.		
	val	The initial value to assign to the coefficient. If you have set low and high limits, and the value is outside these limits, it is set to the nearer limit.		
	lo,hi	If present, these two values set the acceptable range of values for this coefficient. If omitted, or if both values are set to the same value, there is no limit. The value of the coefficient is tested against these limits after each iterative step, and if it exceeds a limit, it is set to that limit.		
Iterate to a solution	Once yo start the values the	bu have set up the problem and given initial values to your coefficients, you can iteration process that will move the coefficients from their starting values to new hat minimise the error (optionally scaled by the s argument).		
	<pre>func FitExp(a[], &err{, maxI%{, &iTer%{, covar[][]}}});</pre>			
	a[]	An array of size at least nCoef [§] that is returned holding the current set of coefficient values. The first amplitude is in a [0], the first exponent in a [1], the second amplitude in a [2], the second exponent in a [3] and so on.		
	err	A real variable returned as the sum over the data points of $(y_{x[i]}-y[i])^2/s[i]^2$ if $s[]$ is used or holding the sum of $(y_{x[i]}-y[i])^2$ if $s[]$ is not used, where $y_{x[i]}$ is the value predicted from the coefficients at the x value $x[i]$.		
	max1%	This is the maximum number of fitting iterations to do before returning from the function. If you omit this value, the function sets 1000. You can set any value from 1 to 10000. If you set more than 10000, the number is limited to 10000.		
	iTer%	An optional integer variable that is returned holding the number of iterations done before the function returned.		
	covar	An optional two dimensional array of size at least $[nCoef][nCoef]$ that is returned holding the covariance matrix when the fit is complete. It is not changed unless the function return value is 0 or -1.		
	Returns 1 if the number of iterations are exhausted and the fit has not conve fit has converged, -1 if the fit is not improving but the result may b negative numbers indicate failure.			
Select coefficients to fit	Normall coefficie until the	y the command fits all the coefficients, but you can use this variant to select the ents to fit. You can use this command at any time after you have set the problem iteration variant returns 0 or -1.		

func FitExp(fit%[]);

fit%[] An array of at least nCoef% integers. If fit%[i] is 0, coefficient i is held constant, otherwise it is fitted. If all elements are 0, all arguments are fitted.

The effect of this command persists until either the iteration variant returns a value less than 1, or you set up a new problem, or you call this variant again.

An example This is a template for using these commands to fit all the coefficients:

	const nData%:=50;	'set number of data elements					
	<pre>var x[nData%], y[nData%];</pre>	'space for our arrays					
	<pre>var coefs[4],err;</pre>	'coefficients and squared error					
		'in here goes code to get the data					
	<pre>FitExp(4, y[], x[]);</pre>	'fit two exponentials (no sigma array)					
	FitExp(0, 1.0, 0.2, 4);	'set first amplitude and limit range					
	FitExp(1, .01, .001, .03);	'set first time constant and range					
	FitExp(2, 2.0, 0.1, 6);	'set second amplitude and limit range					
	FitExp(3, .08, .03, .15);	'set second time constant and range					
	repeat						
	DrawMyData(coets[], x[]	, y[]); 'Some function to show progress					
	until FitExp(coeis[], err,	1) < 1;					
	<pre>DrawMyData(coefs[], x[], y[]); 'Show the final state</pre>						
	<pre>See also:ChanFit(), FitGauss(), FitLinear(), FitNLUser(), FitPoly(),</pre>						
FitGauss()	This command fits multiple Gaussi point. Fitting is by an iterative met	ians to x,y data points with an optional weight for each hod. The input data is fitted to the equation:					
	$x = a \exp(\frac{1}{(x + a)^2/a^2}) + a \exp(\frac{1}{(x + a)^2/a^2})$	$1/(1-2)^2/2^2$					
	$y = \mathbf{a}_0 \exp(-\gamma_2(x - \mathbf{a}_1) / \mathbf{a}_2) + \mathbf{a}_3 \exp(-\gamma_2(x - \mathbf{a}_1) / \mathbf{a}_2)$	$-7_2(x-a_4)/a_5 + \dots$					
	The fitted parameters (coefficient commands to implement this are:	ts) are the \mathbf{a}_i . You can fit up to 3 Gaussians. The					
	Ĩ						
Sat up the problem	The first command gots the number	or of Coursians to fit and the data set to be fitted. You					
Set up the problem	must call this function before you c	call any of the others.					
	<pre>func FitGauss(nCoef%, y[],</pre>	x[]{, s[] s});					
	- Coof [®] The number of coofficien	ts to fit. The legal values are 2, 6 and 0 for one two					

- nCoefs The number of coefficients to fit. The legal values are 3, 6 and 9 for one, two and three Gaussians.
- $y [\] \qquad \text{An array of } y \text{ data values. The length of the array must be at least nCoef8}.$
- $x [] \qquad \text{An array of x data values. The length of the array must be at least nCoef%}.$
- S An optional argument which is either an array with one value for each y data point or a single value for all data points. If the value is the standard deviation (σ) , then the error value returned when you iterate to find the best fit is the χ^2 value and the fit is a chi-squared fit.

If this value is proportional to the error in the y values, then the fit is still a chisquared fit, and the error returned is proportional to the chi-squared value. If you omit this array, the fit is a least-squares error fit, and the error value returned is the sum of the squares of the errors in the y values.

Returns The function returns 0. There is no other return value as all errors stop the script.

The number of data points is set by the smallest of the sizes of the y[], x[] or s[] (if present) arrays. The number of data points must be at least the number of coefficients. If it is not, you will get a fatal error, so check this before calling the function.

Set coefficient values and ranges	This variant of the function sets the initial value of each coefficient and optionally set the range of allowed values. You can call this function at any time after the setup call a before the iterate call (below) has returned 0 indicating that the fit is completed.		
	In this t coefficie each dis they car Gaussia at the fin commar	ype of fitting it is very important that you give reasonable starting values for the ents. In particular, when fitting multiple Gaussians, it is usual that the centre of stribution is easy to determine. If you can set the centres and limit them so that not overlap, the fit usually will proceed without any problems, even for multiple ns. If you do not give starting values, the command will make a simplistic guess tting values. As we expect that you know more about the "right" answer than the nd does, we suggest that you set the values you want.	
	func F	itGauss(coef%, val{, lo, hi});	
	coef%	The coefficient to set. The first coefficient is number 0, the last is $nCoef$ -1.	
	val	The initial value to assign to the coefficient. If you have set low and high limits, and the value is outside these limits, it is set to the nearer limit.	
	lo,hi	If present, these two values set the acceptable range of values for this coefficient. If omitted, or if both values are set to the same value, there is no limit. The value of the coefficient is tested against these limits after each iterative step, and if it exceeds a limit, it is set to that limit.	
	As long problem	g as you make a reasonable estimate of the centre points, there should be no is fitting multiple Gaussians.	
Iterate to a solution	Once yo start the values th	bu have set up the problem and given initial values to your coefficients, you can iteration process that will move the coefficients from their starting values to new hat minimise the error (optionally scaled by the s argument).	
	func F	itGauss(a[], &err{, maxI%{, &iTer%{, covar[][]}}});	
	a[]	An array of size at least $nCoef$ ⁸ that is returned holding the current set of coefficient values. The first amplitude is in a[0], the first centre in a[1], the first sigma in a[2], the second amplitude in a[3] and so on.	
	err	A real variable returned as the sum over the data points of $(y_{x[i]}-y[i])^2/s[i]^2$ if $s[]$ is used or holding the sum of $(y_{x[i]}-y[i])^2$ if $s[]$ is not used, where $y_{x[i]}$ is the value predicted from the coefficients at the x value $x[i]$.	
	max1%	This is the maximum number of fitting iterations to do before returning from the function. If you omit this value, the function sets 1000. You can set any value from 1 to 10000. If you set more than 10000, the number is limited to 10000.	
	iTer%	An optional integer variable that is returned holding the number of iterations done before the function returned.	
	covar	An optional two dimensional array of size at least $[nCoef%][nCoef%]$ that is returned holding the covariance matrix when the fit is complete. It is not changed unless the function return value is 0 or -1.	
	Returns	1 if the number of iterations are exhausted and the fit has not converged, 0 if the fit has converged, -1 if the fit is not improving but the result may be OK. Other negative numbers indicate failure.	
	Remem minimu	ber that even when a minimum is found, there is no guarantee that this is <i>the</i> m. It is the best minimum that this algorithm can find given the starting point.	
Select coefficients to fit	Normall coefficie until the	y the command fits all the coefficients, but you can use this variant to select the ents to fit. You can use this command at any time after you have set the problem iteration variant returns 0 or -1.	
	func F	itGauss(fit%[]);	

fit%[] An array of at least nCoef% integers. If fit%[i] is 0, coefficient i is held constant, otherwise it is fitted. If all elements are 0, all arguments are fitted.

The effect of this command persists until either the iteration variant returns a value less than 1, or you set up a new problem, or you call this variant again.

An example This is a template for using these commands to fit all the coefficients:

```
const nData%:=50;
                             'set number of data elements
                             'space for our arrays
var x[nData%], y[nData%];
var s[nData%];
                             'space for sigma of each point
var coefs[4], err;
                             'coefficients and error squared
                             'in here goes code to get the data
FitGauss(3, y[], x[], s[]); 'fit one gaussian
                             'set amplitude and limit range
FitGauss(0, 1.0, 0.2, 4);
FitGauss(1, 2, 1.5, 2.5);
                             'set centre of the gaussian and range
FitGauss(2, 0.5, 0.3, 1.9); 'set width and limit range
repeat
   DrawMyData(coefs[], x[], y[]); 'Some function to show progress
until FitGauss(coefs[], err, 1) < 1;</pre>
                                    'Show the final state
DrawMyData(coefs[], x[], y[]);
See also:ChanFit(), FitExp(), FitLinear(), FitNLUser(), FitPoly(),
      FitSin()
```

FitLine()

This function calculates the least-squares best-fit line to a set of data points from a time view waveform, RealWave or WaveMark channel or a result view. It fits: $y = \mathbf{m} x + \mathbf{c}$ through the data points (x_i, y_i) so as to minimise the error given by: $\operatorname{Sum}_i(y_i - \mathbf{m} x_i - \mathbf{c})^2$. In this expression, **m** is the gradient of the line and **c** is the y axis intercept when x is 0.

Func FitLine(chan%, start, finish, &grad{, &inter{, &corr}});

- chan% A channel holding suitable data in the current time or result view.
- start The start time for processing in a time view, the start bin in a result view.
- finish The end time for processing in a time view, the end bin in a result view. Data at the finish time, or in the finish bin, is included in the calculation.
- grad This is returned holding the gradient of the best fit line (**m**).
- inter Optional, returned holding the intercept of the line with the y axis (c).
- corr Optional, returned holding correlation coefficient indicating the "goodness of fit" of the line. Values close to 1 or -1 indicate a good fit; values close to 0 indicate a very poor fit. This parameter is often referred to as r in textbooks.

Returns 0 if all was OK, or -1 if there were not at least 2 data points.

The results are in user units, so in a time view with a waveform measured in Volts, the units of the gradient are Volts per second and the units of the intercept are Volts. In a result view, the units are y axis units per x axis unit. They are not y units per bin.

See also:ChanFit(), ChanMeasure(), FitLinear(), FitPoly()

FitLinear()

This command fits $y = \mathbf{a}_0 f_0(x) + \mathbf{a}_1 f_1(x) + \mathbf{a}_2 f_2(x)...$ to a set of (x,y) data points. If you can provide error estimates for each y value, you can use the covariance output from this command to provide confidence limits on the calculated coefficients and you can use the returned χ^2 value to test if the model is likely to fit the data. The command is:

func FitLinear(coef[], y[], x[][]{, s{, covar[][]{, r[]{, mR}}});

- coef[] A real array which sets the number of coefficients to fit and which return the best fit set of coefficients. The array must be between 2 and 10 elements long. The coefficient \mathbf{a}_0 is returned in coef[0], \mathbf{a}_1 in coef[1] and so on.
- y[] A real array of y values.
- x[][] This array specifies the values of the functions f(x) at each data point. If there are *nc* coefficients and *nd* data values, this array must be of size at least [nd] [nc]. Viewed as a rectangular grid with the coefficients running from left to right and the data running from top to bottom, the values you must fill in are:

$f_0(x_0)$	$f_1(x_0)$	$f_2(x_0)$	$f_3(x_0)$		$f_{\text{nc-1}}(x_0)$
$f_0(x_1)$	$f_1(x_1)$	$f_2(x_1)$	$f_3(x_1)$		$f_{\text{nc-1}}(x_1)$
$f_0(x_2)$	$f_1(x_2)$	$f_2(x_2)$	$f_3(x_2)$		$f_{\rm nc-1}(x_2)$
				•••	
$f_0(x_{\rm nd-1})$	$f_1(x_{nd-1})$	$f_2(x_{\rm nd-1})$	$f_3(x_{nd-1})$		$f_{\rm nc-1}(x_{\rm nd-1})$

- s This is either a real array holding the standard deviations of the y[] data points, or a real value holding the standard deviation of all data points. If s is omitted or zero, a least-squares error fit is performed, otherwise a chi-squared fit is done.
- covar An optional two dimensional array of size at least [nc] [nc] (nc is the number of coefficients fitted) that is returned holding the covariance matrix.
- r[] An optional array of size at least [nc] (*nc* is the number of coefficients fitted) that is returned holding diagnostic information about the fit. The less relevant a fitting function f(x) is to the fit, the smaller the value returned. The element of the array that corresponds to the most relevant function is returned as 1.0, smaller numbers indicate less relevance.

If your fitting functions are not independent of each other, several coefficients may have low r values. The solution is to remove one of the functions from the fit, or to set the mR argument to exclude one of the functions, then fit again. If the remaining coefficients become relevant, you have excluded a function that was a linear combination of the others. If the remaining coefficients still are not relevant, you have eliminated a function that did not contribute to the fit.

- mR You can use this optional variable to set the minimum relevance for a function. Functions that have less relevance than this are "edited" out of the fit and their coefficient is returned as 0. If you do not provide this value, the minimum is set to 10⁻¹⁵, which will probably not exclude any values.
- Returns The function returns the chi-square value for the fit if s[] or s is given (and non-zero), or the sum of squares of the errors between the data points and the best fit line if s is omitted or is zero.

The smallest of the sizes of the y[] array (and s[] array, if provided) and the second dimension of x[][] sets the number of data points. It is a fatal error for the number of data points to be less than the number of coefficients.

An example This demonstrates how to fit data to the function y = a*sin(x/10) + b*cos(x/20). The x values vary from 0 to 49 in steps of 1. The function MakeFunc() calculates a trial data set plus noise. We set s to 1.0, so FitLinear() returns the sum of squares of the errors between the fitted and input data. If you run this example, you will notice that the returned value is slightly less than the sum of squares of the added errors.

```
const NCOEF%:=2,NDATA%:=50;' coefficient and data sizes
var x[NDATA%][NCOEF%]; ' array of function information
const noise := 0.01; ' controls how much noise we add
var data[NDATA%], err := 0;' space for our function and errors
' Generate raw data. Fit y = a*sin(x/10)+b*cos(x/20)
var coef[NCOEF%], i%, r; ' coefficients, index, random noise
coef[0]:=1.0; coef[1]:=2; ' set coefficients for generated data
MakeFunc(data[], coef[], x[][]); ' Generate the data, then...
for i%:=0 to NDATA%-1 do ' ...add noise to it
```

```
r := (rand()-0.5)*noise;' the noise to add
data[i%] += r; ' add noise to the data
err := err + r*r; ' accumulate sum of squared noise
   next;
var covar[NCOEF%][NCOEF%]; ' covariance array
                               ' sigma, fitted coefficients
' array for "relevance" values
var sig2, a[NCOEF%];
var rel[NCOEF%];
sig2 := FitLinear(a[], data[], x[][], 1, covar[][], rel[]);
Message("sig^2=%g, err=%g\ncoefs=%g\nrel=%g",sig2,err,a[],rel[]);
halt;
'y is the output array (x values are 0, 1, 2...), a is the array
'of coefficients. Y = a*sin(x/10)+b*cos(x/20)
proc MakeFunc(y[], a[], x[][])
var nd%,v;
                               ' coefficient index, work space
for nd% := 0 to NDATA%-1 do
    v := Sin(nd% / 10.0);
                                ' first funcion
                                ' save the value;
    x[nd%][0] := v;
    y[nd%] := a[0] * v;
                               ' start to build the result
                               ' second function
    v := Cos(nd%/20.0);
                                ' save it
    x[nd%][1] := v;
    y[nd%] += a[1]*v;
                                ' full result
next;
end;
See also:ChanFit(), FitExp(), FitGauss(), FitNLUser(), FitPoly(), FitSin()
```

FitNLUser()	This cor data poi coefficie with res each dat	mmand uses a non-linear fitting algorithm to fit a user-defined function to a set of nts. The function to be fitted is of the form $y = f(x, \mathbf{a}_0, \mathbf{a}_1, \mathbf{a}_2)$ where the \mathbf{a}_i are ents to determine. You must be able to calculate the differential of the function f pect to each of the coefficients. You can optionally supply an array to weight a point. The commands to implement this are:
Set up the problem	The first fit, the 1 point. Y	t command sets the user-defined function, the number of coefficients you want to number of data points and optionally, you can set the weight to give each data ou must call this function before you call any of the others.
	func F	itNLUser(User(ind%, a[], dyda[]), nCoef%, nData%{, s[] s});
	User()	A user-defined function which is called by the fitting routine. The function is passed the current values of the coefficients. It returns the error between the function and the data point identified by ind% and the differentials of the function with respect to each of the coefficients at that point. The return value should be the y data value at the index minus the calculated value of the function at the x value, using the coefficients passed in.
	ind%	The index into the data points at which to evaluate the error and differentials. If there are n data points, ind% runs from 0 to n-1. You can rely on the function being called with the same coefficients as ind% increments from 0 to n-1, which may be useful if you have complex functions of the coefficients to evaluate.
	a	An array of length $nCoef$ holding the current values of the coefficients. The coefficients are refreshed for each call to the user-defined function, so it is not an error to change them; however this is usually not done.
	dyda	An array of length nCoef% which your function should fill in with the values of the partial differential of the function with respect to each of the coefficients. For example, if you were fitting $y = \mathbf{a}_0 \exp(-\mathbf{a}_1 \times \mathbf{x})$ then set dyda[0] = $\delta y/\delta \mathbf{a}_0 = \exp(-\mathbf{a}_1 \times \mathbf{x})$ and dyda[1] = $\delta y/\delta \mathbf{a}_1 = -\mathbf{a}_0 \times \mathbf{a}_1 \times \exp(-\mathbf{a}_1 \times \mathbf{x})$.
	nCoef%	The number of coefficients to fit in the range 1 to 10.

nData%	The number of data points you will be fitting. If s [] is provided as an array, the
	value of nData% used is the smaller of nData% and the length of the s[] array.
	It is a fatal error for the number of data points used to be less than nCoef8.

This argument is optional. It is either an array of weights to be given to each data point in the fit or a single weight to apply to all data points. If this value is the expected standard deviation (σ) of the y value of the data points, then the error value returned is the χ^2 value and the fit is a chi-squared fit. If this value is proportional to the expected error at the data point, then the fit is still a chi-squared fit, and the error returned is proportional to the chi-squared value. If you omit this argument, the fit is a least-squares error fit, and the error value returned is the sum of the squares of the errors in the y values.

Returns The function returns 0. There is no other return value as all errors stop the script.

Unlike the other fitting routines, you will notice that the x and y data values are not passed into the command. Instead, the user-defined function is passed an index to the data values. It is assumed that the data is accessible by the user function.

Due to restrictions in the implementation of the script language, you cannot debug through the user-defined function. If you set a break point in it, or attempt to step into it you will get errors. We recommend that you check the returned values from the user-defined function by calling it from your own script code.

Set coefficient values and ranges This variant of the function sets the initial value of each coefficient and optionally sets the range of allowed values. You can call this function at any time after the setup call and before the iterate call (below) has returned 0 indicating that the fit is completed.

In this type of fitting it is very important that you give reasonable starting values for the coefficients. If you do not give starting values, the command will set them all to zero, which is unlikely to be correct.

func FitNLUser(coef%, val{, lo, hi});

- coef% The coefficient to set. The first coefficient is number 0, the last is nCoef%-1.
- val The initial value to assign to the coefficient. If you have set low and high limits, and the value is outside these limits, it is set to the nearer limit.
- 10, hi If present, these two values set the acceptable range of values for this coefficient. If omitted, or if both values are set to the same value, there is no limit. The value of the coefficient is tested against these limits after each iterative step, and if it exceeds a limit, it is set to that limit.

Iterate to a solution Once you have set up the problem and given initial values to your coefficients, you can start the iteration process that will move the coefficients from their starting values to new values that minimise the error (optionally scaled by the s argument).

func FitNLUser(a[], &err{, maxI%{, &iTer%{, covar[][]}}});

- a[] An array of size at least nCoef% that is returned holding the current set of coefficient values.
- err A real variable returned as the sum over the data points of $(y_{x[i]}-y[i])^2/s[i]^2$ if s[] is used or holding the sum of $(y_{x[i]}-y[i])^2$ if s[] is not used, where $y_{x[i]}$ is the value predicted from the coefficients at the x value x[i].
- max1% This is the maximum number of fitting iterations to do before returning from the function. If you omit this value, the function sets 1000. You can set any value from 1 to 10000. If you set more than 10000, the number is limited to 10000.
- iTer% An optional integer variable that is returned holding the number of iterations done before the function returned.

- covar An optional two dimensional array of size at least [nCoef%] [nCoef%] that is returned holding the covariance matrix when the fit is complete. It is not changed unless the function return value is 0 or -1.
- Returns 1 if the number of iterations are exhausted and the fit has not converged, 0 if the fit has converged, -1 if the fit is not improving but the result may be OK. Other negative numbers indicate failure.

Remember that even when a minimum is found, there is no guarantee that it is *the* minimum. It is the best minimum that this algorithm can find given the starting point.

Select coefficients to fit Normally the command fits all the coefficients, but you can use this variant to select the coefficients to fit. You can use this command at any time after you have set the problem until the iteration variant returns 0 or -1.

func FitNLUser(fit%[]);

fit%[] An array of at least nCoef% integers. If fit%[i] is 0, coefficient i is held constant, otherwise it is fitted. If all elements are 0, all arguments are fitted.

The effect of this command persists until either the iteration variant returns a value less than 1, or you set up a new problem, or you call this variant again.

An example The following is an example of using this set of commands to fit the user-defined function $y = \mathbf{a} + \exp(-\mathbf{b} + \mathbf{x})$. In this example we generate some test data and add to it a random error. There are two coefficients to be fitted (\mathbf{a} and \mathbf{b}).

```
const NDATA%:=100, NCOEF% := 2;
var x[NDATA%],y[NDATA%],i%;
                             ' Generate data
for i% := 0 to NDATA%-1 do
                              ' a:=1, b:=0.05 and add some noise
   x[i%] := i%;
   y[i\%] := exp(-0.05*i\%) + (rand() - 0.5)*0.01;
   next;
FitNLUser(UserFnc, NCOEF%, NDATA%); 'Link user function
FitNLUser(0, 0.5, 0.01, 2); 'Set range of amplitude
FitNLUser(1, 0.01, 0.001, 1); 'Set range of exponent
var coefs[NCOEF%], err, iter%;
i% := FitNLUser(coefs[], err, 100, iter%);
Message("fit=%d, Err=%g, iter=%d, coefs=%g",i%,err,iter%,coefs[]);
halt;
' The user-defined function: y = a * exp(-b*x);
dy/da = exp(-b*x), dy/db = -x * a * exp(-b*x)
func UserFnc(ind%, a[], dyda[])
var xi,yi,r;
xi := x[ind%];
                            ' local copy of x value
                           ' local copy of y value
yi := y[ind%];
dyda[0] := exp(-a[1]*xi); ' differential of y with respect to a
r := dyda[0] * a[0];
                           ' intermediate value
                            ' differential of y with respect to b
dyda[1] := -xi * r;
return yi-r;
end
See also:FitExp(), FitGauss(), FitLinear(), FitPoly(), FitSin()
```

FitPoly()

This command fits $y = \mathbf{a}_0 + \mathbf{a}_1 x + \mathbf{a}_2 x^2 + \mathbf{a}_3 x^3 \dots$ to a set of (x,y) data points. If you can provide error estimates for each y value, you can use the covariance output from this command to provide confidence limits on the calculated coefficients and you can use the returned χ^2 value to test if the model is likely to fit the data. The command is:

func FitPoly(coef[], y[], x[]{, s[]|s{, covar[][]}});

- coef[] A real array that sets the number of coefficients and returns the coefficient values. The array must be between 2 and 10 elements long. The coefficient \mathbf{a}_0 is returned in coef[0], \mathbf{a}_1 in coef[1] and so on.
- y[]A real array of y values. The smaller of the sizes of the x[] and y[] arrays (and s[] array, if provided), sets the number of data points. It is a fatal error for the number of data points to be less than the number of coefficients.
- x[] A real array of x values.
- s This is an optional argument. It is either a real array holding the standard deviations of each of the y[] data points, or it is a real value holding the standard deviation of all of the data points. If the argument is omitted or set to zero, a least-squares error fit is performed, otherwise a chi-squared fit is done.
- covar An optional two dimensional array of size at least [nc] [nc] (nc is the number of coefficients fitted) that is returned holding the covariance matrix.
- Returns The function returns chi-squared if s[] or s is given (and non-zero), otherwise it returns the sum of squares of the errors between the raw and fitted data.

```
An example This example generates test data, adds random noise, then fits a polynomial to the data.
```

```
const NCOEF% := 5;
                          ' number of coefficients
const NDATA%:=50;
                            number of data points
var y[NDATA%], x[NDATA%]; ' space for x and y values for fit
                          ' noise to add
const noise := 1;
                          ' will be sum of squares of added noise
var err := 0.0;
var cf[NCOEF%], i%, r;
cf[0]:=1.0; cf[1]:=-80; cf[2]:=-2.0; cf[3]:=0.5; cf[4]:=-0.009;
MakePoly(cf[],x[],y[]);
                         ' generate ideal data as polynomial
for i%:=0 to NDATA%-1 do ' now add some noise to it
  r := (rand()-0.5)*noise;
                          ' add noise to the data
  y[i%] += r;
                          ' sum of squares of added noise
  err += r*r;
  next;
                          ' a[] will be the fitted coefficients
var sig2, a[NCOEF%];
sig2 := FitPoly(a[], y[], x[]);
Message("sig2=%g, noise=%g\nfitted=%8.4f\nideal =%8.4f",
                                         sig2, err, a[], cf[]);
halt:
'a[] input array of coefficients
'x[] output x co-ordinates, y[] output data values
proc MakePoly(a[], x[], y[])
var i%,j%,xv,s;
for i% := 0 to Len(y[])-1 do
  s := 0.0;
  xv := 1;
   for j% := 0 to NCOEF%-1 do
      s += a[j%]*xv;
     xv *= i%;
     next;
  v[i%] := s;
      x[i%] := i%;
  next:
end;
See also:ChanFit(), FitExp(), FitGauss(), FitLinear(), FitSin()
```

FitSin()	This con for each	mmand fits multiple sinusoids to arrays of x,y data points, with an optional weight point. Fitting is by an iterative method. The input data is fitted to the equation:
	$y = \mathbf{a}_0 \operatorname{si}$	$\ln(\mathbf{a}_1 x + \mathbf{a}_2) + \mathbf{a}_3 \sin(\mathbf{a}_4 x + \mathbf{a}_5) + \dots \{+ \mathbf{a}_{3n}\}$
	The fitt sinusoic function \mathbf{a}_5 , \mathbf{a}_8) a	ed coefficients are the \mathbf{a}_i . Angles are evaluated in radians. You can fit up to 3 ls and an optional offset. Although the function is given in terms of sine is, you can easily convert to cosines by subtracting $\pi/2$ from the phase angle (\mathbf{a}_2 , fter the fit. The commands to implement this are:
Set up the problem	The firs must ca	t command sets the number of sinusoids to fit and the data set to be fitted. You ll this function before you call any of the others.
	func F	<pre>itSin(nCoef%, y[], x[]{, s[] s});</pre>
	nCoef%	The number of coefficients to fit. The only legal values are 3, 6 and 9 for one, two and three sinusoids or 4, 7 and 10 to include an offset as the last coefficient.
	у[]	An array of y data values. The length of the array must be at least $nCoef$.
	x[]	An array of x data values. The length of the array must be at least $nCoef$.
	S	An optional argument which is either an array with one value for each y data point or a single value for all data points. If the value is the standard deviation (σ) , then the error value returned when you iterate to find the best fit is the χ^2 value and the fit is a chi-squared fit.
		If this value is proportional to the error in the y values, then the fit is still a chi- squared fit, and the error returned is proportional to the chi-squared value. If you omit this array, the fit is a least-squares error fit, and the error value returned is the sum of the squares of the errors in the y values.
	Returns	The function returns 0. There is no other return value as all errors stop the script.
	The nur present) it is not,	nber of data points is set by the smallest of the sizes of the $y[], x[]$ or $s[]$ (if arrays. The number of data points must be at least the number of coefficients. If you will get a fatal error, so check this before calling the function.
Set coefficient values and ranges	This van the rang before t	riant of the function sets the initial value of each coefficient and optionally sets ge of allowed values. You can call this function at any time after the setup call and he iterate call (below) has returned 0, indicating that the fit is complete.
	In this t coefficie or have that the give sta we expo suggest	ype of fitting it is very important that you give reasonable starting values for the ents. In particular, when fitting multiple sinusoids you will usually either know, a good idea of the frequencies. You should limit the range of each frequency so y cannot overlap. If you can do this, the fit will proceed quickly. If you do not rting values, the command will make a simplistic guess at the fitting values. As ect that you know more about the "right" answer than the command does, we that you set the values you want.
	func F	<pre>itSin(coef%, val{, lo, hi});</pre>
	coef%	The coefficient to set. The first coefficient is number 0, the last is nCoef%-1.
	val	The initial value to assign to the coefficient. If you have set low and high limits, and the value is outside these limits, it is set to the nearer limit.
	lo,hi	If present, these two values set the acceptable range of values for this coefficient. If omitted, or if both values are set to the same value, there is no limit. The value of the coefficient is tested against these limits after each iterative step, and if it exceeds a limit, it is set to that limit.

Iterate to a solution	1 Once you have set up the problem and given initial values to your coefficients, you can start the iteration process that will move the coefficients from their starting values to new values that minimise the error (optionally scaled by the s argument).		
	func F	itSin(a[], &err{, ma	xI%{, &iTer%{, covar[][]}});
	a[]	An array of size at least coefficient values. The fi the first phase angle in a [t nCoef% that is returned holding the current set of rst amplitude is in a[0], the first frequency in a[1], 2], the second amplitude in a[3] and so on.
	err	A real variable returned a if $s[]$ is used or holding $y_{x[i]}$ is the value predicted	s the sum over the data points of $(y_{x[i]}-y[i])^2/s[i]^2$ g the sum of $(y_{x[i]}-y[i])^2$ if $s[]$ is not used, where d from the coefficients at the x value $x[i]$.
	max1%	This is the maximum num function. If you omit this from 1 to 10000. If you see	ber of fitting iterations to do before returning from the value, the function sets 1000. You can set any value at more than 10000, the number is limited to 10000.
	iTer%	An optional integer variated one before the function r	ble that is returned holding the number of iterations eturned.
	covar	An optional two dimension returned holding the con- changed unless the function	onal array of size at least $[nCoef%][nCoef%]$ that is variance matrix when the fit is complete. It is not on return value is 0 or -1.
	Returns	1 if the number of iteratio fit has converged, -1 if th negative numbers indicate	ns are exhausted and the fit has not converged, 0 if the e fit is not improving but the result may be OK. Other e failure.
	Even wh it is the	nen a minimum is found, the best minimum that this alg	here is no guarantee that this is <i>the</i> minimum, only that orithm can find given the starting point.
Select coefficients to fit	fit Sometimes you may wish to hold some coefficients fixed while you fit others. Normally the command will fit all the coefficients, but you can use this command variant to select the coefficients to fit. You can use this command at any time after you have set the problem until the iteration variant returns 0 or -1.		
	<pre>func FitSin(fit%[]);</pre>		
	fit%[]	An array of at least nCc constant, otherwise it is fi	ef% integers. If fit%[i] is 0, coefficient i is held tted. If all elements are 0, all arguments are fitted.
	The effe than 1, c is likely	ect of this command persis or you set up a new proble that you will know the free	ts until either the iteration variant returns a value less m, or you call this variant again. For a sinusoidal fit it quency to fit, so you may well hold this constant.
An example	The follo frequence	owing is a template for using you which we assume you k	ng this command (assuming you don't want to fit the now).
	const var x[var s[var fi var co var er	<pre>nData%:=50; nData%], y[nData%]; nData%]; t%[3]; efs[4]; r;</pre>	<pre>'set number of data elements 'space for our arrays 'space for sigma of each point 'we want to hold the frequency 'space for coefficients 'will hold error squared 'in here goes code to get the data</pre>
	<pre>FitSin(3, y[], x[], s[]); 'fit one sinusoid</pre>		'fit one sinusoid
	'Note FitSin FitSin FitSin	that we let the phas (0, 1.0, 0.2, 4); (1, .02, .01, .03); (2, 0., 0.3, 1.9);	e take any value 'set amplitude and limit range 'set frequency 'set width and limit range

```
'Now we say that we don't want to fit the frequency
ArrConst(fit%[],1); 'set all elements to 1
```

	<pre>fit%[1] := 0;</pre>
	DrawMyData(coefs[], x[], y[]); 'Some function to show progress until FitSin(coefs[], err, 1) < 1;
	<pre>DrawMyData(coefs[], x[], y[]); 'Show the final state</pre>
	<pre>See also:FitExp(), FitGauss(), FitLinear(), FitNLUser(), FitSin()</pre>
FitValue()	This function returns the value at a particular x axis point of the fitted function set by the last FitData() command.
	<pre>Func FitValue(x);</pre>
	x The x axis value at which to evaluate the current fit. You should be aware that some of the fitting fuctions can overflow the floating point range if you ask for x values beyond the fitted range of the function.
	Returns The value of the fitted function at x. If the result is out of floating point range, the function may return a floating point infinity or a NaN (Not a Number) value or a 0. If there is no fit, the result is always 0.
	<pre>See also:FitCoef(), FitData(), FitExp(), ChanFitValue()</pre>
FocusHandle()	This function returns the view handle of the script-controllable window with the input focus (the active window). Unlike FrontView(), it can return any type of window, for example multimedia and spike shape windows.
	<pre>Func FocusHandle();</pre>
	Returns The handle of a window that the script can manipulate, or 0 if the focus is not in such a window.
	See also:FrontView()
FontGet()	This function gets the name of the font, and its characteristics for the current view. Func FontGet(&name\$, &size, &style%);
	name \$ This string variable is returned holding the name of the font.
	size The real number variable is returned holding the point size of the font.
	style% Returned holding the style: 0=Normal, 1=Italic, 2=Bold, 3=Bold and Italic. Returns The function returns 0 if all was well or a negative error code. If an error occurs, the variables are not changed.
	See also:FontSet()
FontSet()	This function sets the font for the current view. Text-based views (text, sequencer and script) should be set to non-proportionally spaced fonts only. If you set proportional fonts, the caret position in the view will be incorrect and the text will be difficult to edit.
	<pre>Func FontSet(name\$ code%, size, style%);</pre>
	name\$ A string holding the font name to use or you can specify a font code:
	code% This is an alternative method of specifying a font. We recognise these codes:

- 0 The standard system font, whatever that might be
- 1 A non-proportionally spaced font, usually Courier-like
- 2 A proportionally spaced non-serifed font, such as Helvetica or Arial
- 3 A proportionally spaced serifed font, such as Times Roman
- 4 A symbol font
- 5 A decorative font, such as Zapf-Dingbats or TrueType Wingdings

size The point size required. Your system may limit the allowed range.

style% The type style: 0=Normal, 1=Italic, 2=Bold, 3=Bold and Italic. Returns The function returns 0 if the font change succeeded, or a negative error code.

See also:FontGet()

Frac()

Returns the fractional part of a real number or converts a real array to its fractional parts.

Func Frac(x | x[] | x[]]);

x A real number or an array of reals.

Returns For arrays, it returns 0 or a negative error code. If x is not an array it returns a real number equal to x-Trunc(x). Frac(4.7) is 0.7, Frac(-4.7) is -0.7.

See also:Abs(), ATan(), Cos(), Exp(), Ln(), Log(), Max(), Min(), Pow(), Rand(), Round(), Sin(), Sqrt(), Tan(), Trunc()

FrontView()

The front view is the time, result, XY or text-based view that decides which menu appears in the main window. You can use this command to get the front view, or to set it. When a view becomes the front view, it is moved to the front and made the current view. If an invisible or iconised view is made the front view, it is made visible automatically, (equivalent to WindowVisible(1)).

You can also use this command to move other windows (for example the multimedia and Spike shape windows) to the top and make them the current view, but such windows are not returned by this command. Use FocusHandle() to return these windows.

Func FrontView({vh%});

- vh% Either 0 or omitted to return the front view handle, a view handle to be set, or -n, meaning the nth duplicate of the time view associated with the current view.
- Returns 0 if there are no visible views, -1 if the view handle passed is not a valid view handle, otherwise it returns the view handle of the view that was at the front.

See also:FocusHandle(), View(), WindowVisible()

GammaP(), GammaQ()

$$P(a,x) = \int_0^{\infty} e^{-t} t^{a-1} dt / \Gamma(a)$$

This is named GammaP() in the script. We also define the complement of this function, GammaQ(a, x), which is 1.0-GammaP(a, x). From them are obtained the error function, the cumulative Poisson probability function and the Chi-squared probability function.

The error function $erf(x) = 2/\sqrt{\pi} \int_0^x e^{-t^2} dt = \text{GammaP}(0.5, x*x)$

The incomplete gamma function P(a, x) is defined mathematically as:

The cumulative Poisson probability function relates to a Poisson process of random events and is the probability that, given an expected number of events r in a given time

period, the actual number was greater than or equal to n. This turns out to be GammaP(n, r). Also, the probability that there are less than n events is GammaQ(n, r).

The Chi-squared probability function is useful where we are fitting a model to data. Given a fitting function that fits the data with n degrees of freedom (if you have nData data points and nCoef coefficients you usually have nData-nCoef degrees of freedom), and given that the errors in the data points are normally distributed, the probability of a Chi-squared value less than chisq is GammaP(n/2, chisq/2). Similarly, the probability of a chisq value at least as large as chisq is GammaQ(n/2, chisq/2). So, if you know the chi-squared value from a fitting exercise, you can ask "What is the probability of getting this value (or a greater one) given that my model fits the data?" If the probability is very small, it is likely that your model does not fit the data, or your fit has not converged to the correct solution.

Func GammaP(a, x);
Func GammaQ(a, x);

- a This must be positive, it is a fatal error if it is not.
- x This must be positive, it is a fatal error if it is not.
- Returns These functions return the incomplete Gamma function and the complement of the incomplete Gamma function.

See also:LnGamma()

Grid()

This function turns the background grid on and off for the current time, result or XY view. It also returns the state of the grid.

Func Grid({on%});

on \$ Optional, sets the grid state. 0 = off, 1 = on, omit for no change.

Returns The state of the grid at the time of the call, or a negative error code. Changes made by this function do not cause an immediate redraw.

See also:XAxis(), XScroller(), YAxis()

Gutter()

The gutter is the area on the left of a text-based window where bookmarks and script break points appear. This function returns and optionally sets the gutter visible state. If you set a large font size you may wish to hide the gutter.

Func Gutter({show%});

show% Optional, sets the gutter state. 0 = hide, 1 = show, -1 or omitted for no change. Returns The gutter state at the time of the call: 0 = hidden, 1 = visible.

HCursor()

This function returns the position of a horizontal cursor, and optionally sets a new position. You can get and set positions of cursors attached to invisible channels or channels that have no y axis. It is also used with spike shape window to set trigger levels.

Func HCursor(num% {,where {,chan%}});

num% The cursor to use. It is an error to attempt this operation on an unknown cursor.

In a spike shape window, cursor 1 is the low trigger level, 2 is the high trigger, 3 is the low limit (if enabled), 4 is the high limit (if enabled). To set the cursor for trace n (0-3) add 4*n to the cursor number.

	where If this parameter is given it sets the new position of the cursor.
	chan% If this parameter is given, it sets the channel number. In XY or spike shape views you should omit this argument as it is ignored.
	Returns The function returns the position of the cursor at the time of the call.
	See also:HCursorChan(),HCursorDelete(),HCursorLabel(), HCursorLabelPos(),HCursorNew(),HCursorRenumber()
HCursorChan()	This function returns the channel number that a particular cursor is attached to.
	num% The horizontal cursor number
	Returns It returns the channel number that the cursor is attached to, or 0 if this cursor is not attached to any channel or if the channel number is out of the allowed range. XY views return 1 as the channel number.
	See also:HCursor(), HCursorDelete(), HCursorNew(), HCursorRenumber()
HCursorDelete()	This deletes the designated horizontal cursor. It is not an error to delete an unknown cursor; it just has no effect.
	<pre>Func HCursorDelete({num%});</pre>
	num% The number of the cursor to delete. If this is omitted, the highest numbered cursor is deleted. Set -1 to delete all horizontal cursors.
	Returns The number of the deleted cursor or 0 if no cursor was deleted.
	See also:HCursor(),HCursorChan(),HCursorLabel(),HCursorLabelPos(), HCursorNew(),HCursorRenumber()
HCursorExists()	Use this function to determine if a horizontal cursor exists. Func HCursorExists(num%);
	num% The cursor number in the range 1-4.
	Returns 0 if the cursor does not exist, 1 if it does.
	See also:HCursorNew(), CursorExists()
HCursorLabel()	This command sets (or gets) the cursor label style for the current view. Func HCursorLabel({style%{, num%{, form\$}}})
	style% The cursor style. Cursors can be annotated with a position or the cursor number or a user-defined style. The styles are: 0=Neither, 1=Position, 2=Number, 3=Both, 4=User-defined (not Macintosh yet). Unknown styles cause no change.
	num? The cursor number for style 4. 0 means all cursors, 1-4 for a single cursor.
	form\$ The label string for style 4. The string has replaceable parameters %p, and %n for position and number. We also allow %w.dp where w and d are numbers that set the field width and decimal places. You cannot read back a label format string.
	Returns The previous cursor style. If you omit $style$, the style does not change.
	See also:HCursorLabelPos(),HCursorNew(),HCursorRenumber()

HCursorl abelPos()	This lets you set and read the position of the cursor label.		
	Func HCursorLabelPos(num% {,pos});		
	num% The cursor number. Setting a silly number does nothing and returns -1.		
	pos If present, the command sets the position. The position is a percentage of the distance from the left of the cursor at which to position the value. Out of range values are set to the appropriate limit.		
	Returns The cursor position before any change was made.		
	See also:HCursor(),HCursorChan(),HCursorDelete(),HCursorLabel(), HCursorNew(),HCursorRenumber()		
HCursorNew()	This function creates a new horizontal cursor in a time or result view and assigns it to a channel. You can create up to 4 horizontal cursors.		
	<pre>Func HCursorNew(chan% {,where});</pre>		
	chan% A channel for the new cursor. If the channel is hidden, the cursor is not visible. You should use a channel number of 1 for XY and result views.		
	where An optional argument setting the cursor position. If this is omitted, the cursor is placed in the middle of the y axis or at zero if there is no y axis.		
	Returns It returns the horizontal cursor number or 0 if all cursors are in use.		
	See also:HCursor(),HCursorChan(),HCursorDelete(),HCursorLabel(), HCursorLabelPos(),HCursorRenumber()		
HCursorRenumber()	This command renumbers the cursors from bottom to top. There are no arguments.		
	<pre>Func HCursorRenumber();</pre>		
	Returns The number of cursors found in the view.		
	See also:HCursor(), HCursorChan(), HCursorDelete(), HCursorLabel(), HCursorLabelPos(), HCursorNew()		
Help()	The help available depends on the system. The Windows version uses the standard Windows help system. The Macintosh version may provide help through Apple Guide in future releases.		
	<pre>Func Help(topic% topic\$ {,file\$});</pre>		
	topic% A numeric code for the help topic. These codes are assigned by the help system author. Code 0 changes the default help file to file\$.		
	topic\$ A string holding a help topic keyword or phrase to look-up.		
	file\$ If this is omitted, or the string is empty, the standard Spike2 help file is used. If this holds a filename, this filename is used as the help file.		
	Returns 1 if the help topic was found, 0 if it was not found, -1 if no help file was found.		

The Windows SDK has some help-authoring tools, and third-party tools are available.

IIR commands The IIRXXXX() script commands make it easy for you to generate and apply IIR (Infinite Impulse Response) filters to data held in arrays of real numbers. The data values are assumed to be a sampled sequence, spaced at equal intervals. You can create digital filters that are modelled on Butterworth, Bessel, Chebyshev type 1 and Chebyshev type 2 highpass, lowpass, bandstop and bandpass filters. You can also create digital resonators and notch filters. The commands are:

IIRBp()	Bandpass filter	IIRHp()	Highpass filter	IIRNotch()	Notch filter
IIRBs()	Bandstop filter	IIRLp()	Lowpass filter	IIRReson()	Resonator

The algorithms used to create the filters are based on the mkfilter program, written by Tony Fisher of York University. The basic idea is to position the s-plane poles and zeros for a normalised low-pass filter of the desired characteristic and order, then to transform the filter to the desired type.

The theory of IIR filters is beyond the scope of this manual; a classic reference work is *Theory and Application of Digital Signal Processing* by Rabiner and Gold, published in 1975. The IIR filters generated by these commands can be modelled by:

$$y[\mathbf{n}] = \sum_{i=0}^{N} a_i x[\mathbf{n} \cdot i]/G + \sum_{i=1}^{M} b_i y[\mathbf{n} \cdot i]$$

where the x[n] are the sequence of input data values, the y[n] are the sequence of output values, the a_i and the b_i are the filter coefficients (some of which may be zero) and G is the filter gain. Although G could be incorporated into the a_i , for computational reasons we keep it separate. In the filters designed by the IIRXXXX() commands, N=M and is the *order* of the filter for lowpass and highpass designs and is twice the order for bandpass and bandstop designs. The order of these filters determines the sharpness of the filter cut-off.

IIR and FIR filters When compared to FIR filters, IIR filters have advantages:

- They can generate much steeper edges and narrower notches for the same computational effort.
- The filters are causal, which means that the filter output is only affected by current and previous data. If you run a step change through FIR filters you typically see ringing before the step as well as after it.

However, they also have disadvantages:

- IIR filters are prone to stability problems, FIR filters are unconditionally stable. All IIR filters generated by these commands should be stable, but high order band pass and band stop filters with narrow pass or stop bands may have problems.
- They impose a group delay on the data that varies with frequency. This means that they do not preserve the shape of a waveform, in particular, the positions of peaks and troughs will change.

You can remedy the group delay problem by running a filter forwards, then backwards, through the data. However, this makes the filter non-causal, removing one of the advantages of using an IIR filter. The commands allow you to check the impulse, step, frequency and phase response of the filters, and we recommend that you do so before using a generated filter for a critical purpose.

The lowpass, highpass, bandpass and bandstop filters generate digital filters modelled on four types of analogue filter: Butterworth, Bessel, Chebyshev type 1 and Chebyshev type 2. The resulting digital filters are not identical to the analogue filters as the mapping from the analogue to the digital domain distorts the frequency scale. In many cases, this improves the performance of the digital filter over the analogue counterpart.

- **Filter types** You can generate notch and resonator filters plus lowpass, highpass, bandpass and bandstop filters modelled on Butterworth, Bessel and Chebyshev analogue filters. The examples for Butterworth, Bessel and Chebyshev filters show a fifth order lowpass filter with the edge set to 0.2 with inset examples of high-pass, bandpass and bandstop filters.
- **Butterworth** These have a maximally flat pass band, but pay for this by not having the steepest possible transition between the pass band and the stop band. The example shows a low pass fifth order Butterworth filter with a cut-off frequency set to 0.2 of the sample rate. The x axis is frequency, the y axis is the filter gain. Both axes are linear.
 - **Bessel** An analogue Bessel filter has the property that the group delay is maximally flat, which means that it tends to preserve the shape of a signal passed through it. This leads to filters with a gentle cut-off. When digitised, the constant group delay property is compromised; the higher the filter order, the worse the group delay. The example shows a fifth order low pass filter at 0.2 of the sample rate.
- **Chebyshev type 1** Filters of this type are based on Chebyshev polynomials and have the fastest transition between the pass band and the stop band for a given ripple in the pass band and no ripples in the stop band. In the example, the ripple has been set to 3 dB, to match the other examples, though you would normally choose less ripple than this.
- **Chebyshev type 2** Filters of this type are defined by the start of the stop band and the stop band ripple. The filter has the fastest transition between the pass and stop bands given the stop band ripple and no ripple in the pass band. The example shows a fifth order filter with a 40 dB stop band ripple and with the stop band starting at 0.2 of the sample rate.
 - **Notch** Notch filters are defined by a centre frequency and a q factor. q is the width of the stop band at the -3 dB point divided by the centre frequency: the higher the q, the narrower the notch. Notch filters are often used to remove mains hum, but if you do this you will likely need to set notches at the first few odd harmonics of the mains frequency. The example has a centre frequency of 0.2 and a q of 10, so the width the -3 dB point is 0.02 of the sample rate.
 - **Resonator** A resonator is the inverse of a notch. It is defined in terms of a centre frequency and a q factor. q is the width of the pass band at the -3 dB point divided by the centre frequency: the higher the q, the narrower the resonance. Resonators are sometimes used as alternatives to a narrow bandpass filter. The example shows a centre frequency of 0.2 of the sample rate and a q of 10, so the width of the pass band at the -3 dB point is 0.02 of the sample rate.



0.25

Common command variants	There are six IIR commands. The commands are all independent and they each remember the last filter you created. They all support the following variants:			
Get filter information	You can use a command variant of this form to return filter information:			
	<pre>Func IIRxxxx(get%{, arr[]{[]}});</pre>			
	 get% The form of the command using this argument is used to read back information about the last filter created with this command. The argument can be: arr[] is filled in with the impulse response of the filter. The return value is the magnitude of the largest value in arr[]. For example, if the impulse response ranged in values from -0.5 to 0.3, 0.5 would be returned. arr[] is filled in with the step response of the filter. The return value is the magnitude of the largest value in arr[]. arr[] [] is a matrix with 2 columns and r rows. The frequency response is returned as complex numbers in the columns; column 0 holds the real part and column 1 the imaginary part. The first row corresponds to a frequency of 0; the final row corresponds to a frequency of half the sampling rate. The frequencies are spaced as 0.5/(r-1). The return value is the maximum magnitude of the returned frequency response (if you wish to normalise the response curve to a maximum of 1.0). The same as 3 except that the results are returned as the amplitude response in column 0 and the phase response in column 1. If the real and imaginary parts of the response are r and i, column 0 holds sqrt(r*r+i*i) and column 1 holds atan(i, r). The return value is the maximum returned amplitude response. Returns the number of filter coefficients (2*order%+1) to apply to the filter input values and fills in arr[] with these values. These correspond to the <i>a</i>_i in the filter expression. However, we return the values in reverse order to match the <i>a</i>_i. The final value is always -1.0 (corresponding to <i>b</i>₀, which is not used when implementing the filter). Returns the filter gain <i>G</i>. arr[] is ia matrix with 2 columns. The return value is the number of poles in the solute in the poles as complex numbers with the real part in column 0 and the imaginary part in column 1. 			
Apply existing filter	This command applies the current filter to an array of equally spaced data. The variant with a single argument applies the filter forward through the array; with two arguments you have to option of running it backwards. Running the filter forwards introduces a phase shift into the output. Running the filter a second time, but backwards, cancels the phase shift at the expense of a non-causal filter.			
	<pre>Func IIRxxxx(data[]{, flags%});</pre>			
	data An array of data to filter.			
	flags% Optional, taken as 0 if omitted. Set to 1 to run the filter backwards.			
	<pre>See also:FIRMake(), IIRBp(), IIRBs(), IIRHp(), IIRLp(), IIRNotch(), IIRReson()</pre>			

IIRBp()	This fur arrays of	nction creates and applies IIR (Infinite Impulse Response) band pass filters to f data. You can run the filter forwards or backwards through the data.
	Func I Func I Func I	IRBp(data[] 0, lo, hi, order%{, type%{, ripple}}); IRBp(data[]{, flags%}); IRBp(get%{, arr[]{[]}});
	data	An array of data to filter. If there are only 1 or 2 arguments, the last created band pass filter is used. Otherwise, the filter defined by the remaining arguments is used forwards. Replace data with 0 to create a filter without applying it.
	flags%	Optional, taken as 0 if omitted. Add 1 to run the filter backwards.
	lo	The low corner frequency of the band stop filter. This is expressed as a fraction of the sample rate and is limited to the range 0.0002 to 0.4998. For Chebyshev type 2 filters, this is the point at which the attenuation reaches the <code>ripple</code> value, for all other filters this sets the -3 dB point.
	hi	The high corner frequency of the band pass filter. This is expressed as a fraction of the sampling rate and is limited to the range $10+0.0002$ to 0.4998. For Chebyshev type 2 filters, this is the point at which the attenuation reaches the ripple value, for all other filters this sets the -3 dB point.
	order%	The order of the lowpass filter used as the basis of the design, in the range 1 to 10. The order of the filter implemented is $order\%^2$. High orders ($order\% > 7$) and narrow bands may cause inaccuracy in the filter. Narrow means that ($hi-lo$)/sqrt($lo\%hi$) is less than 0.2, for example.
	type%	Set 0 for Butterworth, 1 for Bessel and 2 for Chebyshev type 1, 3 for Chebyshev type 2. The default value is 0 for a Butterworth filter.
	ripple	The desired pass band ripple in dB for Chebyshev type 1 filters (default 3 dB) or the desired minimum cut in the stop bands for Chebyshev type 2 filters (default 40 dB). The value here must be greater than 0 .
	get%	The command variant with this argument returns information about the last filter you created with this command. See the discussion of IIR commands for details.
	arr	This is an option vector or matrix used to return information about the last filter you created with this command. See the discussion of IIR commands for details.
	Returns	All forms of the command return negative numbers for errors. The forms that apply or create filters return 0 for success. The other command forms have their return values included in the description of the get% argument or return 0.
	See also	:FIRMake(), IIRBs(), IIRHp(), IIRLp(), IIRNotch(), IIRReson()

llRBs()

This function creates and applies IIR (Infinite Impulse Response) band stop filters to arrays of data. You can run the filter forwards or backwards through the data.

```
Func IIRBp(data[]|0, lo, hi, order%{, type%{, ripple}});
Func IIRBp(data[]{, flags%});
Func IIRBp(get%{, arr[]{[]}});
```

- data An array of data to filter. If there are only 1 or 2 arguments, the last created band stop filter is used. Otherwise, the filter defined by the remaining arguments is applied forwards. Replace data with 0 to create a filter without applying it.
- flags% Optional, taken as 0 if omitted. Add 1 to run the filter backwards.
- 10 The low corner frequency of the band pass filter. This is expressed as a fraction of the sample rate and is limited to the range 0.0002 to 0.4998. For Chebyshev filters, this is the point at which the attenuation reaches the ripple value, for other filters this sets the -3 dB point.

- hi The high corner frequency of the band pass filter. This is expressed as a fraction of the sampling rate and is limited to the range 10+0.0002 to 0.4998. For Chebyshev type 2 filters, this is the point at which the attenuation reaches the ripple value, for all other filters this sets the -3 dB point.
- order% The order of the lowpass filter used as the basis of the design, in the range 1 to 10. The order of the filter implemented is order%*2. High orders and narrow pass bands may lose numerical accuracy in the filter output.
- type% Set 0 for Butterworth, 1 for Bessel and 2 for Chebyshev type 1, 3 for Chebyshev type 2. The default value is 0 for a Butterworth filter.
- ripple The desired pass band ripple in dB for Chebyshev type 1 filters (default 3 dB) or the desired minimum cut in the stop band for Chebyshev type 2 filters (default 40 dB). The value here must be greater than 0.
- get% The command variant with this argument returns information about the last filter you created with this command. See the discussion of IIR commands for details.
- arr This is an option vector or matrix used to return information about the last filter you created with this command. See the discussion of IIR commands for details.
- Returns All forms of the command return negative numbers for errors. The forms that apply or create filters return 0 for success. The other command forms have their return values included in the description of the get% argument or return 0.

See also:FIRMake(), IIRBp(), IIRHp(), IIRLp(), IIRNotch(), IIRReson()

llRHp()

This function creates and applies IIR (Infinite Impulse Response) high pass filters to arrays of data. You can run the filter forwards or backwards through the data.

Func IIRHp(data[]|0, edge, order%{, type%{, ripple}});
Func IIRHp(data[]{, flags%});
Func IIRHp(get%{, arr[]{[]}});

- data An array of data to filter. If there are only 1 or 2 arguments, the last created high pass filter is used. Otherwise, the filter defined by the remaining arguments is applied forwards. Replace data with 0 to create a filter without applying it.
- flags% Optional, taken as 0 if omitted. Add 1 to run the filter backwards.
- edge The corner frequency of the high pass filter. This is expressed as a fraction of the sample rate and is limited to the range 0.0002 to 0.4998. For Chebyshev filters, this is the point at which the attenuation reaches the ripple value, for other filters this sets the -3 dB point.
- order% The order of the filter in the range 1 to 10.
- type% Set 0 for Butterworth, 1 for Bessel and 2 for Chebyshev type 1, 3 for Chebyshev type 2. The default value is 0 for a Butterworth filter.
- ripple The desired pass band ripple in dB for Chebyshev type 1 filters (default 3 dB) or the desired minimum cut in the stop bands for Chebyshev type 2 filters (default 40 dB). The value here must be greater than 0.
- get% The command variant with this argument returns information about the last filter you created with this command. See the discussion of IIR commands for details.
- arr This is an option vector or matrix used to return information about the last filter you created with this command. See the discussion of IIR commands for details.
- Returns All forms of the command return negative numbers for errors. The forms that apply or create filters return 0 for success. The other command forms have their return values included in the description of the get% argument or return 0.

See also:FIRMake(), IIRBp(), IIRBs(), IIRLp(), IIRNotch(), IIRReson()

llRLp()	This fur arrays of	nction creates and applies IIR (Infinite Impulse Response) low pass filters to f data. You can run the filter forwards or backwards through the data.
	Func I Func I Func I	IRLp(data[] 0, edge, order%{, type%{, ripple}}); IRLp(data[]{, flags%}); IRLp(get%{, arr[]{[]}});
	data	An array of data to filter. If there are only 1 or 2 arguments, the last created low pass filter is used. Otherwise, the filter defined by the remaining arguments is applied forwards. Replace data with 0 to create a filter without applying it.
	flags%	Optional, taken as 0 if omitted. Add 1 to run the filter backwards.
	edge	The corner frequency of the low pass filter. This is expressed as a fraction of the sample rate and is limited to the range 0.0002 to 0.4998. For Chebyshev filters, this is the point at which the attenuation reaches the ripple value, for other filters this sets the -3 dB point.
	order%	The order of the filter in the range 1 to 10.
	type%	Set 0 for Butterworth, 1 for Bessel and 2 for Chebyshev type 1, 3 for Chebyshev type 2. The default value is 0 for a Butterworth filter.
	ripple	The desired pass band ripple in dB for Chebyshev type 1 filters (default 3 dB) or the desired minimum cut in the stop bands for Chebyshev type 2 filters (default 40 dB). The value here must be greater than 0.
	get%	The command variant with this argument returns information about the last filter you created with this command. See the discussion of IIR commands for details.
	arr	This is an option vector or matrix used to return information about the last filter you created with this command. See the discussion of IIR commands for details.
	Returns	All forms of the command return negative numbers for errors. The forms that apply or create filters return 0 for success. The other command forms have their return values included in the description of the get% argument or return 0.
	See also	:FIRMake(), IIRBp(), IIRBs(), IIRHp(), IIRNotch(), IIRReson()

IIRNotch()

This function creates and applies IIR (Infinite Impulse Response) notch filters to arrays of data. You can run the filter forwards or backwards through the data. The gain of the notch filter is zero at the notch frequency and 1 at low and high frequencies.

```
Func IIRNotch(data[]|0, fr, q);
Func IIRNotch(data[]{, flags%});
Func IIRNotch(get%{, arr[]{[]}});
```

- data An array of data to filter. If there are only 1 or 2 arguments, the last created notch filter is used. Otherwise, the filter defined by the remaining arguments is applied forwards. Replace data with 0 to create a filter without applying it.
- flags% Optional, taken as 0 if omitted. Add 1 to run the filter backwards.
- fr The frequency of the notch. This is expressed as a fraction of the sample rate and is limited to the range 0.0002 to 0.4998.
- The desired q factor for the notch filter. If the two -3 dB points either side of the notch are at frequencies Flo and Fhi, q is given by fr/(Fhi-Flo). The higher the q, the narrower the notch. Try 100 as a starting point.
- get% The command variant with this argument returns information about the last filter you created with this command. See the discussion of IIR commands for details.
- arr This is an option vector or matrix used to return information about the last filter you created with this command. See the discussion of IIR commands for details.

Returns All forms of the command return negative numbers for errors. The forms that apply or create filters return 0 for success. The other command forms have their return values included in the description of the get% argument or return 0.

See also:FIRMake(), IIRBp(), IIRBs(), IIRHp(), IIRLp(), IIRReson()

IIRReson()

This function creates and applies IIR (Infinite Impulse Response) resonator filters to arrays of data. You can run the filter forwards or backwards through the data. The gain of the filter is 1 at the resonator frequency and zero at low and high frequencies.

```
Func IIRReson(data[]|0, fr, q);
Func IIRReson(data[]{, flags%});
Func IIRReson(get%{, arr[]{[]}});
```

- data An array of data to filter. If there are only 1 or 2 arguments, the last created band stop filter is used. Otherwise, the filter defined by the remaining arguments is applied forwards. Replace data with 0 to create a filter without applying it.
- flags% Optional, taken as 0 if omitted. Add 1 to run the filter backwards.
- fr The centre frequency of the resonator. This is expressed as a fraction of the sample rate and is limited to the range 0.0002 to 0.4998.
- The desired q factor for the resonator. If the two -3 dB points either side of the resonance are at frequencies Flo and Fhi, q is given by fr/(Fhi-Flo). The higher the q, the narrower the resonance. Try 100 as a starting point.
- get% The command variant with this argument returns information about the last filter you created with this command. See the discussion of IIR commands for details.
- arr This is an option vector or matrix used to return information about the last filter you created with this command. See the discussion of IIR commands for details.
- Returns All forms of the command return negative numbers for errors. The forms that apply or create filters return 0 for success. The other command forms have their return values included in the description of the get% argument or return 0.

See also:FIRMake(), IIRBp(), IIRBs(), IIRHp(), IIRLp(), IIRNotch()

Inkey()

This function is provided for compatibility with the MS-DOS version of Spike2. Do not use it in new scripts. It returns the ASCII code for the key pressed, or -1 if no key was pressed. In some cases Spike2 absorbs keystrokes, for example if you are sampling and the current window is the sampling time window all keystrokes are taken as markers.

Func Inkey();

Returns The key code or -1 if there is no pending key. The codes are:

1-31	Ctrl+ABCXYZ[\]^_	65-95	ABCXYZ[\]^_
32	space	96	`
33-47	!"#\$%&'()*+,/	97-126	abcxyz{ }~
48-64	0123456789:;<=>?@	127	Rubout

See also:Interact(), Keypress(), Toolbar(), ToolbarSet()

Input()	Input() This function reads a number from the user. It opens a window with a message, and displays the initial value of a variable. You can limit the range of the result.		
	<pre>Func Input(text\$, val {,low {,high}});</pre>		
	text\$ A string holding a prompt for the user. If the string includes a vertical bar, the text before the vertical bar is used as the window title.		
	val The initial value to be displayed for editing. If limits are given, and the initial value is outside the limits, it is set to the nearer limit.		
	low An optional low limit for the result. If low>=high, the limits are ignored.		
	high An optional high limit for the result.		
	Returns The value typed in. The function always returns a value. If an out-of-range value is entered, the function warns the user and a correct value must be given. When parsing the input, leading white space is ignored and the number interpretation stops at the first non-numeric character or the end of the string.		
	See also:DlgReal(),DlgInteger(),Input\$(),Inkey()		
Input\$()	This function reads user input into a string variable. It opens a window with a message, and displays a string. You can also limit the range of acceptable characters.		
	<pre>Func Input\$(text\$, edit\${, maxSz%{, legal\$}});</pre>		
	text\$ A string holding a prompt for the user. If the string includes a vertical bar, the text before the vertical bar is used as the window title.		
	edit ^{\$} The starting value for the text to edit.		
	maxSz% Optional, maximum size of the response string.		
	legal\$ An string holding acceptable characters. edit\$ is filtered before display. A hyphen indicates a range of characters. To include a hyphen in the list, place it first or last in the string. Upper and lower case characters are distinct. For upper and lower case characters and integer numbers use: "a-zA-Z0-9".		
	If this string is omitted, all printing characters are allowed, equivalent to " -~" (space to tilde). For simple use, the sequence of printing characters is: space !"#\$%&'()*+,/0123456789:;<=>?@ ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^_` abcdefghijklmnopqrstuvwxyz{ }~		
	The order of extended or accented characters is system and country dependent.		
	Returns The result is the edited string. A blank string is a possible result.		
	See also:DlgString(), Input(), Inkey()		
InStr()	This function searches for a string within another string. This function is case sensitive.		
U	<pre>Func InStr(text\$, find\$ {,index%});</pre>		
	text\$ The string to be searched.		
	find\$ The string to look for.		
	index% If present, the start character index for the search. The first character is index 1.		

Returns The index of the first matched character, or 0 if the string is not found.
Interact()

This function provides a quick and easy way to interact with a user. It displays the interact toolbar at the top of the Spike2 window and pauses the script until a button or a key linked to a button is pressed. Cursors can always be dragged as we assume that they are one of the main ways of interacting with the data. You can limit the user actions when the bar is active.

Func Interact(msg\$, allow% {,help {, lb1\$ {,lb2\$ {,lb3\$...}}}});

- msg\$ A message to display in the message area of the toolbar. The message area competes with the button area. With many buttons, the text may not be visible.
- allow% A number that specifies the actions that the user can and cannot take while interacting with Spike2. 0 allows the user to inspect data and position cursors in a single, unmoveable window. The codes are shown in both decimal and hexadecimal format. The number is the sum of possible activity codes.
 - 1 0x0001 Can change application
 - 2 0x0002 Can change the current window
 - 4 0x0004 Can move and resize windows
 - 8 0x0008 Can use File menu
 - 16 0x0010 Can use Edit menu
 - 32 0x0020 Can use View menu but not ReRun
 - 64 0x0040 Can use Analysis menu
 - 128 0x0080 Can use Cursor menu and add cursors
 - 256 0x0100 Can use Window menu
 - 512 0x0200 Can use Sample menu
 - 1024 0x0400 No changes to y axis
 - 2048 0x0800 No changes to x axis
 - 4096 0x1000 No horizontal cursor channel change
- help This is either the number of a help item (CED internal use) or it is a help context string. This is used to set the help information that is presented when the user presses the F1 key. Set 0 to accept the default help. Set a string as displayed in the Help Index to select a help topic, for example "Cursors: Adding".
- 1b1\$ These label strings create buttons, from right to left, in the tool bar. If no labels are given, one label is displayed with the text "OK". The maximum number of buttons is 17. Buttons can be linked to the keyboard using & and by adding a vertical bar followed by a key code to the end of the label. You can also set a tooltip. The format is "Label|code|tip". See the documentation for label\$ in the ToolbarSet() command for details.
- Returns The number of the button that was pressed. Buttons are numbered in order, so 1b1\$ is button 1, 1b2\$ is button 2 and so on.

If a toolbar created by Toolbar() is present, it is hidden during the Interact() command and restored after Interact() returns.

See also:Toolbar(), ToolbarSet()

Keypress()

This function returns 1 if the lnkey() function would return a character, or 0 if it would not. This function and lnkey() are provided for compatibility with the MS-DOS version of Spike2 and are not recommended for new scripts.

Func Keypress();

Returns 1 if a key is ready to read, 0 if there is no key.

See also:Inkey(), Interact(), Toolbar(), ToolbarSet()

LastTime()	This function finds the first item on a channel before a time. If a marker filter is applied to the channel, only data in the filter is visible. This function is for time views only.		
	Func L	<pre>astTime(chan%, time{,&val code%[]{,data[] data%[] &data\$}});</pre>	
	chan%	The channel number in the view to use.	
	time	The time to search before. Items at the time are ignored. To start a backward search that guarantees to iterate through all items, start at Maxtime(chan%)+1.	
	val	Optional: for waveform channels it returns the waveform value. For event level channels, it is returned 0 if the transition is low to high, and 1 if the transition is high to low. If there is no event it returns the level at time; 0 for low, 1 for high.	
	code%	This optional parameter is only used if the channel is a marker type (marker, RealMark, TextMark, WaveMark). This is an array with at least four elements that is filled in with the marker codes.	
	data	Filled with data from RealMark and WaveMark channels. If there is insufficient data to fill it, unused entries are unchanged. An integer array can be used with WaveMark data to collect a copy of the 16-bit data that holds the waveform. If WaveMark data has multiple traces, use data[points%][traces%] to get real data and data%[points%][traces%] to get the integer data.	
	data\$	A string returned holding the text from a TextMark channel.	
	Returns	The function returns either the time of the next item, or -1 if there are no more items to be found or a negative error code.	
	See also	:ChanData(), MaxTime(), NextTime()	
LCase\$()	This fun	action converts a string into lower case.	
	+ ov + \$	The string to convert	
	Returns	A lower cased version of the original string.	
	See also	:Asc(), Chr\$(), Str\$(), UCase\$(), Val()	
Left\$()	This fun	action returns the first n characters of a string.	
	Func L	eft\$(text\$, n);	
	text\$	A string of text.	
	n	The number of characters to extract.	
	Returns	The first n characters, or all the string if it is less than n characters long.	
	See also	:DelStr\$(),Len(),Mid\$(),Right\$()	
l on()	This fun	ction returns the length of a string or the size of a one dimensional array.	
Len()	Func L	en(text\$):	

arr[] A one dimensional array. It is an error to pass in a two dimensional array.Returns The length of the string or the array, as an integer.

Func Len(arr[]);
text\$ The text string.

You can find out the size of each dimension of a two dimensional array as follows:

```
proc something(arr[][]) 'function passed a 2-d array
var n%; n% := Len(arr[][0]); 'get size of first dimension
var m%; m% := Len(arr[0][]); 'get size of second dimension
See also:Asc(), Chr$(), DelStr$(), InStr(), LCase$(), Left$(), Mid$(),
Print$(), Right$(), Str$(), UCase$(), Val()
```

This function calculates the natural logarithm (inverse of Exp()) of an expression, or replaces the elements of an array with their natural logarithms.

Func Ln(x|x[]|x[][]);

- x A real number or a real array. Zero or negative numbers cause the script to halt with an error unless the argument is an array, when an error code is returned.
- Returns When used with an array, it returns 0 if all was well, or a negative error code. When used with an expression, it returns the natural logarithm of the argument.

See also:Abs(), ATan(), Cos(), Exp(), Frac(), Log(), Max(), Min(), Pow(), Rand(), Round(), Sin(), Sqrt(), Tan(), Trunc()

LnGamma()

Ln()

The gamma function $\Gamma(x)$ has the useful property that $\Gamma(n+1)$ is the same as n! (*n* factorial) for integral values of *n*. However, it increases very rapidly with *x*, reaching floating-point infinity when *x* is 172.62. To avoid this problem, the script returns the natural logarithm of the gamma function. The definition of the gamma function is:

$$\Gamma(a) = \int_0^\infty e^{-t} t^{a-1} dt$$

Func LnGamma(a);

a A positive value. The script stops with a fatal error if this is negative.

Returns The natural logarithm of the Gamma function of a.

See also:GammaP()

Log()

Takes the logarithm to the base 10 of the argument.

Func Log(x|x[]|x[]]);

- x A real number or a real array. Zero or negative numbers cause the script to halt with an error unless the argument is an array, when an error code is returned.
- Returns With an array, this returns 0 if all was well or a negative error code. With an expression, this returns the logarithm of the number to the base 10.
- See also:Abs(), ATan(), Cos(), Exp(), Frac(), Ln(), Max(), Min(), Pow(), Rand(), Round(), Sin(), Sqrt(), Tan(), Trunc()

LogHandle()

This returns the view handle of the log window (which always exists).

Func LogHandle();

See also:Print(), PrintLog(), View()

MarkEdit()	This changes the data stored in a marker at a particular time. You can get the data using LastTime() and NextTime().
	<pre>Func MarkEdit(chan%, time, code%[]{, data\$ data[] data%[]});</pre>
	chan% The marker, TextMark, WaveMark or RealMark channel to edit.
	time The time of the marker (must match exactly).
	code% Array of 4 marker codes (bottom 8 bits used) to replace markers in the channel.
	data The data to replace the data in the marker. If you use integer data with a WaveMark, the bottom 16 bits of each integer replace the data. To update a WaveMark with multiple traces use a two-dimensional array, for example with 32 points and 4 traces use var data[32][4]; to declare the array.
	Returns 0 if a marker was edited, or -1 if no marker exists at time.
	See also:LastTime(),MarkInfo(),MarkMask(),MarkSet(),NextTime()
MarkSet()	This sets the marker codes of a marker, WaveMark, TextMark or RealMark channel in a time range. If a marker filter mask is set, only data that is passed by the mask is changed. Func MarkSet(chan%, sT, eT, code%[]);
	Func MarkSet(chan%, sT, eT, c0%{, c1%{, c2%{, c3%}}});
	chan% The channel to process.
	sT, eT The time range to process.
	code% An array of 4 integers holding the new marker codes in the range 0 to 255 or -1 for a code that is unchanged.
	c0-c3% One to four marker codes as an alternative to code%[]. Use values 0 to 255 to change a code and -1 (or omit the value) to leave a code unchanged.
	Returns The number of markers that were changed.
	See also:LastTime(), MarkEdit(), MarkMask(), NextTime()
MarkMaak()	This function sets the mask for a marker. WaveMark, TextMark or RealMark channel.
MarkMask()	Each data item in one of these channels has four marker codes. Each code has a value from 0 to 255. Each data channel (and duplicated channel) has its own marker filter that determines the visible data items. A marker filter has four masks, one for each of the four marker codes. For each mask, you can specify which codes are wanted. There are two marker filter modes. In the diagrams, a marker data item is represented as a time stamp, four marker codes and data values that depend on the marker type.
Mode 0 (AND)	A marker data item is allowed through the mask if each of the four codes in the data item is present in the corresponding mask. We think of this as <i>and</i> mode because for the filter to pass the marker, marker code 0 must be is mask 0 <i>and</i> marker code 1 in mask 1 <i>mask</i> 2 mask 3 are usually set to accept all codes and masking is used for layer 0.
Mode 1 (OR)	A marker data item is allowed through the mask if any of the four marker codes is present in mask 0. A code 00 is only accepted for the first of the four marker codes. Masks 1 to 3 are ignored. We think of this as or mode because marker code 0 or 1 or 2 or 3 must be present in mask 0 for the filter to pass the marker for display or

analysis. This mode can be used with spike shape data (WaveMark) where two spikes have collided and the marker represents a spike of two different templates.

There are three command variants, the first sets the filter codes, the second the filter mode and the third allows you to read and set the mask:

Func MarkMask(chan%, layer%, set%, code%|code\$ {,code%|code\$...});
Func MarkMask(chan%, mode%);
Func MarkMask(chan%, mask%[]{[]}{, write%});

- chan% The channel number to work on. If this is not a suitable channel, the function does nothing and returns a negative error code.
- layer% The layer of the mask in the range 0 to 3 or -1 for all layers.
- set% 1 to include codes in the mask, 0 to exclude codes or -1 to invert the mask. Inverting a mask changes all included codes to excluded and vice versa.
- code% A number in the range 0 to 255 setting a code to include or exclude. You can specify more than one code at a time. -1 is also allowed, meaning all codes.
- code\$ Each character in the string is converted to its ASCII value, and used as a code.
- mode% The variant with two arguments returns the current mode of the marker filter and optionally sets the mode of matching as 0 or 1 or -1 for no change.
- mask% This is an array that is usually declared as mask%[256][4] or mask%[256] that maps onto the mask. The second dimension sets the layer. A onedimensional array maps onto layer 0. When reading, elements are set to 0 or 1. When writing, zero element clear and non-zero elements set the corresponding item in the mask. If the array size does not match 256 by 4, data is transferred between array items that map onto the mask.
- write% If this optional argument is omitted or zero, the mask% array is filled with data. If it is 1, the mask data is copied to the mask.

Returns 0 if all is OK, or a negative error code.

A common requirement is to allow all markers to be used. This is achieved by:

MarkMask(chan%, -1, 1, -1); 'Set all codes for all layers

To fill or empty or invert a complete layer use:

MarkMask(chan%, layer%, set%, -1); 'Apply set% to entire layer

This example sets the keyboard marker channel mask (channel 31) to show only markers 0 and 1 (start and stop recording markers) and key presses for A, B, C and D:

MarkMask(31, 0); 'set mode 0
MarkMask(31,-1, 1,-1); 'include everything (reset)
MarkMask(31, 0, 0, -1); 'exclude everything in layer 0
MarkMask(31, 0, 1, 0, 1, "ABCD"); 'include the codes we want

You can use this command together with ChanDuplicate() to split a marker channel into several channels based on marker codes.

MarkInfo()	This function is used to get information on the extended marker types (TextMark, WaveMark and RealMark).
	<pre>Func MarkInfo(chan% {,⪯% {,trace%}});</pre>
	chan% The channel to get the information from.
	pre% If present, returned as the number of pre-peak points for WaveMark data.
	trace% If present, returned as the number of traces (electrodes) in the WaveMark data.
	Returns For WaveMark data it returns the number of waveform points, for TextMark data it returns the maximum string length and for RealMark data it returns the number of reals attached to each marker. For all other channel types it returns 0.
	See also:LastTime(),MarkEdit(),MarkMask(),MarkSet(),NextTime()
MATDet()	This calculates the determinant of a matrix (a two dimensional array).
	<pre>Func MATDet(mat[][]);</pre>
	mat A two dimensional array with the same number of rows and columns.
	Returns The determinant of mat or 0.0 if the matrix is singular.
	See also:ArrAdd(),MATInv(),MATMul(),MATTrans()
	This is such a matrix (a true dimensional and a structure due data main at
MATInv()	This inverts a matrix (a two dimensional array) and optionally returns the determinant.
	<pre>Func MATInv(inv[][] {,src[][]{, &det}});</pre>
	inv A two dimensional array to hold the result. If src is omitted, inv is replaced by its own inverse. The number of rows and columns of inv must be the same.
	src If present, the matrix to invert. The numbers of rows and columns of this two dimensional array must be at least as large as inv.
	det If present, returned holding the determinant of the inverted matrix.
	Returns 0 if all was OK, -1 if the matrix was singular or very close to singular.
	See also:ArrAdd(), MATMul()
	This function multiplies matrices and/or vectors. In matrix terms, this evaluates $\mathbf{A} = \mathbf{B}C$
MA I Mul()	where A is an m rows by n columns matrix, B is an m by p matrix and C is a p by n matrix. Vectors of length v are treated as a v by 1 matrix.
	<pre>Proc MATMul(a[]{[]}, b[]{[]}, c[]{[]});</pre>
	a A m by n matrix of reals or a vector of length m (n is 1) to hold the result.

- b A m by p matrix or a vector of length m (p is 1).
- c A p by n matrix or a vector of length p (n must be 1).

If you pass any of a, b or c as a vector, they are treated as a n by 1 matrix, where n is the length of the vector. Use the trans() operator to convert a vector to a 1 by n matrix.

See also:ArrMul(), MATInv()

MATSolve()

This function solves the matrix equation Ax = y for x, given A and y. Both x and y are vectors of length n and A is an n by n matrix.

Func MATSolve(x[], a[][], y[]);

- x A one dimensional real array of length n to hold the result.
- a A two dimensional (n by n) array of reals holding the matrix.
- Y A one dimensional real array of length n.
- Returns The functions returns 0 if all is OK or -1 if a is a singular matrix.

See also:ArrMul(), MATInv()

MATTrans()

This transposes a matrix (a two dimensional array), swapping the rows and columns. This procedure physically moves the data, unlike the trans() or ` operator, which remaps the matrix without moving any data. It is usually much more efficient to use trans().

Proc MATTrans(mat[][]{, src[][]});

- mat A m by n matrix returned holding the transpose of src. If src is omitted, m must be equal to n and the rows and columns of mat are swapped.
- src Optional, a n by m matrix to transpose.

See also:ArrAdd(), MATMul()

Max()

This function returns the index of the maximum value in an array, or the maximum of several real and/or integer variables.

Func Max(arr[]|arr%[]|val1 {,val2 {,val3...}});

arr A real or integer array.

valn A list of real and/or integer values to scan for a maximum.

Returns The maximum value or array index of the maximum.

Maxtime()

In a time view, this returns the maximum time in seconds in the file or in a channel, or the current sample time during sampling. In a result view, it returns the number of bins.

Func MaxTime({chan%});

- chan% Optional channel number for time views, ignored in result views. If present, the function gets the maximum time in the channel ignoring any marker filter. If there is no data in the channel, the return value is -1.
- Returns The value returned is negative if the channel does not exist. If the current view is of the wrong type the script stops with an error.

To find the time of the last item in the marker filter on a channel with a marker filter set:

time := LastTime(chan%, MaxTime(chan%)+1.0);

With a marker filter set, Spike2 has to search the data to find a marker that is in the filter. MaxTime() returns the last data item in the channel regardless of the marker code. The +1.0 is because LastTime() finds data before the search time.

See also:Len(), LastTime(), NextTime(), Seconds()

MeasureChan()

This function adds or changes a measurement channel in an XY view created with MeasureToXY() using the settings previously defined by using MeasureX() and MeasureY(). The XY view must be the current view. This command implements some of the functionality of the XY plot setting dialog.

Func MeasureChan(chan%, name\${, pts%});

- chan% This is 0 to add a new channel or the number of an existing channel to change settings. MeasureToXY() creates an XY view with one channel, so you will usually call this function with chan% set to 1. You can have up to 32 measurement channels in the XY view.
- name\$ This sets the name of the channel and can be up to 9 characters long.
- pts% Sets the maximum number of points for this channel, if omitted or set to zero then all points are used. When a points limit is set and more points are added, the oldest points are deleted.

Returns The channel number these settings were applied to or a negative error code.

See also:CursorActive(), MeasureToXY(), MeasureX(), MeasureY()

MeasureToChan()

This creates a new Event, Marker or RealMark channel with an associated measurement process and cursor 0 iteration method for the current time view. The Process() command adds items to the new channel based on active cursor measurements. Before calling this function, use MeasureX() to set the measurement method to 102 (Time) and set expr1\$ to generate the time stamps of the new items. If you select any other method this command will return an error code. The iteration must produce times in ascending order unless the output is to a memory channel.

If you are creating a RealMark channel, you must call MeasureY() before this function to define the measurement to attach to each data item added to the new channel.

- dest% This is the output channel number or zero for the lowest numbered, unused memory channel. Set 1 to 400 for channels in the data file and 401 to 800 for specific memory channels. It is a fatal error to set a channel that is in use.
- name\$ The output channel name. Channel units, if required, are inferred from chan% and the measurement method.
- type% This sets the output channel type, as for ChanKind(). The allowed types are: 2 or 3 for Event, 5 for a Marker, and 7 for RealMark. If you wish to store data other than event times in the file, you should set type% to 7.
- mode% This is the cursor 0 iteration mode. Modes are the same as in CursorActive() but not all modes can be used. Valid modes are:

4 Peak find	8 Falling threshold	14 Data points	19 Outside levels
5 Trough find	12 Slope peak	16 Expression	20 Inside levels
7 Rising threshold	13 Slope trough	17 Turning point	

- chan% This is the channel searched by the cursor 0 iterator. In expression mode (16), this is ignored and should be set to 0.
- min This is the minimum allowed step for cursor 0 used in all modes except 16.
- exp This is the expression that is evaluated in mode 16.
- Iv This number or string expression sets the threshold level for threshold modes and the peak size for peak and trough modes. It is in the y axis units of channel chan% or y axis units per second for modes 12 and 13. In mode 14 it sets the points as a number and defaults to 1. This argument is ignored and should be 0 or omitted for modes that do not require it.
- hw This sets the hysteresis level in y axis units for threshold modes. This argument is ignored and should be 0 or omitted for modes that do not require it. For backwards compatibility, if width is omitted, it sets the width in seconds for slope measurements.
- flgs% This is the sum of option flags. Add 1 to force a common x axis Add 2 for user checks on the cursor positions. The default value is zero.
- qu\$ This sets the qualification expression for the iteration. If left blank then all iteration positions will be used. If not blank, and it evaluates to non-zero, then the iteration is skipped.
- width This sets the width, in seconds, for slope measurements. Set this 0 or omit it in modes that do not require it.
- 1v2 This number or string expression is used with 1v to set the two threshold levels for cursor iteration modes 19 and 20.

Returns The function result is the number of the created channel.

See also:CursorActive(), MeasureChan(), MeasureX(), MeasureY(), Process()

MeasureToXY()

This creates a new XY view with a measurement process and cursor 0 iteration method for channels in the current time view. It creates one output channel with a default measurement method. Use MeasureX(), MeasureY() and MeasureChan() to edit the method and add channels. Use Process() to generate the plot. The new XY view is the current view and is invisible. Use WindowVisible(1) to make it visible. These commands implement the functionality of the Measurements to XY view dialog.

mode% This is the cursor 0 iteration mode. Modes are the same as in CursorActive() but not all modes can be used. Valid modes are:

4 Peak find	8 Falling threshold	14 Data points	19 Outside levels
5 Trough find	12 Slope peak	16 Expression	20 Inside levels
7 Rising threshold	13 Slope trough	17 Turning point	

- chan% This is the channel searched by the cursor 0 iterator. In expression mode (16), this is ignored and should be set to 0.
- min This is the minimum allowed step for cursor 0 used in all modes except 16.
- exp\$ This is the expression that is evaluated in mode 16.
- 1v This number or string expression sets the threshold level for threshold modes and the peak size for peak and trough modes. It is in the y axis units of channel chan% or y axis units per second for modes 12 and 13. In mode 14 it sets the points as a number and defaults to 1. This argument is ignored and should be 0 or omitted for modes that do not require it.
- hw This sets the hysteresis level in y axis units for threshold modes. This argument is ignored and should be 0 or omitted for modes that do not require it. For backwards compatibility, if width is omitted, it sets the width in seconds for slope measurements.
- flgs% This is the sum of option flags. Add 1 to force a common x axis Add 2 for user checks on the cursor positions. The default value is zero.
- qu\$ This sets the qualification expression for the iteration. If left blank then all iteration positions will be used. If not blank, and it evaluates to non-zero, then the iteration is skipped.
- width Set this value in seconds; use 0 or omit it in modes that do not require it. For slopes it sets the time over which the slope is measured. For peaks and troughs, it sets the maximum peak width (use 0 for no maximum). For threshold modes 7, 8, 19 and 20 it sets the minimum crossing time (Delay in the dialog).
- 1v2 This number or string expression together with 1v sets the two threshold levels for cursor iteration modes 19 and 20.

Returns The function result is an XY view handle or a negative error code.

Arguments passed as strings are not evaluated until data is processed. Invalid strings generate invalid measurements and no data points in the XY view.

Example This code generates a plot of the peak values from channel 1 of the current time view. Peaks must be at least 0.1 seconds apart and the data must fall by at least 1 y axis unit after each peak.

```
var xy%; 'Handle of new xy view
xy%:=MeasureToXY(4, 1, 0.1, 1); 'Peak, chan 1, min step 0.1, amp 1
WindowVisible(1); 'Window is invisible, so show it
MeasureX(102, 0, "Cursor(0)"); 'x = Time, no channel, at cursor 0
MeasureY(100, 1, "Cursor(0)"); 'y = Value of chan 1 at cursor 0
MeasureChan(1, "Peaks", 0); 'Set the title, no point limit
Process(0.0, View(-1).MaxTime(),0,1); 'Process all the data
See also:CursorActive(), MeasureChan(), MeasureX(), MeasureY(), Process()
```

MeasureX(

The MeasureX() and MeasureY() functions set the x and y part of a measurement for a measurement channel. The settings are saved but have no effect until MeasureChan() or MeasureToChan() are used to change or create a channel. This command implements some of the functionality of the XY plot setting dialog. The current view must be the target of the measurements.

Func MeasureX(type%, chan%, expr1|coef% {,expr2 {,width}});
Func MeasureY(type%, chan%, expr1|coef% {,expr2 {,width}});

type% This sets the x or y measurement type. Values less than 100 match those for the ChanMeasure() command (types 5 and 1 are identical for this command). If you are using MeasureX() for use with MeasureToChan(), the only valid type% is 102.

 Area 5 Area (scaled) Mean 6 Curve area Slope 7 Modulus Sum 8 Maximum 	9 Minimum 10 Peak to Peak 11 RMS Amplitude 12 Standard deviation	13 Abs max.14 Peak15 Trough16 RMS error	17 Mean in X 18 SD in X 19 Mean of abs
Values from 100 up are: 100 Value 101 Value difference 102 Time 103 Time difference	104 0-based fit coeffi 105 User entered valu 106 Expression 107 Value ratio	cient 108 Va ue	alue product

- chan% This is the channel number for the measurement. For time, user entered and expression measurements it is ignored and should be set to 0.
- expr1 Either a real value or a string expression that evaluates as the start time for measurements over a time range, as the position for time (102) and value measurements and as the expression used for measurement type 106.
- coef% The zero-based coefficient number for measurement type 104.
- expr2 Either a real value or a string expression that evaluates as the end time for measurements over a time range and as the reference time for position for single-point measurements and differences. Set this to an empty string when width is required and this is not.
- width This is the measurement width for value and value difference measurements. The default value is zero.

Returns The function return value is zero or a negative error code.

See also:CursorActive(), MeasureChan(), MeasureToChan(), MeasureToXY()

MeasureY()

This is identical to MeasureX() and sets the y part of a measurement for a measurement channel. The settings are saved but have no effect until MeasureChan() is used to change or create a channel. See the MeasureX() documentation for details.

Func MeasureY(type%, chan%, expr1\$ {,expr2\$ {,width}});

See also:MeasureX()

MemChan()	This fur time vie	nction creates a new channel in memory and a w. You can have up to 400 memory channels p	ttaches it to the file in the current er file.
	Func M	emChan(type%{, size%{, binsz{, pre%	{, trace%}}});
	type%	The type of channel to create. Codes are:1 Waveform4 Level (Event+-)2 Event (Event-)5 Marker (default)3 Event (Event+)6 WaveMark	7 RealMark 8 TextMark 9 Real wave
	size%	For TextMark, RealMark and WaveMark c characters, reals or waveform points to attach other channel types and should be set to 0.	channels this sets the number of a to each item. It is ignored for all
	binsz	Used for waveform and WaveMark data to sp waveform points. This is rounded to the near resolution. If you set this 0 or negative, the sm	ecify the time interval between the est multiple of the underlying time hallest bin size possible is set.
	pre%	This must be present for WaveMark data to se	t the number of pre-trigger points.
	trace%	Optional, default 1, sets the number of WaveM	Mark traces in range 1 to 4.
	Returns	Either the channel number of the newly create channels, or a negative error code.	ed channel, or 0 if there are no free
	Channel these wi ChanOf wavefor	s created in this way are given default titles, ith the ChanTitle\$(), ChanUnits\$(), Char fset() routines. The following code create m channel) as a memory channel:	units and comments. You can set aComment\$(), ChanScale() and as a copy of channel wChan% (a
	<pre>func C var mc if Cha mc% := if mc% Cha Cha Cha Cha Cha cha endif; return end;</pre>	<pre>opyWave%(wChan%) 'Copy waveform to %; nKind(wChan%)<>1 then return 0 endi MemChan(1,0,BinSize(wChan%)); >0 then nScale(mc%, ChanScale(wChan%)); nOffset(mc%, ChanOffset(wChan%)); nUnits\$(mc%, ChanUnits\$(wChan%)); nTitle\$(mc%, "Copy"); nComment\$(mc%, "Copied from channel Import(mc%, wChan%, 0, MaxTime()); nShow(mc%); mc%; 'Return the new channel num</pre>	a memory channel f; 'Not a waveform! 'Create waveform channel 'Created OK? 'Copy scale 'and offset 'and units 'set our own title "+Str\$(wChan%)); 'Copy data 'display new channel ber
	See also	:ChanNew(), MemDeleteItem(), MemDelete MemImport(), MemSave(), MemSetItem()	eTime(),MemGetItem(),
MemDeleteltem()	This fun the entir	action deletes one or more items from a channel re channel use ChanDelete().	l created by MemChan (). To delete
	Func M	<pre>emDeleteItem(chan% {,index% {,num%}</pre>	});

chan% The channel number of a channel created by MemChan().

- index% The ordinal index into the channel to the item to delete. The first item is numbered 1. If you specify index -1 or omit this argument, all items are deleted.
- num% The number of items to delete from index%. Default value is 1.

Returns The number of items deleted or a negative error code.

MemDeleteTime()

This function deletes one or more items from a memory channel based on a time range. If a marker filter is set, only items in the filter are deleted unless you add 4 to mode%.

Func MemDeleteTime(chan%, mode%, t1 {,t2});

chan% The channel number of a channel created by MemChan().

mode% This sets the items to delete and how to interpret the time range t1 and t2:

- 0 A single item is deleted. t1 is the item time and t2 is a timing tolerance (0 if t2 is omitted). The nearest item to t1 in the time range t1-t2 to t1+t2 is deleted. If there are no items in the range, nothing is deleted.
- 1 Delete all items from time t1-t2 to t1+t2. If t2 is omitted, 0 is used.
- 2 Delete the first item from time t1 to t2. If t2 is omitted it is taken as t1.
- 3 Delete all items from time t1 to t2. If t2 is omitted, it is taken as t1.
- +4 If you add 4 to the mode, any marker filter for the channel is ignored.

t1, t2 Two times, in seconds, that set the time range for items to delete.

Returns The number of items deleted, or a negative error code.

See also:MemChan(), MemDeleteItem(), MemGetItem(), MemSetItem()

	This ast	une information should a manager should it on The items are identified by their
MemGetItem()	ordinal p	position in the channel, not by time, and any mask set for markers is ignored.
	Func M Func M	emGetItem(chan% {,index% {,code%[] {,&data\$ data[]}}); emGetItem(chan%, index%, &wave &wave% wave[] wave%[]{,&n%});
	chan%	The channel number of a channel created by ${\tt MemChan}$ ().
	index%	The ordinal index into the channel to the item required. The first item is numbered 1. If you omit the index, or specify index 0, the function returns the number of items in the channel.
	The rem	aining fields are only allowed if the index is non-zero.
	code%	This is an integer array of at least 4 elements that is returned holding the channel marker codes. If there are no markers, the codes are all set to 0.
	data	This must be a string variable for TextMark data or a real array for RealMark or WaveMark data. It is returned holding the data from the item. If the data type is incorrect for the channel, it is not changed. For an array, the points returned is the smaller of the size of the array and the number of values in the item. You can use a two dimensional array to collect WaveMark data with multiple traces. The second dimension is for the traces: var data[points%][traces%];
	wave	This argument collects waveforms from waveform channels. If a real variable or array is passed, the waveform data is in user units. If an integer variable or array is passed, the data is a copy of the 16-bit integer data used by Spike2 to store waveforms. The maximum number of elements copied is the array size.
		When an array is used, only contiguous data is returned. A gap in the data (when the interval between two points is greater that the sample interval for the channel) terminates the data transfer.
		You may find that ChanData() is easier to use when you want to collect waveform (or event) data based on a time range.
	n%	If the previous argument is an array this optional argument returns the number of data items copied into the array.
	Returns	For index% of 0 or omitted, the function returns the number of items in the channel. If an index is given that is outside the range of items present, the function returns -1. Otherwise it returns the time of the item.
	See also	:ChanData(),MemChan(),MemDeleteItem(),MemDeleteTime(),

MemImport(), MemSave(), MemSetItem()

MemImport()

This function imports data into a channel created by MemChan(). There are some restrictions on the type of data channel that you can import from, depending on the type of the destination channel.

Destination	Source	Restrictions
Waveform	Waveform	Data copied, but must match sample interval
	All others	Not available, see EventToWaveform()
Event	Waveform	Can extract event times
	All others	Times are extracted from the channel.
Level	Waveform	Can extract event times
	All others	Times are extracted from the channel. First time in
		destination is assumed to be a low to high transition.
Marker	Waveform	Can extract event times, coded for peak/trough etc.
	Event, Level	Marker codes all set to 0.
	All others	Marker codes are copied.
TextMark	Waveform	Can extract event times, coded for peak/trough etc.
	Event, Level	Copies times, marker codes set 0, empty strings
	TextMark	Copies all data, strings may be truncated if too long.
	All others	Copies marker information, empty strings
RealMark	Waveform	Can extract event times, coded for peak/trough etc.
	Event, Level	Times copied, marker codes and reals set to 0
	RealMark	Copied, real data truncated or zero padded as needed
	WaveMark	Copied, waveform to reals, padded/truncated
	All others	Marker portion copied, reals filled with 0
WaveMark	Waveform	Special option, event channel marks waveforms
	Event, Level	Copies times, marker codes set to 0, waveform set 0
	WaveMark	Copies all data, waveforms aligned on trigger point
	All others	Copies marker, waveform filled with zeros

You can extract events or markers from waveform data using peak search and level crossing techniques. You can convert waveform data to WaveMark data with a special option that chops sections of waveform data out based on event times on a third channel.

Import compatible channel The first command variant imports data from a compatible channel. All event and marker channels can be copied to each other, but the information transferred is the lowest common denominator of the two channel types. Missing data is padded with zeros. Waveform data is compatible with itself if both channels have the same sampling rate. Func MemImport (chan%, inCh%, start, end); chan% The channel number of a channel created by MemChan(). inCh% The channel number to import data from. start The start time to collect data from. The end time to collect data up to (and including). end Extracting events from The mode%, time, level and code% arguments are used when extracting events from waveforms waveform data. The level crossing modes use linear interpolation between points to find the exact time of the waveform crossing. The peak and trough modes fit a parabola to the three points around the peak or trough to estimate the time more accurately. When extracting events to a WaveMark channel from waveform data, the time saved is the time of the start of the waveform section, not the peak/trough or level crossing time. The saved waveform data starts the number of points before each event set by the pre% parameter to MemChan(). The WaveMark time is adjusted to match the waveform. If both channels do not have the same sampling rate, the copied waveform is set to zero.

Func MemImport(chan%,inCh%,start,end{,mode%,time,level{,code%}});

	mode% time level code%	 The mode of data extraction. The modes are: 0 Extract events based on the time of a peak in the waveform. If the destination is a marker, these events are coded as 2 unless code% is set. 1 Extract events based on the time of a trough in the waveform. If the destination is a marker, the events are coded 3 unless code% is set. 2 Extract events based on times when the waveform rises through level. If the destination is a marker, the events are coded 4 unless code% is set. 3 Extract events based on times when the waveform falls through level. If the destination is a marker, the events are coded 5 unless code% is set. 3 Extract events based on times when the waveform falls through level. If the destination is a marker, the events are coded 5 unless code% is set. The minimum time between detected events. Use this to filter noisy signals. In modes 0 and 1, this is the distance that the waveform must fall after a peak or rise after a trough. In modes 2 and 3 it is the level to cross to detect an event. If present and positive it overrides the codes based on mode% that are applied to the events. The low byte of code% sets the first marker code; the remaining
WaveMark from events and waveform data	The spe the way as set for before e	marker codes are always 0. ecial mode to convert waveform to WaveMark data uses an extra channel to mark reform sections to be extracted. The waveform must have the same sampling rate or the WaveMark channel. The saved waveform data starts the number of points each event set by the pre% parameter to MemChan().
	Func M	<pre>MemImport(chan%, inCh%, start, end, eCh%);</pre>
	eCh%	A channel holding event times to mark the waveform sections to extract. The time saved is the time of the first point in each waveform section. If the channel contains marker codes, these are copied to the memory channel.
	Returns	It returns the number of items added to the channel, or a negative error code.
	See also):MemChan(), MemGetItem(), MemSave(), MemSetItem()
MemSave()	This w window ChanDe	rites a channel created by MemChan() to the data file associated with the current y, making the data permanent. The memory channel is not changed; use elete() to remove it.
	Func M	<pre>lemSave(chan%, dest%{, type%{, query%{, bufSz%}});</pre>
	chan%	A channel created by the MemChan () function
	dest%	The destination channel in the file. This must be in the range 1 to the maximum number of channels allowed in the data file.
	type%	The type of data to save the data as. The type selected must be compatible with the data in the memory channel. Codes are:
		0 Same type (default)3 Event (Evt+)6 WaveMark9 RealWave1 Waveform4 Level (Evt+-)7 RealMark2 Event (Evt-)5 Marker8 TextMark
		The special code -1 means append the memory channel to an existing channel. The new data must occur after the last item in the dest% channel and the dest% channel must be of a compatible type to the memory channel.
	query	If this is not present or zero, and the dest% channel is already in use, the user is queried about overwriting it. If this is non-zero, no query is made.
	bufSz१	If present and greater than 0, this sets the disk block size in bytes for waveform and RealWave data. If omitted (recommended) or zero, a default size is used.
	Returns	The number of items written, or a negative error code.
	See also	<pre>D:ChanDelete(),ChanWriteWave(),MemChan(),MemImport()</pre>

MemSetItem()	This function edits or adds an item in a channel created by MemChan(). The item is identified by its ordinal position in the channel and any mask set for markers is ignored.
	<pre>Func MemSetItem(chan%, index%, time{, code%[]{, data\$ data[]}}); Func MemSetItem(chan%, index%, time, wave wave% wave[] wave%[]);</pre>
	chan% The channel number of a channel created by MemChan().
	index% The index of the item to edit. The first item is number 1. An index of 0 adds a new item to the buffer at a position set by time (which must be positive).
	time The item time, or -1 for no change. If index% is 0, you must supply a time. For a waveform channel, it sets the time of the first data point but if there is already data in the channel, time is adjusted by up to half the sampling interval to be compatible with the sampling interval of the channel and the existing data.
	code% This is an integer array of at least 4 elements that hold the marker codes for the channel. If the channel does not require marker codes, this argument is ignored. If this parameter is omitted for a channel with markers, the codes are set to 0.
	data A string for TextMark data or a real array for RealMark or WaveMark, holding the item data. If the data type is incorrect it is ignored. The number of points or characters set is the smaller of the number passed and the number expected. For RealMark and WaveMark data, if the array is too short, the extra values are unchanged when index% > 0 (editing) and have the value 0 if index% is 0.
	WaveMark and waveform real values are limited to the range -5*scale+offs to 4.99985*scale+offs. A two dimensional array is allowed for WaveMark data: var_data[point.s%][traces%]: passed in as data[][].
	Wave For a waveform you can set one value, or an array. Real values are limited as described above. For integers, the lower 16-bits of the 32-bit integer are copied to the channel (values greater than 32767 or less than -32768 will overflow).
	Returns The function returns the index at which the data was stored. If an index is given that is outside the range of items present, the function returns -1.
	See also:ChanWriteWave(), MemChan(), MemGetItem(), MemImport(), MemSave()
Message()	This function displays a message in a box with an OK button that the user must click to remove the message. Alternatively, the user can press the Enter key.
	<pre>Proc Message(form\$ {,arg1 {,arg2}});</pre>
	form\$ A string that defines the output format as for Print(). If the string includes a vertical bar, the text before the vertical bar is used as the window title.
	arg1, 2 The arguments used to replace %d, %f and %s type formats.
	You can split the message into multiple lines by including n in the forms string. Long messages are truncated.
	<pre>See also:Print(), Input(), Query(), DlgCreate()</pre>
Mid\$()	This function returns a sub-string of a string.
	<pre>Func Mid\$(text\$, index% {,count%});</pre>
	text\$ A text string.
	index% The starting character in the string. The first character is index 1.
	count & The maximum characters to return. If omitted, there is no limit on the number.
	Returns The sub-string. If Index% is larger than Len(text\$), the string is empty.
	<pre>See also:Asc(), Chr\$(), DelStr\$(), InStr(), LCase\$(), Left\$(), Len(), Print\$(), Right\$(), Str\$(), UCase\$(), Val()</pre>

Min()

This function returns the index of the minimum value in an array, or the minimum of several real and/or integer variables.

```
Func Min(arr[]|arr%[]|val1 {,val2 {,val3...}});
```

arr A real or integer array.

valn A list of real and/or integer values to scan for a minimum.

Returns The minimum value or array index of the minimum.

Minmax()

Minmax() finds the minimum and maximum values for result views and time view channels with a y axis, or the minimum and maximum intervals for an event or marker channel drawn as dots or lines. Min() and Max() are preferred in result views.

- chan% The channel number in the time or result view.
- start The start position in time for a time view, in bins for a result view.
- finish The end position in time for a time view, in bins for a result view.
- min The minimum value is returned in this variable.
- max The maximum value is returned in this variable.
- minP The position of the minimum is returned in this variable.
- maxP The position of the maximum is returned in this variable.
- mode% If present, this sets the drawing mode in which to find the minimum and maximum. If mode% is absent or inappropriate, the display mode is used. This parameter is ignored in a result view. The modes in a time view are:
 - 0 The standard mode for the channel.
 - 1 Dots mode for events. The dotSz% argument can be used.
 - 2 Lines mode.
 - 3 Waveform mode. This is the only mode for waveform channels.
 - 4 WaveMark mode.
 - 5 Rate mode. The binSz argument sets the width of each bin.
 - 6 Mean frequency mode. binSz sets the time period.
 - 7 Instantaneous frequency mode. dotSz% can be used.
 - 8 Raster mode. trig% sets the trigger channel, dotSz% is used.
- binSz This sets the width of the rate histogram bins and the smoothing period for mean frequency mode when specifying your own mode.
- trig% The trigger channel for raster displays, level data raster displays are impossible.
- edge% For level data event channels. It sets which edges of the level signal are used for mean frequency, instantaneous frequency and rate modes. The values are:
 - 0 Use both edges (same as omitting the parameter).
 - 1 Use rising edges.
 - 2 Use falling edges

Returns Zero if all was well or a negative error code.

See also:Min(), Max(), XYRange()

MMAudio()

This returns information about the audio content of the current multimedia window.

Func MMAudio({&rate});

rate If present, this is returned as the sample rate per channel, in Hz.

Returns The number of audio channels or 0 if no audio or not a multimedia window.

See also:MMOpen(), MMImage(), MMPosition(), MMVideo()

MMImage()

This copies the image at the current position in the current multimedia view to one or more data arrays. The two-dimensional arrays that collect the image are treated as x[v][h], where v is a vertical co-ordinate and h is a horizontal co-ordinate. x[0][0] is at the bottom left of the image. Use MMVideo() to get the image size. Sections of the array that do not map onto the image are not changed.

```
Func MMImage(rgb%[][]);
Func MMImage(map[][]|mode%{, r[][]{, g[][]{, b[][]}});
Func MMImage(-1, h[][]{, l[][]{, s[][]}});
```

- rgb% If the first argument is a two-dimensional integer array, it is returned with each array element holding one pixel in RGB32 format. Bits 0-7 hold the Red intensity, bits 8-15 hold the Green and bits 16-23 hold the Blue. Bits 24-31 are 0. Any additional arguments are ignored.
- mode% If the first argument is 1, the r,g and b arrays are filled with Red, Green and Blue intensities. If mode% is set to 2, the arrays are filled with Cyan, Magenta and Yellow intensities. The colour components will be in the range 0 to 1.0, inclusive. If mode% is 0, the r array is filled with a monochrome representation of the data with 0 representing black and 1.0 representing white; the g and b arrays are set to 0, if they are present.



If the first argument is -1, the h, 1 and s arrays return HLS (Hue, Lightness and Saturation) data. Hue is in degrees around a colour hexagon with

Blue at 0, Red at 120 and Green at 240 degrees. Lightness represents the amount of white in a colour, from 0 (black) to 1.0 (white). Saturation is a measure of the purity of a colour, from 0 (gray) to 1.0, which is a pure colour.

map If this is present, it should be a 3x3 array that maps the red, green and blue (RGB) in the image into the r, g and b arrays that you provide. The first array index sets the array: 0 for r, 1 for g and 2 for b. The second array index sets the colour: 0 for red, 1 for green and 2 for blue. The array values set how much of each colour should be added into the result.

If R, G and B (all in the range 0.0 to 1.0), represent the colour of a pixel, the r array pixel is set to R*map[0][0]+G[]*map[0][1]+B*map[0][2]. The g and b arrays (if present) are calculated with the first index of map set to 1 and 2.

- r,g,b These arrays, when present are filled with RGB colour information. If map% is used, the values in the arrays depend on the map values. If mode% is used, the values will be in the range 0 to 1.
- h, l, s These arrays are used when the first argument is -1. They are returned holding the hue (0 to 360), lightness (0 to 1.0) and saturation (0 to 1.0).

Returns 1 if the copy succeeded and 0 if it failed.

See also:EditCopy(), MMOpen(), MMPosition(), MMVideo()

MMOpen(

This counts, opens and returns the view handles of multimedia windows associated with the current time or multimedia view. Multimedia files associated with fname.smr are named fname-1.avi, fName-2.avi and so on. If a multimedia window is opened or located, it becomes the current view. If multiple windows open, the last-opened window becomes the current view. This does not access the s2video application window.

Func MMOpen({n%{, mode%}});

- n $\$ Set to 0 to open or count all associated multimedia files. Set to n (greater than 0) for the nth multimedia file. It is not an error to open a file that is already open.
- mode% This optional argument has the default value 0 and is the sum of the following:
 - 1 Make new windows visible. Omitting this does not hide an existing window.
 - 2 Reserved.
 - 4 Do not open any new windows, do not change the current view. Use this to find a view handle or count the multimedia windows.
- Returns If n% is 0 or omitted, the return value is the number of associated multimedia files that are open. If n% is greater than 0, and a file was opened, or was already open, the return value is the view handle.

See also:FileClose(), MMAudio(), MMImage(), MMPosition(), MMVideo()

MMPosition()

This sets and gets the play position and state of the current multimedia window or sets the state of all multimedia windows associated with the current time view. Multimedia windows track any change in the associated time view position. Use this command to change the play position or state independently of the time view. Use Rerun() to replay a time view and all associated multimedia windows together. Multimedia windows play in real time or at the rate set by the last Rerun() to the window.

Func MMPosition({pos{, sPlay%{, &gPlay%}}});

- pos The time, in seconds, in the multimedia stream to move to. Omit this argument or set it negative for no position change.
- sPlay% This argument controls the play state. Omit it or set it to -1 for no change in the play state. Set 0 to stop playing and set 1 to start playing.
- gplay% If this argument is present it is returned holding the play state after the command has run as 0 for not playing and 1 for playing.
- Returns The return value is the play position after any changes made by the command have been made or -1 if an error occurred or no window was found. If the command is applied to a time view, the return value is for the lowest numbered multimedia window found.

See also:MMAudio(), MMImage(), MMOpen(), MMVideo(), Rerun()

MMRate()

This requests the attached s2video applications to change the video frame rate. You cannot run faster than the camera frame rate. The rate is achieved by dropping frames from the video stream; the result may be faster or slower than you request.

Func MMRate(fps);

fps The desired frames per second. Rates of 0 or less set the minimum of the camera rate or the maximum frame rate set in the s2video application.

Returns The number of external windows that this request was sent to.

See also:MMAudio(), MMImage(), MMOpen(), MMPosition(), MMVideo()

MMVideo()	This returns information about the video content of the current multimedia window. Func MMVideo({&hPix%{, &vPix%{, &fps}}}); hPix% If present, returned as the number of horizontal pixels in the image.
	VPix% If present, returned as the number of vertical pixels in the image.
	fps If present, returned as the frames per second value held in the image file.
	Returns 0 if no video is present or this is not a multimedia window or 1 if there is video in the multimedia window.
	See also:MMAudio(), MMImage(), MMOpen(), MMPosition()
MoveBy()	This moves the text caret in a text window relative to the current position by lines and/or a character offset. You can extend or cancel the current selection.
	<pre>Func MoveBy(sel%{, char%{, line%}});</pre>
	sel% With char% present, if sel% is zero, all selections are cleared. If non-zero the selection is extended to the destination of the move. With char% omitted sel%=0 returns the character offset, 1 the line number and 2 the column.
	char* If line* is absent, the new position is obtained by adding char* to the current character offset in the file. You cannot move the caret beyond the existing text.
	line% If present it specifies a line offset. The new line is the current line number plus line% and the new character position is the current character position in the line plus char%. The new line number is limited to the existing text. If the new character position is beyond the start or end of the line it is limited to the line.
	Returns It returns the new position. MoveBy(1,0) returns the current position without changing the selection. See sel% (above) to get line and column numbers.
	See also:MoveTo(), Selection\$()
MoveTo()	This moves the text caret in a text window. You position the caret by lines and/or a character offset. You can extend or cancel the text selection. The first line is number 1. Func MoveTo(sel%, char%{, line%});
	sel% If zero, all selections are cleared. If non-zero the selection is extended to the
	destination of the move and the new current position is the start of the selection.
	char% If line% is absent, this sets the new position in the file. You cannot move the caret beyond the existing text. 0 places the caret at the start of the text.
	line% If present it specifies the new line number. The new line number is limited to the existing text. char% sets the position in this line (and is limited to the line).
	Returns The function returns the new position in the file.
	See also:MoveBy(), Selection\$()
NextTime()	This finds the next item on a channel after a time. If a marker filter is in use, only data that is included in the filter is visible. It is an error to use this function in a result view.
	<pre>Func NextTime(chan%, time{,&val code%[]{,data[] data%[] &data\$}});</pre>
	chan% The channel number in the view to use.
	time The time to start the search after. Items at the time are ignored. To ensure that items at time 0 are found, set the start time of the search to a negative time.

- val Optional: for waveform channels it returns the waveform value. For event level channels, it is returned 0 if the transition is low to high, and 1 if the transition is high to low. If there is no event it returns the level at time; 0 for low, 1 for high.
- code% This optional parameter is only used if the channel is a marker type. This is an array with at least four elements that is filled in with the marker codes.
- data Gets data from RealMark and WaveMark channels. If there is insufficient data, unused entries are unchanged. Integer arrays are for WaveMark channels and return the 16-bit data that holds the waveform. If WaveMark data has multiple traces, use a two dimensional array. The trace number is the second index.
- data\$ A string returned holding the text from a TextMark channel.

Returns The time of the next item, -1 if there are no more items or a negative error code.

See also:ChanData(), LastTime(), MaxTime()

Optimise()

This optimises y axes over an x axis range in time, result or XY views in the same way as the Y Range dialog optimise button. You can optimise a hidden channel. In an XY view, all channels share the same y range, so optimising affects all channels. You can also use this on a spike shape window (with no arguments) to optimise the display.

Proc Optimise({chan%{, start{, finish}}});

- chan% The channel to optimise. We also allow -1 for all, -2 for visible and -3 for all selected channels. If omitted, -1 is used, for all channels.
- start The start of the region to optimise. This is in x axis units for a time or XY view and in bins for a result view. If omitted, this is the start of the window.
- finish The end of the region to optimise. If omitted, this is the end of the window.

See also:ChanOffset(), ChanScale(), YRange(), YLow(), YHigh()

|--|

This reads back the percentages of red, green and blue in a colour in the palette.

Proc PaletteGet(col%, &red, &green, &blue);

- col% The colour index in the palette in the range 0 to 39.
- red The percentage of red in the colour.
- green The percentage of green in the colour.
- blue The percentage of blue in the colour.

See also:ChanColour(), Colour(), PaletteSet(), XYColour()

PaletteSet()	This set cannot model. Black is	his sets the colour of one of the 40 palette colours. Colours 0 to 6 form a grey scale and annot be changed. Colours are specified using the RGB (Red, Green, Blue) colour odel. For example, bright blue is $(0\%, 0\%, 100\%)$. Bright yellow is $(100\%, 100\%, 0\%)$. lack is $(0\%, 0\%, 0\%)$ and white is $(100\%, 100\%, 100\%)$.			
	Proc E	<pre>?aletteSet(col%, red, green, blue {,solid%});</pre>			
	col%	The colour index in the palette in the range 0 to 39. Attempting to change a fixed colour or a non-existent colour has no effect.			
	red	The percentage of red in the colour.			

- green The percentage of green in the colour.
- blue The percentage of blue in the colour.

solid% If present and non-zero, sets the nearest solid colour (all pixels have the same hue in a solid colour). Systems that don't need to do this ignore solid%.

See also:ChanColour(), Colour(), PaletteGet(), XYColour()

PCA()

This command performs Principal Component Analysis on a matrix of data. This can take a long time if the input matrix is large.

Func PCA(flags%, x[][]{, w[]{, v[][]}});

flags% Add the following values to control pre-processing of the input data:

- 1 Subtract the mean value of each row from each row
- 2 Normalise each row to have mean 0.0 and variance 1.0
- 4 Subtract the mean value of each column from each column
- 8 Normalise each column to have mean 0.0 and variance 1.0

You would normally pre-process the rows or the columns, not both. If you set flags for both, the rows are processed first.

- x[][] A m rows by n columns matrix of input data that is replaced by the output data. The first array index is the rows; the second is the columns. There must be at least as many rows as columns ($m \ge n$). If you have insufficient data you can use a square matrix and fill the missing rows with zeros. If you were computing the principal components of spike data, on input, each row would be a spike waveform. On output, each row holds the proportion of each of the n principal components scaled by the w[] array that, when added together, would best (in a least-squares error sense) represent the input data.
- w[] This is an optional array of length at least n that is returned holding the variance of the input data that each component accounts for. The components are ordered such that w[i] $\geq w[i+1]$.
- v[][] This is an optional square matrix of size n by n that is returned holding the n principal components in the rows.
- Returns 0 if the function succeeded, -1 if m < n, -2 if w has less than n elements or v has less than n rows or columns.

You can find an overview of PCA in the *Clustering spikes* chapter of the Spike2 manual. In the terms of that overview, x[][] corresponds with the X matrix on input and the U matrix on output, w[] is the diagonal of W and v[][] is the matrix V.

PlayOffline()

This command plays an area of the current time view through the 1401 DACs or through the sound card in the computer. If you close the view during sampling, output will cease. You cannot use the 1401 for sampling and PlayOffline() simultaneously. You cannot use the sound card if it is already in use. If you choose the sound card, the default wave out device is used.

Unless the waveform to replay is small (total points less than 32000), it is copied to the output device in chunks, as required. This operation happens in the background, but you must make time for it by allowing Spike2 to idle, or by checking the play position.

There are two versions of the command; the first starts output and the second reports the play position and allows you to stop the output early.

func PlayOffline(cSpc,dacs%,sTime,eTime{,rep%{,scale{,flag%}}}); func PlayOffline({what%{, &rep%}});

cSpc A channel specification for up to 4 channels to be played. These can be any mix of waveform, WaveMark or RealWave channels. The channels can have different sample rates; the rate of the first channel is taken as the rate of all

channels and the data for the other channels is matched to the first by linear interpolation. Gaps in the data are output as zeros.

dacs% This is either an integer array of 1401 DAC channel numbers (the first DAC is numbered 0), or a single integer, being the number of the first DAC channel to use. If a single integer is used, and more than 1 channel is in use, the second, third and fourth channels use DACs (dacs%+1) mod 4, (dacs%+2) mod 4, (dacs%+3) mod 4.

If the DAC number for the first channel is negative, output is to the system wave out device (usually a sound card). The maximum number of channels allowed depends on your system, but will usually be at least 2.

- sTime The start time in the current time view. This, together with eTime defines the data region to play.
- eTime The end time of the region in the current time view to output.
- rep% In the setup call, this sets the number of times to play the output. If omitted, the value 1 is used. To play for as many repeats as possible use the value 0. The total number of points to play is limited to 2 to the power 32 (about 4 billion).

In the status call, this returns the number of repeats completed.

- scale You can scale the replay rate in the range 0.25 to 4.0 (quarter speed to four times the channel rate). If you omit this argument, 1 is used, for output at the rate of the first channel.
- flag% Set this to 1 to position cursor 0 at the current replay position. The cursor is displayed if it is hidden. When the play finishes, the cursor is hidden if it was originally hidden.
- what Set this to -1 to stop output. Omit or set to 0 to report the position only. Both forms of the command return the position at the time of the call or -1 if there was no play in progress.
- Returns The setup call returns 0 if all was OK, else -1 if there is a problem with the number of channels or DAC channels, -2 for an internal setup problem, -3 if the output rate is too fast, -4 if the rate is too fast to replay in chunks, -5 if memory was exhausted, -6 for a problem with the sound card, -7 for a stupid points count. The status call returns the replay position, in seconds, or -1 if there is no replay in progress, or it has finished.

See also:DlgShow(), Interact(), PlayWaveAdd(), Toolbar(), Yield()

PlayWaveAdd()

This command adds a new area to the on-line play wave list. When you create a data file, Spike2 reserves 1401 memory and transfers any stored waveform to it. The area is set to play once and is not linked to any other area. The area replay speed factor is set to 1.0 and the wave is set to non-triggered. You must use this command before you use FileNew() to create the sampling window.

There are three command variants. The first adds a wave from the current time view or from a Spike2 data file, the second adds a wave from a data array, and the third reserves space without setting any data.

func PlayWaveAdd(key\$, lb\$, dac%, sT, eT, wch%{, mem% {,path\$}}); func PlayWaveAdd(key\$, lb\$, dac%, rate, data{%}[]{[]}); func PlayWaveAdd(key\$, lb\$, dac%, rate, size%);

- key\$ The first character of key\$ identifies this wave and triggers the wave playing in SampleKey(). It is an error to use Chr\$(0) or an empty string as a key.
- 1b\$ The label for the play wave control bar button that will play this wave, and record the character code as a keyboard marker. Labels can be up to 7 characters

long. If you include & as a character, it will not appear on the button, but the next character will be underlined and can be used as a keyboard short-cut.

- dac% Either a single DAC channel number or an array of DAC channel numbers. These are the outputs that will be used to play the data. The channel numbers must be in the range 0 to 3, and if more than one channel is specified, the channel numbers must be different.
- sT, eT The start and end times of the data in either the current time view, or in the file identified by the path\$ variable to be used as a source of output data.
- wch% Either a single channel, or an array of channels to use as a data source for playing. There must be one source channel for each output channel set by the dac% variable. The channels can be either waveform or WaveMark data. The sample rate is taken from the sampling rate of the first channel in the list. If subsequent channels in the list have different rates, data is interpolated.
- mem% If present, and non-zero, the data is converted to a memory image and becomes independent of the data file. Otherwise, Spike2 stores the file name and extracts the data as required for sampling.
- path\$ This optional argument sets the name of the file to extract data from. If absent, the current view (which must be a time view) is used as the data source. The file must exist and hold suitable data channels.
- rate When the data does not come from a file, this value sets the sample rate for each channel in Hz. Spike2 will get as close to this rate as it can.
- data This is either an integer or a real array. If there is more than one DAC channel to play, the array must have the same number of rows as there are DACs, that is if there are three DACs, the array must be equivalent to var data[n][3]; where n is the number of data points. If this is an integer array, the bottom 16 bits of each element is played through the DACs. If the data is in a real array, we assume that the full range of the DACs is ± 5 Volts and that the data is the required output value in Volts.
- size% This is the number of data points per channel to reserve for this area. The data values are not specified and you must use PlayWaveCopy() to transfer data to the 1401 for playing after sampling has started.

Returns The memory bytes in the 1401 used to hold this data, or a negative error code.

This command does not transfer the data to the 1401, that happens when the start sampling command is given. Spike2 always keeps a minimum memory area for data sampling, so there is a limit to size of the waveforms that can be copied to the 1401. This size limit is not known until sampling starts. There is also a limit on the size of a waveform that can be stored in the list, which is 32,000,000 bytes (was 2,000,000 bytes before version 5.09).

If you wish to link areas together or set the number of times the area is to be repeated or change the area speed factor, use one of the other PlayWave...() commands. The SampleClear() command removes all stored waveforms.

To start a waveform playing from the script you can use the SampleKey() function with the same key as set for the area or the output sequencer WAVEGO instruction.

```
See also:FileNew(), PlayWaveChans(), PlayWaveCycles(), PlayWaveDelete(),
        PlayWaveEnable(), PlayWaveInfo$(), PlayWaveLabel$(),
        PlayWaveLink$(), PlayWaveRate(), PlayWaveSpeed(),
        PlayWaveTrigger(), SampleClear()
```

PlayWaveChans()	This function lets you read back or set the DAC channels assigned to a particular play wave area. You cannot change the DAC assignments after sampling has started.						
	func P	layWaveChans(key\${, ch%[] {, set%}});					
	key\$	The first character of the string identifies the play wave area.					
	ch%	An optional integer array used to collect or set the DAC channel numbers. The array size must match the number of channels.					
	set%	If omitted or 0, the ch% array is filled in with the DAC channels used. If non-zero, the DAC channels are changed to the list defined by the ch% argument.					
	Returns	The command returns the number of DACs in the area or a negative error code.					
	See also	<pre>also:PlayWaveAdd(), PlayWaveInfo\$(), PlayWaveLabel\$(), PlayWaveLink\$(), PlayWaveRate()</pre>					
PlayWaveCopy()	This condone at	mmand is used to update a play wave data area in the 1401 memory. This can be any time that SampleStatus () returns 0 or 2, even while a wave is playing.					
	func P	layWaveCopy(key\$, data{%}[]{[]}{, offs%});					
	key\$	The first character of the string identifies the area to be updated.					
	data	This is either an integer or a real array. If there is more than one DAC channel to play, the array must have the same number of columns as there are DACs, that is if there are three DACs, the array must be equivalent to var data[n][3]; where n is the number of data points. If this is an integer array, the bottom 16 bits of each element is played through the DACs. If the data is in a real array, we assume that the full range of the DACs is ± 5 Volts and that the data is the required output value in Volts. The data in the array is copied to 1401 memory. It is an error for the array size to be larger than the memory area in the 1401.					
	offs%	The destination offset within the data area, in data points per channel. If the size of the data is such that the copy operation would extend beyond the end of the target area, the extra data is copied to the start of the area. The first offset is 0. Use PlayWaveStatus\$() to find the next offset to be written to the DACs.					
	Returns	The function returns 0 if all went well or if you call when there is no sampling or sampling is stopping. It returns a negative error code if there is a problem or you have requested sampling to start and it hasn't yet started.					
	See also:PlayWaveCycles(), PlayWaveLink\$(), PlayWaveSpeed(),						
		<pre>PlayWaveStatus\$(), PlayWaveStop(), SampleStatus()</pre>					
PlayWaveCycles()	This fur particula next pla	action gets and sets the number of times to play a waveform area associated with a ar key. If this is used on-line, it will also change the number of repeats for the y of the waveform.					
	func P	layWaveCycles(key\$, {new%});					

- key The first character of the string identifies the play wave area.
- new% If present, this sets the number of cycles to play. The value 0 sets a very large number of cycles.

Returns The number of cycles set at the time of the call.

PlayWaveDelete()	This function deletes one or more play wave areas from the sampling configuration. If you do this while sampling is in progress it will not change the waves loaded to the 1401.					
	<pre>func PlayWaveDelete({keys\$});</pre>					
	keys\$ If this is omitted, all play wave areas are deleted. If it is present, all areas with a key that is in this string are deleted. Case is significant for play areas; "abc" and "ABC" are not the same.					
	Returns The number of areas deleted or a negative error code.					
	See also:PlayWaveAdd(), SampleClear()					
PlayWayeEnable()	This function reports on the enabled state of a play wave area and optionally enables and					
T layWaveEllable()	disables it. Enabled areas are setup in the 1401 when sampling starts. Changes made with this function once sampling has started will not change the waves loaded to the 1401.					
	<pre>func PlayWaveEnable(keys\$, {set%});</pre>					
	key\$ The first character of the string identifies the play wave area.					
	set% If present and zero, the area is disabled. If non-zero, the area is enabled.					
	Returns The enabled state (1=enabled, 0=disabled) of the area at the time of the call, or a negative error code.					
	See also:PlayWaveAdd(), PlayWaveInfo\$(), PlayWaveStatus\$()					
PlayWaveInfo\$()	This function returns the list of keys associated with play wave areas, or gets the type of a particular area and the name of any associated data file.					
	<pre>func PlayWaveInfo\$({key\$, &size%, &type%});</pre>					
	key\$ If omitted, the function returns the list of keys associated with the play wave areas. If present, information about a specific area is returned in the remaining arguments and the function returns any file name associated with the area.					
	size% If present, this is returned as the number of data points per channel that will be played out for this channel.					
	type% This returns the type of the area as: 0= unused, 1 = data is taken from a data file, 2 = area has a memory image of data (and may also have an associated data file), 3 = area is reserved by the script but there is no associated data.					
	Returns If there is no key\$ value, then a string is returned holding one key character for each area. If there is a key, then the function returns the name of any data file associated with the area, in which case type% will be returned as 1 or 2.					
	See also:PlayWaveAdd(),PlayWaveChans(),PlayWaveCycles(), PlayWaveEnable(),PlayWaveLabel\$(),PlayWaveLink\$(), PlayWaveRate(),PlayWaveSpeed(),PlayWaveStatus\$()					
PlayWaveLabel\$()	This function returns and/or changes the label associated with each play area. The label can be up to 7 characters long and is used to label the buttons that appear on the Play waveform control bar. If you include & as a character, it does not appear in the label and the next character is underlined and can be used as a short cut to the button when the control bar is the current window.					
	<pre>func PlayWaveLabel\$(key\$, {new\$});</pre>					
	key\$ The first character of the string identifies the play wave area.					

new\$ If this is present, the label for the area is changed to the string. If the string is more than 7 characters long, only the first 7 are used.

Returns The label at the time of the function call.

See also:PlayWaveAdd()

PlayWaveLink\$(

You can link play wave areas together. Linked areas must have the same number of output channels and the same output channel list. The sample rate used is the sample rate for the area that is played first. You can change links during replay if SampleStatus() is 0 or 2. Both the area to link from and the area to link to must exist at the time of the call.

func PlayWaveLink\$(key\${, to\$});

- key\$ The first character of the string identifies the play wave area.
- to\$ If present, the first character of this string sets the area to link to. Use Chr\$(0) to cancel the link from the area set by key\$.
- Returns The key character of the area that was linked at the time of the call or an empty string if no area was linked or there was an error.

See also:PlayWaveAdd(), PlayWaveStatus\$(), SampleStatus()

PlayWaveRate()

This function gets or sets the base play rate for a play wave area. This is the standard play rate that can be changed by PlayWaveSpeed(). Changes to the rate made after sampling starts have no effect on the output; use PlayWaveSpeed() for on-line changes.

func PlayWaveRate(key\$, {new});

- key\$ The first character of the string identifies the play wave area.
- new If present, this is the new play rate for the area, in samples per second. You can set any value you like (not 0 or negative) and Spike2 will get as close as it can with the available hardware.

Returns The rate for the channel at the time of the function call.

See also:PlayWaveAdd(), PlayWaveLink\$(), PlayWaveSpeed(), SampleKey()

PlayWaveSpeed()

You can alter the sample rate for a play wave area by a factor of 0.25 to 4.0 with this command. Spike2 may not be able to play at the rate you request; it will set the closest rate it can. The H_Z argument returns the achieved rate. On-line changes are allowed.

func PlayWaveSpeed(key\${, new{, wait%{, Hz}}});

- key\$ The first character of the string identifies the play wave area. If this area is playing, or an area that this area links to, the rate will change during playing.
- new If present, this is the new speed factor for the area, in the range 0.25 to 4.0. Spike2 gets as close to this speed factor as it can with the available hardware.
- wait% If present and non-zero, any on-line speed change is postponed until the end of the current cycle and will happen within a few milliseconds of the cycle end.
- Hz If present and a sampling document is open, it returns the real replay rate in Hz.
- Returns The speed factor for the area at the time of the function call or 0 if there is no area defined by the key.

See also:PlayWaveAdd(), PlayWaveLink\$(), PlayWaveRate(), SampleKey()

PlayWaveStatus\$()	This function returns information about waveform output during sampling.					
	<pre>func PlayWaveStatus\$({&pos%{, &cyc%}});</pre>					
	If present, this returns the next position, in terms of the number of points per channel in the area, to write to the DACs. The first position is 0. The DACs hold one data point ready for output on the next clock tick, so pos% is the next but one index to output. If you play an area with a triggered start but do not trigger it, pos% is 1, not 0. This is because the DACs are already holding index 0, ready to output when the trigger arrives.					
	cyc% If present, this integer variable is returned holding the number of cycles left to play (including the current cycle).					
	Returns The key code of the area that is playing or waiting for a trigger, or an empty string if no area is playing or sampling is not active.					
	See also:PlayWaveAdd(), PlayWaveCycles(), PlayWaveLink\$(), PlayWaveRate(), PlayWaveSpeed(), PlayWaveStop(), SampleKey()					
PlayWaveStop()	when the current cycle finishes.					
	<pre>func PlayWaveStop(cEnd%);</pre>					
	cEnd% If present and non-zero, the current cycle for the playing area will be the last cycle for that area, otherwise output will stop immediately.					
	Returns 1 for OK, 0 if not playing or a negative error code.					
	See also:PlayWaveAdd(), PlayWaveCycles(), PlayWaveStatus\$(), SampleKey()					
PlayWaveTrigger()	This function reports and optionally changes the trigger state of a play wave area. If an area is triggered, a play request prepares the area but output does not start until a trigger signal is received by the 1401. This trigger is the E3 front panel input for the 1401 <i>plus</i> and is the Trigger input for the micro1401 and Power1401 unless it is routed to the rear panel by the Edit Preferences menu.					
	<pre>func PlayWaveTrigger(keys\$,{set%});</pre>					
	key\$ The first character of the string identifies the play wave area.					
	set% If present this sets the triggered state. 0 = not triggered and non-zero = triggered.					
	Returns The trigger state (1=triggered, 0=not triggered) of the area matching the first character in keys\$ at the time of the call, or a negative error code.					
	See also:PlayWaveAdd(), PlayWaveEnable(), PlayWaveStatus\$()					
	This function misses to the neuron of If the coloradation underflower the neurltic O					
Pow()	Function raises x to the power of y. If the calculation under nows, the result is 0. Func Pow(x x[] x[]], y):					
	Δ real number or a real array to be raised to the power of π					
	The exponent If y is negative must be integral					
	y The exponent. If x is negative, y must be integral. Returns If x is an array it returns 0 or a negative error code. If y is a number, it returns y					
	to the power of y unless an error is detected, when the script halts.					
	See also: Abs (), ATan (), Cos (), Exp (), Frac (), Ln (), Log (), Max (), Min (), Rand (),					

Round(),Sin(),Sqrt(),Tan(),Trunc()

Print()

This command prints to the current view at the text caret. If the first argument is a string (not an array), it is used as format information for the remaining arguments. If the first argument is an array or not a string or if there are more arguments than format specifiers, Spike2 prints the arguments without a format specifier in a standard format and adds a new line character at the end. If you provide a format string and you require a new line character at the output, include n at the end of the format string.

Func Print(form\$|arg0 {,arg1 {,arg2...}});

- form\$ A string that specifies how to treat the following arguments. The string contains two types of characters: ordinary text that is copied to the output unchanged and format specifiers that convert the following arguments to text. Format specifiers start with % and end with one of the letters d, x, c, s, f, e or g in upper or lower case. For a % in the output, use %% in the format string.
- arg1, 2 The arguments used to replace %c, %d, %e, %f, %g, %s and %x type formats.

Returns 0 or a negative error code. Fields that cannot be printed are filled with asterisks.

Format specifiers The full format specifier is: %{flags}{width}{.precision}format

- flags The flags are optional and can be placed in any order. They are single characters that modify the format specification as follows:
 - Specifies that the converted argument is left justified in the output field.
 - + Valid for numbers, and specifies that positive numbers have a + sign.
 - *space* If the first character of a field is not a sign, a space is added.
 - 0 For numbers, causes the output to be padded on the left to the field width with 0.
 - # For x format, 0x is prefixed to non-zero arguments. For e, f and g formats, the output always has a decimal point. For g formats, trailing zeros are not removed.
- width If this is omitted, the output field will be as wide as is required to express the argument. If present, it is a number that sets the minimum output field width. If the output is narrower than this, the field is padded on the left (on the right if the flag was used) to this width with spaces (zeros if the 0 flag was used). The maximum width for numbers is 100.
- precision This number sets the maximum number of characters to be printed for a string, the number of digits after the decimal point for e and f formats, the number of significant figures for g format and the minimum number of digits for d format (leading zeros are added if required). It is ignored for c format. There is no limit to the size of a string. Numeric fields have a maximum precision value of 100.
 - format The format character determines how the argument is converted into text. Both upper and lower cased version of the format character can be given. If the formatting contains alphabetic characters (for example the e in an exponent, or hexadecimal digits a-f), if the formatting character is given in upper case the output becomes upper case too (e+23 and 0x23ab become E+23 and 0x23AB). The formats are:
 - c The argument is printed as a single character. If the argument is a numeric type, it is converted to an integer, then the low byte of the integer (this is equivalent to integer mod 256) is converted to the equivalent ASCII character. You can use this to insert control codes into the output. If the argument is a string, the first character of the string is output. The following example prints two tab characters, the first using the standard tab escape, the second with the ASCII code for tab (8): Print("\t%c", 8);
 - d The argument must be a numeric type and is printed as a decimal integer with no decimal point. If a string is passed as an argument the field is filled with asterisks. The following prints "23,0002":

Print("%4d,%.4d", 23, 2.3);

The argument must be a numeric type; otherwise the field is filled with asterisks. The argument is printed as $\{-\}$ m.dddddde \pm xx $\{x\}$ where the number of d's is set by the precision (which defaults to 6). A precision of 0 suppresses the decimal point unless the # flag is used. The exponent has at least 2 digits (in some implementations of Spike2 there may always be 3 digits, others use 2 digits unless 3 are required). The following prints "2.300000e+01,2.3E+00": Print("%4e,%.1E", 23, 2.3); The argument must be a numeric type; otherwise the field is filled with asterisks. The f argument is printed as {-}mmm.ddd with the number of d's set by the precision (which defaults to 6) and the number of m's set by the size of the number. A precision of 0 suppresses the decimal point unless the # flag is used. The following prints "+23.000000,0002.3": Print("%+f,%06.1f", 23, 2.3); The argument must be a numeric type; otherwise the field is filled with asterisks. α This uses e format if the exponent is less than -4 or greater than or equal to the precision, otherwise f format is used. Trailing zeros and a trailing decimal point are

> not printed unless the # flag is used. The following prints "2.3e-06,2.300000": Print ("%g, %#g", 0.0000023, 2.3);

- s The argument must be a string; otherwise the field is filled with asterisks.
- x The argument must be a numeric type and is printed as a hexadecimal integer with no leading 0x unless the # flag is used. The following prints "1f, 0x001F": Print("%x, %#.4X", 31, 31);

Arrays in the argument list The d, e, f, g, s and x formats support arrays. One dimensional arrays have elements separated by commas; two dimensional arrays use commas for columns and new lines for rows. If there is a format string, the matching format specifier is applied to all elements.

See also:Message(), ToolbarText(), Print\$(), PrintLog()

Print\$()

This command prints formatted output into a string. The syntax is identical to the Print() command, but the function returns the generated output as a string.

Func Print\$(form\$|arg0 {,arg1 {,arg2...}});

- form\$ An optional string with formatting information. See Print() for a description.
- arg1, 2 The data to format into a string.
- Returns It returns the string that is the result of the formatting operation. Fields that cannot be printed are filled with asterisks.

See also:Asc(), Chr\$(), DelStr\$(), LCase\$(), Left\$(), Len(), Mid\$(), Print(), PrintLog(), Right\$(), Str\$(), UCase\$(), Val()

PrintLog()

This command prints to the log window. The syntax is identical to the Print() command, except that the output always goes to the end of the log window.

Func PrintLog(form\$|arg0 {,arg1 {,arg2...}});

form\$ An optional string with formatting information. See Print() for a description. arg1,2 The data to print.

Returns 0 or a negative error code. Fields that cannot be printed are filled with asterisks. See also:Print(), Print\$(), Message()

Process()

Processes the current view. For result and XY views, the time view to process data from must be open. The times for start and end of processing are times in the time view. Use View(-1). Cursor(1) to refer to time view times from result and XY views.

- The time to start processing from, in seconds. Negative times are treated as zero. Times greater than Maxtime() cause no processing. In triggered modes with no trigger channel, this sets the trigger time and eTime is ignored. In Phase histograms with no Cycle channel, this sets the start of a cycle.
- eTime The end time for processing. In Phase histogram mode with no Cycle channel, this sets the end time of a single cycle.
- clear% If present, and non-zero, the result view bins are cleared before the results of the analysis are added to the result view and Sweeps () result is reset.
- opt% If present and non-zero, the result view is optimised after processing the data.
- dest% Used when processing to a channel from MeasureToChan(). If not used it should be set to 0. This identifies the destination channel for processing. If omitted or 0, all suitable channels are processed. It is not an error if this doesn't match a destination channel, but no processing happens.
- gate% If present, this is the channel number of an event or marker-based channel in the associated time view to use as a gate. Both len and pre must be present.
- len The time to process for each gate event in the time sTime to eTime.
- pre How far before each gate event to start processing.
- mCd% If present and gate% is a marker or derived type, it holds the marker code to use as the gate. Set this to -1 or omit it to use the current channel marker filter.
- Returns This returns the number of items processed. This is the number of intervals considered for an INTH (even if they didn't fall in the histogram), the number of sweeps for sweep-based analysis, the number of data blocks for SetPower(). In the case of an error, a negative error code is returned.

See also:MeasureToChan(), MeasureToXY(), ProcessAll(), SetAverage(), SetEvtCrl(), SetINTH(), SetPhase(), SetPSTH(), SetPower(), SetResult(), SetWaveCrl(), Sweeps()

ProcessAll()

This function runs all processes that use the current time view as a data source. It is as if a Process() command were used for each target view using the current processing settings. These are the settings you would see if you opened the Process Settings dialog for the target. If you have already used a Process() command on a target, this sets the current processing settings.

Func ProcessAll(sTime, eTime);

- The time to start processing from, in seconds. Negative times are treated as zero. Times greater than Maxtime() cause no processing. In triggered modes with no trigger channel, this sets the trigger time and eTime is ignored. In Phase histograms with no Cycle channel, this sets the start of a cycle.
- eTime The end time for processing. In Phase histogram mode with no Cycle channel, this sets the end time of a single cycle.

Returns Zero if no errors or a negative error code.

```
See also:MeasureToChan(), Process(), MeasureToXY(), SetAverage(),
        SetEvtCrl(), SetINTH(), SetPhase(), SetPSTH(), SetPower(),
        SetResult(), SetWaveCrl(), Sweeps()
```

ProcessAuto()	This is equivalent to the Process Dialog for a New file in Automatic mode. The current view must be the view where the results appear. The processing parameters set by this command are used when this view is given a chance to update.						
	Func P	<pre>rocessAuto(delay, mode%{, opt%{, last{, leeway{, dest%}}});</pre>					
	delay	The minimum time between updates, in seconds.					
	mode%	mode% 0=accumulate all data. 1=clear the result, process the most recent last se					
	opt%	If present and non-zero, the result view is optimised after each process.					
	last	The length of time to process in mode 1, in seconds; ignored in mode 0.					
	leeway	How close cursor $0\ \text{can}\ \text{be}\ \text{to}\ \text{the}\ \text{file}\ \text{end}\ \text{for}\ XY\ \text{views}\ \text{and}\ \text{MeasureToChan}\ ()$.					
	dest%	The channel for MeasureToChan(), ignored for other process types. If omitted or 0, all suitable channels are processed. It is not an error if this doesn't match a destination channel, but no processing happens.					
	Returns	0 or a negative error code.					
	See also	:MeasureToChan(),MeasureToXY(),ProcessTriggered(),Process()					
ProcessTriggered()	This is e view me comman	equivalent to the Process Dialog used with a New file in Gated mode. The current ust be the view where the results appear. The processing parameters set by this and are used when this view is given a chance to update.					
	<pre>Func ProcessTriggered(len,pre,gate%{,clr%{,opt%,{mCd%{,dest%}}});</pre>						
	len	The length of data to process around each gate event, in seconds.					
	pre	The pre-event time, in seconds.					
	gate%	The channel holding gate events or markers.					
	clr%	If present, and non-zero, the result view bins are cleared before the results of the analysis are added to the result view and $Sweeps()$ result is reset.					
	opt%	If present and non-zero, the result view is optimised after processing the data.					
	mCd%	If present and gate% is a marker or derived type, it holds the marker code to use as the gate. Set this to -1 or omit it to use the current channel marker filter.					
	dest%	The channel for MeasureToChan(), ignored for other process types. If omitted or 0, all suitable channels are processed. It is not an error if this doesn't match a destination channel, but no processing happens.					
	Returns 0 or a negative error code.						
	See also	:MeasureToChan(),MeasureToXY(),ProcessTriggered(),Process()					

Profile()

Spike2 saves information in the HKEY_CURRENT_USER\Software\CED\Spike2 section of the system registry. The registry is a tree of keys with values attached to each key. If you think of the registry as a filing system, the keys are folders and the values are files. Keys and values are identified by case-insensitive text strings. The Profile() command manipulates the keys and values within the Spike2 section of the registry.

You can also view and edit the registry with the regedt32 program, which is part of your system. Select Run from the start menu and type regedt32 then click OK. Please read the regedt32 help information before using this program. It is a very powerful tool; careless use can severely damage your system.

Do not write vast quantities of data into the registry; it is a system resource and should be treated with respect. If you must save a lot of data, write it to a text or binary file and save the file name in the registry. If you think that you may have messed up the Spike2 section

	of the registry, use regedt32 to locate the Spike2 section and delete it. The next you run Spike2 the section will be restored; you will lose any preferences you had se					
	<pre>Proc Profile(key\${, name\${, val%{, &read%}}}); Proc Profile(key\${, name\${, val\${, &read\$}}});</pre>					
	key\$ This sets the Spik "My bit may not marks; a key that	This sets the key to use in the Spike2 registry section. Set key\$ empty to the Spike2 key. You can use nested keys separated by a backslash, for exam "My bit\\stuff" to use the key stuff inside the key My bit. The key na may not start with a backslash. Remember to use two backslashes inside qu marks; a single backslash is an escape character. It is never an error to refer key that does not exist; the system creates missing keys for you.				
	name\$ This strin name, th	ig identifies the data in the key to read or write. If you set an empty s refers to the (default) data item for the key set by key .				
	val This is e registry.	ither a string or an integer. If read is omitted, val is written to the if read is present, val is returned if the registry item does not exist.				
	read If presen value hel	t, it must have the same type as val. This is a variable that is set to the d in the registry. If name\$ is not in the registry, read is set to val.				
	 Profile() can be used with 1 to 4 arguments. It has a diffe The key identified by key\$ is deleted. All sub-keys ar key and sub-keys are also deleted. Nothing is done if key The value identified by name\$ in the key key\$ is deleted The value identified by name\$ in the key key\$ is set to The value identified by name\$ in the key key\$ is return 					
	The following scr var path\$, cor Profile("My d. Profile("My d. Profile("My d. Profile("My d.	<pre>pt example collects values at the start, then saves them at the end: .nt%; .ta", "path", "c:\\work", path\$); 'get initial path .ta", "count", 0, count%); 'and initial count</pre>				
Registry use by Spike2	HKEY_CURRENT_	JSER\Software\CED\Spike2 holds the following keys:				
BarList	This key holds the list of scripts to load into the script bar when Spike starts.					
Editor	This key holds the editor settings for scripts, output sequences and general text editing					
PageSetup	This key holds the margins in units of 0.01 mm for printing data views, and the margin in mm and header and footer text for text-based views.					
Preferences	The values in this Edit menu Prefer immediately. The	key are mainly set by the Edit menu preferences. If you change any ences value in this key, Spike2 will use the changed information values are all integers except the file path, which is a string:				
	E3 trigger at rear Enhanced Metafil Enter debug on er Event ports at rear Fast online scroll Line thickness coo	 0=Sample and PlayWave trigger on front panel, 1=on rear 0=Windows metafile, 1=enhanced metafile for clipboard. 0=Do not enter debug, 1=enter debug 0=Events 0&1 on front panel, 1=on rear panel. 0=Normal, 1=faster (inaccurate) scrolling algorithm Bits 0-3 = Axis code, bits 4-7 = Data code. The codes 0-15 map onto the 16 values in the drop down list. Bit 7=1 to use lines not rectangles to draw axes 				
	Low channels at t Metafile NoComp Metafile Scale New file path No Save Prompt	 0=Standard display shows low channel at bottom, 1=at top. 0=Compress when multiple points per x pixel, 1=no compress. 0-11 selects from the list of allowed scale factors. New data file directory or blank for current folder. 0=Prompt to save derived views, 1=no prompt. 				

Save modified scripts Ten volt 1401 Use colour	0=Do not save, 1=save 0=5 Volt 1401, 1=10 Volt 1401, 2=same as last 1401. 0=Use Black and White, 1=Use Colour (from the View menu).					
The keys with names starting "Bars-" are used by system code to restore dockabl toolbars. You can delete them all safely; any other change is likely to crash Spike2.						
This key holds the list of recently used files that appear at the bottom of the file menu.						
This key holds the information to recover data from interrupted sampling sessions.						
This is where the evaluate bar saves the last few evaluated lines. The colour palette is also saved here. Set "No MThread Sampling" to 1 to force single threaded sampling.						
The Tip of the Day dialog uses this key to remember the last tip position.						
Spike2 uses this key to detect when a new version of the program is run for the first time.						
In Windows NT derive menu About Spike2 d dialog to read more abo Minimum working set Maximum working set	d systems, this key holds the desired working set sizes. The Help ialog displays the current working set sizes. Click Help in this out using these registry values. Minimum size in kB (units of 1024 bytes), default is 800 Maximum size in kB, default is 4000 (4 MB)					
	Save modified scripts Ten volt 1401 Use colour The keys with names toolbars. You can delet This key holds the list of This key holds the infor This is where the eval also saved here. Set "No The <i>Tip of the Day</i> dial Spike2 uses this key to In Windows NT derive menu About Spike2 d dialog to read more abo Minimum working set Maximum working set					

See also:ViewUseColour()

ProgKill()

This function terminates a program started using ProgRun(). This can be dangerous, as it will terminate a program without giving it the opportunity to save data.

Func ProgKill(pHdl%);

pHdl% A program handle returned by ProgRun ().

Returns Zero or a negative error code.

See also:ProgRun(), ProgStatus()

ProgRun()

This function runs a program using command line arguments as if from a command prompt. Use ProgRun() to test the program status, ProgKill() to terminate it. The program inherits its environment variables from Spike2, see System\$() for details.

Func ProgRun(cmd\$ {,code% {,xLow, yLow, xHigh, yHigh}});

- cmd\$ The command string as typed at a command prompt. To run shell command x
 use "cmd /c x" in Windows NT and XP, "command /c x" in Windows 9x.
 To run another copy of Spike2, use "sonview.exe /M {filename}".
- code% If present, this sets the initial application window state: 0=Hidden, 1=Normal, 2=Iconised, 3=maximised. Some programs set their own window state so this may not work. The next 4 arguments set the Normal window position:
- xLow Position of the left window edge as a percentage of the screen width.
- yLow Position of the top window edge as a percentage of the screen height.
- xHigh The right hand edge as a percentage of the screen width.
- yHigh The bottom edge position as a percentage of the screen height.
- Returns A program handle or a negative error code. ProgStatus() releases resources associated with the handle when it detects that the program has terminated.

See also:FileCopy(), FileDelete(), ProgKill(), ProgStatus(), System\$()

ProgStatus()	This fun resources	ction tests if a program started with ProgRun() is still running. If it is not, s associated with the program handle are released.					
	Func ProgStatus(pHdl%);						
	pHdl%	The program handle returned by ProgRun().					
	Returns	1=program is running, 0=terminated, resources released, handle now invalid. A negative error code (-1525) means that the handle is invalid.					
	See also:	ProgKill(), ProgRun()					
Query()	This fun	ction is used to ask the user a Yes/No question. It opens a window with a					
Query()	message	and two buttons. The window is removed when a button is pressed.					
	Func Qu	<pre>nery(text\$, {,Yes\$ {,No\$}});</pre>					
	text\$	This string forms the text in the window. If the string includes a vertical bar, the text before the vertical bar is used as the window title.					
	Yes\$ This sets the text for the first button. If this argument is omitted, "Yes"						
	No $\$$ This sets the text for the second button. If this is omitted, "No" is used.						
	Returns	Returns 1 if the user selects Yes or presses Enter, 0 if the user selects the No button.					
	See also:	<pre>Print(),Input(),Message(),DlgCreate()</pre>					
Rand()	This returned	arns pseudo-random numbers with a uniform density in a set range. The values are $\mathbf{R} * \text{scl} + \text{off}$ where \mathbf{R} is in the range 0 up to but not including 1. Spike?					
	initialises a repeata	s the generator with a random seed based on the time. You must set the seed for ble sequence. The sequence is independent of RandExp() and RandNorm().					
	Func Ra Func Ra Func Ra	<pre>and(seed); and({scl, off}); and(arr[]{[]}{, scl{, off}});</pre>					
	seed	If present, this is a seed for the generator in the range 0 to 1. If seed is outside this range, the fractional part of the number is used. If you use 0.0 as the seed, the generator is initialised with a seed based on the system time.					
	arr	This 1 or 2 dimensional real or integer array is filled with random numbers. If an integer array is used, the random number is truncated to an integer.					
	scl	This scales the random number. If omitted, it has the value 1.					
	off	This offsets the random number. If omitted it has the value 0.					
	Returns	If the first argument is not an array, the return value is a random number in the range off up to off+scl. If a seed is given, a random number is still returned. If the first argument is an array, the return value is 0.0.					

See also:RandExp(), RandNorm()

RandExp()

This function returns pseudo-random numbers with an exponential density, suitable for generating Poisson statistics. The values returned are $\mathbf{R} * \text{mean} + \text{off}$ where \mathbf{R} is a random number with the density function $p(x) = e^{x}$. When you start Spike2, the generator is initialised with a random seed based on the time. For repeatable sequences, you must set a seed. The sequence is independent of Rand() and RandNorm().

Func RandExp(seed);
Func RandExp({mean, off});
Func RandExp(arr[]{[]}{, mean{, off}});

seed	If present, this is a seed for the generator in the range 0 to 1. If seed is outside
	this range, the fractional part of the number is used. If you use 0.0 as the seed,
	the generator is initialised with a seed based on the time.

- arr This 1 or 2 dimensional real or integer array is filled with random numbers. If an integer array is used, the random number is truncated to an integer.
- mean This scales the random number. If omitted, it has the value 1.
- off This offsets the random number. If omitted it has the value 0.
- Returns If the first argument is not an array, the return value is a random number. If a seed is given, a random number is still returned. If the first argument is an array, the return value is 0.0.

The following example fills an array with event times with a mean interval t:

<pre>RandExp(arr[], t);</pre>	'Fill	arr	with	ever	nt	intervals
<pre>ArrIntgl(arr[]);</pre>	'conv	ert	interv	vals	to	times

See also:Rand(), RandNorm()

RandNorm(

This function returns pseudo-random numbers with a normal density. The values returned are $\mathbf{R} * \text{scl} + \text{off}$ where \mathbf{R} is a random number with a normal probability density function $p(x) = (2\pi)^{-\frac{1}{2}}e^{x^2/2}$; this has a mean of 0 and a variance of 1. When you start Spike2, the generator is initialised with a random seed based on the time. For a repeatable sequence, you must set a seed. The sequence is independent of Rand() and RandExp().

```
Func RandNorm(seed);
Func RandNorm({scl, off});
Func RandNorm(arr[]{[]}{, scl{, off}});
```

- seed If present, this is a seed for the generator in the range 0 to 1. If seed is outside this range, the fractional part of the number is used. If you use 0.0 as the seed, the generator is initialised with a seed based on the time.
- arr This 1 or 2 dimensional real or integer array is filled with random numbers. If an integer array is used, the random number is truncated to an integer.
- scl This scales the random number. If omitted, it has the value 1.
- off This offsets the random number. If omitted it has the value 0.
- Returns If the first argument is not an array, the return value is a random number. If a seed is given, a random number is still returned. If the first argument is an array, the return value is 0.0.

See also:Rand(), RandExp()

RasterAux()

This is now deprecated; use RasterSort() or RasterSymbol() instead. It is here for version 4 compatibility, but will be removed in a future release. This function returns and optionally sets some of the raster auxiliary values for a channel in the current result view.

Values 0 and 1 are used to sort the sweeps and can be selected in addition to time order by the DrawMode() command. Unset values read back as zero.

Values 2 to 5 are the times, in seconds relative to the start of the file, of four markers (circle, cross, square and triangle) to display in the sweep. The markers are drawn in the colours set for WaveMark codes 1 to 4. Markers appear if their position relative to the sweep trigger lies within the sweep. Unset values read back as a negative time.

Func RasterAux(chan%, sweep%, num% {,new});

chan% The channel in the result view. The channel must have raster data enabled.
sweep% The sweep number in the range 1 to the number of sweeps in the channel.

num% The auxiliary value to return and optionally change in the range 0 to 5.

new If present, this changes the auxiliary value for the sweep.

Returns The auxiliary value at the time of the call.

Replace RasterAux(c%,s%,n%{,new}) with RasterSort(c%,s%,n%{,new}) if n% is 0 or 1, and with RasterSymbol(c%,s%,n%-2{,new}) if n% is 2 to 5.

See also:DrawMode(), RasterGet(), RasterSet(), RasterSort(), RasterSymbol(), SetEvtCrl(), SetPSTH()

RasterGet()

This returns a sweep of raster data for a result view channel for which raster data has been enabled. Using this when the current view is not a result view causes a fatal error.

Func RasterGet(chan%{,sweep%{,&sT{,data{ $\%}[]$ {,&tick{,&eT}}}});

- chan% The channel number in the result view. If this is the only argument, the return value is the number of sweeps in the channel.
- sweep% The sweep number in the result view in the range 1 to the number of sweeps. If this argument is present, the return value is the number of times in the sweep.
- This is returned holding the time in seconds of the start of the sweep in the original data. This is not the sweep trigger time (the time that corresponds to the x axis 0 in the result view). The trigger time is sT-BinToX(0).
- data This optional array is filled with sweep event times. A real array returns times in seconds; an integer array returns times in units of the microseconds per time used when the view was created. The number of times copied into the array is the lesser of the number of times in the sweep and the length of the array.
- tick This optional real value is the number of seconds per time units used when times are returned as an integer array. That is, if you multiply each integer time returned in data%[], you would get the time in seconds. If you created the result view using SetResult() it is the tick value you passed in, otherwise it is the BinSize() value from the associated time view.
- This optional real value is returned holding the end time of the sweep, in seconds. This is usually only of interest for data created with SetPhase().

Returns The count of sweeps in the channel, times in the sweep or a negative error code.

See also:RasterSet(), RasterSort(), SetEvtCrl(), SetPhase(), SetPSTH()

RasterSet()

This function sets the raster data for a sweep for a channel of a result view. You can replace the data for an existing sweep, or add a new sweep. It is a fatal script error to call this function when the current view is not a result view.

Func RasterSet(chan%, sweep%, sT {,data[]|data%[] {,eT}});

- chan% The channel number in the result view.
- sweep% This is either in the range 1 to Sweeps () to replace the existing raster data for a sweep, or it can be 0 to add a new sweep. If you add a new sweep, the new sweep is added to all channels. The remaining channels have a sweep added that has the same start time, and no data. If your result view has multiple channels, only use a sweep% value of 0 for one of the channels.
- This sets the time of the start of the sweep, in seconds. This is not the trigger time of the sweep (the time that corresponds to the x axis 0 in the result view). The trigger time is sT-BinToX(0).

- data This is a real or an integer array holding the events times for the sweep. The size of the array sets the number of items. If this is a real array, the times are in seconds. If this is an integer array, the times are in the underlying tick units (see RasterGet() and SetResult() for details). Times outside the time range start to start+MaxTime()*BinSize() are ignored.
- eT This optional value is the end time of the sweep. If you supply this, all times in the data array are mapped into the time period sT to eT, regardless of the x axis scaling set for the result view. This argument is usually only used when working with or emulating a Phase histogram.

Returns The sweep number that received the data or a negative error code.

See also:RasterGet(), RasterSort(), SetEvtCrl(), SetPhase(), SetPSTH()

RasterSort()

Each raster sweep has 4 values that can be used to sort displayed rasters. The DrawMode() command selects between time order and one of these values. This function returns and optionally sets these values for the current result view.

Func RasterSort(chan%, sweep%, num% {,new});

- chan% The channel in the result view. The channel must have raster data enabled.
- sweep% The sweep number in the range 1 to the number of sweeps in the channel.

num% The sort value to return and optionally change in the range 1 to 4.

new If present, this changes sort value num% for the sweep.

Returns The sort value at the time of the call. Unset values read back as zero.

See also:DrawMode(), RasterGet(), RasterSet(), RasterSymbol(), SetEvtCrl(), SetPSTH()

RasterSymbol()

Each raster sweep has 8 values that can be displayed as symbols. This function returns and optionally sets these values for the current result view. The values are times, in seconds relative to the start of the file. The 8 markers are: circle, cross, square, up triangle, plus, diamond, down triangle, filled square. The markers are drawn in the colours set for WaveMark codes 1 to 8. Markers appear if their position relative to the sweep trigger lies within the sweep. Unset values read back as a negative time.

Func RasterSymbol(chan%, sweep%, num% {,new});

- chan% The channel in the result view. The channel must have raster data enabled.
- sweep% The sweep number in the range 1 to the number of sweeps in the channel.

num% The symbol number in the range 1 to 8 to get or set.

- new If present, this sets the symbol time for the sweep. Set -1 to cancel the symbol.
- Returns The symbol time at the time of the call. Unset symbols have a negative time.

See also:DrawMode(), RasterGet(), RasterSet(), SetEvtCrl(), SetPSTH()

Read()

This function reads the next line from the current text view or external text file and converts the text into variables. The read starts at the beginning of the line containing the text cursor. The text cursor moves to the start of the next line after the read.

Func Read({&var1 {, &var2 {,&var3 ...}}});

- varn Arguments must be variables. They can be any type. One dimensional arrays are allowed. The variable type determines how to convert the string data. In a successful call, each variable matches a field in the string, and the value of the variable changes to the value found in the field. A call to Read() with no arguments skips a line.
- Returns The function returns the number of fields in the text string that were successfully extracted and returned in variables, or a negative error code. Attempts to read past the end of the file produce the end of file error code.

It is not an error to run out of data before all the variables have been updated. If this is a possibility you must check that the number of items returned matches the number you expected. An array of length n is treated as n individual values.

The source string is expected to hold data values as real numbers, integer numbers (decimal or hexadecimal introduced by 0x) and strings. Strings can be delimited by quote marks, for example "This is a string", or they can be just text. If a string is not delimited, it is deemed to run to the end of the source string, so no other items can follow it. You can use ReadSetup() to change the characters that delimit a string and also to define hard separator characters within non-delimited strings. String delimiters are not returned as part of the string.

Normally, the fields in the source string are separated by white space (tabs and spaces) and commas. Space characters are "soft" separators. You can have any number of spaces between fields. Tabs and commas are treated as "hard" separators. Two consecutive hard separators (with or without intervening soft separators), imply a blank field. You can use ReadSetup() to redefine the soft and hard separators. When reading a field, the following rules are followed:

- 1. Soft separator (space) characters are skipped over
- 2. If the field is a string and the next character is a delimiter, it is skipped.
- 3. Characters that are legal for the destination variable are extracted until a non-legal character or a separator or a required string delimiter or end of data is found. The characters read are converted into the variable type. If an error occurs in the translation, the function returns the error. Blank fields assigned to numbers are treated as 0. Blank fields assigned to strings produce empty strings.
- 4. Characters are skipped until a separator character is found or end of data. If a soft separator is found, it and any further soft separators are skipped. If the next character is a hard separator it is also skipped.
- 5. If there are no more variables or no more data, the process stops, else back to step 1.
- **Example** The following example shows a source line, followed by a Read() function, then the assignment statements that would be equivalent to the Read():

"This is text" , 2 3 4,, 4.56 Text too 3 4 5 The source line n := Read(fred\$, jim[1:2], sam, dick%, tom%, sally\$, a, b, c); n := 7; This is equivalent to the result fred\$:= "This is text"; jim[1] := 2; jim[2] := 3; sam := 4; dick% := 0; tom% := 4; sally\$:= "Text too 3 4 5" a, b and c are not changed

See also:EditCopy(), FileOpen(), ReadSetup(), ReadStr(), Selection\$()

ReadSetup()	This sets the separators and delimiters used by Read() and ReadStr() to convert text into numbers and strings. You can also set string delimiters and set a string separator.					
	Proc R	eadSetup({hard\${, soft\${, sDel\${, eDel\${, sSep\$}}}});				
	hard\$	The characters to use as hard separators between all fields. If this is omitted or the string is empty, the standard hard separators of comma and tab are used.				
	soft\$	The characters to use as soft separators. If this is omitted, the space character is set as a soft separator. If $soft$ is empty, no soft separators are used.				
	sDel\$	The characters that delimit the start of a string. If omitted, a double quote is used. If empty, no delimiter is set. Delimiters are not returned in the string.				
	eDel\$	If empty, no delimit is set. If sDel\$ and eDel\$ are the same length, only the end delimiter character that matches the start delimiter position is used. For example, to delimit strings with <text> or 'text' set sDel\$ to "<'" and eDel\$ to ">'". You can repeat a character to force different lengths.</text>				
	sSep\$	The list of hard separator characters for strings that have no start delimiter. For example, setting " " lets you read one two three into three separate strings.				
	See also	<pre>:Read(),ReadStr(),Val()</pre>				
	TT1 · C					
ReadStr()	This fun	ction extracts data fields from a string and converts them into variables.				
	FUNC R	<pre>unc ReadStr(text\$, &var1 {, &var2 {, &var3}});</pre>				
	text\$	The string used as a source of data.				
	be one dimensional arrays. The type of each variable determines function tries to extract data from the string. See Read () for details.					
	Returns	Returns The function returns the number of fields in the text string that were successful extracted and returned in variables, or a negative error code.				
	It is not an error to run out of data before all the variables have been updated. If this is a possibility you must check the returned value. If an array is passed in, it is treated as though it was the number of individual values held in the array.					
	See also	See also:Read(),ReadSetup(),Val()				
ReRun()	This fur menu R	action controls the rerun of the current time view and is equivalent to the View eRun command. You cannot use this on a file that is being sampled.				
	Func R	eRun({run%{, sTime{, eTime{, scale}}});				
	run%	If >0, the current time view starts to rerun unless it is already rerunning. If zero, rerun is cancelled for the time view. Omit to leave the state unchanged. Use -1 to return the current sTime, -2 to return eTime and -3 to return scale. The sTime, eTime and scale arguments are used when run% > 0.				
	sTime	Sets the start time for the rerun. If omitted, 0 is used.				
	eTime	Sets the end time of the rerun. If omitted $MaxTime()$ is used.				
	scale	Sets the time scale for the rerun. If omitted, 1.0 is used. A value of 2 reruns twice as fast. Values from 0.01 to 100.0 are allowed.				
	Returns	If run% is positive or omitted the command returns the state at the time of the call: 0=not rerunning, 1=rerunning, -1= rerunning is not allowed. Negative values of run% return the current rerun settings.				

Right\$()	<pre>This function returns the rightmost n characters of a string. Func Right\$(text\$, n); text\$ A string of text. n The number of characters to return. Returns The last n characters of the string, or all the string if it is less than n characters. See also:Asc(), Chr\$(), DelStr\$(), InStr(), LCase\$(), Left\$(), Len(), Mid\$(),</pre>
Round()	<pre>Rounds a real number or an array of reals to the nearest whole number. Func Round(x x[] x[][]); x A real number or an array of reals. Returns If x is an array it returns 0. Otherwise it returns a real number with no fractional part that is the nearest to the original number. See also:Frac(), Max(), Min(), Trunc()</pre>
SampleAbort()	<pre>This cancels sampling a file or file sequence and closes any associated data, result and cursor views without saving them. If a result view derived from the data has been saved, the saved file remains. It is equivalent to the Abort button in the sampling control panel. Func SampleAbort(); Returns 0 if sampling was aborted, or a negative error code. See also:SampleKey(), SampleReset(), SampleStart(), SampleStop(), SampleStatus(), SampleText(), SampleWrite()</pre>
SampleAutoComment()	<pre>This gets or sets file auto-commenting state as set in the sampling configuration dialog. Func SampleAutoComment({yes%}); yes% If present a non-zero value turns on automatic prompting for file comments when sampling ends and zero turns it off. If absent, no change is made. Returns the automatic commenting flag at the time of the function call, 0=off, 1= on. See also:FileComment\$(), SampleAutoFile(), SampleAutoName\$()</pre>
SampleAutoCommit()	<pre>This gets or sets the automatic file commit period as set in the sampling configuration dialog. There is an overhead associated with the commit operation; do not set needlessly short time periods. A period of 0 means no file commit. Func SampleAutoCommit({every%}); every% If present this sets the automatic file commit interval in seconds. The value is</pre>

SampleAutoFile()	Func SampleAutoFile({yes%});			
	<pre>gets or sets the automatic filing state as set in the sampling configuration dialog. SampleAutoFile({yes%});</pre>			
	Returns the automatic filing state at the time of the function call $0=$ off $1=$ on			
	See also:FileSave(), FileSaveAs(), FilePathSet(), SampleAutoName\$(), SampleAutoComment(), SampleRepeats()			
SampleAutoName\$()	This gets or sets the template for file auto-naming as in the sampling configuration dialog. The path is set by FilePathSet(). The name must be set before you open the file for sampling for the automatic name to be used.			
	<pre>Func SampleAutoName\$({name\$});</pre>			
	name\$ If present, this is the template for file auto-naming. An empty string turns off auto-naming. The maximum name length is 20 characters. See the sampling configuration documentation for details on the template string.			
	turns The auto-naming template at the time of the function call.			
	See also:FilePathSet(),SampleAutoFile(),SampleAutoComment()			
SampleBar()	This gives you access to the Sample toolbar. The format of strings passed by this command is the button label (up to 8 characters), followed by a vertical bar, followed b the full path name to a sampling configuration file, including the .s2c file extension followed by a vertical bar then, a comment to display when the mouse pointer is over the button. If you call the command with no arguments it returns the number of buttons in the toolbar.			
	tooldar. Func SampleBar({n% {, &get\$}}); Func SampleBar(set\$);			
	n% If set to -1, get\$ must be omitted, all buttons are cleared and the function returns 0. When set to the number of a button (the first button is 0), get\$ is as described above. In this case, the function returns -1 if the button does not exist, 0 if it exists and is the last button, and 1 if higher-numbered buttons exist.			
	set\$ The string passed in should have the format described above. The function returns the new number of buttons or -1 if all buttons are already used.			
	Returns See the descriptions above. Negative return values indicate an error.			
	For example, the following code clears the script bar and sets two buttons:			
	SampleBar(-1); 'clear all buttons SampleBar("Fast C:\\Spike3\\Fast.s2c Fast 4 channel sampling"); SampleBar("Faster C:\\Spike3\\FastXX.s2c Very fast sampling");			
	See also:App()			
	This sets the waveform or WaveMark channel calibration. It changes the units scale and			
SampleCalibrate()	offset fields only. See the Sampling Configuration dialog for a description of these fields.			
	Func SampleCalibrate(chan%, units\$, scale, offset);			
	change. The channel number of a waveform or WaveMark channel (1-100).			

units\$ The units to use. If the string is longer than 5 characters only the first 5 are used.

- scale The scale factor for the channel.
- offset The offset for the channel.

Returns 0 if all was well, or a negative error code.

See also:SampleTimePerAdc(), SampleTitle\$(), SampleUsPerTime(), SampleWaveform(), SampleWaveMark()

SampleChannels()

This function returns and optionally sets the maximum number of channels in the data file. The next time you sample, the data file will have space for this number of channels. Spike2 version 5 allows up to 256 channels. Version 4 reads files with up to 100 channels. Version 3 can only read files with 32 channels. The maximum number of channels that you can sample is 100.

Func SampleChannels({nChan%});

nChan[&] If present, this sets the number of channels in the next data file created for sampling in the range 32-256. Numbers outside this range are set to the nearer limit. The maximum channel number for sampling is the smaller of the maximum channels in the file or 100.

Returns The number of channels set at the time of the call.

See also:FileNew(), FileSaveAs(), SampleClear(),

SampleClear()

This sets the sampling configuration to a standard state. The data file is set to 32 channels (all Off except the keyboard), sampling mode is Continuous, microseconds per time is 10, time per ADC is 10, sequence files are disconnected, *Stop sampling when*... fields are disabled and the automatic file name template is cleared. All sample rate optimising is disabled (equivalent to SampleOptimise(0,0,0)).

```
Proc SampleClear();
```

```
See also:SampleLimitSize(), SampleLimitTime(), SampleOptimise(),
SampleTimePerAdc(), SampleUsPerTime()
```

SampleComment\$()

This function gets and sets the comment attached to a channel in the Sampling Configuration dialog.

Func SampleComment\$(chan% {,new\$});

chan% The channel number in the window (1 to 100).

new\$ If present, the new comment. If the comment is too long, it is truncated.

Returns The original comment, or an empty string if the channel number is illegal.

See also:SampleDigMark(), SampleEvent(), SampleTextMark(), SampleTitle\$(), SampleWaveform(), SampleWaveMark()

SampleDigMark()

This adds the digital marker channel to the sampling configuration.

Func SampleDigMark(rate)

rate The expected sustained rate for digital markers on this channel in Hz. Returns 0 if all went well, or a negative error code.

See also:SampleEvent(), SampleTextMark(), SampleTitle\$()

SampleEvent()	This function sets a channel to sample event data.						
Sample Vent()	<pre>Func SampleEvent(chan%, port%, type%, rate);</pre>						
	chan% The channel number in the file to use for this data (1-100).						
	port% The event port number.						
	type%The type of the event channel: 0=Events on a falling edge (Event-), 1=Events on a rising edge (Event +), 2=Events on both edges (Level)rateThe expected maximum sustained event rate on the channel in Hz.						
	Returns 0 if all went well, or a negative error code.						
	<pre>See also:SampleClear(),SampleDigMark(),SampleTitle\$()</pre>						
SampleHandle()	Returns view handles linked to sampling. This can be used to position, show and hide the						
	output sequencer and sampling control panels. App () also returns control panel nandles.						
	Func SampleHandle (which%);						
	which Selects which view handle to return: 0 Sampling time view 1 Sampling control panel 2 Sequencer control panel						
	Returns The view handle or 0 if the view does not exist						
	See also:App(), View(), ViewList(), Window(), WindowVisible()						
SampleKey()	Adds an event to the keyboard marker channel of a sampling file as if you had typed it, including triggering the output sequencer and arbitrary waveform output. There was no return value before version 5.04.						
	<pre>Func SampleKey(key\$);</pre>						
	key\$ The first character of the string is added to the keyboard marker channel.						
	Returns The time stamp of the added marker in seconds or -1 if no file is sampling.						
	See also:PlayWaveAdd(),SampleText(),SampleWrite()						
SampleLimitSize()	This corresponds to <i>Starting and stopping: Stop at file size</i> in the Automation dialog.						
	Func SampleLimitSize({size});						
	size The size limit for the output file, in KB. A positive value sets the size and enables the limit. A negative value sets the limit to the positive value of size, but disables the limit. A zero value, or omitting the argument, means no change.						
	Returns The limit before the call. If the limit is disabled, the size is returned negated.						
	See also:SampleClear(),SampleLimitTime(),SampleMode(),SampleRepeats()						
Samplel imitTime()	This corresponds to Starting and stopping: Stop at time in the Automation dialog.						
	<pre>Func SampleLimitTime({time});</pre>						
	time The time in seconds to set as a limit. A positive time sets the limit and enables it. A negative time sets the limit to the positive time, but disables the limit. A value of zero, or omitting the argument, leaves the time limit unchanged.						
	Returns The limit before the call. If the limit is disabled, the time is returned negated.						
	See also :SampleClear(), SampleLimitSize(), SampleMode(), SampleRepeats()						

SampleMode()

This function sets and gets the sampling mode in the Sampling configuration dialog. To set Triggered mode you should use SampleTrigger() as mode 3 is maintained here for backwards compatibility with previous versions of Spike2.

Func SampleMode(mode%{, time, trig%|every});

- mode% This argument determines the action of the command:
 - 0 Returns the current mode as 1, 2, or 3 for Continuous, Timed or Triggered. Any additional arguments are ignored.
 - 1 Sets Continuous recording mode.
 - 2 Sets Timed recording mode. All three arguments are required.
 - 3 Sets Triggered recording mode. All three arguments are required.
 - -1 Returns the For field value, valid in Timed mode.
 - -2 Returns the value of the Every field, valid for Timed sampling mode.
 - -3 Returns the Trigger field value, valid in Triggered sampling mode.
- time In modes 2 and 3 it sets the For field of the Sampling mode group, in seconds.
- every In mode 2 it sets the Every field for timed sampling, in seconds.
- trig[®] The trigger channel number for mode 3. This channel must exist in the sampling configuration and can be of any type except a waveform or WaveMark channel.
- Returns For modes 0-3 it returns the current sampling mode as 1-3, or a negative error code. In mode -1 to -3 it returns the information described for mode[®], above.

If you set triggered mode, triggers 2 to 4 are disabled and the command is equivalent to calling SampleTrigger(1, trig%, -1, 0.0, time, -1).

See also:SampleClear(), SampleTrigger()

SampleOptimise()

This function sets and gets the sample rate optimising settings in the Resolution tab of the Sampling configuration dialog. Spike2 can match the requested sample rates for waveform and WaveMark channels by changing the microseconds per time unit, time units per ADC convert and by making waveform channels into Quick or Slow channels.

The available ADC sample rate is divided by giving one share to each WaveMark channel, one to each Quick channel and one to the entire Slow channel group. The actual rates for each waveform channel (Quick or Slow) can then be further down-sampled by taking one point in n. For example, with 4 Slow channels, the fastest rate for a Slow channel is one quarter the fastest rate for a Quick channel.

The rates are optimised each time you change any aspect of the sampling configuration that affects the sample rate. With full optimise set and slow sampling rates it can take an appreciable time for Spike2 to search all possible combinations for the best sample rates. For this reason, SampleClear() sets no optimise and version 3 compatibility. We recommend that you use this command after all the other sample setup commands so that you pay the time penalty for optimising once.

Func SampleOptimise(opt%, group%, type% {,usLo% {,usHi%}});
Func SampleOptimise(get%);

- opt% This sets the optimise method. Values are:
 - 0 No optimise of microseconds per time or time per ADC.
 - 1 Partial, optimise time per ADC, no optimise of microseconds per time.
 - 2 Full, optimise both time per ADC and microseconds per time.

group% This controls the use of Quick channels. Values are:

- 0 Version 3 compatible; no Quick channels and channel divides up to 65535.
- 1 Group channels with the same ideal rate so they get the same actual rate.
- 2 Optimise for the least error in rate; channels with the same ideal rate may get different actual rates.

CAbc	Set the type of 1401 to optimise for:							
	0 Works with all supported 1401s	0 Works with all supported 1401s except a 1401 <i>plus</i> with an old ADC.						
	1 1401 <i>plus</i> with an old ADC.	4 Micro1401 mk II						
	2 Power1401	5 Power1401 625						
	3 1401plus or a micro1401	6 micro1401 (not mk II)						
usLc	If present, it sets the low limit, in microseconds, for optimising microseconds per time unit (when opt% is 2). Values outside the range 1 to 1000 are ignored.							
usHi	If present, it sets the high limit, in microseconds, for optimising microseconds per time unit. Values outside the range 1 to 1000 are ignored.							
get%	Use the function with one argument 0 to 4 read back the current values of	t to read back the current settings. The values of opt%, group%, type%, usLo% and usHi%.						

- Returns The setting version returns the current optimise method. When called with one argument, the return values are as documented for get%.
- See also:SampleClear(), SampleLimitSize(), SampleLimitTime(), SampleTimePerAdc(), SampleUsPerTime()

SampleRepeats()

This gets and optionally sets the Sampling configuration Automation tab Repeats field, which is used to sample a sequence of files with automatically incrementing names. You can use it online to return the number of files left to sample. You can choose to trigger each file in the sequence, or just the first with SampleStartTrigger().

Func SampleRepeats({files%});

- files% The number of files to sample; 0 and 1 both sample a single file. Set -1 to return the number of files that remain to be sampled (including the current file). If you omit this argument, the number of files is not changed.
- Returns The number of repeats set at the time of the call or the number of files that remain to be sampled.

You must also set a file name with SampleAutoName\$(), set automatic file saving with SampleAutoFile(1) and set a time or size limit with SampleLimitTime() or SampleLimitSize(), or load a suitable sampling configuration with these values set.

To make life easier for the script programmer, each file in the sequence has the same view handle, which is the handle of the first file you create with FileNew(). As each file reaches the time or size limit, it is saved, closed and replaced with the next file. The view handle always references a file, but its file name and sampling status will vary. The last sampled file in the sequence is not closed. You can use SampleStop() or SampleAbort() to end a sequence early and SampleReset() to restart a file.

```
See also:SampleAutoFile(), SampleAutoName$(), SampleLimitSize(),
SampleLimitTime(), SampleStartTrigger(), SampleStop()
```

SampleReset()

This abandons sampling, deletes any data that has been written to disk, and returns to a state as if FileNew() had just been used to create a new data file. If it is used when sampling a sequence of files, the current file is restarted and the sequence continues.

Func SampleReset();

Returns 0 if the reset operation completed without a problem, or a negative error code.

See also:SampleAbort(),SampleKey(),SampleStart(),SampleStop(),
 SampleStatus(),SampleText(),SampleWrite()

SampleSeqCtrl()	This sets and gets options that control the use of the output sequence. Currently there is only one option.			
	<pre>Func SampleSeqCtrl(opt%{, new%});</pre>			
	opt% There is currently only one option: 1 = get or set the sequencer jump control. This is the same as the Sequencer jumps controlled by setting in the Sequencer tab of the Sampling configuration dialog.			
	new% The new value for the control option. For the jump control: 0 = keyboard, control panel and script, 1 = control panel and script, 2 = script only.			
	Returns If you are setting a value or this is used at an inappropriate time, the function returns 0. If you are reading a value, the function returns the value.			
	See also:SampleKey(),SampleSeqStep(),SampleSequencer\$(),SampleStart()			
SampleSeqStep()	This returns the current sequencer step or -1 if not sampling. If no sequence is running the result is usually 0 (but this is not guaranteed).			
	<pre>Func SampleSeqStep();</pre>			
	Returns The current sequence step number, or -1 if not sampling.			
	<pre>See also:SampleKey(), SampleSequencer\$(), SampleSeqVar()</pre>			
SampleSeqTable()	If there is a sampling document with an output sequence, you can use this function to find the size of any table set in the sequence by the TABSZ directive. You can also use this to transfer data between an integer array and the table.			
	<pre>Func SampleSeqTable({tab%[]{, offs%{, get%}});</pre>			
	tab%[] An integer array holding items to transfer to the 1401 sequencer table or to hold items read back from the table. The array size sets the maximum item count			
	offs% This sets the index into the sequencer table to start the transfer. The first index in the table is 0. If this value is negative or greater than or equal to the sequencer table size, no data is transferred. If omitted, the value 0 is used.			
	get% Set 0 or omit this argument to transfer data to the sequencer table, set to 1 to transfer data from the sequencer table.			
	Returns If you call this with no arguments, the return value is the size of the sequence table. Otherwise, the returned value is the number of items transferred between the sequencer table and the array. A negative error code is also possible, fo example -1 if there is no sampling document.			
	See also:SampleKey(),SampleSequencer\$(),SampleSeqVar()			
SampleSequencer()	You can use this function to set the sequencer file to attach to the Sampling configuration. Use SampleSequencer\$() to get the name of the current sequencer file.			
	<pre>Func SampleSequencer(new\$);</pre>			
	new\$ The name of the sequence file. Pass an empty string to set no sequencer file.			
	Returns It returns 0 if all was well, or a negative error code.			
	See also:SampleKey(), SampleSequencer\$(), SampleSeqVar()			

SampleSequencer\$()	This function returns the name of the sequencer file that is currently attached to the sampling configuration. Use SampleSequencer() to set the file.				
	<pre>Func SampleSequencer\$();</pre>				
	Returns It returns the current sequencer file name, or an empty string if there is no file. The returned name includes the full path.				
	See also:SampleKey(),SampleSequencer(),SampleSeqVar()				
SampleSeqVar()	This is used during sampling with an output sequence, to get or set the value of an output sequencer variable. Values set before the sampling window exists are ignored. From 4.04, values set before SampleStart() set the initial variable values. Previously they took effect after sampling started.				
	<pre>Func SampleSeqVar(sVar%{, new%});</pre>				
	sVar% The sequencer variable to set or read, in the range 1 to 64.				
	new% The new value for the output sequencer variable. If present, the value of the variable is updated. Omit to return the variable value. A common error when setting variables for the DAC instruction is to set a value 65536 times too small.				
	Returns If you are setting a value, or this is used at an inappropriate time, the function returns 0. If you are reading a value, the function returns the value.				
	See also:SampleKey(), SampleSeqStep(), SampleSequencer\$(), SampleStart()				
SampleStart()	This function can be used after FileNew() has created a new time view based on the current Sampling configuration. It starts sampling immediately, or on a trigger. Func SampleStart({trig%});				
	<pre>trig% 0 = start sampling immediately (default), 1 = wait for a Trigger input signal (1401plus Event 3 input), -1 = use the sample configuration trigger setting.</pre>				
	See also: Sample Depart (). Sample Departs (). Sample Depart (). Sample Stop ().				
	SampleStartTrigger(), SampleStatus(), SampleWrite()				
SampleStartTrigger()	This gets and optionally sets the trigger option used when sampling starts. It is equivalent to <i>Starting and stopping: Triggering</i> in the Sampling configuration Automation tab.				
	<pre>Func SampleStartTrigger({trig%});</pre>				
	trig% Omit this argument for no change. The following values can be used:				
	 -1 Use whatever state the user previously set in the sampling control panel 0 Sampling is not triggered 1 Sampling starts on a trigger; repeated files are not triggered 2 Sampling starts on a trigger; repeated files are triggered 				
	Returns The trigger state (as for the trig% argument) at the time of the call.				
	See also:SampleAbort(), SampleRepeats(), SampleReset(), SampleStop(),				

SampleStart(), SampleStatus(), SampleWrite()

	This function enquires about the state of any sampling			
SampleStatus()	Funce SampleStatus():			
	Potume A code indicating the compline state on 1 if there is no compliant			
	 A time view is ready to sample, but it has not been told to start yet Sampling is waiting for an Event 3 trigger Sampling is now in progress Sampling is gtorning (the code changes to 1 when it has storned) 			
	See also: Comple Theast (). Comple Decet (). Comple Chart (). Comple Chart ().			
	See also.sampleAbort(), sampleReset(), samplestart(), samplestop()			
SampleStop()	This function stops sampling in progress, and is equivalent to the Stop button of the floating command window. The function does not return until sampling has stopped. If used in a file sequence, the sequence is cancelled and the current file is saved.			
	<pre>Func SampleStop();</pre>			
	Returns 0 if sampling stopped correctly or a negative error code.			
	<pre>See also:SampleAbort(), SampleKey(), SampleReset(), SampleStart(), SampleStatus(), SampleText(), SampleWrite()</pre>			
SampleText()	This function adds text to a text marker channel during sampling. Text marker channels are created with the SampleTextMark() command.			
	<pre>'unc SampleText(text\$ {,time {,code%[]}});</pre>			
	text\$ The text string to attach to the text marker. If the string is longer than the maximum length set for the channel, extra characters are ignored.			
	Time The time for the text marker. If this argument is omitted, negative, less than the time of the last text marker, greater than the current sampling time or this is the triggered sampling mode trigger channel, then the current sampling time is used			
	code% The first four elements of this array set the marker codes stored with the text string. If this argument is omitted the codes are set to 0. Codes are limited to the range 0-255. Only the lower 8 bits of codes outside this range are stored.			
	Returns 0 if all was OK, or a negative error code.			
	<pre>See also:SampleAbort(), SampleKey(), SampleReset(), SampleStart(), SampleStop(), SampleStatus(), SampleTextMark()</pre>			
SampleTextMark()	This function sets channel 30 as a text marker channel. Each event on a text marker channel holds a time, marker codes and a text string. You can add text markers to this channel using the SampleText() command and from a serial line.			
	Func SampleTextMark(max%{, port%{, term\$ {, baud%{, bits%{, par%{, stop%{, hsk%}}}}});			
	max% This sets the maximum number of characters that can be attached to each text marker in the range 1 to 200. If 0 is passed, the channel is deleted from the list.			
	port [®] Serial port, in the range 1 to 9, to use to read on-line TextMark data. If omitted, no serial data is read. When reading, characters with codes less than 32 are ignored unless they are the terminator character.			
	term ^{\$} Optional terminating character for serial line read. If omitted, "\r" is used (carriage return, character code 13).			

baud%	This sets the serial line Baud rate (number of bits per second). The maximum character transfer rate is of order one-tenth this figure. All standard Baud rates from 50 to 115200 are supported. If you do not supply a Baud rate, 9600 is used.			
bits%	The number of data bits to encode a character (7 or 8). If omitted, 8 is set.			
par%	Set this to 0 for no parity check, 1 for odd parity or 2 for even parity. If you do not specify this argument, no parity is set.			
stop%	This sets the number of stop bits as 1 or 2. If omitted, 1 stop bit is set.			
hsk%	This sets the handshake mode, sometimes called flow control. 0= no handshake, 1=hardware handshake, 2=XON/XOFF protocol. If omitted, 0 is set.			
Returns 0 or a negative error code.				
<pre>See also:SampleText(), SampleComment\$(), SampleClear()</pre>				

SampleTimePerAdc()

This sets and gets the number of time units set by SampleUsPerTime() for each ADC conversion. The product of the time per ADC and the microseconds per unit time must be at least 3 for Power1401 and 6 for 1401*plus* and micro1401. Lower values will cause sampling to fail. If the optimise method is set by SampleOptimise() is not 0, this call will have no effect as time per ADC will be recalculated.

Func SampleTimePerAdc({new%});

new% The number of clock ticks per conversion in the range 1 to 32767. If this is omitted the value is not changed. Illegal values stop the script with a fatal error.

Returns The value of the time per ADC convert at the time of the call.

See also:SampleOptimise(), SampleUsPerTime(), SampleWaveform()

SampleTitle\$()

This gets and sets the title attached to a channel in the Sampling Configuration dialog.

Func SampleTitle\$(chan% {,new\$});

chan% The channel number.

new\$ If present, the new title. If the title is too long, it is truncated.

Returns The title at the time of the call, or an empty string for illegal channel numbers.

See also:SampleCalibrate(), SampleComment\$(), SampleClear()

SampleTrigger()

There are four triggers that can be set to determine which data for a channel is written to disk. Each trigger is associated with a list of channels. The first trigger is also set by SampleMode() for backwards compatibility. You cannot change the triggering after sampling has started. Channels not associated with a trigger are written continuously.

Func SampleTrigger(trig%, scr%, code%, relS, relE{, cSpc});
Func SampleTrigger(trig%, get%);

- trig% A trigger number in the range 1 to 4. Setting a trigger sets triggered sampling mode. This is cancelled by SampleClear() or by calling SampleMode().
- A source channel for triggers. The channel is not checked here, so you can put in anything that could be a channel number. However, it is checked when you start sampling and if this is not an event-based channel the trigger is ignored.
- code% If the src% channel is Marker-based, you can choose to only trigger when the first marker code matches this value (in the range 0 to 255). Set -1 to accept all markers. If the source is not a marker, code% is ignored.

relS	The start of the triggered area relative to the trigger time, in seconds. This can be
	negative or positive, but must be less than relE. There are limits on how far
	into the past a trigger can go; there is a warning message when sampling starts if
	the triggering requests are impossible.

- relE The end of the triggered area relative to the trigger time. This must be greater than rels. Spike2 saves data in buffers so you may get more data saved to disk that you requested.
- cSpc A channel specification for the channels to trigger. If you omit this argument all channels are selected for triggering. You cannot use the standard channel specification values -2 or -3.
- get % You can use the second command variant to return trigger values. Get values: -1 -2 -3 -4 -5

scr% code% relS relE CSpc, -1=all, 1-n for one, 0=multiple

If more than one channel is selected, the channel specification returns 0.

Returns The first variant returns 0 or an error code. The second variant returns the requested value.

See also:SampleClear(), SampleMode()

SampleUsPerTime() This gets and optionally sets the basic time unit used for sampling. If the optimise method is set to full by SampleOptimise(2,...), changes to the basic time units have no effect as Spike2 chooses the best value to optimise the waveform sample rates.

Func SampleUsPerTime({new});

new If present, this sets the basic time unit in the range 2-1000 microseconds. Out of range values cause a fatal script error. If you sample with a Power1401 or Micro1401 mk II, the value is rounded to the nearest 0.1 microseconds. For all older 1401 types the value is rounded to the nearest microsecond.

Returns The current value of microseconds per time unit.

See also:SampleLimitTime(), SampleOptimise(), SampleTimePerAdc()

SampleWaveform()

This function adds a waveform channel to the list of channels required. If the channel is already in use, it is replaced. The units, scale and offset fields are set to the standard Spike2 defaults (input in Volts). The title and comment of the channel are not changed.

Func SampleWaveform(chan%, port%, ideal);

- chan% The channel number to use for the new channel in the range 1 to 100.
- port% An unused (by a waveform channel) 1401 waveform port in the range 0-31.
- ideal The ideal sampling rate that you would like for the port in Hz. Remember that you may not get this rate. Spike2 will set the nearest rate it can.

Returns 0 if all went well, or a negative error code.

See also:SampleOptimise(), SampleTimePerAdc(), SampleUsPerTime()

SampleWaveMark()

This adds a WaveMark channel to the sample list. The units, scale and offset fields are set to the Spike2 defaults (input in Volts). The channel title and comment are not changed. See the Sampling Configuration dialog description for sample rate details.

Func SampleWaveMark(chan%, port%, rate,size%,pre%{,ideal{,nTr%}});

chan%	The channe	l number to	use for the	new channel	in the range 1	to 100.
-------	------------	-------------	-------------	-------------	----------------	---------

port%	The	1401	waveform	port in	the range	0-3	1.
-------	-----	------	----------	---------	-----------	-----	----

- rate The estimated maximum sustained spike rate, used to allocate buffer space.
- size% The number of waveform samples to save for each WaveMark This must be an even number in the range 10 to 126.
- pre% The number of points before the peak/trough of each spike, range 0 to size%-1;
- ideal If present, this sets the ideal sample rate, in samples per second, for all WaveMark channels. This must be in the range 1.0 to 500000.0 Hz.
- nTr% If present, this sets the number of traces to sample as 1, 2 or 4. The default is 1.

Returns 0 if all went well, or a negative error code.

See also:SampleOptimise(), SampleTimePerAdc(), SampleUsPerTime()

SampleWrite()

This controls writing data to the file during sampling. You can enable and disable writing of some or all channels, and obtain the sampling state of channels and channel groups. Func SampleWrite(write% {, chan%|chan%[]});

write% This determines the action taken by the command:

- -1 Report on the state of the channel (or channels) given by the next argument.
- 0 Disable writing to disk (pause) the channels given by the next argument
- 1 Enable writing to disk (un-pause) the channels given by the next argument
- chan% This sets the channels to operate on. If omitted, all channels are used. If this is an integer, it is the channel number. If it is an array, index 0 holds the number of channels following; the remaining elements hold a list of channel numbers.
- ReturnsThe state of writing for the channel or channels:0Disabled1Enabled2Some channels in list are enabled

See also:SampleAbort(), SampleKey(), SampleReset(), SampleStart(), SampleStop(), SampleStatus(), SampleText()

ScriptBar()

This controls the Script toolbar. Call the command with no arguments to return the number of toolbar buttons. The first button is numbered 0.

Func ScriptBar({nBut%{, &get\$}});
Func ScriptBar(set\$);

- nBut% Set -1 and omit get\$ to clear all buttons and return 0. Otherwise it is a button number and returns -1 if the button does not exist, 0 if it is the last button, and 1 if higher-numbered buttons exist. get\$ returns the information as for set\$.
- set\$ This holds up to 8 characters of button label, a vertical bar, the path to the script file including .s2s, a vertical bar and a pop-up comment. The function returns the new number of buttons or -1 if all buttons are already used.

Returns See the descriptions above. Negative return values indicate an error.

For example, the following code clears the script bar and sets a button: ScriptBar(-1); 'clear all buttons ScriptBar("ToolMake|C:\\Scripts\\ToolMake.s2s|Build a toolbar");

See also:App()

ScriptRun() This sets the name of a script to nu when the current script terminates. You can pass information to the new script using disk files or by using the Profile() command. You can call this function as often as you like; only the last use has any effect. Proc ScriptRun (name%(, flags%)); name% The script file. If you supply a plath relative to the current folder or a full plath to the script file. If you supply a relative path, it must still be valid at the end of the current script. Set name% to "" to cancel running a script. flags% Optional flags that control the new script. If omitted, 0 is used. The only flag defined now is 1 = run new script even if the current script and in an error. If the file you name does not exist when Spike2 tries to run it, nothing happens. If the nominated script is not already loaded, Spike2 will load it, run it and unload it. Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the time: If you want the time into sampling, use the Maxtime() function. Func Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the time: If you want the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(), MaxTIme(), 7ime\$(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCout(), EditFaste(), MoveBy(), MoveTo() See also:Ed		
 Proc ScriptRun(name\${, flags%}); name\$ The script file to run. You can supply a path relative to the current folder or a full path to the script file. If you supply a relative path, it must still be valid at the end of the current script. Set name5 to "" to cancel running a script. flags\$ Optional flags that control the new script If omitted, 0 is used. The only flag defined now is 1 = run new script even if the current script ends in an error. If the file you name does not exist when Spike2 tries to run it, nothing happens. If the nominated script is not already loaded, Spike2 will load it, run it and unload it. See also:App(), Profile() Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the timer. If you want the time into sampling, use the Maxtime() function. Func Seconds((set)); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(), MaxTime(), Time\$(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func \$election\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:Date\$(), DateCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(); port\$ The serial port oclose a defined for SerialOpen(). Returns O an anguive error code See also:SerialOpen(), SerialPered(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial p	ScriptRun()	This sets the name of a script to run when the current script terminates. You can pass information to the new script using disk files or by using the Profile() command. You can call this function as often as you like; only the last use has any effect.
 name5 The script file to run. You can supply a path relative path, it must still be valid at the end of the current script. Stame5 0 "" to cancel running a script. flags% Optional flags that control the new script. If omitted, 0 is used. The only flag defined now is 1 = run new script even if the current script. State to run it, nothingpens. If the file you name does not exist when Spike2 tries to run it, nothingpens. If the nominated script is not already loaded, Spike2 will load it, run it and unload it. See also:App(), Profile() Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the timer. If you want the time into sampling, use the Maxtime() function. Funce Seconds((set)); set the time. If You want the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(), MaxTime(), Time\$(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), BditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a soft and for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialTead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waining for serial data. There is an internal buffor of 1024 charactery nev or thot is filled wainty or serial data. There is an in		<pre>Proc ScriptRun(name\${, flags%});</pre>
<pre>flags% Optional flags that control the new script. If omitted, 0 is used. The only flag defined now is 1 = run new script even if the current script ends in an error. If the file you name does not exist when Spike2 will load it, run it and unload it. See also:App(), Profile()</pre> Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the timer. If you want the time into sampling, use the Maxtime() function. Func Seconds((set)); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:DateS(), MaxTime(), TimeS(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is imappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); ports The serial port oclose as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount()		name\$ The script file to run. You can supply a path relative to the current folder or a full path to the script file. If you supply a relative path, it must still be valid at the end of the current script. Set name\$ to "" to cancel running a script.
If the file you name does not exist when Spike2 tries to run it, nothing happens. If the nominated script is not already loaded, Spike2 will load it, run it and unload it. See also:App(), Profile() This returns the time in seconds. This is used for relative time measurements. You can also set the timer. If you want the time into sampling, use the Maxtime() function. Func Seconds((set)); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(), MaxTime(), Time\$(), TimeDate() This function returns the text in the current view that is currently selected. Func Selection\$() Returns The current text selection, If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose() The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(). SerialWrite(), SerialRead(), SerialCount() See also:SerialOpen(). SerialCount() This counts the characters or items buffered in a serial po		flags% Optional flags that control the new script. If omitted, 0 is used. The only flag defined now is 1 = run new script even if the current script ends in an error.
See also:App(), Profile() Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the timer. If you want the time into sampling, use the Maxtime() function. Func Seconds((set)); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(),MaxTime(),Time\$(),TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(),EditCut(),EditPaste(),MoveBy(),MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port\$); port\$ The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a scrial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use to the is an internal buffer of 1024 characters per port that is filled when you use.		If the file you name does not exist when Spike2 tries to run it, nothing happens. If the nominated script is not already loaded, Spike2 will load it, run it and unload it.
Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the time: If you want the time into sampling, use the Maxtime() function. Func Seconds((set)); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(),MaxTime(),Time\$(),TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(),EditCut(),EditPaste(),MoveBy(),MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port\$); port% The serial port to closes as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. They is an internal buffer of 1024 characters per port that is filled when you use.		See also:App(), Profile()
Seconds() This returns the time in seconds. This is used for relative time measurements. You can also set the time. If you want the time into sampling, use the Maxtime() function. Func Seconds((set)); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(),MaxTime(),Time\$(),TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(),EditCut(),EditPaste(),MoveBy(),MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port\$); port The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen().		
<pre>Func Seconds({set}); set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(),MaxTime(),Time\$(),TimeDate() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(),EditCut(),EditPaste(),MoveBy(),MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose() Returns 0 or a negative error code See also:SerialOpen(),SerialWrite(),SerialRead(),SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters pro tot that is filled when you mere serial count is filled when you mere that is filled when you mere there is an internal buffer of 1024 characters pro the tasks while waiting for serial data. There is an internal buffer of 1024 characters pro thet is filled when you mere there is fully and the prove the protein the start of the above you mere there is fully and the prove the p</pre>	Seconds()	This returns the time in seconds. This is used for relative time measurements. You can also set the timer. If you want the time into sampling, use the Maxtime() function.
<pre>set If present, this sets the time in seconds. The time is 0 when Spike2 starts. Returns The time in seconds. This is the value before any new time is set. See also:Date\$(), MaxTime(), Time\$(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters pro trut is filled when you use</pre>		<pre>Func Seconds({set});</pre>
Returns The time in seconds. This is the value before any new time is set. See also:Date\$(), MaxTime(), Time\$(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). There is an internal buffer of 1024 characters near nor that is filled when you use		set If present, this sets the time in seconds. The time is 0 when Spike2 starts.
See also:Date\$(), MaxTime(), Time\$(), TimeDate() Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port\$); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per nort that is filled when you use		Returns The time in seconds. This is the value before any new time is set.
Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		<pre>See also:Date\$(),MaxTime(),Time\$(),TimeDate()</pre>
Selection\$() This function returns the text in the current view that is currently selected. Func Selection\$(); Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Let his to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters pre nor that is filled when you use		
Returns The current text selection. If there is no text selected, or if the view is inappropriate for this action, an empty string is returned. See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use	Selection\$()	This function returns the text in the current view that is currently selected.
See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		Returns The current text selection. If there is no text selected, or if the view is
SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). This counts the characters or items buffered in a serial port opened by SerialOpen(). This counts the characters or items buffered in a serial port opened by SerialOpen(). SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen().		See also:EditCopy(), EditCut(), EditPaste(), MoveBy(), MoveTo()
SerialClose() This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it. Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		
<pre>Func SerialClose(port%); port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use</pre>	SerialClose()	This function closes a serial port opened by SerialOpen(). Closing a port releases memory and system resources. Ports are automatically closed when a script ends, however it is good practice to close a port when your script has finished with it.
<pre>port% The serial port to close as defined for SerialOpen(). Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use</pre>		Func SerialClose(port%);
Returns 0 or a negative error code See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		port% The serial port to close as defined for SerialOpen().
See also:SerialOpen(), SerialWrite(), SerialRead(), SerialCount() SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		Returns 0 or a negative error code
SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		See also:SerialOpen(),SerialWrite(),SerialRead(),SerialCount()
SerialCount() This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use		
SerialCount(). The size of this buffer limits the number of characters that this function can tell you about. To avoid character loss when you are not using a serial line handshake do not buffer up more than a few hundred characters with SerialCount().	SerialCount()	This counts the characters or items buffered in a serial port opened by SerialOpen(). Use this to detect input so your script can do other tasks while waiting for serial data. There is an internal buffer of 1024 characters per port that is filled when you use SerialCount(). The size of this buffer limits the number of characters that this function can tell you about. To avoid character loss when you are not using a serial line handshake do not buffer up more than a few hundred characters with SorialCount().

Func SerialCount(port% {,term\$});

- port% The serial port to use as defined for SerialOpen().
- term\$ An optional string holding the character(s) that terminate an input item.
- Returns If term\$ is absent or empty, this returns the number of characters that could be read. If term\$ is set, this returns the number of complete items that end with term\$ that could be read.

See also:SerialOpen(), SerialWrite(), SerialRead(), SerialClose()

SerialOpen()

This function opens a serial port and configures it for use by the other serial line functions. It is not an error to call SerialOpen more than once on the same port. The serial routines use the host operating system serial line support. Consult your system documentation for information on serial line connections and Baud rate limits.

Func SerialOpen(port%{, baud%{, bits%{, par%{, stop%{, hsk%}}}});

- port[®] The serial port to use, in the range 1 to 9. The number of ports depends on the computer. Two ports (1 and 2) are common on both PC and Macintosh systems.
- baud% This sets the serial line Baud rate (number of bits per second). The maximum character transfer rate is of order one-tenth this figure. All standard rates from 50 to 115200 Baud are supported. If you omit baud%, 9600 is used.
- bits% The number of data bits to encode a character. Windows supports 4 to 8 bits, the Macintosh supports 7 or 8. If bits% is omitted, 8 is set. Standard values are 7 or 8 data bits. If you set 7 data bits, character codes from 0 to 127 can be read. If you set 8 data bits, codes from 0 to 255 are possible.
- pars Set this to 0 for no parity check, 1 for odd parity or 2 for even parity. If you do not specify this argument, no parity is set.
- stop% This sets the number of stop bits as 1 or 2. If omitted, 1 stop bit is set. If you specify 5 data bits, a request for 2 stop bits results in 1.5 stop bits being used.
- hsk% This sets the handshake mode, sometimes called "flow control". 0 sets no handshake, 1 sets a hardware handshake, 2 sets XON/XOFF protocol.

Returns 0 or a negative error code.

See also:SerialWrite(), SerialRead(), SerialCount(), SerialClose()

SerialRead()

This function reads characters, a string, an array of strings, or binary data from a nominated serial port that was previously opened with SerialOpen(). Binary data can include character code 0; string data never includes character 0.

Func SerialRead(port%, &in\$|in\$[]|&in%|in%[]{,term\${, max%}});

- port% The serial port to read from as defined for SerialOpen().
- in\$ A single string or an array of strings to fill with characters. To use an array of strings you must set a terminator or all input goes to the first string in the array.
- in% A single integer (term\$ and max% are ignored) or an array of integers (term\$ and max% can be used) to read binary data. Each integer can hold one character, coded as 0 up to 255. The function returns the number of characters returned.
- term\$ If this is an empty string or omitted, all characters read are input to the string, integer array or to the first string in the string array and the number of characters read can be limited by max%. The function returns the number of characters read.

If term\$ is not empty, the contents are used to separate data items in the input stream. Only complete items are returned and the terminator is not included. For example, set the terminator to "\n" if lines end in line feed, or to "\r\n" if

input lines end with carriage return then line feed. If in\$ is a string, one item at most is returned. If in\$[] is an array, one item is returned per array element. The function returns the number of items read unless in is an integer (in% or in\$[]) when it returns the number of characters returned.

max% If present, it sets the maximum number of characters to read into each string or into the integer array. If a terminator is set, but not found after this many characters, the function breaks the input at this point as if a terminator had been found. There is a maximum limit set by the integer array size, the size of the buffers used by Spike2 to process data and by the size of the system buffers used outside Spike2. This is typically 1024 characters.

Returns The function returns the number of characters or items read or a negative error code. If there is nothing to read, it waits 1 second for characters to arrive before timing out and returning 0. Use SerialCount() to test for items to read to avoid a time out.

See also:SerialOpen(), SerialWrite(), SerialCount(), SerialClose()

SerialWrite()

This writes one or more strings or binary data to a serial port opened by SerialOpen().

Func SerialWrite(port%, out\$|out\$[]|out%|out%[]{, term\$});

- port% The serial port to write to as defined for SerialOpen().
- out \$ A single string or an array of strings to write to the output. The return value is the number of strings written. If a time-out occurs, the function returns the number of complete strings sent before the time-out.
- out% A single integer or an integer array to write as binary. One value is written per integer. The output written depends on the number of data bits set for the port; 7-bit data writes as out% band 127, 8-bit data writes as out% band 255. The return value is 1 if the transfer succeeded.
- term\$ If present, it is written to the output port after the contents of out%, out%[] or out\$ or after each string in out\$[].
- Returns As defined above or a negative error code. If the output system becomes full, the function waits for one second before timing out.

See also:SerialOpen(), SerialRead(), SerialCount(), SerialClose()

SetXXXX commands These commands create result views attached to the current time view or to the time window associated with the current result view. Apart from SetResult(), which creates a result view that is not dependent on a time view, these routines correspond with the Analysis menu New Result view commands. They do not update the display (use Draw() or DrawAll()), nor do they perform any analysis (see the Process() command).

The commands return a view handle, or a negative error code. Errors include: Bad channel number, illegal number of bins, out of memory, illegal bin size. The new view is made the current view. However, it is created invisibly and must be made visible with WindowVisible(1) before it will appear when drawn.

Most commands accept a channel list specifier, written as cSpc. It can be an integer channel number, or -1 for all, -2 for visible or -3 for selected channels. It can also be a channel specification string, such as "1..4, 6, 10" or an integer array where the first element is the channel count, the remaining elements are the channel numbers. Whatever channels are specified, only channels of the correct type for the command are used. If a channel is duplicated, the first occurrence is used. It is an error if the resulting list is empty. Channels are created in the result view in the order they appear in the specification. The first channel in the list generates result view channel 1, the second generates channel 2, and so on.

SetAverage()

This command creates a result view to hold the sum or average of sweeps of waveform data. The Process() command does the analysis. Sweeps() reports the number of sweeps accumulated. Omitted optional arguments have the value 0. The function is:

Func SetAverage(cSpc, bins% {,offset {,trig% {,flags%}}});

- CSpc A channel list specifier defining channels to average. The channels must be waveform or WaveMark channels in the current time view, or in the time view associated with the current result view. Invalid channels are removed from the list. The resulting channel list must hold at least one valid channel. All the channels must have the same sampling rate; any channels that do not match the sample rate of the first channel in the list are ignored.
- bins% The number of bins required. There must be at least one bin. The maximum is limited by available memory. The bin width is the waveform sampling interval.
- offset This sets the pre-trigger time to show in the result. If omitted, 0 is used.
- trig% The channel number to use as a trigger for each sweep. If this is omitted, or set to 0, then each call to Process () takes the start time as the trigger time.
- flags% This is the sum of flag values. Add 1 to display the mean data and not the sum. Add 4 to enable error bars. If flags% is omitted, the value 0 is used.

Returns The function returns a handle for the new view, or a negative error code.

See also:BinError(), ChanData(), DrawMode(), Process(), Sweeps()

SetEvtCrl()

This creates a result view of event correlation histograms and optional raster displays. The Process() command does the analysis. Sweeps() reports the number of triggers processed. You can take measurements from an auxiliary channel from each sweep and use this value to sort rasters, display a symbol and/or discard sweeps. See the Event correlation in the Analysis menu for details. See RasterSort() and RasterSymbol() for a description of sort values and symbol times.

- cSpc A channel list specifier of event or marker channels in the current time view, or in the time view associated with the current result view. Invalid channels are ignored. The resulting channel list must hold at least one valid channel.
- bins% The number of bins in the histogram. There must be at least 1 bin.
- binsz The width of each bin in seconds. This is converted into underlying time units.
- offset This sets the pre-trigger time to show in the histogram.
- trig% The channel number to use as a trigger for each sweep. If this is omitted, or set to 0, then each call to Process () takes the start time as the trigger time.
- flags% This is the sum of: 1 to scale the result as spikes per second, 2 to enable raster displays, 8 for backwards aCh% event search, 16 to exclude rather than include values that lie in the Mn to Mx range. If omitted, 0 is used.
- aCh% Auxiliary measurement channel. The measured value sets raster sort value 1. For waveform or RealWave channels, this is the waveform value at the trigger time. For all other channel types, it is the latency of the first event on the channel before/after the trigger and it also sets the time of symbol 1. If you omit aCh%, there is no auxiliary channel.
- Mn, Mx If present, these arguments set the range of measured values from aCh% that control if a sweep is included or excluded.

Returns The function returns a handle for the new view, or a negative error code.

For an auto-correlation, set trig% the same as chan%. In this case, we do not count the result of correlating each event with itself. If you must count self-correlations, add Sweeps() to the bin holding the time shift of zero.

See also:SetEvtCrlShift(), DrawMode(), Process(), RasterSet(), Sweeps()

SetEvtCrlShift()

This can be used when the current result view was created by SetEvtCrl(). It sets a time shift for the next Process() command to produce shuffled correlations by shifting one channel with respect to the other without changing the result view time axis.

Proc SetEvtCrlShift(shift);

shift This value can be positive or negative. It causes the data on the trig% channel to be correlated with data from channel chan% at a time shift later. The shift is set back to zero after the Process() command on this result view.

Consider two channels of events (to be correlated) and a channel marking stimuli:



If the stimulus channel holds events at some constant interval int, by shifting all events on one of the two channels by int, the only correlation between the two channels left is the correlation due to the stimulus. Thus by forming the difference between the unshifted correlation and the shifted correlation, you can obtain the correlation between a and b with the effect of the stimulus removed.

See also:SetEvtCrl(), Process()

SetINTH()	This function creates a result window to hold an interval histogram. The Process() command does the analysis. Each interval processed increments the value returned by Sweeps(), even intervals that are too short or too long to contribute to the histogram.				
	Func S	etINTH(cSpc, bins%, binsz{, minInt});			
	cSpc	A channel list specifier of event or marker channels in the current time view, or in the time view associated with the current result view. Invalid channels are ignored. The resulting channel list must hold at least one valid channel.			
	bins%	The number of bins required. There must be at least one bin. The maximum is limited by available memory.			
	binsz	The width of each bin in seconds. This is converted into underlying time units.			
	minInt	This sets the start of the first analysis bin, in seconds. Shorter intervals are counted for $Sweeps()$, but not included in the histogram. If omitted, 0 is used.			
	Returns	The function returns a handle for the new view, or a negative error code.			
	See also	:Process(),Sweeps()			
SetPhase()	This fun The Pro analysed value to Analysis for a des	action creates a result view to hold phase histograms with optional raster displays. bccess() command does the analysis. Sweeps() reports the number of cycles d. You can measure a value for each sweep from an auxiliary channel and use this b sort rasters, display a symbol and discard sweeps. See the description of the s menu Phase histogram for details. See RasterSort() and RasterSymbol() scription of sort values and symbol times.			
	<pre>Func SetPhase(cSpc, bins%{, minCyc{, maxCyc{, cycle%{, flags%</pre>				
	cSpc	A channel list specifier of event or marker channels in the current time view, or in the time view associated with the current result view. Invalid channels are ignored. The resulting channel list must hold at least one valid channel.			
	bins%	The number of bins required. There must be at least one bin. The maximum is limited by available memory. The bins are each 360/bins% wide.			
	minCyc	The minimum cycle time to use in seconds. If two consecutive cycle markers are closer than this time, the cycle is ignored, and is not counted by $Sweeps()$. If this is omitted, a value of 0.0 is used.			
	maxCyc	The maximum cycle time to use in seconds. If two consecutive cycle markers are further apart than this, the cycle is ignored, and is not counted by Sweeps(). If this is omitted, there is no limit on cycle size.			
	cycle%	The channel number to use as the cycle start and end marker. If this is omitted, or set to 0, then each call to Process() takes the start time and end time as marking a single cycle.			

- flags% This is the sum of: 2 to enable raster displays, 8 for backwards aCh% event search, 16 to exclude rather than include sweeps with measured values that lie in the Mn to Mx range. If omitted, 0 is used.
- aCh% Auxiliary measurement channel. The measured value sets raster sort value 1. For waveform or RealWave channels, this is the waveform value at the cycle start. For all other channel types, it is the latency of the first event on the channel before/after the cycle start and it also sets the symbol 1 time. If you omit aCh%, there is no auxiliary channel.
- Mn, Mx If present, these arguments set the range of measured values from aCh% that control if a sweep is included or excluded.

Returns The function returns a handle for the new view, or a negative error code.

See also:DrawMode(), SetEvtCrl(), Process(), RasterSet(), Sweeps()

SetPower()

This function creates a result view to hold a power spectrum. The Process () command does the analysis. The function is as follows:

Func SetPower(cSpc, fftsz% {,wnd%});

- cSpc A channel list specifier of waveform channels in the current time view, or in the time view associated with the current result view. Invalid channels are removed from the list. The resulting channel list must hold at least one valid channel. All the channels must have the same sampling rate; any channels that do not match the sample rate of the first channel in the list are ignored.
- fftsz% The size of the transform used in the FFT. This must be a power of 2 in the range 16 to 16384. The result view has half this number of bins. The width of each bin is the sampling rate of the channel divided by fftsz%. Each block of fftsz% data points processed increments the value for Sweeps().
- wnd% The window to use. 0 = none, 1 = Hanning, 2 = Hamming. Values 3 to 9 set Kaiser windows with -30 dB to -90 dB sideband ripple in steps of 10 dB. If this is omitted a Hanning window is applied.

Returns The function returns a handle for the new view, or a negative error code.

See also:ArrFFT(), Process(), Sweeps()

SetPSTH()

This creates a result view to hold a peri-stimulus time histogram with an optional raster display. The Process() command does the analysis. Sweeps() returns the number of triggers processed. You can measure a value for each sweep from an auxiliary channel and use this to sort rasters, display a symbol and discard sweeps. See the description of the PSTH in the Analysis menu for details. See RasterSort() and RasterSymbol() for a description of sort values and symbol times.

- cSpc A channel list specifier of event or marker channels in the current time view, or in the time view associated with the current result view. Invalid channels are ignored. The resulting channel list must hold at least one valid channel.
- bins% The number of bins required. There must be at least one bin.
- binsz The width of each bin in seconds. This is converted into underlying time units.
- offset This sets the pre-trigger time, in seconds. If omitted, 0.0 is used.
- trig% The channel number to use as a trigger for each sweep. If this is omitted, or set to 0, then each call to Process() takes the start time as the trigger time.
- flags% This is the sum of: 1 to scale the result as spikes per second, 2 to enable raster displays, 8 for backwards aCh% event search, 16 to exclude rather than include values that lie in the Mn to Mx range. If omitted, 0 is used.
- aCh% Auxiliary measurement channel. The measured value sets raster sort value 1. For waveform or RealWave channels, this is the waveform value at the trigger time. For all other channel types, the value is the latency of the first event on the channel before/after the trigger and it also sets the symbol 1 time. If you omit aCh%, there is no auxiliary channel.
- Mn, Mx If present, these arguments set the range of measured values from aCh% that control if a sweep is included or excluded.

Returns The function returns a handle for the new view, or a negative error code.

See also:DrawMode(), SetEvtCrl(), Process(), RasterAux(), Sweeps()

SetResult()	This fur with no	This function creates a result view of user-defined type, attached to no time view and with no implied Process(). Sweeps() will return 0 unless you set the sweep count.					
	Func S	Func SetResult({chans%,} bins%, binsz, offset, title\$, xU\$ {,yU\$ {,xT\$ {,yT\$ {,flags% {,tick}}}});					
	chans%	chans% This sets number of channels in the new result view. Omit for 1 channel. See the View() documentation for access to multiple channel result view data.					
	bins%	The number of bins in the view.					
	binsz	The width of each bin. Bins should have a positive non-zero width.					
	offset The x axis value at the start of the first bin.						
	title\$	The window title.					
	xU\$	The x axis units.					
	YU\$	Optional, y axis units, blank if omitted.					
	xT\$	Optional, x axis title (otherwise blank).					
	yT\$	Optional, y axis title (otherwise blank).					
	flags%	Add 2 to enable raster data, add 4 to enable error bars. The default value is 0.					
	tick	This is required if you enable raster data. Result views store raster data as 32-bit integers. tick is the time resolution for this data. When working with a time view, set tick to BinSize(), the time resolution of the time view. The maximum time of a raster event is 2147483647*tick seconds.					
	Datuma	The function returns a handle for the new view, or a negative error and					

Returns The function returns a handle for the new view, or a negative error code.

See also:BinError(), RasterGet(), RasterSet(), Sweeps(), View()

SetWaveCrl()

This creates a result view to hold the waveform correlation between waveform channels sampled at the same rate. The Process() command does the analysis. Sweeps() returns the number of data points used on the reference channel.

Func SetWaveCrl(cSpc, ref%, bins% {,offset});

- cSpc A channel list specifier of waveform channels in the current time view, or in the time view associated with the current result view. Invalid channels are removed from the list. There must be at least one valid channel. Channels that do not match the sample rate of the first channel in the list are ignored.
- ref? A reference waveform channel that is correlated with all the waveform channels.
- bins% The number of bins in the correlation. The bin width is the channel sampling interval. This calculation is slow compared to the other SetXXXX commands. The time taken is proportional to the length of the area processed times bins%.
- offset The time to show before the zero time shift, in seconds. If omitted, 0 is used.
- Returns The function returns a handle for the new view, or a negative error code.

See also:SetWaveCrlDC(), Process(), Sweeps()

SetWaveCrIDC()

To use this function the current view must be a waveform correlation result view. It sets whether the DC levels of the two signals is removed before the correlation or not.

Func SetWaveCrlDC({useDC%})

useDC% If present, a zero value removes the DC, a non-zero value (default) includes it.

Returns The useDC% value before the call. It returns a negative error if this view is a result view but not a waveform correlation.

If you change the setting the view is drawn at the next opportunity.

See also:SetWaveCrl(), Process()

Sin()

This calculates the sine of an angle in radians, or converts an array of angles into sines.

Func Sin(x|x[]|x[][]);

- The angle, expressed in radians, or a real array of angles. The best accuracy of the result is obtained when the angle is in the range -2π to 2π .
- Returns When the argument is an array, the function replaces the array with the sines and returns 0 or a negative error code. When the argument is not an array the function returns the sine of the angle.

See also:ATan(), Cos(), Cosh(), Sinh(), Tan(), Tanh()

Sinh()

This calculates the hyperbolic sine of a value or of an array of values.

Func Sinh(x | x[] | x[]];

- x The value, or an array of real values.
- Returns When the argument is an array, the function replaces the array elements with their hyperbolic sines and returns 0. When the argument is not an array the function returns the hyperbolic sine of the argument.

See also:ATan(), Cos(), Cosh(), Sinh(), Tan(), Tanh()

Sound()

This has two variants. The first plays a short "beep" in Windows 9x and a tone of set pitch and duration in later systems using the system speaker. The second plays a .WAV file or system sound if your system has multimedia support.

Func Sound(freq%, dur{, midi%});Tone output (not using sound card)Func Sound(name\${, flags%};Multimedia sound card output

- freq% If midi% is 0 or omitted this holds the sound frequency in Hz. If midi% is nonzero this is a MIDI value in the range 1-127. A MIDI value of 60 is middle C, 61 is C# and so on. Add or subtract 12 to change the note by one octave.
- dur The sound duration, in seconds. The script stops during output.
- midi% If present and non-zero, freq% is interpreted as a MIDI value.
- name\$ The name of .wav file or of a system sound. You can supply the full path to the file or just a file name and the system will search for the file in the current directory, the Windows directory, the Windows system directory, directories listed in the PATH environmental variable and the list of directories mapped in a network. If no file extension is given, .wav is assumed. The file must be short enough to fit in available physical memory, so this function is suitable for files of a few seconds duration only.

A blank name halts sound output. If name\$ is any of the following (case is important), a standard system sound plays:

"s*"	Asterisk	"ss"	Start	"SE"	Exit	"SH"	Hand
"sw"	Welcome	"s?"	Query	"SD"	Default	"s!"	Exclamation

flags% This optional argument controls how the data is played. It is the sum of:

0x0001 1 Play asynchronously (start output and return). Without this flag, Spike2 does nothing until replay ends.

0x0002	2	Silence when sound not found. Normally Sound() plays the
		system default sound if the nominated sound cannot be found.
0x0008	8	Loop sound until stopped by another $Sound()$ command. You
		must also supply the asynchronous flag if you use loop mode.
0x0010	16	Don't stop a playing sound. Normally, unless the "No wait" flag
		is set, each command cancels any playing sound.
0x2000	8192	No wait if sound is already playing. Sound ("", 0x2010) returns
		1 if a previous asynchronous sound is finished and 0 if not.
10 1	• •	

If you don't supply this argument, the flag value is set to 0x2000.

Returns The tone output returns 0 or a negative error code. The multimedia output returns non-zero if the function succeeded and zero if it failed.

See also:Speak()

Speak()

If your system supports text to speech, this command allows you to convert a text string into speech. Currently we provide no facilities to setup voices or to route the sound output; you must do this from the Speech applet in the control panel.

Func	<pre>Speak(text\${,</pre>	opt%});	Convert text to speech
Func	Speak({what%{,	val}});	Speech status and control

text\$ A string holding the text to output, for example "Sampling has started".

- opt% This optional argument (default value 1) controls the text conversion and output method. It is the sum of the following flags:
 - 1 Speak asynchronously. Without this flag the command waits until speech output is over before returning.
 - 2 Cancel any pending speech output.
 - 4 Speak punctuation marks in the text.
 - 8 Process embedded SAPI XML; the www.microsoft.com/speech web site has more information on this advanced topic. For example: Speak("Emphasis on <EMPH>this</EMPH> word", 8);
 - 16 Reset to the standard voice settings before speaking.
- what% An optional variable, taken as 0 if it is omitted:
 - 0 Returns 1 if speech is playing and 0 if it is not.
 - 1 Wait for up to val seconds (default value 3.0) for output to end. The return value is 0 if playing is finished, 1 if it continues after val seconds.
 - 2 Returns the current speech speed in the range -10 to 10; 0 is the standard speed. If val is present, it sets the new speed.
 - 3 Returns the current speech volume in the range 0 to 100, 100 is the standard volume. If val is present, it sets the new volume.
- val An optional argument used when what% is greater than 0.
- Returns If there is no speech support available, or a system error occurs, the command returns -1. Otherwise the first command variant returns 0 if all is well and the second variant returns the values listed for what%.

To use TTS (text to speech), you need a suitable sound card and the Microsoft SAPI software support. Windows XP has this software included with the operating system. You can get text to speech support as a download for other versions of Windows. In August 2004, the speech support was available as SpeechSDK51MSM.exe from the web page www.microsoft.com/speech/download/sdk51/ but this location may change.

See also:Sound()

Sqrt()

Forms the square root of a real number or an array of real numbers. Negative numbers halt the script with an error when \times is not an array. With an array, negative numbers are set to 0 and an error is returned.

Func Sqrt(x|x[]|x[][]);

x A real number or a real array to replace with an array of square roots.

Returns With an array, this returns 0 if all was well, or a negative error code. With an expression, it returns the square root of the expression.

See also:Abs(), ATan(), Cos(), Exp(), Frac(), Ln(), Log(), Max(), Min(), Pow(), Rand(), Round(), Sin(), Tan(), Trunc()

SMOpen()

This function gets the handle and optionally opens the spike monitor window. The current view must be a time view or a spike monitor window. Close the spike monitor window with FileClose(). Use the SMControl() command to control the window. You can use the Window() and WindowVisible() commands to size, show and hide the window.

Func SMOpen({type%{, mode%});

type% This argument is assumed to be zero if omitted. Allowed values are:

- 0 The return value is the window handle of the spike monitor window associated with the current view, or 0 if there is no open window.
- 1 Open the spike monitor window. Return its handle or 0 if we failed.
- mode% If this is zero or omitted, the dialog is created invisibly. Set 1 to make it visible. This is ignored unless type% is 1.

Returns The return value is the view handle of the spike monitor window or 0.

See also:FileClose(), SMControl(), ViewLink(), Window(), WindowVisible()

SMControl()

This function sets and reads the toolbar and edit field states of the current spike monitor window. See the documentation of the spike monitor window for a full description.

Func SMControl(item%{, new});

item% This identifies the item:

0	Draw mode	0=3D, 1=2D, 2=2D separated
1	Fade to background	0=no fade, 1=fade
2	Colour last spike only	0=colour all, 1=colour last, the rest are grey
3	Use thick lines	0=use thin lines, 1=use thick lines (slow)
4	Timed (smooth) update	0=update on new data, 1=update on timer
5	At end mode	0=use cursor 0, 1=at end of file
6	Lock y axes	0=track data size, 1=fix at current size
7	Show duplicate channels	0=Show originals only, 1=show duplicates
8	Maximum spikes	In the range 1 to 40
9	Time range	In the range 0.003 to MaxTime()
10	Display rectangles	0=no, 1=show front and back rectangles
11	Front rectangle x offset	In the range 0 to 0.9
12	Front rectangle scale	In the range 0.1 to 1.0
13	Back rectangle x offset	In the range 0 to 0.9
14	Back rectangle scale	In the range 0.1 to 1.0
16		

new If present, this sets the state of the item. See the table for acceptable values.

Returns The item state at the time of the call or -1 if the item does not exist.

See also:SSChan(), SSClassify(), SSOpen(), SSParam(), SSRun()

SSButton()	This fun window.	This function sets and reads the toolbar and checkbox states of the current spike shape window. See the documentation of the spike shape dialogs for a description of the items.		
	Func SS	unc SSButton({item%{, new%}});		
	item%	If or 0 1 2 3 4 5 6 7	nitted, all items are reset to 0. Otherwise Overdraw spikes Show template boundaries Show non-matching spikes Play sound for each spike Scroll time view to track cursor 0 At end mode (online) Circular replay Make templates	e this identifies the item: 0=no overdraw, 1=overdraw 0=hide, 1=show 0=hide, 1=show 0=no sound, 1=sound 0=no track, 1=track 0=not at end, 1=at end 0=stop at end, 1=circular 0=no, 1=Yes
	new%	8 9 10 11 12 If pr butt	Build by code Number of horizontal cursors Set which marker code to use Overlay drawing of traces Set size of displayed templates resent, this sets the state of the item. Us on or check the box) and 0 to unselect it	0=by shape, 1=by code 0=2, 1=4 0-3 for layers 0-3 0=no overlay, 1=overlay 0=small, 1=medium, 2=large se 1 to select the feature (depress the t. Item 10 supports values 0-3.

Returns The item state at the time of the call or -1 if the item does not exist or no item%.

See also:SSChan(), SSClassify(), SSOpen(), SSParam(), SSRun()

SSChan()

This function sets and gets the current channel in the current spike shape dialog.

Func SSChan({ch%});

ch% If omitted or 0, the return value is the current channel. If greater than 0, ch% sets the channel and the return value is the ch% (or the original channel if channel ch% is not suitable). If the channel changes, the old channel settings are saved.

If ch% is -1, the return value is the number of interleaved traces in the dialog. In the New Stereotrode or New Tetrode dialogs, this is 2 or 4. In the Edit WaveMark dialog it is the number of interleaved traces for the current channel.

If ch% is -2, this forces the channel configuration to be saved; the return value is 0. Use this before FileClose() as that does not save the current settings.

In a New Stereotrode or New Tetrode dialog ch% is a zero based channel index into the list of channels to combine and the return value is the channel number.

Returns The return value depends on the ch% argument.

See also:SSButton(), SSClassify(), SSOpen(), SSParam(), SSRun()

SSClassify()

This function is equivalent to the New Channel button in the New WaveMark dialog and to the Reclassify button in the Edit WaveMark dialog. The arguments from new% onward are ignored in the Edit WaveMark dialog.

Func SSClassify({speed%{, new%{, type%{, noQu%}}});

speed% If this is omitted, the function returns 1 if the dialog is currently reclassifying or creating a new channel and 0 if it is not. If present, use the values:

- -1 To cancel any ongoing reclassification or writing to a new channel.
- 0 Start reclassifying or writing a new channel with display updates. This takes place in idle time, so you must provide some with Yield() or by using Toolbar(), Interact() or DlgShow().

- 1 Reclassify or write a new channel in fast mode, with no display updates. The operation will complete before control returns to the script.
- new? The channel number to create in the New WaveMark dialog.
- type% The channel type to create. You can use types: 2 (Event-), 3(Event+), 4(Level), 5(Marker) or 6 (WaveMark) only. If omitted, type 6 (WaveMark) is used.
- noQu[®] If the channel set by new[®] is already in use and this argument is omitted or 0, the user is asked if they want to overwrite it. Set non-zero for no query.
- Returns If speed% is omitted, the return value is 1 if reclassification or writing a new channel with display updates is in progress, otherwise 0.

See also:SSButton(), SSChan(), SSOpen(), SSParam(), SSRun()

SSOpen()

This function opens and returns information about spike shape dialogs. The current view must be a time view or a spike shape dialog. Close this dialog with FileClose().

Func SSOpen({type%{, mode%});

- type% This argument is assumed to be zero if omitted. Allowed values are:
 - -1 The function returns the type of any open spike shape dialog associated with the current time view or 0 if there is no open dialog. Returned values are 1=Edit WaveMark, 2=New WaveMark, 3=New n-trode.
 - 0 The return value is the window handle of any open spike shape dialog associated with the current view, or 0 if there is no open dialog.
 - 1 Open the Edit WaveMark dialog and return its handle or 0.
 - 2 Open the New WaveMark dialog and return its handle or 0.
 - 3 Open the New n-trode dialog with selected channels. Return its handle or 0.
- mode% If this is zero or omitted, the dialog is created invisibly. Set 1 to make it visible. This is ignored unless type% is greater than 0.

Returns The return value depends on the type% argument.

See also:SSButton(), SSChan(), SSClassify(), SSParam(), SSRun(), ViewLink()

SSParam()

This function sets the template parameters for the current spike shape dialog and reads back the state of individual parameters. This is equivalent to the template parameters dialog; see that dialog for a full description of the arguments. The command has two variants. In the first you supply all the arguments (use the value -1 to leave an argument unchanged). In the second variant you can get or set the value of a single item.

Func SSParam(new%,wide,rare%,AMax,PMin,flg%,mode%,nAT%,int%,tc%);
Func SSParam(item%{, value});

- new% If new% is greater than 0, this sets number of spikes for a new template.
- wide The width of a new template as a percentage of the template amplitude.
- rare% Templates with spikes rarer than 1 in rate% are discounted.
- AMax Maximum amplitude change for tracking or 0 for no amplitude tracking.
- PMin Minimum percentage of point in the template for a possible match.
- flg% This is the sum of: 1 = use minimum percentage only when building templates, 2 = Remove DC value from spike data before matching.
- mode% The mode for adding spikes to the template: 0=All, 1=Autofix, 2=track.
- nAT% The number of spikes for Autofix and track modes.
- int% The interpolation method to use: 0=linear, 1=parabolic, 2=cubic spline.

tc%	The high-pass time constant to apply	y to the c	lata expressed	as 2 to the power tc%
	data samples with tc% in the range	1 to 30. U	Jse 31 for no l	nigh pass filter.

- item% Used in the second command variant. Set -2 to copy the current settings to all channels. Set -1 to reset the values to a standard state. Set 0 to close the parameters dialog if it is open. Otherwise set 1 to 10 to select one of the arguments new% (1) to tc% (10) and the return value will be the value of that item at the time of the call.
- value If present, this sets the new value of the parameter when item% is greater than 0.
- Returns The first command variant returns 0. The second variant, with item% greater than 0, returns the item value at the time of the call and item% less than zero returns 0. For item%=0, the return value is 1 if the parameter dialog was open.

See also:SSButton(), SSChan(), SSClassify(), SSOpen(), SSRun()

SSRun()

This function gets and sets the run state of the current view, which must be a spike shape window. Use Cursor(0) to get the current run position.

Func SSRun({run%});

run% If omitted, no change is made, otherwise negative values are backwards and positive are forwards. 0 means stop, 1=step, 2=run as fast as possible, 3= run in real time, 4-= run at 91 Hz, 5=31 Hz and so on up to 15=1 Hz.

Returns The run state at the time of the call.

See also:SSButton(), SSChan(), SSClassify(), SSOpen(), SSParam()

SSTempDelete()

This function deletes one or more templates for the current channel of the current view, which must be a spike shape window. If this command is used with no arguments, or with both arguments set to -1, the entire template system is cleared, equivalent to clicking the clear all templates button in the dialog.

Func SSTempDelete({n%{, code%}});

- n% The zero-based index of a template or -1 for all templates. Templates that match this index and the code% argument are deleted. If this is omitted, -1 is used.
- code% If this is omitted or set to -1, all templates codes are deleted. Otherwise only templates selected by n% that match code% are deleted.

Returns The number of deleted templates or -1 if n% is not a template index.

```
See also:SSButton(), SSChan(), SSClassify(), SSOpen(), SSParam(), SSRun(),
SSTempGet(), SSTempInfo(), SSTempSet(), SSTempSizeGet(),
SSTempSizeSet()
```

SSTempGet()

This function gets information about the current set of templates for the current channel of the current view, which must be a spike shape window. It can also return the current raw data displayed in the template window.

Func SSTempGet(n%{, temp[]{[]}{, what%{, &count%}});

- n% This is the zero-based index of the template to return information from or to -1 to get the currently displayed waveform. Use SSTempSizeGet() to get the number of displayed points and the number of points in the template.
- temp This optional array is filled with template data. If there is insufficient data to fill it, unused entries are unchanged. An integer or real array can be used. If the

template has multiple traces, use temp[points%][traces%] to get real data and temp%[points%][traces%] to get integer data. See ChanScale() for an explanation of the use of integer and real data for waveform values.

- what% Omit this or set it to 0 to return the mean template waveform. It is ignored if n% is -1. 1 = return the template upper limit, 2 = return the lower limit, 3 = return the template width. The width is half the distance between the upper and lower template limits.
- count% This optional integer variable is returned set to the number of events that were accumulated into this template. It is ignored if n% is -1.
- Returns If n% is positive, the return value is the template code or -1 if the template does not exist. If n% is -1, in Edit WaveMark dialogs, the return value is the currently selected sort code of the spike or -1 if you are on-line and the current data did not trigger (background data). In Create WaveMark dialogs the return value is 1 for triggered data and 0 for background data.
- See also:ChanScale(), SSButton(), SSChan(), SSClassify(), SSOpen(), SSParam(), SSRun(), SSTempDelete(), SSTempInfo(), SSTempSet(), SSTempSizeGet(), SSTempSizeSet()

SSTempInfo()

This function gets and sets template information for the current channel of the current view, which must be a spike shape window.

Func SSTempInfo({n%, item%{, val%{, noLim%}});

- If this is omitted, the return value is the number of templates. If present, it is the zero-based index of a template or -1 for all templates (if appropriate). When returning a value, if -1 is used, the value for template 0 is returned.
- item% This selects the information to get or set. Items 0-2 operate on all templates. Items 3-6 operate on visible templates only (n% in the range 0 to 19).
 - 0 Get the index of the first template, starting at n%, that has a code of val%.
 - 1 Get the code for template n%. It sets the code if val% is greater than 0.
 - 2 Get the count of spikes for template n%.
 - 3 Get the event counter. If val% is positive, it sets the event count.
 - 4 Get the lock state of template n%. val%=0 to clear the lock and >0 to set it.
 - 5 Change the template width by val% pixels. Set val% to 1 and -1 to mimic the increase and decrease width buttons. noLim% can be used.
 - 6 Change the template width by val% channel units. noLim% can be used.
- val% Use this argument when setting values and to locate codes when item% is 0.
- noLim% Set non-zero with item codes 5 and 6 to remove the template width limits that are usually enforced to stop the template becoming too wide or too narrow. Each point of a template has a mean value and a width. It also has a minimum width, which is initialised to the original width when the template was created. The width is allowed to increase to up to 4 times the minimum width, but not to become less that the original. If you allow the width to override these limits, the minimum width is adjusted, as required.
- Returns The value selected by item% at the time of the call. The return value is 0 for items 5 and 6. The return value is -1 if a template is not found.
- See also:SSButton(), SSChan(), SSClassify(), SSOpen(), SSParam(), SSRun(), SSTempDelete(), SSTempGet(), SSTempSizeGet(), SSTempSizeSet()

SSTempSet()

This function creates a template or adds the current spike in the window to a nominated template in the current channel of the current view, which must be a spike shape view. The first command variant creates a new template from user-supplied data.

The second variant adds the currently displayed data in the window to a nominated template or creates a new template from it. In the Edit WaveMark dialog this also sets the classification code of the displayed data. If this code is not in the current Marker Filter for the source channel, the displayed data changes to the next available spike. It can also renumber and sort the displayed templates and to update online templates in the 1401 to match the templates in the dialog.

```
Func SSTempSet(temp[]{[]}, code%, wide{, count%});
Func SSTempSet({n%{, code%}});
```

- temp Used when creating a new template from arrays of data. This is a one or two dimensional array holding the template shape. An integer or real array can be used. If the template has multiple traces, use temp[points%][traces%] for real data and temp%[points%][traces%] for integer data. See ChanScale() for an explanation of the use of integer and real data for waveform values. If the array is too short or has too few traces, missing values are set to 0.
- code% This is used when creating a new template. Set it to 0 (or omit it in the second command variant) to use the lowest unused code or set it to the code to use for the new template in the range 1-255.
- count% This optional argument has a default value of 1, and sets the number of spikes that contributed to the template.
- wide A number or an array holding the template width, which is half the distance between the upper and lower template limits. If you supply an array, the shape of this array should match the temp array.

If you supply a number, all points are given the same width. This is in user units if temp is a real array and is in channel units if temp is an integer array.

- n% Optional. If omitted, the currently displayed waveform is added to the best-fit template. Set it to the zero-based index of an existing template to add the current waveform to the template. Set to -1 to add a new template, in which case code% determines the code of the new template. Set -2 to renumber the existing templates, in which case code% sets the lowest numbered code used. Set -3 to sort the templates into ascending code order. Set -4 to update online templates during sampling.
- Returns The zero-based index of the template that was added or modified or -1 if there was a problem (no matching template, unknown template index or no current waveform). The second variant with n% less than -1 returns 0.
- See also:SSTempDelete(), SSTempGet(), SSTempInfo(), SSTempSizeGet(), SSTempSizeSet()

SSTempSizeGet()

This function gets the template size and the displayed data area for the current channel. The current view must be a spike shape window. All arguments are integer variables that are set to the template and display sizes.

Func SSTempSizeGet({&start%{, &show%{, &pre%}}});

start% If present, set to the point offset to the template start in the raw data display area.

show% If present, set to the number of data points in the raw data display area.

pre% If present, set to the number of pre-trigger points in the raw data display area.

Returns the number of data points in the template.

See also:SSTempDelete(), SSTempGet(), SSTempInfo(), SSTempSet(), SSTempSizeSet()

SSTempSizeSet()	This sets the template size and the displayed data area. The template must lie within the displayed data area with at least 2 data points before and after the template. If you set a			
	template area that exceeds these limits, the command attempts to increase the display area (if this is possible). Use YRange() and Optimise() to set the displayed y range,			
	<pre>Proc SSTempSize(start%{, n%{, show%{, pre%}}});</pre>			
	start% This is the offset to the beginning of the template from the start of the displayed area. It must be at least 2. Set -1 to leave the start point unchanged.			
	n% This sets the number of data points in the template. To keep the current number of points, either omit this argument or set it to -1.			
	show% The requested number of points to display in the main data window; it has no effect in the Edit WaveMark dialog. Omit show% or use -1 for no change. If odd, show% is reduced by 1. It is also adjusted if it is too large or too small.			
	pre%This sets the number of pre-trigger points to display in the main data window in the range 0 to show%-1. This has no effect in the Edit WaveMark dialog. Omit pre% or set -1 for no change. Out of range values are set to the appropriate limit.			
	<pre>See also:SSClassify(), SSTempDelete(), SSTempGet(), SSTempInfo(), SSTempSet(), SSTempSizeGet(), YRange()</pre>			
Str\$()	This converts a number to a string.			
	Func Stra (x {, wratha {, Srgra}});			
	x A number to be converted.			
	sigf% Optional number of significant figures (default is 6).			
	Returns A string holding a representation of the number.			
	<pre>See also:Asc(), Chr\$(), DelStr\$(), InStr(), LCase\$(), Left\$(), Len(), Mid\$(), Print\$(), Print(), Right\$(), UCase\$(), Val()</pre>			
Sweeps()	This function returns the number of items accumulated in the result view.			
	<pre>Func Sweeps({set%});</pre>			
	set% If present, this sets the number of sweeps held by the view. You might use it to sum waveform averages in a result view so that the sweep count is correct.			
	Returns The value returned depends on the type of the result view (see the SetXXXX commands). It is an error to use this in any type of view other than a result view.			
	<pre>See also:SetAverage(), SetEvtCrl(), SetINTH(), SetPhase(), SetPower(), SetPSTH(), SetWaveCrl()</pre>			

System\$()

This returns the operating system name and accesses Spike2 environment variables. The environment holds a list of strings of the form "name=value". You can get or set the value associated with name. You can also read all the strings into a string array.

Func System\$({var\$ {,value\$}});
Func System\$(list\$[] {,&n%});

var\$ If present, this is the name of an environment variable (case insensitive).

value\$ If present, the new value. An empty string deletes the environment variable.

- list\$ An array of strings to fill with environment strings of the form "name=value".
- n% An optional integer that is returned holding the number of elements copied.
- Returns With no arguments, it returns: "Windows SS build n" where SS is the operating system and n is the build number. Otherwise it returns the value of the environment variable identified by var\$ or an empty string.

Each process has its own environment. A process started with ProgRun () inherits a copy of the Spike2 environment, so you can pass it information. However, you cannot see environment changes made by the new process. Here are some examples of use:

```
var list$[200], value$, n%, i%;
PrintLog("%s\n", System$()); 'Print OS name
System$("fred","good"); 'Assign the value "good" to fred
PrintLog("%s\n", System$("fred")); 'get value of fred
System$("fred",""); 'Delete fred from the environment
System$(list$[], n%); 'Print all environment strings
for i%:=0 to n%-1 do PrintLog("%s\n",list$[i%]) next;
```

See also:ProgRun(), System()

System()	This function returns the operating system version as a number and gets information about desktop screens. Use the $App()$ command to get the version number of Spike2.			
	Func System({get%{, scr%{, sz%[]}}});			
	If omitted or 0, the return value is the operating system revision times 100: 351=NT 3.51, 400=95 and NT 4, 410=98, 490=Me, 500=NT 2000, 501=XP. If 1, the function returns information about installed desktop monitors.			
	Scr% Set 0 to return the number of desktop monitors; sz%[] gets the pixel co- ordinates of the desktop. Set to n (>0) to get the pixel co-ordinates of screen n; returns 1 for the primary monitor, 0 if not and -1 if it does not exist.			
	sz%[] Optional array of at least 4 elements to return pixel positions. Elements 0 and hold the top left x and y, 2 and 3 hold bottom right x and y.			
	<pre>See also:App(), System\$(), Window(), WindowVisible()</pre>			
Tan()	This calculates the tangent of an angle in radians or converts an array of angles into tangents. Tangents of odd multiples of $\pi/2$ are infinite, so cause computational overflow. There are 2π radians in 360°. π is approximately3.14159265359 (4.0*ATan(1)).			
	Func Tan(x x[] x[][]);			
	The angle, expressed in radians, or a real array of angles. The best accuracy of the result is obtained when the angle is in the range -2π to 2π .			
	Returns For an array, it returns a negative error code (for overflow) or 0. When the argument is not an array the function returns the tangent of the angle.			
	See also:ATan(), Cos(), Ln(), Log(), Pow(), Sin(), Sqrt()			
	This calculates the human aligner and a solution on an array of solution			
Tanh()	This calculates the hyperbolic tangent of a value of an array of values.			
	Func Tanh(x x[] x[][]);			
	x I ne value or an array of real values.			
	Returns For an array, it returns 0. Otherwise it returns the hyperbolic tangent of x .			
	See also: ATan(), Cos(), Cosh(), Ln(), Log(), Pow(), Sin(), Sinh(), Sgrt()			

Time\$()	This function returns the current system time of day as a string. If no arguments are supplied, the returned string shows hours, minutes and seconds in a format determined by the operating system settings. To obtain the time as numbers, use TimeDate(). To obtain relative time and fractions of a second, use Seconds().							
	<pre>Func Time\$({tBase%, {show%, {amPm%, {sep\$}}});</pre>							
	tBase	Specifies the time base to show the time in. You can choose between 24 hour or 12 hour clock mode. If this argument is omitted, 0 is used.						
		 Operating system settings 2 12 hour format 24 hour format 						
	show%	Specifies the time fields to show. Add the values of the required options together and use that as the argument. If this argument is omitted or a value of 0 is used, 7 (1+2+4) is used for 24 hour format and 15 (1+2+4+8) for 12 hour format.						
		1Show hours4Show seconds2Show minutes8Remove leading zeros from hours						
	amPm%	This sets the position of the "AM" or "PM" string in 12 hour format and has no effect in 24 hour format. If omitted, a value of zero is used. The string that gets printed ("AM" or "PM") is specified by the operating system.						
		 Operating system settings 2 Show to the left of the time Show to the right of the time 3 Hide the "AM" or "PM" string 						
	sep\$	This string appears between adjacent time fields. If sep = ":" then the time will appear as 12:04:45. If an empty string is entered or sep is omitted, the operating system settings are used.						
	See also	<pre>D:Date\$(), FileTime\$(), Seconds(), TimeDate()</pre>						

TimeDate()

This returns the time and date in seconds, minutes, hours, days, months, and years plus the day of the week. You can use separate variables for each field or an integer array. To get the data or time as a string, use Date\$() or Time\$(). To measure relative times, or times to a fraction of a second, see the Seconds() command. To get the current sampling time, see MaxTime().

```
Proc TimeDate(&s%, {&m%, {&h%, {&d%, {&mon%, {&y%, {&wDay%}}}}});
Proc TimeDate(now%[])
```

- s% If s% is the only argument, it is set to the number of seconds since midnight. Otherwise it is set to the number of seconds since the start of the current minute.
- m%If this is the last argument, it is set to the number of minutes since midnight.Otherwise it is set to the number of full minutes since the start of the hour.
- h% If present, the number of hours since Midnight is returned in this variable.
- d% If present, the day of the month is returned as an integer in the range 1 to 31.
- mon % If present, the month number is returned as an integer in the range 1 to 12.
- y[%] If present, the year number is returned here. It will be an integer such as 2002.
- wDay% If present, the day of the week will be returned here as 0=Monday to 6=Sunday.
- now%[] If an array is the first and only argument, the first seven elements are filled with time and date data. The array can be less than seven elements long. Element 0 is set to the seconds, 1 to the minutes, 2 to the hours, and so on.

See also:Date\$(), MaxTime(), Seconds(), Time\$()

The toolbar The toolbar is at the top of the screen, below the menu. The bar has a message area and can hold buttons that are used in the Interact() and Toolbar() commands. When you start a script, the toolbar is invisible and contains no buttons. When a script stops running, the toolbar becomes invisible (if it was visible).

You can define up to 40	Spike2 - DEMO.SMR	
visible buttons, numbered	<u>File E</u> dit <u>V</u> iew <u>A</u> nalysis <u>C</u> urso	r <u>S</u> ample Script <u>W</u> indow <u>H</u> elp
from 1, in your toolbar.	Message area	<u>Z</u> ero <u>F</u> ilter data Cancel Ok
There is an invisible button		

0, which sets a function that is called while the toolbar waits for a button to be pressed.

Buttons can be linked to the keyboard. However, any keys linked to buttons belong to the toolbar when it is active and waiting for a button press. If you link the A key to a button, each time you press A, the button is pressed, even if the text caret is in a text window.

See also:Interact(), Toolbar(), ToolbarClear(), ToolbarEnable(), ToolbarSet(), ToolbarText(), ToolbarVisible()

Toolbar()

This function displays the toolbar and waits for the user to click a button or press a linked key. If button 0 is defined with an associated function, that function is called repeatedly while no button is pressed. If no buttons are defined or enabled, or if all buttons become undefined or disabled, the toolbar state is illegal and an error is returned. If the toolbar was not visible, it becomes visible when this command is given.

When the Toolbar is active, Spike2 ignores the "escape" key (command-dot on the Macintosh, Esc on the PC) unless an "escape button" has been set by ToolbarSet().

Func Toolbar(text\$, allow% {,help%|help\$});

- text\$ A message to display in the message area of the toolbar. The message area competes with the button area. With many buttons, the text may not be visible.
- allow% A code that defines what the user can do (apart from pressing toolbar buttons). See Interact() for the allowed values.
- help This is either the number of a help item (CED internal use) or it is a help context string. This is used to set the help information that is presented when the user presses the F1 key. Set 0 to accept the default help. Set a string as displayed in the Help Index to select a help topic, for example "Cursors: Adding".
- Returns The function returns the number of the button that was pressed to leave the toolbar, or a negative code returned by an associated function.

The buttons are displayed in order of their item number. Undefined items leave a gap between the buttons. This effect can be used to group related buttons together.

```
See also:Interact(),ToolbarClear(),ToolbarEnable(),ToolbarSet(),
ToolbarText(),ToolbarVisible()
```

ToolbarClear()

This function removes some or all of the toolbar buttons. If you delete all buttons, the Toolbar() function inserts an OK button, so you can get out of the Toolbar() function. Use ToolbarText("") to clear the toolbar message.

Proc ToolbarClear({item%});

item% If present, this is the button to clear. Buttons are numbered from 0. If omitted, all buttons are cleared. If the toolbar is visible, changes are shown immediately.

See also:Interact(), Toolbar(), ToolbarSet(), ToolbarText()
This function enables and disables toolbar buttons, and reports on the state of a button. Enabling an undefined button has no effect. If you disable all the buttons and then use the Toolbar() function, or if you disable all the buttons in a function linked to the toolbar, and there is no idle function set, a single OK button is displayed.
<pre>Func ToolbarEnable(item% {,state%});</pre>
<pre>item% The number of the button or -1 for all buttons. You must enable and disable button 0 with ToolbarSet() and ToolbarClear().</pre>
state% If present this sets the button state. 0 disables a button, 1 enables it.
Returns The function returns the state of the button prior to the call, as 0 for disabled and 1 for enabled. If all buttons were selected the function returns 0. If an undefined button, or button 0 is selected, the function returns -1.
<pre>See also:Interact(),Toolbar(),ToolbarClear(),ToolbarSet(), ToolbarText(),ToolbarVisible()</pre>

ToolbarSet()

This function adds a button to the toolbar and optionally associates a function with it. When a button is added, it is added in the enabled state.

Proc ToolbarSet(item%, label\$ {,func ff%()});

item% The button number in the range 1 to 40 to add or replace or 0 to set or clear a function that is called repeatedly while the toolbar waits for a button press.

You can set an "escape" key as described in Toolbar(), by negating item%. For example ToolbarSet(-2, "Quit"); sets button 2 as the escape key.

label\$ The button label plus optional key code and tooltip as "Label|code|tip". Labels compete for space with each other; use tooltips for lengthy explanations. The label is ignored for button 0. Tooltips can be up to 79 characters long. To use a tooltip with no code use "Label||A tooltip with no code field".

To link a key to a button, place & before a character in the label or add a vertical bar and a key code in hexadecimal (e.g. 0x30), octal (e.g. 060) or decimal (e.g. 48) to the end of the label. Characters set by & are case insensitive. For example "a&Maze" generates the label **a**<u>M</u>aze and responds to m or M; the label "F1:G0|0x70" generates the label F1:G0 and responds to the F1 key. Useful key codes include (nk = numeric keypad):

0x08 Backspace	0x09 Tab	0x0d Enter	0x1b Escape
0x20 Spacebar	0x21 Page up	0x22 Page down	0x23 End
0x24 Home	0x25 Left arrow	0x26 Up arrow	0x27 Right arrow
0x28 Down arrow	0x2e Del	0x30-0x39 0-9	0x41-0x5a A-Z
0x60-0x69 nk 0-9	0x6a nk *	0x6b nk +	0x6c nk separator
0x6d nk -	0x62 nk .	0x6f nk /	0x70-0x87 F1-F24

Use of other keys codes or use of & before characters other than a-z, A-Z or 0-9 may cause unpredictable and undesirable effects.

Beware: When the toolbar is active, it owns all keys linked to it. If A is linked, you cannot type a or A into a text window with the toolbar active.

ff%() This is the name of a function with no arguments. The name with no brackets is given, for example ToolbarSet(1, "Go", DoIt%); where Func DoIt%() is defined somewhere in the script. When the Toolbar() function is used and the user clicks on the button, the linked function runs. If the item% 0 function is set, that function runs while no button is pressed. The function return value controls the action of Toolbar() after a button is pressed.

If it returns 0, the Toolbar() function returns to the caller, passing back the button number. If it returns a negative number, the Toolbar() call returns the

negative number. If it returns a number greater than 0, the Toolbar() function does not return, but waits for the next button. An item 0 function must return a value greater than 0, otherwise Toolbar() will return immediately.

If this argument is omitted, there is no function linked to the button. When the user clicks on the button, the Toolbar() function returns the button number.

See also:Asc(), Interact(), Toolbar(), ToolbarClear(), ToolbarEnable(), ToolbarText(), ToolbarVisible()

ToolbarText()

ToolbarVisible()

This replaces any message in the toolbar, and makes the toolbar visible if it is invisible. It can be used to give a progress report on the state of a script that takes a while to run.

Proc ToolbarText(msg\$);

msg\$ A string to be displayed in the message area of the toolbar.

```
See also:Interact(), Toolbar(), ToolbarClear(), ToolbarEnable(),
ToolbarSet(), ToolbarVisible()
```

This function reports on the visibility of the toolbar, and can also show and hide it. You cannot hide the toolbar if the Toolbar() function is in use.

Func ToolbarVisible({show%});

- show% If present and non-zero, the toolbar is made visible. If this is zero, and the Toolbar() function is not active, the toolbar is made invisible.
- Returns The state of the toolbar at the time of the call. The state is returned as 2 if the toolbar is active, 1 if it is visible but inactive and 0 if it is invisible.
- See also:Interact(), Toolbar(), ToolbarClear(), ToolbarEnable(), ToolbarSet(), ToolbarText()

Trunc()

Removes the fractional part of a real number or array. To truncate a real number to an integer, assign the real to the integer. ArrConst() copies a real array to an integer array.

Func Trunc(x|x[]);

- x A real number or a real array.
- Returns 0 or a negative error code for an array. For a number it returns the value with the fractional part removed. Trunc (4.7) is 4.0; Trunc (-4.7) is -4.0.

See also:Abs(), ATan(), Cos(), Exp(), Frac(), Ln(), Log(), Max(), Min(), Pow(), Rand(), Round(), Sin(), Sqrt(), Tan(), Trunc()

UCase\$()

This function converts a string into upper case. The upper-case operation may be system dependent. Some systems may provide localised uppercasing, others may only provide the minimum translation of the ASCII characters a-z to A-Z.

Func UCase\$(text\$);

text\$ The string to convert.

Returns An upper case version of the original string.

See also:Asc(), Chr\$(), DelStr\$(), InStr(), LCase\$(), Left\$(), Len(), Mid\$(), Print\$(), ReadStr(), Right\$(), Str\$(), Val()

Val()

This converts a string to a number. The converter allows the same number format as the script compiler and leading white space is ignored.

```
Func Val(text${, &nCh%});
```

text\$ The string to convert to a number. The expected format is:
 {white space}{-}{digits}{.digits}{e|E{+|-}digits}

nCh% If present, it is set to the number of characters used to construct the number.

Returns It returns the extracted number, or zero if no number was present.

View(),	V	iew(().X() and
			Viev	v().[]

The View() function sets the current view and returns the last view handle. A view handle is a positive integer > 0. Changing the current view does not change the focus or bring the view to the front, use FrontView() to do that.

Func View({vh%});

- vh% An integer argument being:
 - >0 A valid view handle of a view that is to be made the current view. Use ViewKind() to test for a valid handle.
 - 0 (or omitted) no change of the current view is required.
 - <0 If vh% is -n, this selects the nth duplicate of the time view that is associated with the current view. The current view can be a time, result or XY view. If the current view is a time view, this is equivalent to Dup(n). View(-1) in a result or XY view returns the time view from which it was created.
- Returns the view handle of the view that was current at the time of the call. If an argument is passed in which is not valid, the script stops with an error.
- **View().x()** The View().x() construction overrides the current view for the evaluation of x(). For example, View(vh%).Draw(1,2) draws view number vh%. It is an error if the selected view does not exist or if the function closes the original view, and the script stops.

View(vh%).x()

vh% A view handle of an existing view, 0 for the current view, or -n for the nth duplicate of the time view associated with the current view.

The equivalent code to View (vh%).x() is:

var temp%; temp% := View(vh%); 'Save the current view x(); 'call the user-defined or built-in function View(temp%); 'restore the original view

This means that View (vh%).FileClose() causes an error if vh% is the current view.

View().[] The View().[] construction overrides the current view to give you access to the data arrays that form result view channels.

View(vh%{,ch%}).[]

- vh% A view handle of an existing result view or 0 for the current result view.
- ch% An optional channel number in the result view. If omitted, channel 1 is assumed.

For example:

ArrConst(View(0,2).[6:20], 0); 'zero 20 elements of channel 2
View(rv%).[4] := 3; 'set fifth element of channel 1

See also:App(), Dup(), FileClose(), FrontView(), SampleHandle(), ViewFind(), ViewKind(), ViewList()

ViewColour()	This function gets and sets the colours of time, result and XY view items, overriding the application-wide colours set by Colour(). Currently you can set the background colour.			
	<pre>Func ViewColour(item%{, col%});</pre>			
	item% The colour item to get or set; 0=background			
	col% If present, the new colour index for the item. There are 40 colours in the palette with indices 0 to 39. Use -1 to revert to the application colour for the item.			
	Returns The palette colour index at the time of the call or -1 if no view colour is set.			
	See also:ChanColour(), Colour(), PaletteGet(), PaletteSet(), XYColour()			
ViewFind()	This function searches for a window with a given title and returns its view handle.			
Ű	<pre>Func ViewFind(title\$);</pre>			
	title\$ A string holding the view title to search for.			
	Returns The view handle of a view with a matching title, or 0 if no view matches.			
	See also:WindowTitle\$(),View(),ViewList()			
ViewKind()	 This function returns the type of the current view or a view identified by a view hand Types 5-7 are reserved. Type 11 windows include the sampling configuration window and control panel, the sequencer control panel and multimedia windows. Func ViewKind({vh%}); vh% An optional view handle. If omitted, the type of the current view is returned Use -n in a result or XY view to find the nth duplicate of the time view from the type of the type of the time view from the type of the type of the time view from the type of the type of the time view from the type of the type of the type of the time view from type of the type of type of the type of the type of the type of the type of t			
	Returns. The type of the current view or the view identified by yb View types are:			
	-2 Invalid handle1 Text view4 Result view10 Application window-1 Unknown type2 Output sequence8 External text file11 Other types0 Time view3 Script9 External binary file12 XY view			
	See also:ChanKind(),SampleHandle(),ViewList()			
ViewLink()	This function returns the view handle of the view that owns the current window. For example, you can use this to get the time view that owns a multimedia, spike shape or spike monitor window or the time view that created a result or XY view. This is slightly different from $\forall iew(-n)$, which finds the n^{th} duplicate of the time view linked to the current time, result or XY view.			
	FUNC VIEWLINK();			

Returns The handle of the linked view, or 0 if there is no such view.

See also:App(), SampleHandle(), View(), ViewKind(), ViewList()

ViewList()	This function fills an integer array with a list of view handles. It never returns the view handle of the running script; use the $App()$ command to get this.
	<pre>Func ViewList({list%[] {,types%}});</pre>
	list% An integer array that is returned holding view handles. The first element of the array (element 0) is filled with the number of handles returned. If the array is too small to hold the full number, the number that will fit are returned.
	types% The types of view to include. This is a code that can be used to filter the view handles. The filter is formed by adding the types from the list below. If this is omitted or if no types are specified for inclusion, all view handles are returned.
	1 Time views8 Script views512 External binary4096 XY views2 Text views16 Result views1024 Application view4 Sequencer256 External text2048 Other view types
	 You can also exclude views otherwise included by adding: 8192 Exclude views not directly related to the current view 16384 Exclude visible windows 32768 Exclude hidden windows 65536 Exclude duplicates
	Returns The number of windows that match types%.
	The following example prints all the window titles into the log view:
	<pre>var list%[100],i%; 'Assume no more than 99 views ViewList(list%[]); 'Get a list of all the views for i%:=1 to list%[0] do 'Iterate round all the views PrintLog(view(list%[i%]).WindowTitle\$()+"\n"); next;</pre>
	<pre>See also:App(), SampleHandle(), ViewKind()</pre>
ViewStandard()	This sets the current time, result or XY view to a standard state by making all channels and axes visible in their standard drawing mode, axis range and colour. All channels are given standard spacing and are ungrouped and the channels are sorted into the numerical order set by the Edit menu Preferences. In a time view, duplicate channels are deleted, triggered mode is disabled and any channel processing is removed. In an XY view the key is hidden and the axes are optimised. It has no effect on other view types. Proc ViewStandard();
	See also:ChanOrder(),ChanWeight(),DrawMode(),XYDrawMode()
ViewTrigger()	This controls the triggered drawing mode available for time views; it is a fatal error to use this command in a different view type. This command is the equivalent of the View menu Display Trigger command.
	Func ViewTrigger(chan%, pre, hold, xZero% {,cur0%});Set and enableFunc ViewTrigger({mode%});Enable/disable and get information
	chan% The event, marker, WaveMark, TextMark or RealMark trigger channel. Set this to 0 for a paged display on-line or during rerun.
	pre The pre-trigger time to display, in seconds.
	hold The minimum hold time for on-line displays, the minimum time to the next/previous trigger event for off-line use.
	xZero% If this is set to 1, the x axis zero point (for display only) is moved to the trigger time. All time measurements are still relative to the start of the file.

- cur0% This optional argument controls the cursor 0 action each time the view is triggered. 0=no action, 1=move cursor 0 without causing the active cursors to iterate, 2=move cursor 0 and update any active cursor positions. If this argument is omitted, the value 0 is used.
- mode[®] This command version gets information, moves to the next and previous trigger, and enables and disables the trigger. With no arguments, the command returns the enabled or disabled state as 1 or 0. If there is a single argument, it can be:
 - 3 Move to the next trigger and return the trigger time or -1 if no trigger is found or if the trigger is not enabled.
 - 2 Move to the previous trigger and return the trigger time or -1 if no trigger is found or if the trigger is not enabled.
 - 1 Enable the View trigger with the current settings.
 - 0 Disable the trigger.
 - -1 Return the trigger channel.
 - -2 Return the pre trigger time.
 - -3 Return the hold time.
 - -4 Return the xZero% state.
 - -5 Return the cur0% state.
- Returns When the command is used with a negative argument, the command returns the requested information, of -1 if the information is not available. All other use returns the enabled/disabled state as 1 or 0, or a negative error code.

See also:DrawMode(), ViewStandard(), XYDrawMode()

ViewUseColour()

This function can be used to force the display to use black and white only, or to use the colours set in the Colour dialog or by the Colour () command.

Func ViewUseColour({use%});

use% If present, a value of 0 forces Spike2 to display all windows in black and white. Any other value allows the use of colour. If omitted, no change is made.

Returns The current state as 1 if colour is in use, 0 if black and white is used.

See also:ChanColour(), Colour(), XYColour()

VirtualChan()

This command controls virtual channels in the current time view. Virtual channels are formed by combining existing waveform and RealWave channels using an algebraic expression. This command has the functionality of the Virtual Channel dialog. The first command variant creates and modifies a virtual channel, the second reports the state of a virtual channel.

Func VirtualChan(chan%, expr\${, match%{, binsz{, align}});
Func VirtualChan(chan%, get%{, &expr\$});

- chan% In the first command variant, set this to 0 to create a new virtual channel and return the channel number. The remaining arguments set the initial channel settings, otherwise default values are used. If not 0, this is the number of an existing virtual channel to modify or from which to read back the settings.
- expr\$ In the first command variant, this is a string expression defining the output. See the Analysis chapter of the Spike2 manual for details. In the second command variant, this is an optional string variable that is returned holding the current expression for the channel.

- match% This is the number of an existing waveform, RealWave or WaveMark channel to match for sample interval and alignment. If 0, the sample interval and alignment are set by the binsz and align arguments. If negative, no change is made.
- binsz If match% is 0, this sets the sample interval of the channel. Values of binsz less than the file time resolution are ignored, as are values less than or equal to 0.
- align If match% is 0, this sets the channel alignment. Values less than 0 are ignored.
- get% In the second command variant, get% determines the returned value. 0 = the state of the expression parsing, 1 = match%, 2 = binsz, 3 = align.
- Returns When creating a new channel, the return value is the new channel number or a negative error code. When modifying an existing channel, the return value is a negative code if the match%, binsz or align arguments are illegal, a positive code if the expression contains an error and 0 if there was no error. The second variant returns the requested information or a negative error code. If get% is 0, the return value is 0 if the expression is acceptable and a positive code if not.

See also:MemChan()

Window()

This sets the position and size of the current view. Normally, positions are percentages from the top left-hand corner of the application window size. You can also set positions relative to a monitor. This can also be used to position, dock and float dockable toolbars.

- Proc Window(xLow, yLow{, xHigh{, yHigh{, scr%{, rel%}}});
- xLow Position of the left hand edge of the window. When docking a dockable toolbar, the xLow and xHigh values correspond to the position of the top left corner of the window when dropped with the mouse.
- yLow Position of the top edge of the window.
- xHigh If present, the right hand edge. If omitted the previous width is maintained. If the window is made too small, the minimum allowed size is used. If the current view is a dockable control and yHigh is 0, values less than 1 or greater than 4 float the window at (xLow, yLow), otherwise xHigh sets the docking state:
 - 1 Docked to the left window edge 3 Docked to the right edge
 - 2 Docked to the top window edge 4 Docked to the bottom edge
- yHigh If present, the bottom edge position. If omitted the previous height is maintained. If the window is made too small, the minimum allowed size is used. If the Window is dockable and yHigh is 0, this command sets the docked state of the window (see xHigh). Otherwise the window is floated with the nearest allowed width that is no more than xHigh-xLow. If xHigh-xLow is 0 or negative YHigh sets the height of the dockable window.
- SCT% Optional screen selector for views, dialogs and the application window. If omitted or -1, positions are relative to the application window. Otherwise, 0 selects the entire desktop rectangle and greater values select a particular screen rectangle (but see rel%). See System() for more screen information.
- rel% Ignored unless scr% is greater than 0. Set 0 or omit for positions relative to the intersection of the selected monitor and the application window, 1 for positions relative to the monitor. If there is no intersection, there is no position change. When positioning the application window, rel% is treated as if it were 1.

Examples:

view(App()).Window(0,0,100,100, 0); 'use all desktop view(App()).Window(0,0,100,100,2); 'use all second monitor

See also:App(), System(), WindowDuplicate(), WindowGetPos(), WindowSize(), WindowTitle\$(), WindowVisible()

	This duplicates the surrant time window, greating a new window that has all the actings
WindowDuplicate()	of the current window. It does not duplicate channels; these are shared with the existing window. The new window becomes the current view and is created invisibly.
	<pre>Func WindowDuplicate();</pre>
	Returns The view handle of the new window, a negative error code or 0 if there are no free duplicates. There is a limit of 64 duplicates per window.
	See also:Window(),WindowGetPos(),WindowTitle\$(),WindowVisible()
	This gets the window position of the surrent view with respect to the employtion window.
WindowGetPos()	Positions are measured from the top left-hand corner as a percentage.
	Proc WindowGetPos(&xLow, &yLow ,&xHigh ,&yHigh);
	xLow A real variable that is set to the position of the left hand edge of the window.
	yLow A real variable that is set to the position of the top edge of the window.
	xHigh A real variable that is set to the position of the right hand edge of the window or
	that returns a docking code for a docked control bar if yHigh is returned as 0.
	 Docked to the left window edge Docked to the top window edge Docked to the bottom edge
	yHigh A real variable that is set to the position of the bottom edge of the window or to 0 if the window is docked.
	See also:Window(),WindowDuplicate(),WindowSize(),WindowVisible()
WindowSize()	This resizes the current view without changing the top left-hand corner position. Setting a
	maximum allowed sets the appropriate limit. There are no errors from this function.
	<pre>Proc WindowSize(width, height);</pre>
	width The width of the window as a percentage of the available area.
	height The height of the window as a percentage of the available area.
	You can also use this to resize the application window or a dockable control bar when it is floating; use App() to get the handles. For control bars, if width is greater than zero, its sets the width, otherwise height is used to set the height. If the control bar can be resized, it will use the width or the height and will calculate the other dimension itself.
	<pre>See also:App(),Window(),WindowDuplicate(),WindowGetPos(), WindowTitle\$(),WindowVisible()</pre>
WindowTitle\$()	This function gets and sets (where allowed) the current window title. Most windows can return a title. If you change a title, dependent window titles change, for example, cursor windows belonging to time views track the title of the time view.
	<pre>Func WindowTitle\$({new\$});</pre>
	new\$ If present, this sets the new window title. Window titles must follow any system rules for length or content. Illegal titles (for example titles containing control characters) are mangled or ignored at the discretion of the system.
	Returns The window title as it was prior to this call.
	See also:Window(),WindowDuplicate(),WindowGetPos(),WindowSize()

WindowVisible()	This function is used to get and set the visible state of the current window. This function can also be used on the application window, however the effect will vary with the system and on some, there may be no effect at all.
	<pre>Func WindowVisible({code%});</pre>
	 code% If present, this sets the window state. The possible states are: 0 Hidden, the window becomes invisible. A hidden window can be sent data, sized and so on, the result is just not visible. 1 Normal, the window assumes its last normal size and position and is made visible if it was invisible or iconised.
	 Iconised under Windows, hidden on the Macintosh. An iconised window can be sent data, sized and so on; the result is not visible. This will dock a dockable window at its last docked position. Maximise, make it as large as possible or float a dockable window. Application window only; extend over all available desktop monitors. Returns The window state prior to this call.
	<pre>See also:FrontView(),Window(),WindowDuplicate(),WindowGetPos(), WindowSize(),WindowTitle\$()</pre>
XAxis()	This turns on and off the x axis of the current view and returns the state of the x axis. Func XAxis({on%});
	on% Set the axis state. If omitted, no change is made. 0=Hide axis, 1=Show axis.
	Returns The axis state at the time of the call (0 or 1, as above) or a negative error code. It is an error to use this function on a view that has no concept of an x axis.
	See also:XAxisMode(),XAxisStyle(),XHigh(),XLow(),XRange()
XAxisMode()	This function controls what is drawn in an x axis. Func XAxisMode({mode%});
	 mode% Optional argument that controls how the axis is displayed. If omitted, no change is made. Possible values are the sum of the following: 1 Hide all the title information. 2 Hide all the unit information. 4 Hide small ticks on the x axis. Small ticks are hidden if big ticks are hidden. 8 Hide numbers on the x axis. Numbers are hidden if big ticks are hidden. 16 Hide the big ticks and the horizontal line that joins them. 32 Scale bar axis. If selected, add 4 to remove the end caps.
	Returns The x axis mode value at the time of the call or a negative error code.
XAxisStyle()	This function controls the x axis mode (seconds, hours minutes and seconds, time of day) for a time view, and the major and minor tick spacing for all views that have an x axis. If you set values that would cause illegible or unintelligible axes, they are stored but not used unless the axis range or scaling changes to make the values useful.
	<pre>Func XAxisStyle({style%{, nTick%{, major}}});</pre>
	style% In a time view, set I for an axis in seconds, 2 for hours minutes and seconds and 3 for time of day. Omit style% or use 0 to leave the style unchanged. A value

of -1 returns the number of minor divisions set or 0 for automatic. A value of -2 returns the major tick spacing or 0 for automatic spacing.

- nTick% The number of minor tick subdivisions or 0 for automatic spacing. Omit nTick% or set it to -1 for no change.
- major If present, values greater than 0 sets the major tick spacing. Values less than or equal to 0 select automatic major tick spacing.
- Returns If style% is positive or omitted it returns the style at the time of the call. See the description of style% for negative values.

See also:XAxis(), XAxisMode(), XHigh(), XLow(), XRange(), YAxisStyle()

XHigh()

This returns the x axis value at the right hand side of the current time, result or XY view or the end of the time range to process for a spike shape view. Use with other views causes a fatal script error.

Func XHigh()

Returns In a time view, the result is in seconds. In a result view, the result is in bins and can be fractional. In an XY view the result is in x axis units.

This example pages through a time or result view from the current position to the end. while (xHigh() < MaxTime()) do Draw(XHigh()) wend;

See also:Draw(), XRange(), BinToX(), XToBin(), XLow()

XLow()

This returns the x axis value at the left hand side of the current time, result or XY view or the start of the time range to process for a spike shape view.

Func XLow()

Returns In a time view, the result is in seconds. In a result view, the result is in bins and can be fractional. It is a fatal error to use this in an inappropriate view.

For example, this code pages a time or result view from the current position to the start of the file. In an XY view it moves the view to the left if the current position is past 0.

while XLow()>0 do Draw(XLow()-(XHigh()-XLow())) wend;

See also:Draw(), XRange(), BinToX(), XToBin(), XHigh()

XRange()

This sets the x axis range to display in a time result or XY view in x axis units (not bins for a result view). Unlike Draw(), the view does not update immediately; updates wait for the next Draw(), DrawAll(), Yield() or some interactive activity. It also sets the time range to process in a spike shape window.

Proc XRange(from {,to});

- from The left hand edge of the view in x axis units (seconds for a time view). You can set from to -1 in a time view that is sampling or re-running to make the view scroll automatically to show the most recent data at the right-hand edge.
- to The right hand edge of the view. If omitted, the view stays the same width.

Values are limited to the axis range for time and result views; there is no limit in an XY view. Without to, it preserves the width, adjusting from if required. If the resulting width is less than the minimum allowed, no change is made.

See also:Draw(), XLow(), XHigh()

XScroller()	This function gets and optionally sets the visibility of the x axis scroll bar and controls. Func XScroller({show%}); show% If present, 0 hides the scroll bar and buttons, non-zero shows it.
	Returns 0 if the scroll bar was hidden, 1 if it was visible.
XTitle\$()	and returns an empty string. The window updates with a new title at the next opportunity.
	<pre>Func XTitle\$({new\$});</pre>
	New\$ If present, this sets the new x axis title in a result view.
	Returns The x axis title at the time of the call.
	See also:ChanTitle\$()
	In a result view, this converts between hin numbers and x axis units. In a time view, it
XTOBIN()	converts between time in seconds and the underlying Spike2 time units.
	<pre>Func XToBin(x);</pre>
	x An x axis value. If it exceeds the x axis range it is limited to the nearer end.
	Returns In a result view it returns the bin position that corresponds to x. In a time view, it converts seconds to the underlying time units used for the file.
	See also:BinToX()
XUnits\$()	This function gets the units of the x axis. You can also set the units in a result or XY view. The window will update with the new units at the next opportunity.
	New $\$$ If present, this sets the new x axis units in a result view.
	Returns The x axis units at the time of the call.
	See also:ChanUnits\$()
	This adds data points to an XV view channel. If the axes are set to automatic expanding
XYAddData()	mode by XYDrawMode(), they will change when you add a new data point that is out of the current axis range. If the channel is set to a fixed size (see XYSize()), adding new points causes older points to be deleted once the channel is full. The first form of the command allows unrestricted x and y positions. The second form was added at version 5.04 and is for data that is equally spaced in the x direction.
	<pre>Func XYAddData(chan%, x x[] x%[], y y[] y%[]); Func XYAddData(chan%, y[], xInc{, xOff});</pre>
	chan% A channel number in the current XY view. The first channel is number 1.
	x The x co-ordinate(s) of the added data point(s). In the first form of the

- command, both x and y must be either single variables or arrays. If they are arrays, the number of data points added is equal to the size of the smaller array.
- The y co-ordinate(s) of the added data point(s). In the second form of the У command, this is an array of equally spaced data in x.
- Sets the x spacing between the y data points in the second form of the command. xInc

omitted, the first position is set to 0.

Returns The number of data points which have been added successfully.

xOff

Sets the x position of the first data point in the second form of the command. If

	<pre>See also:XYColour(), XYDelete(), XYDrawMode(), XYJoin(), XYKey(), XYRange(), XYSetChan(), XYSize(), XYSort()</pre>
XYColour()	This gets or sets a channel colour in the current XY view. The default colour is black. Func XYColour(chan% {, col%}); chan% A channel number in the current XY view. The first channel is number 1. col% The index of the colour in the colour palette. There are 40 colours in the palette; their indexes are numbered from 0 to 39. If omitted, there is no colour change. Returns The colour index in the colour palette at time of call or a negative error code. See also:Colour(), XYDrawMode(), XYJoin(), XYKey(), XYSetChan()
XYCount()	This gets the number of data points in a channel in the current XY view. To find the maximum number of data points, see the XYSize() command. Func XYCount(chan%); chan% A channel number in the current XY view. The first channel is number 1. Returns The number of data points in the channel or a negative error code. See also:XYAddData(), XYDelete(), XYJoin(), XYGetData(), XYInCircle(), XYInRect(), XYRange(), XYSetChan(), XYSize()
XYDelete()	<pre>This command deletes a range of data points or all data points from one channel of the current XY view. Use ChanDelete() to delete the entire channel. Func XYDelete(chan% {, first% {,last%}}); chan% A channel number in the current XY view. The first channel is number 1. first% The zero-based index of the first point to delete. Omit to delete all points. last% The zero-based index of the last data point to delete. If omitted, data points from first% to the last point in the channel are deleted. If last% is less than first% no data points are deleted.</pre> Returns The function returns the number of deleted data points. The index number of a data point depends on the current sorting method of the channel set by XYSort(). For different sorting methods, a data point may have different index numbers. The data points in a channel have continuous index numbers. When a point has been deleted the remaining points re-index themselves automatically. See also:ChanDelete(),XYAddData(),XYCount(),XYSetChan(),XYSize()

XYDrawMode()

This gets and sets the drawing and automatic axis expansion modes of a channel in the current XY view. It is an error to use this command with any other view type.

Func XYDrawMode(chan%, which% {,new%});

1

- A channel number in the current XY view. The first channel is number 1. This is chan% ignored when which is 5, as all XY channels share the same axes. -1 can also be used, meaning all channels.
- which% The drawing parameter to get or set in the range 1 to 5. When setting parameters, the new value is held in the new% argument. The values are:
 - Get or set the data point draw mode. The drawing modes are:
 - 0 dots (default) 1 boxes 2 plus signs +
 - 3 crosses x 4 circles (NT only)
 - 5 triangles 6 diamonds 7 horizontal line 8 vertical line
 - 2 Get or set the size of the data points in points (units of approximately 0.353) mm). The sizes allowed are 0 to 100; 0 is invisible. The default size is 5.
 - 3 Get or set the line style. If the line thickness is greater than 1 all lines are drawn as style 0. Styles are:
 - 0 solid (default) 1 dotted 2 dashed
 - 4 Get or set the line thickness in points. Thickness values range from 0 (invisible) to 10. The default is 1.
 - 5 Get or set automatic axis range mode. This applies to the entire view, so the chan% argument is ignored. Values are:
 - 0 The axes do not change automatically when new data points are added.
 - When new data points are added that lie outside the current x or y axis range, the data and axes screen area update at the next opportunity to display all the data.
- New draw mode or axis expanding mode. If omitted, no change is made. new%
- Returns The value of the relevant channel draw mode or axis expanding mode at the time of the call or a negative error code.

See also:XYAddData(), XYColour(), XYCount(), XYDelete(), XYGetData(), XYInCircle(), XYInRect(), XYJoin(), XYKey(), XYRange(), XYSetChan(),XYSize(),XYSort()

XYGetData()

This gets data points between two indices from a channel in the current XY view. It is an error to use this command with any other view type.

Func XYGetData(chan%, &x|x[], &y|y[] {,first% {,last%}});

- chan% A channel number in the current XY view. The first channel is number 1.
- The returned x co-ordinate(s) of data point(s). When arrays are used, either both x | x [] x and y must be arrays or neither can be. The smaller of the two arrays sets the maximum number of data points that can be returned.
- y | y[] The returned y co-ordinate(s) of data point(s).
- first% The zero-based index of the first data point to return. If omitted, 0 is used.
- last & The zero-based index of the last data point returned, used when x and y are arrays. If omitted or greater than or equal to the number of data points, the final data point is the last one in the channel. If last% is less than first%, no data points are returned.
- Returns The number of data points copied. If the x or y arrays are not big enough to hold all the data points from first% to last%, the return value is the array size. If xor y are not arrays, if a data point with index first% exists, 1 is returned.

The index number of a data point depends on the current sorting method (see XYSort()).

See also:XYAddData(), XYCount(), XYDelete(), XYInCircle(), XYInRect(), XYJoin(), XYRange(), XYKey(), XYSetChan(), XYSize(), XYSort()

XYInCircle()	This gets the number of data points inside a circle defined by xc, yc, and r in the current XY view. A general point (x, y) is considered to be inside the circle if:
	$(x-xc)^{2} + (y-yc)^{2} <= r^{2}$
	Points lying on the circumference are considered inside, but owing to floating-point rounding effects they may be indeterminate.
	<pre>Func XYInCircle(chan%, xc, yc, r);</pre>
	chan% A channel number in the current XY view. The first channel is number 1.
	xc, yc These are the x and y co-ordinates of the centre of the circle.
	r This is the radius of the circle. r must be ≥ 0 .
	Returns The number of data points inside the circle or a negative error code.
	See also:XYAddData(),XYInRect(),XYJoin(),XYRange()
XYInRect()	This function returns the number of data points in a channel of the current XY view that lie inside a rectangle. To be inside, a data point must lie between the low rectangle coordinate up to, but not including the high coordinate. This is so that two rectangles with a common edge will not both count a data point on the boundary.
	<pre>Func XYInRect(chan%, xl, yl, xh, yh);</pre>
	chan% A channel number in the current XY view.
	x1, xh The x co-ordinates of the left and right hand edges of the rectangle. xh must be greater than or equal to x1.
	yl, yh The y co-ordinates of the bottom and top edges of the rectangle. yh must be greater than or equal to yl.
	Returns The number of data points inside the rectangle or a negative error code.
	See also:XYInCircle(),XYGetData(),XYJoin(),XYRange()
XYJoin()	This function gets or sets the data point joining method of a channel in the current XY view. Data points can be separated, joined by lines, or joined by lines with the last point connected to the first point (making a closed loop).
	<pre>Func XYJoin(chan% {,join%});</pre>

- chan% A channel number in the current XY view. The first channel is number 1. -1 is also allowed, meaning all channels.
- join% If present, this is the new joining method of the channel. If this is omitted, no change is made. The data point joining methods are:
 - 0 Not joined by lines (this is the default joining method)
 - 1 Joined by lines. The line styles are set by XYDrawMode().
 - 2 Joined by lines and the last data point is connected to the first data point to form a closed loop.

Returns The joining method at the time of the call or a negative error code.

See also:XYColour(),XYDrawMode(),XYKey(),XYSetChan(),XYSort()

XYKey()

This gets or sets the display mode and positions of the channel key for the current XY view. The key displays channel titles (set by ChanTitle\$()) and drawing symbols of all the visible channels. It can be positioned anywhere within the data area. The key can be framed or unframed, transparent or opaque and visible or invisible.

Func XYKey(which%, {new});

which% This determines which property of the key we are interested in. Properties are:

- 1 Visibility of the key. 0 if the key is hidden (default), 1 if it is visible.
- 2 Background state. 0 for opaque (default), 1 for transparent.
- 3 Draw border. 0 for no border, 1 to draw a border (default)
- 4 Key left hand edge x position. It is measured from the left-hand edge of the x axis and is a percentage of the drawn x axis width in the range 0 to 100. The default value is 0.
- 5 Key top edge y position. It is measured from the top of the XY view as a percentage of the drawn y axis height in the range 0 to 100. The default is 0.

new If present it changes the selected property. If it is omitted, no change is made.

Returns The value selected by which% at the time of call, or a negative error code.

See also:ChanTitle\$(), XYColour(), XYDrawMode(), XYJoin(), XYSetChan()

XYRange()

This function gets the range of data values of a channel or channels in the current XY view. This is equivalent to the smallest rectangle that encloses the points.

Func XYRange(chan%, &xLow, &yLow, &xHigh, &yHigh);

- chan% A channel number in the current XY view or -1 for all channels or -2 for all visible channels. The first channel is number 1.
- xLow A variable returned with the smallest x value found in the channel(s).
- yLow A variable returned with the smallest y value found in the channel(s).
- xHigh A variable returned with the biggest x value found in the channel(s).
- yHigh A variable returned with the biggest y value found in the channel(s).

Returns 0 if there are no data points, or the channel does not exist, 1 if values found.

See also:XYAddData(),XYCount(),XYDrawMode(),XYGetData(),XYInCircle(), XYInRect(),XYSetChan(),XYSize()

XYSetChan()

This function creates a new channel or modifies an existing channel in the current XY view. It is an error to use this function if the current view is not an XY view. This function can modify all properties of an existing channel without calling the XYSize(), XYSort(), XYJoin() and XYColour() commands individually.

Func XYSetChan(chan% {,size% {,sort% {,join% {,col%}}});

- chan% A channel number in the current XY view. If chan% is 0, a new channel is created. Each XY view can have maximum of 256 channels, numbered 1 to 256. Spike2 creates the first channel automatically when you open a new XY view with FileNew(). If chan% is not 0, it must be the channel number of an existing channel to modify or -1 to modify all channels.
- size% This sets the number of data points in the channel and how and if the number of data points can extend. The only limits on the number of data points are the available memory and the time taken to draw the view.

A value of zero (the default) sets no limit on the number of points and the size of the channel expands as required to hold data added to it.

If a negative size is given, for example -n, this limits the number of points in the channel to n. If more than n points are added, the oldest points are deleted to make space for the newer points. If you set a negative size for an existing channel that is smaller than the points in the channel, points are deleted.

If a positive value is set, for example n, this allocates storage space for n data points, but the storage will grow as required to accommodate further points. Using a positive number rather than 0 can save time if you know in advance the likely number of data points as it costs time to grow the data storage.

- sort% This sets the sorting method of data points. The sorting method is only important if the points are joined. If they are not joined, it is much more efficient to leave them unsorted as sorting a large list of points takes time. The sort methods are:
 - 0 Unsorted (default). Data is drawn and sorted in the order that it was added. The most recently added point has the highest sort index.
 - 1 Sorted by x value. The index runs from points with the most negative x value to points with the most positive x value.
 - 2 Sorted by y value. The index runs from points with the most negative y value to points with the most positive y value.

If this is omitted, the default value 0 is used for a new channel. For an existing channel, there is no change in sorting method.

- join% If present, this is the new joining method of the channel. If this is omitted, no change is made to an existing channel; a new channel is given mode 0. The data point joining methods are:
 - 0 Not joined by lines (this is the default joining method)
 - 1 Joined by lines. The line styles are set by XYDrawMode().
 - 2 Joined by lines and also connect the first and last data points to form a loop.
- col% If present, this sets the index of the colour in the colour palette to use for this channel. There are 40 colours with indices 0 to 39. If omitted, the colour of an existing channel is not changed. The default colour for a new channel is the colour that a user has chosen for an ADC channel in a time window.
- Returns The highest channel number that was affected or a negative error code. When you create a channel, the value returned is the new channel number.

See also:XYAddData(), XYColour(), XYDelete(), XYDrawMode(), XYJoin(), XYKey(), XYRange(), XYSize(), XYSort()

XYSize()

This function gets and sets the limits on the number of data points of a channel in the current XY view. Channels can be set to have a fixed size, or to expand as more data is added. The only limit on the number of data points is the available memory and the time taken to draw data.

Func XYSize(chan% {,size});

- chan% A channel number in the current XY view. The first channel is number 1.
- size% This sets the number of data points in the channel and how and if the number of data points can extend. A value of zero sets no limit on the number of points and the size of the channel expands as required to hold data added to it.

If a negative size is given, for example -n, this limits the number of points in the channel to n. If more than n points are added, the oldest points are deleted to make space for the newer points. If you set a negative size for an existing channel that is smaller than the points in the channel, points are deleted.

If a positive value is set, for example n, this allocates storage space for n data points, but the storage will grow as required to accommodate further points.

Using a positive number rather than 0 can save time if you know in advance the likely number of data points as it costs time to grow the data storage.

If this is omitted, there is no change to the size.

Returns If the number of points for the channel is fixed at n points, the function returns -n. Otherwise, the function returns the maximum number of points that could be stored in the channel without allocating additional storage space.

See also:XYAddData(),XYColour(),XYCount(),XYDelete(),XYDrawMode(), XYGetData(),XYInCircle(),XYInRect(),XYJoin(),XYKey(), XYRange(),XYSetChan(),XYSort()

XYSort()

In the current XY view, gets or sets the sorting method of the channel.

Func XYSort(chan% {,sort%});

- chan% the channel number in the current XY view. The first channel is number 1.
- sort% This sets the sorting method of data points. The sorting method is only important if the points are joined. If they are not joined, it is much more efficient to leave them unsorted as sorting a large list of points takes time. The sort methods are:
 - 0 Unsorted (default). Data is drawn and sorted in the order that it was added. The most recently added point has the highest sort index.
 - 1 Sorted by x value. The index runs from points with the most negative x value to points with the most positive x value.
 - 2 Sorted by y value. The index runs from points with the most negative y value to points with the most positive y value.

If this is omitted, there is no change in sorting method.

Returns The function returns the sorting method at time of call or a negative error code.

See also:XYAddData(),XYColour(),XYCount(),XYDelete(),XYDrawMode(), XYGetData(),XYInCircle(),XYInRect(),XYJoin(),XYKey(), XYRange(),XYSetChan(),XYSize()

YAxis()

This function is used to turn the y axes on and off, in the current view and to find the state of the y axes in a view.

Func YAxis({on%});

- on% Optional argument that sets the state of the axes. If omitted, no change is made. Possible values are:
 - 0 Hide all y axes in the view.
 - 1 Show all y axes in the view.

Returns The state of the y axes at the time of the call (0 or 1) or a negative error code.

See also:ChanOffset(), ChanScale(), Grid(), Optimise(), XAxis(), XScroller(), YAxisMode(), YHigh(), YLow(), YRange()

YAxisLock()

This function locks and unlocks the axes of grouped channels and reports on the locked state of grouped channels. If you lock a group, the grouped channels keep their own axis ranges, but display using the axis of the first channel in the group. The YRange(), YHigh() and YLow() commands operate on the information stored with a channel. To change the displayed range of grouped and locked channels, you must use YRange() on the first channel in a group.

Func YAxisLock(chan%{, opt%{, vOffs}});

- chan% A channel that is in the group that you wish to address.
- opt% If present, values of 1 and 0 set and unset the locked state. A value of -1 returns the visual offset per channel for the group. If omitted, no change is made.
- voffs If present, this sets the y axis display offset to apply between channels in the group. The nth channel has a visual offset of (n-1) *offs.
- Returns The current locked state of the group unless opt% was -1, when it returns the y axis visual offset per channel for the group.

See also:ChanOrder(), YAxisMode(), YHigh(), YLow(), YRange()

YAxisMode()

Func YAxisMode({mode%});

respect to the data.

mode[®] Optional argument that controls how the axis is displayed. If omitted, no change is made. Possible values are the sum of the following values:

This function controls what is drawn in a y axis and where the y axis is placed with

- 1 Hide all the title information.
- 2 Hide all the unit information.
- 4 Hide y axis small ticks. They are also hidden when big ticks are hidden.
- 8 Hide y axis numbers. They are also hidden when big ticks are hidden.
- 16 Hide the big ticks and the vertical line that joins them.
- 32 Scale bar axis. If selected add 4 to remove the end caps.
- 4096 Place the y axis on the right of the data

Returns The state of the y axis mode at the time of the call or a negative error code.

See also:ChanNumbers(), YAxis(), YAxisStyle(), YHigh(), YLow(), YRange()

YAxisStyle()

This function controls the y axis major and minor tick spacing. If you set values that would cause illegible or unintelligible axes, they are stored but not used unless the axis range or scaling changes to make the values useful.

Func YAxisStyle(cSpc, opt%{, major});

- CSpc A channel specifier or -1 for all, -2 for visible and -3 for selected channels. When multiple channels are specified, returned values are for the first channel.
- opt% Values greater than 0 set the number of subdivisions between major ticks. 0 sets automatic small tick calculation. Use -1 for no change. Values less than -1 return information, but do not change the axis style
- major If present and opt% is greater than -2, values greater than 0 sets the major tick spacing. Values less than or equal to 0 select automatic major tick spacing.
- Returns If opt% is -2 this returns the current number of forced subdivisions or 0 if they are not forced. If opt% is -3 this returns the current major tick spacing if forced or 0 if not forced. Otherwise the return value is 0 or a negative error code.

See also:YAxis(),YAxisMode(),YHigh(),YLow(),YRange(),XAxisStyle()

Yield()

This function suspends script operation for a user defined time and allows Spike2 to idle. During the idle time, invalid screen areas update, you can interact with the program and Spike2 has the opportunity to do housekeeping such as releasing temporary system resources. If your script runs for long periods without using Interact() or Toolbar(), adding an occasional Yield() can make it feel more responsive and stop the system marking Spike2 as "not responding".

Func Yield({wait{, allow%}});

- wait The minimum time to wait, in seconds. If omitted or 0, Spike2 gives the system one idle cycle. If greater than 0 and there is still waiting time left after an idle cycle completes, other processes are given the opportunity to run between idle cycles. If set negative, there is no idle cycle but the allow% argument is applied.
- allow% This defines what the user can do during the wait period. See Interact() for the allowed values. The allow% value is cancelled after the command unless wait is 0 or negative.

Returns The function returns 1. We may add more return codes in future versions.

See also:Interact(), MaxTime(), Seconds(), Toolbar(), YieldSystem()

YieldSystem()

To share the system CPU(s) among competing processes, the operating system allocates time slices of around 10 milliseconds based on process priorities and recent process activity. A process can surrender a time slice if it has nothing to do. A typical application spends most of its time waiting for user input, which appears as messages in the application message queue; it will surrender a time slice if the message queue is empty unless it has a task to work on. Spike2 normally surrenders time slices, but if you run a script, it runs for the full time slice unless it is in Yield(), Interact(), or you use ToolBar() or DlgShow() without an idle routine.

The YieldSystem() command causes Spike2 to surrender the current time slice and suspends the user interface and script thread for a user-defined time or until a new message arrives in the Spike2 input queue. It has no effect on sampling, which runs in a separate thread. Unlike Yield(), it does not allow Spike2 to idle.

Proc YieldSystem({wait});

wait The time to wait, in seconds, before resuming the thread. Values are rounded to the nearest millisecond. Values greater than 10 are treated as 10 seconds; values less than -10 are treated as -10 seconds.

For wait values greater than 0, the wait is ended by unserviced messages; keyboard and mouse activity, timers for screen updates and the like cause messages. If wait is 0 or omitted, the current time slice is surrendered, but if Spike2 is the highest priority task it will be re-scheduled immediately. Negative values suspend Spike2 for -wait seconds regardless of messages.

YieldSystem() with wait values greater than 0 returns immediately if there are messages in the input queue. Unless you allow Spike2 to idle, either with a Yield() call or with Toolbar() or DlgShow(), there will always be pending messages, so it will have no effect. If you have a script loop that causes 100% CPU usage, inserting:

Yield();YieldSystem(0.001);

will give other processes a chance to run. Increasing the wait time up to 0.05 will further reduce the CPU usage. Larger values have little additional effect due to timer messages

ending the wait early. To give as much system time as possible to other tasks without allowing Spike2 to idle, you can use:

```
YieldSystem(-0.001);
```

In this case, increasing the time to -10.0 will have an effect; Spike2 will feel completely unresponsive until the time period has elapsed.

See also:Interact(), DlgShow(), Seconds(), Toolbar(), Yield()

YHigh()

This function returns the upper y axis limit for the y axis of a channel in a time, result or XY view. The return value ignores the grouped or locked state of the channel.

Func YHigh(chan%);

chan% The channel number (use 0 for an XY view).

Returns The value at the appropriate end of the axis.

See also:YAxisLock(), YRange(), ChanScale(), ChanOffset(), Optimise()

YLow()

This function returns the lower y axis limit for the y axis of a channel in a time, result or XY view. The return value ignores the grouped or locked state of the channel.

Func YLow(chan%);

chan% The channel number. Use 0 for an XY view.

Returns The value at the appropriate end of the axis.

See also:YAxisLock(), YRange(), ChanScale(), ChanOffset(), Optimise()

YRange()

This sets the displayed y axis range for a channel. It has no effect if the current channel display mode has no y axis. If the y range changes, the display updates at the next opportunity. If you omit low and high, Spike2 displays the "full" range of the channel, equivalent to the Show All button in the Y Axis dialog.

Proc YRange(chan%{, low, high});

- chan% The channel number in a time or result view, 0 in an XY or spike shape view. You can also use -1 for all, -2 for visible and -3 for selected channels.
- low The value for the bottom of the y axis. If omitted, and the channel type has a known range, Spike2 sets low and high to suitable limits. For example, for a waveform channel displayed as a waveform the limits are those set by the 16-bit nature of the data. For a sonogram channel, the limits are 0 and half the sampling rate. For an event channel the limits are 0 and the maximum sustained event rate set in the sampling configuration dialog.
- high The value for the top of the y axis. If low and high are the same, or too close, the range is not changed. It is an error to supply low and omit high.

See also:YAxisLock(), YHigh(), YLow(), ChanScale(), Optimise()

Introduction This section is for Spike2 for DOS users who need to move scripts to the Spike2 for Windows environment. If you do not need to translate old scripts you can skip this information. We provide a translator that will do most of the hard work required for the changes, but it cannot cope with all possible scripts and it cannot translate keywords for which no equivalent exists.

To use the translator, open the DOS script as a text file and select the Script menu Convert DOS Script option. This opens a new script window for the output (plus any notes about changes you may need to make) and sends a summary of the notes to the Log window. The translator marks the output with notes like:

DrawMode(-1, mode+1); 'Original was "DRAWMODE mode"
'*** Note 7: No exact equivalent of DRAWMODE
'*** The effect of the new DRAWMODE may not be exactly the same. We have
'*** selected all channels and set a mode, please check it.

Because the mode is held in a variable the translator has no way of knowing what this will be at run time or what type of view will be active. It has made the best guess it can at the required translation, but is warning that you should inspect the result to be sure that it is right. If the same warning is required later in the script, the message is just:

'*** Note 7: No exact equivalent of DRAWMODE

Please don't expect the result to be perfect. It is our experience that the more effort went into the DOS script to make it look good on screen, the less successful the result.

What will translate Spike2 for DOS supports some 146 keywords and 14 operators. Because many of the keywords use numeric codes to make one keyword perform several different functions, there are some 200 different commands to translate. The vast majority translate exactly without any problem. All the maths functions, use of variables, strings and arrays and PROCS translate without any problem, as do expressions using these items.

What will not translate Spike2 for DOS has built into it a structure of 6 views that can hold time or result views plus text and graphic output streams. It is the mapping of this architecture to a windowed environment that causes most of the problems in translating a script.

We cannot translate concepts that have no equivalent in the windowed environment. For example, overdrawing of data has no equivalent (yet) so does not translate. Scripts that contain syntax errors cannot translate! It is quite common for the translator to find errors in scripts, usually in IF or DOCASE statements in branches that were never executed. If the translator complains about your script, correct the errors and try again.

The keywords that do not translate at all come into the following categories:

- Keywords that make little difference to the functioning of a script or are there because of lack of memory in DOS or are intrinsic to the DOS system or are related to debugging: DUMPM, FREEMEM, NUMVAR, NUMARR, NUMSTR, SHOWVARS, PARAMSTR. The lack of these is usually not a problem.
- Keywords associated with user-defined drawing of lines and/or boxes in views and overdrawing of data. These are the ones that will cause most problems. We will address these in future releases: BOX, DRAWTO, DRAWR, MOVER, OVERDRAW, SETAXES. You may be able to work around this using XY views.
- Keywords associated with running multiple scripts from one script to compensate for the lack of space (20,000 characters) allowed in the DOS scripts. From version 3, the only limit to script size built into the program is that a script must be less than 32,000

lines long. We have not yet decided on the best way to produce large multi-file scripts. Missing keywords are: EXECUTE, EXECUTED, GLOBAL. Note that you can have multiple scripts loaded at a time, so some of the need for these keywords is reduced. From version 4.08, the ScriptRun() command lets you set a script to run after the current script ends.

- Keywords associated with plotting and to some extent printing. Plotting is translated to printing as much as possible. There is no exact equivalent of export in HPGL or PIC format, but there is the far superior (in most cases) option of copy to the clipboard. PLOTTO is not supported at all.
- Finally there is the SETPLACE command. Please let us know if you use it.

Everything else either translates completely, or translates with some worries about compatibility that are marked with notes in the output.

issues

General translation This section describes how the translator copes with problems in the translation, and also explains the extra script variables and functions declared to emulate the original Spike2 for DOS environment.

- In the DOS version := is treated as an operator, and has a value. In the := windowed version, := is not an operator. Instead it indicates an assignment statement. Use of := as an operator is flagged as an error, for example: a:=b:=3.0; 'This is legal in the DOS version b:=3.0; a := b; 'You must edit the original to this The following is also not translatable (and is often a mistake): if a := b ... 'Allowed in the DOS version a := b; if a... 'You must edit the original to this Very few scripts make use of the assignment operator in this way.
- +-*/... Other operators (+-*/% and comparisons) differ because in the new script operators have an order of precedence. This means that 4+5*6 is 34 in the new script, not 54 as in the DOS script. Because of this all expressions are analysed and rewritten to insert and remove brackets as required.
- Codes Many DOS script commands use numeric codes to indicate the operation, for example STRFUNC. The new script avoids the use of codes as much as possible to make the commands readable, and because the length of a command name does not slow down the script at run time. For example compare i% := InStr(text\$, find\$); with STRFUNC 4 text\$ find\$ I. If you use a variable in place of a code we cannot translate.
- This keyword is not translated. If you intend to clear out text you could try Clear using SelectAll();EditClear().
- The commands associated with drawing (BOX, DRAWTO, DRAWR, MOVER and to Drawing some extent MOVETO) have no equivalent.
- ESCAPE This keyword has no real equivalent and is translated as 0. You can get exactly the same effect by adding an extra Cancel button to an Interact bar, testing the result from the DlgShow() command or by setting an escape button in the Toolbar() command. However, such transformations are beyond the capabilities of the script converter, so you must do them yourself.

Execute This is not translated. Executed is translated as 0.

- File The new file command does not automatically support the use of numbers as part of a name. We add Func File() to emulate this and sort out the views.
- FRAME This is ignored in ON/OFF and treated as 1 in expressions.
- Global This is translated as var.
- Help This is translated, but most unlikely to be of any use unless you write your own help files.
- names The new script language defines many more commands than the DOS version. These commands have descriptive names that may well clash with names your script used. If a name clash is detected we add QQ to the end of the original name. The translator also preserves your capitalisation of all user-defined names (though the system is not case-sensitive).
- NewEvent The memory channel functions have more functionality than the NEWEVENT buffer. The variable neqq% holds the memory channel number used to emulate both NEWEVENT and EVENT. The memory channels behave just like any other channel, so you must be in the correct time view to use them.
- Normal If your script uses this command we add Proc Normal (v%) to emulate it.
- Numxxx The NUMARR, NUMVAR and NUMSTR keywords are ignored.
- Overdraw This keyword is not supported. It is translated as 0 in expressions. You may be able to emulate overdrawing with an XY view. The overdraw WaveMark mode will overdraw Spike shapes for you.
- PrintTo This keyword does not exist, so we add our own to emulate it. We also add the variable pqq% to hold the view handle of the view to which printing to a file should be sent. Any printing that previously went to the screen is now routed to the Log window.
- Print When printing text, you must now specify new lines in the output with the \n escape sequence. The translator adds this to the end of format strings unless it detects that output is to a device that is specified by an expression that starts with a minus sign. This will work in *almost* all cases, but you should check. The %d format now applies to integer output, use %f for reals. The translator attempts to do this for you by looking for possible numeric format expressions in all literal strings. Please check all strings in your script that contain % to make sure that the translator hasn't corrupted them.
- Plot The difference between printing and plotting is blurred in the windowed environment. Plotting calls are translated to printing calls and the plotter control calls are not supported. To export an image you can copy to the clipboard or save the view as a picture.
- Process There is a Process() command in the new script, but unlike the DOS script where the current view must be a time view to use Process, in the new script you must be in the result view (this is because there can be multiple result views simultaneously attached to a time view). To get round this we use ProcessQQ() to emulate the old situation.

- ScreenL The ScreenL and ScreenR keywords translate to XLow() and XHigh(). However, users can scroll result views to fractional bin positions. In the windowed version, do not assume that the value returned by XLow() or XHigh() in a result view is integral.
- Setxxx If you create result views, we call Proc CreateView(v%, vh%), where v% is the old view number (1 to 6) and vh% is the handle for the new view. This takes care of linking result views and time views together (in DOS you could have only one result view linked to a time view).
- SetAxes This is not translated. You can draw your own images using the XY views, but the mapping between the DOS drawing commands and the Windows ones is not simple and it is probably better done by hand.
- TITLE This is ignored in ON/OFF and treated as 1 in expressions.
- View To emulate the 6 views we add Proc ViewQQ() and a two dimensional array vqq%[7][2] to hold the view handles of the 6 views and the handle of any result view attached to a time view for ProcessQQ(). We also add a variable view% to hold the current view number (1-6). DOS scripts all assume that there is a time window, so if your script uses views, a call to ViewQQ(0) is made at the start to initialise the system. If the current view is already a time view it is used, otherwise you must name a file.

After translating Once Spike2 has translated your script, and you have checked that it runs we would suggest that you consider reworking it so that it no longer depends on the old architecture of 6 views. This usually means declaring variables at the top of your script to hold handles to the views you will use and removing the vqq%[][] array, the view% variable and Proc ViewQQ().

You should also examine the code for places where the script converter has been verbose. For example, DRAW becomes FrontView(view()); Draw(); to make sure that a drawn view is visible. Very often you can delete FrontView(view()); as this is only needed if the current view is hidden or covered by other windows.

Introduction

XY views have a wide variety of uses, from displaying user-defined graphs to drawing pictures. XY views have the following features:

- One x axis and one y axis shared between all data channels in the XY view, so all the channels share the same space, so you can overdraw one channel with another.
- Up to 256 data channels allowed in the view.
- Each data channel is a list of (x,y) data points. The number of data points in a channel is limited only by available memory and drawing time. However, you can limit the number of data points on a channel, in which case new data points replace the oldest data points.
- The data points can be drawn with markers at each data point. The range of marker styles currently includes: dots, boxes, plus signs, crosses, circles (Windows NT only), triangles, diamonds, horizontal lines and vertical lines. The size of the markers can also be set, and they can be made invisible.
- The data points can be joined with solid, dotted or dashed lines, and the line thickness can be varied. You can also choose to join the last point in a channel to the first point to make a loop.
- You can sort the order of the data points in a channel by x, by y or by order of insertion in the channel. This is only important if the data points are joined.
- The colour of the lines and markers can be chosen. If no colour is set, the same colour as for a waveform channel in a time view is set.

There are several example scripts included with Spike2 that illustrate some of the uses of XY views. You will find them in the scripts folder:

- FFTWater This draw a "waterfall" type display showing the variation in the power spectrum of a waveform channel with time.
- DE A script that uses a genetic algorithm to fit a double exponential.
- HRaster This shows a raster display and a PSTH of event data with the raster drawn above the PSTH.
- Clock An analogue clock for your Toolbar idle routine.
- PingPong A very silly example indeed.

Creating an XY view Although you can create an XY view from the File New menu command, the usual way to generate XY views is with the Analysis menu Measurements command or with the script language. The script language can generate an XY view for general use with the FileNew() command and as the target for measurements with a Process() command with the MeasureToXY() command. See the documentation for MeasureToXY() for an example.

The following assumes that you have some familiarity with the script language. An XY view always has at least one data channel, so when you create a view, you also create a channel. The following script code shows you how to make an XY view:

If you want to add additional channels you can do this using the XYSetChan() command. You can also use this command to set a channel to a particular state. The following sets channel 1 (the first channel) to show data points joined by lines with no limit on the number of data points, drawn in the standard colour:

XYSetChan(1, 0, 0, 1); 'chan 1, no size limit, no sort, joined XYDrawMode(1, 2, 0); 'set a marker size of 0 (invisible) To add data points to a channel you use the XYAddData() command. You can add single points, or pass an array of x and y co-ordinates. The following code adds three points to draw a triangle:

XYAddData(1, 0, 0); 'add a point to channel 1 at (0, 0) XYAddData(1, 1, 0); 'add a point at (1,0) XYAddData(1, 0.5, 1); 'add a point at (0.5, 1)

You will notice that the result of this draws only two sides of a triangle. We could complete the figure by adding an additional data point at (0,0), but it is just as easy to change the line joining mode to "Looped", and the figure is completed for you:

'set looped mode

XYJoin(1,2);

Overdrawing data

You can use the XY view to overdraw data. For example, suppose we have a waveform channel and another channel of trigger markers and we want to superimpose the first 100 data points after each trigger. The following example does this. To make it a little easier to see, we have displaced each triggered data sweep slightly to the right and up. We have also turned off the x and y axes.



The code required to produce this image is quite short and can be easily adapted to display other types of data. This example expects to run with the demo.smr data file in the Spike2\Data folder created when you installed the program. To run with other files you would need to adjust the channel numbers.

The first line declares variables to remember the time and XY window handles (identifiers) and a variable to hold the number of points read. The second line checks that the current view is a time view and gives up if it is not. The third line remembers the time view handle.

var tv%, wh%, np%; 'time view, XY view, points read if ViewKind() then Message("Not a time view");Halt endif; tv% := View(); 'save view handle

The next three lines declare an array to hold 100 data points, a variable to hold the trigger time (and we set it to -1 so we see the first event, even if it is at time 0), and we call the MakeWindow%() function (see below) that creates the XY window for us. We set the first argument to the time interval between data points in out time view so that the x spacing of the data points is correct. The second and third arguments to this function set the distance each "slice" of data is moved to the right and up.

var wave[100]; 'space for 100 waveform data points var t := -1; 'start time for trigger search wh% := MakeWindow%(BinSize(1), 0.01, 0.5); 'start waterfall view

To generate the data we run round a loop (repeat..until (t<0) or (t>20);) that finds the time of the next trigger on channel 2, reads in 100 data points into the array wave[] and then as long as the data read correctly, adds the data into the XY window.

```
repeat
View(tv%); 'Make time view current
t := NextTime(2, t); 'Get next trigger time from channel 2
np% := ChanData(1, wave[], t, MaxTime()); 'get channel 1 data
if (np% > 0) then AddSlice(wave[:np%]) endif; 'add to picture
until (t < 0) or (t>20); 'until end or we reach 20 seconds
```

The final task is to make the XY window visible and halt.

```
View(wh%).WindowVisible(1); 'make XY window visible
halt;
```

The remaining code can be copied from the FFTWater.s2s script. The variables with names starting WF are used to remember the initial waterfall settings. There are two functions in this code. MakeWindow%() prepares for the waterfall display by creating the XY window and storing the information needed to position the channels. AddSlice() takes an array of data points and adds it to the display as an additional channel.

```
var WFxInc,WFyInc,WFbinSz,WFSlices%,WFvh%;
'xBinSz Width of each bin (sample interval)
'xInc
       Add to each slice x co-ords to give waterfall effect
'vInc Add to each slice y co-ords to give waterfall effect
Func MakeWindow%(xBinSz, xInc, yInc)
WFSlices% := 0;
                          'no slices yet
WFxInc := xInc;
                          'save x increment per slice
WFyInc := yInc;
                          'save y increment per slice
                          'save data point separation
WFbinSz:= xBinSz;
WFvh% := FileNew(12);
                          'create a new XY window (hidden)
return WFvh%;
                          'return the XY window handle
end:
                           'Add data to the waterfall
Func AddSlice(y[])
                           'select the waterfall view
View(WFvh%);
var ch%:=1;
                           'true if this is the first channel
                          'if first channel no need to create
if WFSlices% = 0 then
  XYSetChan(1,-Len(y[]),0,1);
                                      'set original channel
else
  ch% := XYSetChan(0, -Len(y[]), 0, 1); 'create new channel
  if (ch% <= 0) then return ch% endif; 'No more channels
endif;
WFSlices% := ch%;
                           'number of slices
XYDrawMode(ch%,2,0);
                           'Hide the markers (set size of 0)
var x[Len(y[])];
                           'space for x values, same size as y[]
                           'generate x axis values, set the same
ArrConst(x[], WFbinSz);
x[0]:=(ch%-1)*WFxInc;
                          'set x offset as first value
                           'form the x positions
ArrIntgl(x[]);
ArrAdd(y[],(ch%-1)*WFyInc); 'add the y offset to the y array
                           'add the (x,y) data points
XYAddData(ch%, x[], y[]);
                           'return >0 means all OK
return 1;
end;
```

Introduction	It frequently happens that you have a set of data values $(x_1,y_1), (x_2,y_2) \dots (x_n,y_n)$ that you wish to test against a theoretical model $y = f(x, \mathbf{a}_0, \mathbf{a}_1, \mathbf{a}_2)$ where the \mathbf{a}_i are coefficients that are to be set to constant values which give the <i>best fit</i> of the model to the data values.
	For example, if we were looking at the extension of a spring (y) as it is loaded by weights (x) , we might wish to fit the straight line $y = \mathbf{a}_0 + \mathbf{a}_1 x$ to some measured data points so that we could measure a weight by the extension it caused. A careful experimenter might also wish to know what the probable error was in \mathbf{a}_0 and \mathbf{a}_1 so that the probable error in any weight deduced from an extension would be known. An even more cautious experimenter might want to know if the straight-line formula was likely to model the measured data.
	To avoid repeating definitions throughout the remainder of this chapter the following will be taken as defined. We apologise to the statisticians who may read the following and shudder.
mean	Given a set of <i>n</i> values y_j , the mean is $\sum_j y_j / n$ (the symbol \sum_j means form the sum over all indices _j of the expression that follows).
variance	If the mean of a set of <i>n</i> data values y_j is y_m , then the variance σ^2 (sigma squared) of this set of values is:
	$\sigma^{2} = \sum_{j} (y_{j} - y_{m})^{2} / n$ $\sigma^{2} = \sum_{j} (y_{j} - y_{m})^{2} / (n - 1)$ if y_{m} is known independently of the data values y_{j} if y_{m} is calculated from the data values y_{j}
	For a data set of any reasonable size, the use of $n-1$ or n in the denominator should make little difference.
standard deviation	The standard deviation σ (sigma) of a data set is the square root of the variance. Both the variance and the standard deviation are used as measures of the width of the distribution.
Normal distribution	If you measure a data value in any real system, there is always some error in the measurement. Once you have made a (very) large number of measurements, you can form a curve showing the probability of getting any particular value. One would hope that this error distribution would show a large peak at the "correct" value of the measurement and the width of this distribution would show the spread of likely errors.
	There is a particular error distribution that often occurs, called the Normal distribution. If a set of measurements is normally distributed, with a mean y_m and standard deviation σ , then the probability of measuring any particular value y is proportional to:
	$P(y) \propto \exp(-\frac{1}{2}(y-y_{\rm m})^2/\sigma^2)$
	It is for this distribution of errors that we have the well-known result that 68% of the values lie within one standard deviation of the mean, that 95% lie within two standard

deviations and that 99.7% lie within three standard deviations. Of course, if the error

distribution is not normal, these results do not apply.

Chi-squared The fitting routines given here define *best fit* as the values of \mathbf{a}_i (the coefficients) that minimise the chi-squared value (χ^2), defined as the sum over the measured data points of the square of the difference between the measured and predicted values divided by the variance of the data point:

 $\chi^2 = \Sigma_j \left(\left(\mathbf{y}_j - f(\mathbf{x}_j, \mathbf{a}_i) \right)^2 / \sigma_j^2 \right)$

where (x_i, y_i) is a data point and σ_i^2 is the variance of the measured data at that point.

If the sigma of each data point is unknown, then the fitting routines can be used to minimise $\Sigma_j (y_j - f(x_j, \mathbf{a}_i))^2$ which produces the same result as a chi-squared fit would produce if the variance of the errors at all the data points was the same. This is commonly called least-squares fitting (meaning that the fit minimises the sum of squares of the errors between the fitted function and the data).

Chi-squared fitting is also a maximum likelihood fit if the errors in the data points are normally distributed. This means that as well as minimising the chi-squared value, the fit also selects the most probable set of coefficients that model your data. If your data measurement errors are not normally distributed you can still use this method, but the fit is not maximum likelihood.

If your errors are normally distributed and if you know the variance(s) of the data points, you can form good estimates of the variance of the fitted coefficients, and you can also test if the function you have fitted is likely to model the data.

If your errors are normally distributed but you do not know the variance of the errors at the data points, you can make an estimate of the variance of the errors (based on the assumption that the variance is the same for them all and that the model does fit the data), by fitting your model and calculating the variance from the errors between the best fit and the data. Having done this, you cannot then use this variance to test if the fit is likely to model the data.

Residuals Once your fit is completed, it is a good idea to look at the graph of the errors between your original data and the fitted data (the residuals or residual errors). If your errors are normally distributed and are independent, you would expect this graph to be more or less horizontal with no obvious trends. If this is not the case, you should consider if the correct model function has been selected, or if the fitting function has found the true minimum.

Linear fit Non-linear fit

A linear fit is one in which the theoretical model $y = f(x, \mathbf{a}_0, \mathbf{a}_1, \mathbf{a}_2...)$ can be expressed as $y = \mathbf{a}_0 f_0(x) + \mathbf{a}_1 f_1(x) + \mathbf{a}_2 f_2(x) \dots$ for example $y = \mathbf{a}_0 + \mathbf{a}_1 x + \mathbf{a}_2 x^2$. Linear fits are relatively quick as they are done in one step. Usually, the only thing that can cause a problem is if the functions $f_i(x)$ are not linearly independent. The methods we use can usually detect this problem, and can still give a useful result.

A non-linear fit means all other cases, for example, $y = \mathbf{a}_0 \exp(-\mathbf{a}_1 x) + \mathbf{a}_2$. We solve these types of problem by making an initial guess at the coefficients (and ideally providing a range of values that the result is known to lie in) and then improving this guess. This process repeats until some criterion is met. Each repeat is called an *iteration*, so we call this an iterative process.

Covariance array	Several of the script fitting routines return a covariance array. If you have n coefficients, this array is of size n by n and is diagonally symmetric. If the errors in the original data points are normally distributed, the diagonal elements of this array are the variances in the values of the fitted coefficients. The remaining elements are the co-variances of pairs of the fitting parameters and can be used to estimate errors in derived values that depend on the product of two of the coefficients. If the errors are not normally distributed, the further away from normal the errors are, the less useful is the covariance array as a direct indication of the variances in the fitted coefficients.

For example, in the case of the linear fit $y = \mathbf{a}_0 + \mathbf{a}_1 x + \mathbf{a}_2 x^2$ you might collect your three coefficients in the array coef[], and the covariance in the array covar[][]. In this case, the \mathbf{a}_0 value is returned in coef[0] and its variance in covar[0][0], the \mathbf{a}_1 value is returned in coef[1] and its variance in covar[1][1], and the \mathbf{a}_2 value is returned in coef[2] and its variance in covar[2][2].

Because the array is diagonally symmetric, covar[i][j] is equal to covar[j][i] and the off-diagonal elements are the expected variance in the product of pairs of the coefficients, so covar[1][2] is the variance of a_1a_2 .

If you have not supplied the standard deviations of the errors in the data points, the covariance array is calculated on the assumption that all data points have a standard deviation of 1.0, and the covariance array is incorrectly scaled. In this case, if inspection of the residuals leads you to the conclusion that the function does indeed fit the data and that the errors are more or less the same for all values and not too far from normally distributed, then you can scale the covariance array to the correct values by multiplying all the elements of the array by the sum of squares of the errors between the data and the fitted values divided by the number of data points. If there are nD% data points and the sum of squares of the errors is errSq, then use ArrMul(cover[][], errSq/nD%); to rescale the covariance.

What does the covariance mean? Having fitted our data, we would like some idea of how the errors in the original data feed through to uncertainties in the values of the coefficients. The best way to do this is to obtain many sets of (x,y) data and fit our coefficient to each set. Then we can inspect the values of the coefficients and obtain a mean and standard deviation for each coefficient. However, this is very time consuming.

If the errors in the data are normally distributed (or not too far from this ideal case) and known, then the covariance array gives you some useful information. The square root of the covariance for a particular coefficient is the expected standard deviation in that value (given that the remaining coefficients remain fixed at optimum values). In script language terms, the standard deviation of coef[i] is sqrt(covar[i][i]).

In this case you would expect the coefficient to be within one standard deviation of the "correct" result 68% of the time, within 2 standard deviations 95% of the time and within 3 standard deviations 99.7% of the time.

Testing the fit If the errors in the original data are normally distributed and known (not calculated from the fit), and you know the χ^2 value for the fitted data, you can ask the question, "Given the known errors in the original data, how likely is it that you would get a value of χ^2 at least this large given that the data is correctly modelled by the fitting function plus normally distributed noise?" The answer is (at least in terms of the script language) that the probability is: GammaQ((nData% - nCoef%)/2.0, chiSq/2.0); where nData% is the number of data points to be fitted, nCoef% is the number of coefficients that were fitted and chiSq is the χ^2 value for the fit. GammaQ() is the incomplete Gamma function.

If you want to follow this result up in a statistical textbook, you should look up *chi-squared distribution for n degrees of freedom*. In our case, we have nData%-nCoef% degrees of freedom.

If the fit is reasonable, you should expect a probability value between 0.1 and 1 (but be a bit suspicious if you always get values close to 1.0, as you may have overestimated the errors in the data). If the wrong function has been fitted or if the fit is poor you usually get a very small probability.

Intermediate values (0.0001 to 0.1) may indicate that the errors in the original data were actually larger than you thought, or they may indicate that the model does not explain all the variation in the data.

Fitting routines In addition to the interactive Analysis menu Fit Data command, the following fitting related functions are available in the script language:

- **FitPoly()** This fits a polynomial of the form $y = \mathbf{a}_0 + \mathbf{a}_1 x + \mathbf{a}_2 x^2 + \mathbf{a}_3 x^3 \dots$
- FitLine() This fits a straight line to data in a time or result view.
- **FitLinear()** This fits $y = \mathbf{a}_0 f_0(x) + \mathbf{a}_1 f_1(x) + \mathbf{a}_2 f_2(x)...$
- **FitExp()** This fits multiple exponentials of the form $y = \mathbf{a}_0 \exp(-x/\mathbf{a}_1) + \mathbf{a}_2 \exp(-x/\mathbf{a}_3)$...
- **FitGauss()** This fits Gaussians of the form $y = \mathbf{a}_0 \exp(-\frac{1}{2}(x-\mathbf{a}_1)^2/\mathbf{a}_2^2) + \mathbf{a}_3 \exp(-\frac{1}{2}(x-\mathbf{a}_4)^2/\mathbf{a}_5^2) + \dots$
 - **FitSin()** This fits multiple sinusoids of the form $y = \mathbf{a}_0 \sin(\mathbf{a}_1 x + \mathbf{a}_2) + \mathbf{a}_3 \sin(\mathbf{a}_4 x + \mathbf{a}_5) + \dots$
- **FitNLUser()** This fits a user-defined function of the form $y = f(x, \mathbf{a}_0, \mathbf{a}_1, \mathbf{a}_2...)$ to a set of data points where the \mathbf{a}_i are constant coefficients to be determined. You must be able to calculate the differential of the function f with respect to each of the \mathbf{a}_i coefficients. This is the most general fitting routine, and also the slowest and most complex to use.

GammaP(), GammaQ() Use these functions to calculate the error function erf(x), the cumulative Poisson probability function and the Chi-squared probability function. These are very useful when considering probabilities of fits associated with the normally distributed data.

LnGamma() This calculates the natural logarithm of the Gamma function. It can be used to calculate the log of large factorials and is useful when working with probability distributions.

-Symbols-

 \land string literal escape character, 3-2 \$ string variable designator, 3-2 & reference parameter designator, 3-10 & in Windows dialog prompts, 5-48 {} optional syntax, 1-4 vertical bar, 1-4 = equality test operator, 3-6 + arithmetic add, 3-5 += add and assign, 3-4 += Append string, 3-4 + Concatenate strings, 3-6 := assignment, **3-4**, 3-6 () round brackets, 3-5 [] Array subscripts, 3-5 [start:size] array syntax, 3-3 : array range designator, 3-3 ; Statement separator, 3-4 < <= > >= arithmetic comparison, 3-5 < <= > >= String comparisons, 3-6 / arithmetic divide, 3-5 /= divide and assign, 3-4 = arithmetic equality, 3-5 = String equality, 3-6 % Integer variable designator, 3-2 * multiply, 3-5 *= multiply and assign, 3-4 <> arithmetic inequality, 3-5 <> String inequality, 3-6 ' Comment designator, 3-4 -= subtract and assign, 3-4 - Unary minus and subtract, 3-5

...

1902

Get revision, 5-35 script support, 5-31

—A—

1

Abort sampling, 5-149 Abs(), 5-1 Absolute value of expression or array, 5-1 Active cursors CursorActive(), 5-38 CursorActiveGet(), 5-39 CursorSearch(), 5-40 CursorValid(), 5-41 Add cursor, 5-40 Add to array, 5-1 Amplitude of power spectrum, 5-5 Amplitude of waveform, 5-27, 5-29 Analysis Command synopsis, 4-7 Create arbitrary result view, 5-168 Event correlation, 5-164 Interval histogram, 5-166 Number of items accumulated, 5-177 Peri-stimulus time histogram (PSTH), 5-167 Phase histogram, 5-166

Process all linked views, 5-139 Process time view data, 5-139 Shuffled event correlation, 5-165 Waveform average and accumulate, 5-164 Waveform correlation, 5-168 Waveform correlation DC removal, 5-168 and operator, 3-5 App(), 5-1 Application close, 5-73 Application window handle, 5-1 Arbitrary waveform output, 4-10 Arc tangent function, 5-10 Argument lists, 3-10 ArrAdd(), 5-1 Array arithmetic, 4-8 Arrays, 3-3 Absolute value, 5-1 Add constant or array, 5-1 Arc tangent function, 5-10 Command summary, 4-8 Copy, 5-2 Cosine of array, 5-37 Cubic splines, 5-8 Declaring, 3-4 Difference of two arrays, 5-8, 5-9 Differences between elements, 5-2 Division, 5-2, 5-3 Dot product, 5-3 Exponential function, 5-64 FFT analysis, 5-4 Fill with channel data, 5-17 Fill with channel list, 5-21 FIR filter. 5-6 Fractional part of real number, 5-97 Gain and phase, 5-4 Hyperbolic cosine of array, 5-37 Hyperbolic sine, 5-169 Hyperbolic tangent, 5-178 Integrate, 5-7 Inverse FFT, 5-5 Length of array, 5-110 Logarithm to base 10, 5-111 Logarithm to base e, 5-111 Maximum value, 5-115 Mean and standard deviation, 5-9 Minimum value, 5-125 Multiplication, 5-7 Natural logarithm, 5-111 Negate, 5-7 Passing to functions, 3-3 Power function, 5-136 Power spectrum, 5-4 Resample array, 5-8 Result view as an array, 3-4 Set to constant, 5-2 Sine of array elements, 5-169 Smoothing and filtering, 5-6 Sorting, 5-7 Square root of array elements, 5-171 Subtract array from value or array, 5-8 Subtract value or array from array, 5-9 Sum of product, 5-3 Sum of values, 5-9 Syntax, 3-3

Tangent of the array elements, 5-178 Total of array elements, 5-9 Truncate real array elements, 5-182 Use of [a:b] syntax, 3-3 ArrConst(), 5-2 ArrDiff(), 5-2 ArrDiv(), **5-2** ArrDivR(), 5-3 ArrDot(), 5-3 ArrFFT(), 5-4 ArrFilt(), 5-6 ArrIntgl(), 5-7 ArrMul(), 5-7 ArrSort(), 5-7 ArrSpline(), 5-8 ArrSub(), 5-8 ArrSubR(), 5-9 ArrSum(), 5-9 Asc(), 5-9 ASCII code of character, 5-9 ASCII to string conversion, 5-30 Assignment, 3-4, 3-6 ATan(), 5-10 Auto-correlation of waveform, 5-168 Automatic file save at sample end, 5-150 Automatic file saving folder, 5-72 Automatic processing during sampling, 5-140 Automatic scrolling, 5-190 Average waveform data, 5-164 Axis controls Drawing colour, 5-31 XY view data tracking, 5-192 Axis grid display control, 5-98

<u>–B–</u>

band binary and operator, 3-5 Beep or tone output, 5-169 Bin access in result view, 3-4 Bin number to x axis unit conversion, 5-11 Binary file commands, 4-5 Binary files Close, 5-66 Little or big endian, 5-12 Move current position, 5-12 Open, 5-71 Read data, 5-11, 5-12 Write data, 5-14, 5-15 BinError(), 5-10 Binsize(), 5-11 BinToX(), 5-11 Bitmap output To clipboard, 5-61 To file, 5-74 Black and White display, 5-186 Boltzmann sigmoid, 5-83 bor binary logical or operator, 3-5 BOX DOS script command, 6-1 BRead(), 5-11 BReadSize(), 5-12 Break points, 2-1 Clear all, 2-1 Breaking out of a script, 2-2

BRWEndian(), **5-12** BSeek(), **5-12** Burst analysis Burst statistics, 5-14 Creating bursts from events, 5-13 Revising bursts, 5-13 BurstMake(), 5-13 BurstRevise(), 5-13 BurstStats(), 5-14 Buttons. see toolbar and Interact() BWrite(), **5-14** BWriteSize(), **5-15** bxor binary exclusive or operator, 3-5

—C—

Calibrate waveform, 5-15, 5-23, 5-27 Call stack in debug, 2-4 case. see docase CED 1902 signal conditioner, 5-31 Chan\$(), 5-15 ChanCalibrate(), 5-15 ChanColour(), 5-16 ChanComment\$(), 5-16 ChanData(), 5-17 ChanDelete(), 5-17 ChanDuplicate(), 5-17 ChanFit(), 5-18 ChanFitCoef(), 5-19 ChanFitShow(), 5-20 ChanFitValue(), 5-20 ChanHide(), 5-20 ChanKind(), 5-20 ChanList(), 5-21 ChanMeasure(), 5-21 Channel Add new XY view channel, 5-195 Attach horizontal cursor, 5-99 Comment, 5-16 Copy data from XY view, 5-193 Copy into an array, 5-17 Create, 5-22 Delete, 5-17 Duplicate, 5-17 Event count, 5-37 Find duplicate, 5-60 Generate channel list, 5-21 Get or set title. 5-28 Get or set units, 5-28 Get selected state, 5-28 Get visible state, 5-29 Groups, 5-23 Hide, 5-20 Maximum time of item in channel, 5-115 Measure region, 5-21 Memory based channels, 5-120 Minimum, maximum and positions, 5-125 Modify XY view channel settings, 5-195 Number in data file, 5-151 Number show and hide, 5-23 Order, 5-23 Selecting, 5-28 Selecting in a dialog, 5-51 Set colour, 5-16 Show, 5-28

Time of next item, 5-128 Time of previous item, 5-110 Type of a channel, 5-20 Value at given position, 5-29 Vertical space, 5-30 Weight, 5-30 Write wave data, 5-30 Channel commands, 4-2 Channel number Drawing colour, 5-31 Show and hide, 5-23 Channel specifier, 3-12, 5-164 ChanNew(), 5-22 ChanNumbers(), 5-23 ChanOffset(), 5-23 ChanOrder(), 5-23 ChanPort(), 5-24 ChanProcessAdd(), 5-24 ChanProcessArg(), 5-24 ChanProcessClear(), 5-25 ChanProcessInfo(), 5-25 ChanSave(), 5-25 ChanScale(), 5-27 ChanSearch(), 5-27 ChanSelect(), 5-28 ChanShow(), 5-28 ChanTitle\$(), 5-28 ChanUnits\$(), 5-28 ChanValue(), 5-29 ChanVisible(), 5-29 ChanWeight(), 5-30 ChanWriteWave(), 5-30 Character code, 5-30 Character code (ASCII), 5-9 Check box in a dialog, 5-52 Chi-squared probability function, 5-98 Chi-squared value, 8-2 Chr\$(), 5-30 CLEAR DOS script command, 6-2 Clipboard Copy and Cut data, 5-61 Cut current selection to the clipboard, 5-Paste, 5-61 Close Spike2 application, 5-73 Close window, 5-66 Coefficients, 5-6 Colon, array range designator, 3-3 Colour palette, 5-129 Colour(), 5-31 Colours of screen items, 5-31 Data channels, 5-16 Force Black and White, 5-186 View colours, 5-184 XY view, 5-192 Column number in text view, 5-128 Command line, 5-142 Comment File comment automatic prompting, 5-149 Get and set channel comment, 5-16 Get and set file comment, 5-67 In script language, 3-4 Commit data file, 5-149 Comparison operators, 3-5

Compile script, 2-1 CondFeature (), 5-31 CondFilter(), 5-32 CondFilterList(), 5-32 CondFilterType(), 5-33 CondGain(), 5-33 CondGainList(), 5-33 CondGet(), 5-34 Conditioner commands, 5-31 CondOffset(), 5-34 CondOffsetLimit(), 5-35 CondRevision\$(), 5-35 CondSet(), 5-35 CondSourceList(), 5-36 CondType(), 5-36 Configuration files Load from script, 5-71 Save from script, 5-74 Suppress use of, 5-70 Connections Port for WaveMark data, 5-159 Setting port for event channel, 5-152 Setting port for waveform channel, 5-159 Constant delarations, 3-5 Convert A number to a string, 5-177 A string to a number, 5-183 A string to upper case, 5-182 Between channel types with MemImport(), 5-122 Event to waveform, 5-62, 5-186 Number to a character, 5-30 Parse string into variables, 5-148 Parse string setup, 5-148 Result view bin to x axis units, 5-11 Seconds to Spike2 time units, 5-191 String to lower case, 5-110 X axis units to bins, 5-191 Convert foreign file format, 5-67 Converting between data types, 3-2 Copy Array or result view to another, 5-2 External file, 5-68 The current selection to the clipboard, 5-61 Waveform to result view, 5-164 Correlation of events, 5-164 Cos(), 5-37 Cosh(), 5-37 Cosine of expression, 5-37 Count points in XY circle, 5-194 points in XY rectangle, 5-194 Count events in time range, 5-37 Count(), 5-37 Covariance array, 8-3 Create new permanent channel, 5-123 Create new result view, 5-164 Create result view, 5-168 CreateView() DOS emulation function, 6-4 Creating new channels In memory, 5-120 On disk, 5-22

Cross-correlation of events, 5-164 Cross-correlation of waveforms, 5-168 cSpc channel specifier, 3-12, 5-164 Cumulative Poisson probability function, 5-97 Curly brackets, 1-4 Current directory, 5-72 Current view, 1-2, 2-3, 5-183 Cursor commands, 4-4 Cursor(), 5-37 CursorActive(), 5-38 CursorActiveGet(), 5-39 CursorDelete(), 5-39 CursorExists(), 5-39 CursorLabel(), 5-39 CursorLabelPos(), 5-40 CursorNew(), 5-40 CursorRenumber(), 5-40 Cursors Channel for horizontal cursor, 5-99 Create new cursor, 5-40 Delete, 5-39 Delete horizontal cursor, 5-99 Drawing colour, 5-31 Get and set horizontal position, 5-98 Horizontal cursor label position, 5-100 Horizontal cursor label style, 5-99 Label position, 5-40 Label style, 5-39 New horizontal cursor, 5-100 Position. 5-37 Renumber cursors, 5-40 Renumber horizontal cursors, 5-100 Set number, 5-41 User-defined style, 5-39, 5-99 CursorSearch(), 5-40 CursorSet(), 5-41 CursorValid(), 5-41 CursorVisible(), 5-41 Curve fitting Exponentials, 5-82, 5-84 Gaussian, 5-83, 5-86 Linear, 5-88 Non-linear, 5-90 Polynomial, 5-83, 5-92 Sigmoid, 5-83 Sinusoid, 5-83 Sinusoids, 5-94 Cut The current selection to the clipboard, 5-61 CyberAmp Get revision, 5-35 script support, 5-31

—D–

DAC output during sampling, 4-10 Data types, 3-1 Compatibility, 3-2 Date as a string, 5-42 as numbers, 5-179 data file creation, 5-68 Date\$(), **5-42** DC removal in waveform correlation, 5-168 Debug Call stack, 2-4 Enter debug on error, 2-3 Enter debugger, 2-2, 5-42 Globals window, 2-4 Locals window, 2-4 Operations, 4-11 Script commands, 4-11 Debug script, 2-2 Debug(), 2-2, 5-42 DebugList(), 5-43 DebugOpts(), 5-43 Delay in script. See Yield() Delete Channel, 5-17 Cursor, 5-39 File, 5-68 Horizontal cursor, 5-99 Memory channel, 5-120 Memory channel time range, 5-121 Selection, 5-60 Substring, 5-43 XY view data, 5-192 DelStr\$(), 5-43 Determinant of a matrix, 5-114 diag() operator, 3-3 Diagonal of a matrix, 3-3 Dialogs Buttons, 5-50 Check box, 5-52 Create new dialog, 5-52 Dialog units, 5-48 Display and collect responses, 5-55 Enable and disable items, 5-52 Get item value, 5-56 Get time or x value, 5-57 Group boxes, 5-53 Integer number input, 5-53 Real number input, 5-55 Selecting a channel, 5-51 Selecting one value from a list, 5-54 Set label, 5-54, 5-56 Show or hide items, 5-57 Simple format, 5-47 Text string input, 5-55 User actions and call-backs, 5-49 User-defined, 5-47 Differences between array elements, 5-2 Digital filter bank Apply from script, 5-76 Create digital filter from script, 5-78 Filter name from script, 5-79 Force filter calculation, 5-77 Get filter bank information, 5-78 Sampling frequency range, 5-79 Set attenuation of filter bank filter, 5-76 Set comment in filter bank, 5-78 Digital filtering, 4-9 Directory for files, 5-72 DiscrimChanGet(), 5-44 DiscrimChanSet(), 5-45 DiscrimClear(), 5-46 DiscrimLevel(), 5-46

DiscrimMode(), 5-46 DiscrimMonitor(), 5-47 DiscrimTimeOut(), 5-47 Display channel, 5-28 Divide operator with integer arguments, 3-5 Division of arrays, 5-2, 5-3 DlgAllow(), 5-49 DlgButton(), 5-50 DlgChan(), 5-51 DlgCheck(), 5-52 DlgCreate(), 5-52 DlgEnable(), 5-52 DlgGroup(), 5-53 DlgInteger(), 5-53 DlgLabel(), 5-54 DlgList(), 5-54 DlgReal(), 5-55 DlgShow(), 5-55 DlgString(), 5-55 DlgText(), 5-56 DlgValue(), 5-56 DlgVisible(), 5-57 DlgXValue(), 5-57 docase, 3-7 Dockable toolbars, 5-187 Change floating size, 5-188 Change state, 5-189 Get docked state, 5-188 Get handle, 5-1 DOS command line, 5-142 Dot product of arrays, 5-3 Draw a view, 5-57 Draw(), 5-57 DrawAll(), 5-58 Drawing modes Join XY view points, 5-194 XY view, 5-192 DrawMode(), 5-58 DRAWR DOS script command, 6-1 DRAWTO DOS script command, 6-1 DUMPM DOS script command, 6-1 Dup(), 5-60 DupChan(), 5-60 Duplicate Channel, 5-17 Channel, find, 5-60 View, 5-188 View, get handle, 5-60

-E—

e, Mathematical constant, 3-6
Edit marker codes, 5-112
Edit memory channel, 5-124
Edit toolbar handle, 5-1
EditClear(), 5-60
EditCopy(), 5-61
EditCut(), 5-61
Editing script commands, 4-5
EditPaste(), 5-61
EditSelectAll(), 5-61
else. see if and docase. see if and docase

end, 3-10 endcase. see docase endif.see if Environment variables, 5-177 Environmental functions, 4-11 Erase file, 5-68 Error bars DrawMode(), SD, SEM, 5-59 Script access, 5-10 SetAverage(), 5-164 SetResult(), 5-168 Error codes as text, 5-62 Error function, erf(x), 5-97 Error\$(), 5-62 Errors Enter debug on error, 2-3 Esc key in a script, 2-2 Escape character backslash, 3-2 ESCAPE DOS script command, 6-2 Eval(), 5-62 Evaluate argument, 5-62 Event correlation, 5-164 Event data Burst formation, 5-13 Convert to waveform, 5-62 Drawing colour, 5-31 Read times into an array, 5-17 Event discriminator Script commands, 4-10 EVENT DOS script command, 6-3 Events in time range, 5-37 EventToWaveform(), 5-62 EXECUTE DOS script command, 6-2 Execute program, 5-142 EXECUTED DOS script command, 6-2 Exit from Spike2, 5-73 Exp(), 5-64 Exponential fitting, 5-82, 5-84 Exponential function, 5-64 Export data file, 5-74 ExportChanFormat(), 5-64 ExportChanList(), 5-65 ExportRectFormat(), 5-65 ExportTextFormat(), 5-66 Expressions, 3-5 External files, 5-71 External program Kill, 5-142 Run, 5-142 Status of, 5-143 Extract fields from a string, 5-148

—F—

Feature search, 5-27 Active cursors, 5-38 Start cursor search, 5-40 Test search, 5-41 FFT analysis Of arrays, 5-4 Of waveform data, 5-167 File Apply resource, 5-66 Automatic file commit, 5-149

Index-4

Automatic save at sample end, 5-150 Comment, 5-67 Comment automatic prompting, 5-149 Convert foreign format, 5-67 Copy external file, 5-68 Delete list of files, 5-68 Global resource file, 5-68 List of files, 5-69 Name, 5-69 Name for automatic filing, 5-150 New window, 5-70 Open window or external file, 5-71 Parent and child directories/folders, 5-69 Save resource, 5-75 FILE DOS script command, 6-3 File system commands, 4-5 FileApplyResource(), 5-66 FileClose(), 5-66 FileComment\$(), 5-67 FileConvert\$(), 5-67 FileCopy(), 5-68 FileDate\$(), 5-68 FileDelete(), 5-68 FileGlobalResource (), 5-68 FileList(), 5-69 FileName\$(), 5-69 FileNew(), 5-70 XY view, 7-1 FileOpen(), 5-71 FilePath\$(), 5-72 FilePathSet(). 5-72 FilePrint(), 5-72 FilePrintScreen(), 5-73 FilePrintVisible(), 5-73 FileQuit(), 5-73 FileSave(), 5-73 FileSaveAs(), 5-74 FileSaveResource(), 5-75 FileTime\$(), 5-75 FileTimeDate(), 5-75 FiltApply(), 5-76 FiltAtten(), 5-76 FiltCalc(), 5-77 FiltComment\$(), 5-78 FiltCreate(), 5-78 Filter bank. See Digital filter bank Filter coefficients, 5-6 Filter marker data, 5-112 Filtering of arrays, 5-6 FiltInfo(), 5-78 FiltName\$(), 5-79 FiltRange(), 5-79 Find a view, 5-184 Find feature, 5-27 Active cursors, 5-38 Start cursor search, 5-40 Test search, 5-41 Find func or proc, 2-1 FIR filter Apply from script, 5-76 Frequency response, 5-81 Make coefficients, 5-79 Make coefficients (simplified), 5-80 Script functions, 4-9

FIR filter of array, 5-6 FIRMake(), 5-79 FIRQuick(), 5-80 FIRResponse(), 5-81 FitCoef(), 5-81 FitData(), 5-82 FitExp(), 5-84 FitGauss(), 5-86 FitLine(), 5-88 FitLinear(), 5-88 FitNLUser(), 5-90 FitPoly(), 5-92 FitSin(), 5-94 Fitting routines, 8-4 FitValue(), 5-96 Flow of control statements, 3-7 FocusHandle(), 5-96 Folder for files, 5-72 Font Get font characteristics, 5-96 Set font, 5-96 FontGet(), 5-96 FontSet(), 5-96 for, 3-8 Formatted text output, 5-137, 5-138 Frac(), 5-97 Fractional part of real number or array, 5-97 FRAME DOS script command, 6-3 FREEMEM DOS script command, 6-1 FrontView(), 5-97 func, 3-10 Function argument lists, 3-10 Functions and procedures, 3-10 Functions as arguments, 3-12

—G—

Gain, 5-27 GammaP(), 5-97 Gaussian fitting, 5-83, 5-86 GLOBAL DOS script command, 6-3 GLOBAL DOS script command, 6-2 Global resource files, 5-68 Globals window, 2-4 Grid colour, 5-31 Grid(), **5-98** Groups of channels, 5-23 Gutter(), **5-98**

—H—

Halt, **3-9** HCursor(), **5-98** HCursorChan(), **5-99** HCursorDelete(), **5-99** HCursorExists(), **5-99** HCursorLabel(), **5-99** HCursorLabelPos(), **5-100** HCursorNew(), **5-100** HCursorRenumber(), **5-100** HELP DOS script command, 6-3 Help(), **5-100** Hexadecimal number format, 3-2, 5-138
Hide a channel, 5-20 Hide window, 5-189 Hide X axis scroll bar, 5-191 Hide y axis, 5-197 Horizontal cursor. see Cursor Horizontal cursor commands, 4-4 Host operating system, 5-177, 5-178 Hyperbolic cosine, 5-37 Hyperbolic tangent, 5-178

—I—

if statement, 3-7 IIR filter Overview, 5-101 Script commands, 4-9 IIRBp(), 5-104 IIRBs(), 5-104 IIRHp(), 5-105 IIRLp(), 5-106 IIRNotch(), 5-106 IIRReson(), 5-107 Import data into memory channel, 5-122 Import foreign data file, 5-67 Impulse response, 5-6 Indent text, 2-2 Inkey(), 5-107 Input a single number, 5-108 Input a string, 5-108 Input from the keyboard, 5-107 Input\$(), 5-108 Input(), 5-108 Insert data into memory channel, 5-124 Instantaneous frequency, 5-58, 5-62 Drawing colour, 5-31 InStr(), 5-108 Integer data type, 3-2 Integrate array, 5-7 Interact(), 5-109 Interval histogram (INTH), 5-166 Inverse FFT, 5-5 Invert matrix, 5-114

-J-

Join data points in XY view, 5-194

Key window, 5-28, 5-185 Control of, 5-195 Keyboard input, 5-107 Interact(), 5-109 ToolbarSet(), 5-181 Keyboard test, 5-109 Keypress(), **5-109** Keywords, 3-1 Kind of channel, 5-20

_T _

Label position of cursor, 5-40 Label style of cursor, 5-39 LastTime(), **5-110** LCase\$(), **5-110** Least-squares linear fit, 5-88 Left\$(), 5-110 Legal characters in string input, 5-108 Len(), 3-10, 5-110 Length of array or string, 5-110 Line length in script, 3-1 Line number in text window, 5-128 Linear fit, 8-2 Linear fitting, 5-88 Linear least-squares fit, 5-88 List of channels. 5-21 List of files, 5-69 List of views, 5-185 literal string delimiter, 3-2 Ln(), 5-111 LnGamma(), 5-111 Locals window, 2-4 Log amplitude of the power spectrum in dB, 5-5 Log window, 5-111 Log(), 5-111 Logarithm to base 10, 5-111 Logarithm to base e, 5-111 LogHandle(), 5-111 Lower case version of a string, 5-110

-M---

MarkEdit(), 5-112 Marker codes, 5-112 Edit codes, 5-112 Marker data Set codes, 5-112 MarkInfo(), 5-114 MarkMask(), 5-112 MarkSet(), 5-112 Mask for marker data, 5-112 MATDet(), 5-114 Mathematical constants, 3-6 Mathematical functions, 4-7 MATInv(), 5-114 MATMul(), 5-114 Matrix, 3-3 Determinant of, 5-114 Diagonal of, 3-3 Inverse of, 5-114 Multiplication, 5-114 Solve linear equations, 5-115 Transpose of, 3-3, 5-115 Matrix arithmetic, 4-8 MATSolve(), 5-115 MATTrans(), 5-115 Max(), 5-115 Maximum and minimum of XY channel, 5-195 Maximum of channel or result view, 5-125 Maximum value, 5-115 Maxtime(), 5-115 mean, 8-1 Mean frequency, 5-58, 5-62 Drawing colour, 5-31 Mean of array, 5-9 MeasureChan(), 5-116 Measurements MeasureChan(), 5-116

MeasureToXY(), 5-118 MeasureX(), 5-119 MeasureY(), 5-119 MeasureToChan(), 5-117 MeasureToXY(), 5-118 MeasureX(), 5-119 MemChan(), 5-120 MemDeleteItem(), 5-120 MemDeleteTime(), 5-121 MemGetItem(), 5-121 MemImport(), 5-122 Memory channels Add or edit data, 5-124 Creating, 5-120 Delete, 5-120 Delete time range, 5-121 Get item data, 5-121 Import data, 5-122 Write to file, 5-123 MemSave(), 5-123 MemSetItem(), 5-124 Message(), 5-124 Metafile output To clipboard, 5-61 Mid\$(), 5-124 Min(), 5-125 Minimum of channel or result view, 5-125 Minimum value, 5-125 Minmax(), 5-125 MMAudio(), 5-126 MMImage(), 5-126 MMOpen(), 5-127 MMPosition(), 5-127 MMRate(), 5-127 MMVideo(), 5-128 mod remainder operator, 3-5 Monochrome display, 5-186 Move in text and cursor windows, 5-128 MoveBy(), 5-128 MOVER DOS script command, 6-1 MoveTo(), 5-128 Multimedia sound output, 5-169 Speech, 5-170 Multimedia window Audio information, 5-126 Open window, 5-127 RGB data access, 5-126 Set and get play position, 5-127 Multimedia window Set and get play state, 5-127 Multiple monitor support, 4-12, 5-178, 5-189 Multiplication arrays, 5-7 Matrices, 5-114

---N----Name format, 3-1 Natural logarithm, 5-111 Negate array, 5-7 neqq% DOS emulation variable, 6-3 New cursor, 5-40 New file, 5-70 New horizontal cursor, 5-100 New result view, 5-164 NEWEVENT DOS script command, 6-3 next, **3-8** Non-linear fit, 8-2 Non-linear fitting, 5-90 Normal distribution, 8-1 NORMAL DOS script command, 6-3 Normal settings for a view, 5-185 not logical operator, 3-5 NUMARR DOS script command, 6-1 Numeric input, 5-108 NUMSTR DOS script command, 6-1 NUMVAR DOS script command, 6-1

-0-Offset waveform, 5-23 Open file, 5-71 Operating system, 5-177, 5-178 Operators, 3-5 Order of precedence, 3-5 Precedence Spike2 for DOS, 6-2 Optimise the display, 5-129 Optimise(), 5-129 or operator, 3-5 Order of channels, 5-23 Oscilloscope style triggered display, 5-185 Outdent text, 2-2 Output sequencer Get current step, 5-155 Get file name, 5-156 Keyboard link control, 5-155 Set file name, 5-155 Set variable, 5-156 Table access, 5-155 Overdraw data, 7-2 OVERDRAW DOS script command, 6-1 Overwrite wave data, 5-30

-P-

Palette for colour, 5-129 PaletteGet(), 5-129 PaletteSet(), 5-129 PARAMSTR DOS script command, 6-1 Passing arguments by reference, 3-10 by value, 3-10 functions and procedures, 3-12 Paste from clipboard, 5-61 Path for file operations, 5-72 PCA(), 5-130 Peak to Peak value, 5-21 Peri-stimulus time histogram (PSTH), 5-167 Phase display (FFT) in result view, 5-4 Phase histogram, 5-166 Phase of power spectrum, 5-5 pi, Mathematical constant, 3-6, 5-178 Play wave bar handle, 5-1 Play wave during output, 4-10 PlayOffline(), 5-130 PlayWaveAdd(), 5-131 PlayWaveChans(), 5-133

PlayWaveCopy(), 5-133 PlayWaveCycles(), 5-133 PlayWaveDelete(), 5-134 PlayWaveEnable(), 5-134 PlayWaveInfo\$(), 5-134 PlayWaveLabel\$(), 5-134 PlayWaveLink\$(), 5-135 PlayWaveRate(), 5-135 PlayWaveSpeed(), 5-135 PlayWaveStatus\$(), 5-136 PlayWaveStop(), 5-136 PlayWaveTrigger(), 5-136 PLOTTO DOS script command, 6-2 Polynomial fitting, 5-83, 5-92 Port For WaveMark data, 5-159 Setting for event channel, 5-152 Setting for waveform channel, 5-159 Position of cursor, 5-37 Position of window, 5-187, 5-188 Pow(), 5-136 Power function, 5-136 Power spectra Of arrays, 5-4 Of waveform channels, 5-167 pqq% DOS emulation variable, 6-3 Precedence of operators, 3-5 Spike2 for DOS, 6-2 Preferences Script access, 5-140 Principal Component Analysis, 5-130 Print All views on screen, 5-73 Formatted text output, 5-137, 5-138 Print visible region, 5-73 Range of data, 5-72 To log window, 5-138 To string, 5-138 Print\$(), 5-138 Print(), 5-137 PrintLog(), 5-138 PrintTo() DOS emulation function, 6-3 proc, 3-10 Procedures as arguments, 3-12 PROCESS DOS script command, 6-3 Process time view data to result view, 5-139 Process(), 5-139 ProcessAll(), 5-139 ProcessAuto(), 5-140 ProcessQQ() DOS emulation function, 6-3 ProcessTriggered(), 5-140 Profile(), 5-140 ProgKill(), 5-142 ProgRun(), 5-142 ProgStatus(), 5-143

-0-

QQ suffix for DOS emulation, 6-3 Query(), **5-143** Quit Spike2, 5-73

—R— Radians, 5-10, 5-37, 5-169 Convert to degrees, 5-178 Rand(), 5-143 RandExp(), 5-143 RandNorm(), 5-144 Random number generator, 5-143 Exponential distribution, 5-143 Normal distribution, 5-144 Range of data points in XY view, 5-195 Raster drawing mode, 5-58 Drawing colour, 5-31 in a result view, 5-164, 5-166, 5-167 RasterAux(), 5-144 RasterGet(), 5-145 RasterSet(), 5-145 RasterSort(), 5-146 RasterSymbol(), 5-146 Rate drawing mode Drawing colour, 5-31 Set from script, 5-58 Read binary data, 5-11, 5-12 Read channel into an array, 5-17 Read text file Input from a text file into variable(s), 5-147 Open file, 5-71 Read(), 5-147 ReadSetup(), 5-148 ReadStr(), 5-148 Real data type, 3-1 RealMark data Get information, 5-114 RealWave data Write to channel, 5-30 Reciprocal of array, 5-3 Reference parameters, 3-10 Registry access, 5-140 Renumber cursors, 5-40 Renumber horizontal cursors, 5-100 repeat. 3-8 ReRun a file, 5-148 ReRun(), 5-148 Residuals, 8-2 Resource information suppression, 5-71 Result view Access to contents, 3-4 Array access, 3-3, 3-4 Bin width, 5-11 Convert x axis units to bin number, 5-191 Create user-defined view, 5-168 Drawing colours, 5-31 Error bars, 5-10 Get or set units, 5-28 Minimum, maximum and positions, 5-125 Number of items accumulated, 5-177 Open file from script, 5-71 Printing, 5-72 Process data, 5-139 Script commands, 4-3 Sum of bins, 5-37 Value at given x axis position, 5-29 return, 3-10 Right\$(), 5-149

Rightmost characters from a string, 5-149 RMS Amplitude, 5-21 Calibration method, 5-16 Channel process option, 5-24 Trend plot, 5-119 Round a real to nearest whole number, 5-149 Round(), **5-149** RS232 script commands, 4-12 Run external program, 5-142 Run script, 2-1

-S---

Sample bar handle, 5-1 Sample control panel handle, 5-1 Sample toolbar Control from script language, 5-150 SampleAbort(), 5-149 SampleAutoComment(), 5-149 SampleAutoCommit(), 5-149 SampleAutoFile(), 5-150 SampleAutoName\$(), 5-150 SampleBar(), 5-150 SampleCalibrate(), 5-150 SampleChannels(), 5-151 SampleClear(), 5-151 SampleComment\$(), 5-151 SampleDigMark(), 5-151 SampleEvent(), 5-152 SampleHandle(), 5-152 SampleKey(), 5-152 SampleLimitSize(), 5-152 SampleLimitTime(), 5-152 SampleMode(), 5-153 SampleOptimise(), 5-153 SampleRepeats(), 5-154 SampleReset(), 5-154 SampleSeqCtrl(), 5-155 SampleSeqStep (), 5-155 SampleSeqTable (), 5-155 SampleSequencer\$(), 5-156 SampleSequencer(), 5-155 SampleSeqVar(), 5-156 SampleStart(), 5-156 SampleStartTrigger(), 5-156 SampleStatus(), 5-157 SampleStop(), 5-157 SampleText(), 5-157 SampleTextMark(), 5-157 SampleTimePerAdc(), 5-158 SampleTitle\$(), 5-158 SampleTrigger(), 5-158 SampleUsPerTime(), 5-159 SampleWaveform(), 5-159 SampleWaveMark(), 5-159 SampleWrite(), 5-160 Sampling Automatic processing, 5-140 Naming data file from a script, 5-70, 5-74 Runtime control functions, 4-9 Setting where data is stored during sampling, 5-72 Triggered start, 5-156 View handles, 5-152, 5-154 Sampling configuration

Channel comment, 5-151 Event channel, 5-152 Functions, 4-9 Limit file size, 5-152 Limit sample time, 5-152 Number of channels, 5-151 Number of repeats, 5-154 Optimise rate settings, 5-153 Reset configuration, 5-151 Sample mode, 5-153 Save. 5-74 Triggered start, 5-156 Scale waveform, 5-27 Scope of variables and user-defined functions, 3-11 Screen dump to printer, 5-73 SCREENL and SCREENH DOS script command, 6-4 Script Call stack, 2-4 Clear all break points, 2-1 Compile, 2-1 Debug, 2-2 Enter debug on error, 2-3 Find func or proc, 2-1 Handle of running script, 5-1 Run, 2-1 Set and clear break points, 2-1 Split text window, 2-1 Script Bar Control from script language, 5-160 handle, 5-1 ScriptBar(), 5-160 ScriptRun(), 5-161 Scroll bar, show and hide, 5-191 Scroll display, 5-57 Search data Active cursors, 5-38 For feature, 5-27 Start cursor search, 5-40 Test search, 5-41 Search for view, 5-184 Seconds(), 5-161 Select a channel, 5-28 Select all copyable items, 5-61 Selection\$(), 5-161 Semicolon, statement separator, 3-4 Sequencer control panel handle, 5-1 Serial line script commands, 4-12 Serial number Read, 5-1 SerialClose(), 5-161 SerialCount(), 5-161 SerialOpen(), 5-162 SerialRead(), 5-162 SerialWrite(), 5-163 SetAverage(), 5-164 SETAXES DOS script command, 6-1 SetEvtCrl(), 5-164 SetEvtCrlShift(), 5-165 SetINTH(), 5-166 SetPhase(), 5-166 SETPLACE DOS script command, 6-2 SetPower(), 5-167

SetPSTH(), 5-167 SetResult(), 5-168 SetWaveCrl(), 5-168 SetWaveCrlDC(), 5-168 Show channel (list), 5-28 Show window, 5-189 Show X axis scroll bar, 5-191 Show y axis, 5-197 SHOWVARS DOS script command, 6-1 Shuffled event correlation, 5-165 Sigmoid fitting, 5-83 Signal conditioner, 5-31 Get and set gain, 5-33 Get and set offset, 5-34 Get and set special features, 5-31 Get list of gains, 5-33 Get list of sources, 5-36 Get offset range, 5-35 Get revision, 5-35 Get type, 5-36 List filter frequencies, 5-32 List filter types, 5-33 Low-pass and high-pass filters, 5-32 Read all port settings, 5-34 Script commands, 4-10 Set all parameters, 5-35 Sin(), 5-169 Sine of an angle in radians, 5-169 Single step a script, 2-3 Sinh(), 5-169 Sinusoid fitting, 5-83 Sinusoidal fitting, 5-94 Size of window, 5-188 SMControl(), 5-171 Smoothed frequency, 5-62 Smoothing of arrays, 5-6 SMOpen(), 5-171 Solid colour, 5-129 Solve linear equations, 5-115 Sonogram drawing mode, 5-58 Sort arrays, 5-7 Sort channels, 5-23 Sound output, 5-169 Sound(), 5-169 Spawn program, 5-142 Speak(), 5-170 Speech output, 5-170 Spike monitor window control state, 5-171 Open and get handle, 5-171 Spike shape window Button states, 5-172 Create new spike channel, 5-172 Create, merge or replace template, 5-176 Delete templates, 5-174 Display range, 5-200 Duplicate channels, 5-17 Get and set template information, 5-175 Get template and display size, 5-176 Get template data, 5-174 Open dialog, 5-173 Optimise display, 5-129 Parameters, 5-173 Reclassify channel, 5-172

Run state, 5-174 Script function list, 4-11 Set channel, 5-172 Set template and display size, 5-177 Set trigger levels, 5-98 Splitter control, 2-1 Sqrt(), 5-171 Square root, 5-171 SSButton(), 5-172 SSChan(), 5-172 SSClassify(), 5-172 SSOpen(), 5-173 SSParam(), 5-173 SSRun(), 5-174 SSTempDelete(), 5-174 SSTempGet(), 5-174 SSTempInfo(), 5-175 SSTempSet(), 5-176 SSTempSizeGet(), 5-176 SSTempSizeSet(), 5-177 Standard deviation, 5-21, 8-1 Burst statistics, 5-14 Error bars, 5-10, 5-59 of array, 5-9 Standard error of the mean Error bars, 5-10, 5-59 Standard settings for a view, 5-185 Statements, 3-4 Status bar Handle, 5-1 Show and hide, 5-1 Str\$(), 5-177 Straight line fit, 5-88 STRFUNC DOS script command, 6-2 String functions, 4-8 String input, 5-108 Strings, 3-2 ASCII code, 5-9 Comparison operators, 3-6 Conversions, 4-8 Convert a number to a string, 5-177 Convert ASCII to string, 5-30 Convert to a number, 5-183 Convert to lower case, 5-110 Convert to upper case, 5-182 Currently selected text, 5-161 Delete substring, 5-43 Extract fields from, 5-148 Extract fields setup, 5-148 Extract middle of a string, 5-124 Find string within another string, 5-108 Get rightmost characters, 5-149 Leftmost characters of string, 5-110 Length of a string, 5-110 Printing into, 5-138 Read from binary file, 5-12 Read string from user, 5-108 Reading using a dialog, 5-55 Specifying legal characters in input, 5-108 Write to binary file, 5-15 Substring of a string, 5-124 Subtraction of arrays and values, 5-8, 5-9 Sum of array, 5-9 Sum of array product, 5-3

Sum of bin contents, 5-37 Sweeps(), 5-177 Syntax colouring, 2-2 System\$(), 5-177 System(), 5-178

T

Tan(), 5-178 Tangent of an angle in radians, 5-178 Tanh(), 5-178 Text copy, 5-61 Text file script commands, 4-5 Text to speech, 5-170 Text view Get column number, 5-128 Get line number, 5-128 Move absolute, 5-128 Move relative, 5-128 Move to line number, 5-128 TextMark data Add marker from script on-line, 5-157 Add to memory channel, 5-124 Create channel for sampling, 5-157 Edit existing mark, 5-112 Get information, 5-114 Read string, 5-110, 5-129 Time data file creation, 5-75 Into sampling, 5-115 Maximum time in a file, 5-115 Of next item on a channel, 5-128 Of previous item on a channel, 5-110 Resolution for channels, 5-11 Time of day as a string, 5-179 as numbers, 5-179 Time shift, 5-6 Time view Apply resource file, 5-66 Background colour, 5-31 Convert seconds to Spike2 time units, 5-191 Copy data to array, 5-17 Count of events, 5-37 Printing, 5-72 Process data, 5-139 Save resource file, 5-75 Script commands, 4-3 Time of next item, 5-128 Time of previous item, 5-110 Value at given time, 5-29 Time\$(), 5-179 TimeDate(), 5-179 Title of window, 5-188 Title string for channel, 5-28 Toolbar, 5-180 Add buttons, 5-181 Change text, 5-182 Clear all buttons, 5-180 Enable and disable buttons, 5-181 Show and hide, 5-182 System toolbar, 5-1 User interaction, 5-180 window handle, 5-1

Toolbar(), 5-180 ToolbarClear(), 5-180 ToolbarEnable(), 5-181 ToolbarSet(), 5-181 ToolbarText(), 5-182 ToolbarVisible(), 5-182 Tooltips DlgButton, 5-50 In user-defined dialog, 5-48 Interact, 5-109 Toolbar. 5-181 Trace through a script, 2-3 trans() operator, 3-3 Translating DOS scripts, 6-1 Transpose of matrix, 5-115 Trend plot MeasureChan(), 5-116 MeasureToXY(), 5-118 MeasureX(), 5-119 MeasureY(), 5-119 Triggered start of sampling, 5-156 Triggered time view, 5-185 Trunc(), 5-182 Truncate real number, 5-182 Type compatibility, 3-2 Types of data, 3-1

—U—

UCase\$(), 5-182 Units for waveform or WaveMark channel, 5-28 until, 3-8 Update all views, 5-58 Update invalid regions in a view, 5-57 Upper case a string, 5-182 User interaction Ask user a Yes/No question, 5-143 Command summary, 4-6 Dialogs, 5-47 Input single key, 5-107 Let user interact with data, 5-109 Message in pop-up window, 5-124 Print formatted text, 5-137, 5-138 Read a number in a pop-up window, 5-108 Read a string in a pop-up window, 5-108 Test for key available to read. 5-109 The toolbar, 5-180 User-defined functions and procedures, 3-10

__V__

Val(), **5-183** Value parameters, 3-10 var keyword, **3-4** Variable Inspecting value, 2-4 Names, 3-1 Types, 3-1 Variable declarations, 3-4 variance, 8-1 vector, 3-3 Vertical bar notation, 1-4 Vertical cursor commands, 4-4 Vertical space for channels, 5-30 View handle, 1-2 Close view, 5-66 Create result view, 5-164 Duplicate views, 5-60 Find from view title, 5-184 for Log window, 5-111 for new result view, 5-164, 5-166, 5-167, 5-168 for new view, 5-70 for opened view. 5-71 for sampling windows, 5-152 for system windows, 5-1 Get and set, 5-183 Get linked views, 5-184 Get list of views, 5-185 Get type of view, 5-184 Override current view, 5-183 Sampling windows, 5-152 Set colour for view, 5-184 Set or get current view, 5-183 Update the view, 5-57 View manipulation functions, 4-1 View(), 5-183 View().x(), 5-183 ViewColour(), 5-184 ViewFind(), 5-184 ViewKind(), 5-184 ViewLink(), 5-184 ViewList(), 5-185 ViewQQ() DOS emulation function, 6-4 ViewStandard(), 5-185 ViewTrigger(), 5-185 ViewUseColour(), 5-186 VirtualChan(), 5-186 Visible state of a window, 5-189 Visible state of channel, 5-29 vqq%[7][2] DOS emulation array, 6-4

—W—

Wait in script. See Yield() Waterfall script example, 7-3 WAVE file output, 5-169 Waveform data Amplitude, 5-29 Average, accumulate or copy to result view, 5-164 Copy into an array, 5-17 Correlation, 5-168 Create from events, 5-62 Drawing colour, 5-31 Mean level, 5-37 Offset, 5-23 Output during sampling, 4-10 Power spectrum, 5-167 Sampling interval, 5-11 Sampling setup, 5-159 Scaling, 5-27 Units, 5-28 Write to channel, 5-30 Waveform output Add waveform to list, 5-131 Change and get DAC list, 5-133

Change cycles, 5-133 Control bar buttons, 5-134 Delete area, 5-134 Enable area, 5-134 Get area information, 5-134 Link and unlink areas, 5-135 Output rate variation, 5-135 Play offline, 5-130 Sample rate, 5-135 Status during output, 5-136 Stop output, 5-136 Trigger state, 5-136 Update waveform on-line, 5-133 WaveMark data Drawing colour, 5-31 Get information, 5-114 Offset, 5-23 Overdraw mode, 5-59 Sampling interval, 5-11 Scaling, 5-27 Setup for sampling, 5-159 Show setup dialog, 5-70 Traces, 5-22, 5-110, 5-112, 5-114, 5-120, 5-121, 5-124, 5-129 Units, 5-28 Weighting of channel space, 5-30 wend, 3-8 while, 3-8 Window data for FFT, 5-4 Window(), 5-187 WindowDuplicate(), 5-188 WindowGetPos(), 5-188 Windows Close window, 5-66 Current view, 5-183 Duplication, 5-188 Get linked view, 5-184 Get list of view handles, 5-185 Help, 5-100 Manipulation functions, 4-1 Number input with prompt, 5-108 Pop-up message window, 5-124 Position, 5-187, 5-188 Query user in pop-up window, 5-143 Show and hide, 5-189 Size, 5-188 Standard settings, 5-185 String input with prompt, 5-108 Title, 5-188 View handle, 5-183 WindowSize(), 5-188 WindowTitle\$(), 5-188 WindowVisible(), 5-189 Working Set, 5-142 Write binary data, 5-14, 5-15 Write memory channel to disk, 5-123

—X—

X axis Bin number conversions, 5-11 Display set region, 5-57 Drawing mode, 5-189 Drawing style, 5-189 Increment per bin in result view, 5-11

Range, 5-190 seconds, hms and time of day, 5-189 Show and hide features, 5-189 Show and hide scroll bar, 5-191 Tick spacing, 5-189 Title, 5-191 Units, 5-191 Value at given position, 5-29 XAxis(), 5-189 XAxisMode(), 5-189 XAxisStyle(), 5-189 XHigh(), 5-190 XLow(), 5-190 xor exclusive logical or, 3-5 XRange(), 5-190 XScroller(), 5-191 XTitle\$(), 5-191 XToBin(), 5-191 XUnits\$(), 5-191 XY view Add data, 5-191 Automatic axis expansion, 5-192 Background colour, 5-31 Channel list, 5-21 Create a new channel, 5-195 Create new XY view, 5-70 Creating from script, 7-1 Data joining method, 5-194 Data range, 5-195 Delete data points, 5-192 Drawing styles, 5-192 Get data points, 5-193 Get or set title, 5-28 Get or set units, 5-28 Limit points per channel, 5-196 Modify all channel settings, 5-195 Open file from script, 5-71 Overdraw data, 7-2 Points inside a circle, 5-194 Points inside a rectangle, 5-194 Script commands, 4-4 Set channel colour, 5-192 Set Key properties, 5-195 Sort points, 5-197 XYAddData(), 5-191 XYColour(), 5-192 XYCount(), 5-192 XYDelete(), 5-192 XYDrawMode(), 5-192 XYGetData(), 5-193 XYInCircle(), 5-194 XYInRect(), 5-194 XYJoin(), 5-194 XYKey(), 5-195 XYRange(), 5-195 XYSetChan(), 5-195 XYSize(), 5-196 XYSort(), 5-197

Y axis Channel number show and hide, 5-23 Drawing mode, 5-198 Drawing style, 5-198

-Y---

Get current limits, 5-200 Lock channels, 5-198 Range optimising, 5-129 Right and left, 5-198 Set limits, 5-200 Show all, 5-200 Show and hide, 5-197 Show and hide features, 5-198 Tick spacing, 5-198 YAxis(), **5-197** YAxisLock(), **5-198** YAxisMode(), **5-198** YAxisStyle(), **5-198** YHigh(), **5-200** Yield time to the system, 5-199 Yield(), **5-199** YieldSystem(), **5-199** YLow(), **5-200** YRange(), **5-200**