Signal

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
Measurement plots

The trend and measurement plot functions in Signal generate graphs of measurements taken from the recorded data both on-line and off-line. For each frame selected, measurements are taken from individual or multiple events and plotted as X and Y coordinates. Typical examples would be the generation of an I/V curve or measurements from multiple action potentials.

- Identify positions for measurements either manually or by setting 'active cursors' to seek for data features
- Extensive range of available measurements including peak heights, times, fit coefficients, areas, means, slopes, RMS amplitudes and standard deviations
- Plot up to 32 independent sets of measurements with optional log axes
- Easy export of results to spreadsheet

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- Signal's output sequencer provides advanced control functions including fast access to incoming data. This can be used in many ways, for example, to generate a stimulus ramp and halt ramping immediately upon detection of an action potential.
- Signal's built-in script language provides the user with the ability to customise the system for particular requirements. This is used for automation and to add specific functions not available through the menus, for example non-standard leak subtraction routines. Script control can also enhance standard functions including the option of higher-order curve fits.

CED

Cambridge Electronic Design Limited

Technical Centre, 139 Cambridge Road, Milton, Cambridge CB24 6AZ, UK. Tel: (01223) 420186

Email: info@ced.co.uk Europe & International Tel: [44] (0)1223 420186 USA and Canada Toll free: 1-800-345-7794

Distributors in: Australia, Austria, China, France, Germany, Israel, Italy, Japan, Switzerland & Turkey

www.ced.co.uk